GE Healthcare Coils

GE Signa[®] 1.5T 8-Channel CTL Coil

Service Manual

2377426-5
Rev 6
© 2008, General Electric Company,
All Rights Reserved

Revision History

Date	Rev	ECO	Change Description
08/03	1	N/A	First Issue
10/03	2	N/A	Revised Section 2 and Section 3-2-3, Removed Coil Configuration and Section 3-2-4
02/06	3	N/A	Correct Anterior Bridge FRU GE part number in Table 6-1, Page 19, from 2377426-8 to 2418267
01/07	4	N/A	Changed label image to reflect new label with WEEE symbol and add additional information to section 3 for the SNR scan. Added table locator FRU part number.
11/07	5	N/A	Added Preamp to table 6.1 Added, "Please see the locations of the preamps shown in Figure 9 below. The wrist strap MUST be worn all times when handling Preamps" on page 15. Added figure 9 to page 28. Re-numbered figures. Applied new service manual template including new Damage in Transportation guideline and language policy. Added note "The wrist strap MUST be worn all times when handling Coils" to section 5-1. Added "See section 7" to Section 2-3. Added "see section 3" to section 7-1
8/19/08	6	2008217	Revision History: Added ECO column. Section 1-3: changed Op manual p/n. Table 3-2-1: changed phantom kit p/n. Section 3-5: removed reference to preamps in Symptom #2 table, removed Figure 9. Section 6-1: removed preamp from table 6-1. Table 1-1: Updated p/n's. Section 3-2-3 step 4 revised.

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

File a report with

- Call 1-800-548-3366 and use option 8
- Contact your local service coordinator for more information on this process.

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WARNING

(EN)

- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR, OR PATIENT FROM ELECTRIC SHOCK, OR FROM MECHANICAL OR OTHER HAZARDS.

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TABLE OF CONTENTS

SECTION 1 – INTRODUCTION	14
1-1 PRODUCT IDENTIFICATION AND SHIPPING LIST	14
1-2 COMPATIBILITY	15
1-3 RELATED DOCUMENTATION	15
1-4 ENVIRONMENTAL REQUIREMENTS	15
1-5 THEORY OF OPERATION	16
SECTION 2 – SETUP AND CALIBRATION	20
2-1 COIL INSTALLATION	20
2-1-1 SPECIAL INSTALL NOTES	20
2-1-2 CONFIGURATION FOR 11X SOFTWARE	20
2-1-3 CONFIGURATION FOR 10X SOFTWARE	20
2-2 INSTALLATION FUNCTIONAL CHECKS	21
2-3 PERIODIC QUALITY ASSURANCE CHECK	21
SECTION 3 – FUNCTIONAL CHECKS	22
3-1 SCANNER VERIFICATION	22
3-2 COIL IMAGING PERFORMANCE VERIFICATION	22
3-2-1 TOOLS REQUIRED	22
3-2-2 EXPLANATION OF PROCEDURE	22
3-2-3 SNR SETUP AND TEST EXECUTION FOR 11.0 M3A AND GREATER SOFTWARE VERSIONS	22
3-3 EXTERNAL CABLE CHECK	25
3-4 PIN DIODES CHECK	26
-5 TROUBLESHOOTING TIPS	27
3-6 COIL ELEMENT INFORMATION	29
SECTION 4 – MAINTENANCE	30
4-1 COIL CARE	30
4-2 SPECIAL CARE REQUIREMENTS	30
SECTION 5 – REPLACEMENT	31
5-1 EXTERNAL CABLE REPLACEMENT	31
SECTION 6 – RENEWAL PARTS	33
6-1 FIELD REPLACEABLE UNITS	33
6-2 OTHER REPLACEABLE ACCESSORIES	33
SECTION 7 – APPENDIX	34
7-1 SNR DATA SHEET	34

SECTION 1 – INTRODUCTION

1-1 Product Identification and Shipping List

To identify the 8-Channel CTL Coil, refer to the coil label (see *Figure 1*).

Figure 1: Coil label.

Figure 2 shows a picture of the 8-Channel CTL Coil.

Figure 2: The 8-Channel CTL Coil.

SHIPPING LIST – TABLE 1-1

Box #	Part Name	GE Part #	USAI Part #	Qty
1	Coil with Coil ID Cable	N/A	U1-100341	1
3	Patient Comfort Pad Kit	2377426-6	U1-150500	1
2	Phantom Kit	2422887	N/A	1
1	Adapter Block	N/A	U1-150038	1
1	Operator's Guide Packet	2422287-500	N/A	1
1	Service Manual	2377426-5	N/A	1

1-2 Compatibility

The 8-Channel CTL coil is compatible with the Excite 1.5T system and requires 10.0 M4 or greater software. 11.0 M4 software is recommended.

The automated QA Tool is only available for software 11.0 M3 and greater.

1-3 Related Documentation

Operator's Guide, GE Part Number 2422287.

Signa Excite Service Methods CD, GE Part Number 2333500-3.

1-4 Environmental Requirements**Storage Requirements**

Coil should be stored in the scanner room.

Dimensions

43.3" L x 16.1" W x 13.6" H (110.0 cm x 41.0 cm x 34.5 cm)

Weight

18 lbs. (8.1 kg)

1-5 Theory of Operation

The block diagram of the Premier 9000 8-Channel CTL Coil is shown in *Figure 3*. *Figure 4* shows the output diagram, the coil to connector interface is shown in *Figure 5* and *Figure 6* shows the block diagram of the circuit layout. The coil is designed for obtaining diagnostic images of the whole human spine, specifically: cervical, thoracic and lumbar. The coil consists of a main coil base for imaging the entire spine and a removable anterior neck. The coil is a receive-only, phased array 12 element, 8-channel RF coil design. The 12 elements are arranged primarily as a saddle and loop coil configuration. Depending on the anatomy imaged, the user is able to select different imaging modes to selectively image the spine anatomy. The modes are enabled by the clinical user at the console.

The 8-Channel CTL Coil is actively decoupled from the RF transmit coil during transmit by means of a RF choking circuit. The choking circuit elements are switched on and off by pin diodes or by small-signal diodes. During body coil transmission, the pin diodes are turned on, as shorted circuits, with a forward bias DC source provided by the MR system (active decoupling). The small-signal diodes are turned on by the induced RF voltage in the chokes, coupled from the transmit field (passive decoupling). When either the pin diodes or the small-signal diodes are turned on, each of these RF choking elements becomes very high in impedance (typically above 2 kilo-ohm), compared with the other circuit elements (less than 70 ohm). These high impedance elements in the coil segregate the coil circuitry into several isolated electrical segments, preventing any current flow in the coil circuit.


Figure 3: Block diagram of the coil.


Figure 4: Output diagram.


Figure 5: Block diagram of the coil to connector interface.


Figure 6: Block diagram of the circuit layout.

SECTION 2 – SETUP AND CALIBRATION

2-1 Coil Installation

2-1-1 Special Install Notes

None.

2-1-2 Configuration for 11x Software

The system will automatically recognize the coil using the Coil ID feature. No configuration should be required.

2-1-3 Configuration for 10x Software

The 8-Channel CTL coil requires 10.0 M4 or greater software. 11.0 M4 software is recommended.

Updated coil configurations in the system are required for software versions prior to 11.0 M4. Updated coil configurations are recommended for 11.0 M4 and greater software revisions.

To update the coil configuration files:

1. Using the configuration manager, remove all currently installed (if any) configurations for the 8-Channel CTL coil.
2. Update the coil configurations to the latest version (Refer to Section 9 of “Host Software Install Release 10.X”, direction OW2SCA8G, as found on the Service Methods CD, GE Part Number 2333500-3).
3. Install the new configurations listed in Table 2-1-3.

CONFIGURATION MODES – TABLE 2-1-3

8CTLTOP	8CTLMID	8CTLBOT	8CTL123	8CTL234	8CTL345	8CTL456	8CTL12	8CTL23	8CTL34	8CTL45	8CTL56
---------	---------	---------	---------	---------	---------	---------	--------	--------	--------	--------	--------

2-2 Installation Functional Checks

1. From the Scan Desktop, start new scan by selecting [**New Pt**]; set **Patient ID** to “geservice” and **Patient Weight** to “111” pounds. Click [**Patient Position**] to open protocols window.
2. Plug the coil into the patient table.
3. At the console, verify the Coil has been properly identified by the system: correct name in **Coil** field and correct picture on the screen.
4. Perform Section 3-2 Coil Imaging Performance Verification.

2-3 Periodic Quality Assurance Check

On a periodic basis, such as during planned maintenance, perform the quality assurance checks as outlined below to ensure the coil is operating properly.

1. Check external cable for cracks or cuts.
2. Perform Section 3-2 Coil Imaging Performance Verification and record data values in Data Sheet (See section 7).

SECTION 3 – FUNCTIONAL CHECKS

3-1 Scanner Verification

Perform system level Signal to Noise Check. Refer to Service Methods CD; System Level Procedures; Functional Checks; Signal to Noise Check.

3-2 Coil Imaging Performance Verification

3-2-1 Tools Required

TOOLS REQUIRED – TABLE 3-2-1

Description	GE Part #	Qty
Phantom Kit	2422887	1

3-2-2 Explanation of Procedure

The image quality check uses two consecutive acquisitions for the same slice. The signal is measured at the center of the phantom with a ROI of about 60% of the phantom size. The noise is measured on the subtracted image of the two consecutive signal images using the same ROI. SNR is calculated by dividing the Signal mean by the noise standard deviation. All measurements are performed on composite images. Interim images are acquired for reference only.

3-2-3 SNR Setup and Test Execution for 11.0 M3A and Greater Software Versions

1. Remove any other surface coils (if present) from the cradle.
2. Place the coil on the patient table. Position coil so that the top table locators on both sides of the coil fit into the second notch from the magnet. Position the phantom set on the coil without the pad. (see *Figure 7*).

Figure 7: Positioning of the coil and the phantoms on the coil.


3. Landmark the coil at the indicated line on the C-Spine Bridge and advance to scan.
4. For 12.x or earlier systems: From the Common Service Desktop, select **Image Quality**, and then **Coil SNR**. Then select **Click here to start this tool**.

For 14.x or later systems: From the Common Service Desktop, select **Image Quality**, and then **MCQA Tool**. Then select **Click here to start this tool**. For additional information regarding MCQA, refer to the 14.x or later Service CD at the following location: troubleshooting/Image Quality/Multicoil Quality Assurance.

5. Select the coil type you want to run, enter the serial number and select Start Scan. Note: Additional documentation can be accessed by selecting **File, Help, [coil type]**.
6. Follow the reference example below as a visual aid.

- A. From the Common Service Desktop Manager, select the **Image Quality Tab**.

- B. Select **Coil SNR Test**.

- C. Select **Click here to start this tool**. The Multicoil SNR Tool Control displays.

The screenshot shows the 'Multicoil SNR Tool Control' application window. The 'Scan Panel' includes the following fields and controls:

- Coil Name:** A text field containing '8 Channel HR Brain'.
- Coil serial #:** A text field containing 'none'.
- Passes:** A text field containing '1'.
- Buttons:** 'Start Scan' and 'Abort Scan' buttons.
- View Archive Data:** A section containing a 'View Data' button.
- Results:** A large empty area for displaying scan results, currently showing a diagonal 'Example' watermark.

Callout circles are placed as follows:


- E:** Points to the 'Start Scan' button.
- F:** Points to the 'Coil Name' field.
- D:** Points to the 'Coil Name' field.

Multicoil SNR Tool Control

- D. The name of the coil should populate the **Coil Name** field. If desired, select the coil you are testing from the dropdown list.
- E. Enter the **Coil Serial #** found on the coil.
- F. Click **Start Scan** to begin the QA scan.

Note: The scan may take several minutes to complete. For example, an 8-channel coil may take approximately 8 to 10 minutes to complete.

- G. The results of the scan display in the Results window with a Passed or Failed indication.

Multicoil Results (Example)

3-3 External Cable Check

Diode Check from the Connector:

Step 1 – Using a multi-meter set to the **diode check function**, check the coil diodes in the following manner.

Step 2 – Orient the connector as shown in *Figure 8*.

Step 3 – Follow Table 3-3 for the correct pins and lead locations to verify the function of the diodes.

EXTERNAL CABLE EXPECTED READINGS – TABLE 3-3

PIN Number	Positive Lead	Negative Lead	DMM Reading
1	B12	A12	0.95V +/- 0.2V
2	B14	A14	0.95V +/- 0.2V
3	B16	A16	0.95V +/- 0.2V
4	B18	A18	0.95V +/- 0.2V
5	B20	A20	0.95V +/- 0.2V
6	B22	A22	0.95V +/- 0.2V


Figure 8: 60 pin Bendix connector.

Note: This connector is system side. Coil side has rows A and B reversed.

3-4 PIN Diode Check

Same as External Cable Check.

3-5 Troubleshooting Tips

Symptom #1: The system reports a coil fault during prescan or does not recognize the coil connection to the system when selected in the software.

Probable Cause	Suggested Actions	Resolution
The coil connector has become disconnected from the system interface.	Check to make sure the coil connector is fully engaged.	Engage connector and try the scan again. Make sure the indicator light on the patient table turns green with the proper coil selected on the console.
One or more PIN diodes have failed.	From inside the coil, check the PIN diodes using a DVM. Use diode mode for forward-bias measurements and resistance mode for reverse-bias measurements. Note that there should be no reverse-bias leakage. Bad PIN diodes are usually shorted or exhibit reverse-bias leakage current.	Replace the coil.
There is a DC short or open somewhere within the DC path inside the coil.	Check the coil and associated DC paths with a DVM, in both diode check and resistance modes.	If a short or open is observed, check the output cable by itself (see below). If the cable and PIN diodes are okay, then replace the coil.
The output cable has a short or an open.	Disconnect the cable at the coil and check each coaxial line with a DVM (resistance). Try moving and twisting the cable as you watch the meter to find any intermittent connections. Inspect the center pins of the SMB and Bendix connections for wear and proper engagement into the mating connection.	Replace the cable.
There is a problem with the DC bias from the system interface or the coil ID circuitry.	Check to ensure the correct DC bias is being provided to the coil during transmit and receive per GEMS System Service Manual. Try a similar 2-channel coil to see if the same problem occurs.	Correct MRI system problem.
There is a problem with the coil ID; <i>"The prescribed RF coil and detected RF coil do not match."</i>	Check online instructions by accessing the service browser at Diagnostics; Peripherals; SRI Functional Tests; Coil ID	Start the diagnostic by selecting SRI Coil ID and Run . The diagnostic will display the contents of the coil ID chip whenever the connector is inserted. Follow the troubleshooting instructions if the correct ID is not displayed. When the troubleshooting is complete and the ID is still not displayed, replace the cable.

Symptom #2: The coil does not pass SNR tests or exhibits poor image quality on patient scans.

Probable Cause	Suggested Actions	Resolution
One or more of the coil channels has a high noise level.	Look at the uncombined images to determine which channel has the problem. Compare the noise standard deviation measurements between channels and system performance logs; they should be approximately the same.	Replace the coil.
One or more of the coil channels has low signal.	Look at the uncombined images to determine which channel has the problem. Compare signal mean measurements between channels and system performance logs.	Replace the coil.
Excessive ghosting is causing the noise std. dev. measurements to be artificially high.	Window and level the images down to look at the background for signs of ghosting. Try padding the phantom in the holder. Try running the scan without the phantom holder to see if the ghosting diminishes. Try turning off the cold heads to see if the ghosting is vibration-related.	If padding or changing the phantom position can minimize ghosting, then the problem is caused by excessive vibration elsewhere in the system. If the ghosting is not symptomatic of phantom positioning or padding, then replace the coil.
There is a problem with the DC bias from the system interface or the coil ID circuitry.	Check to ensure the correct DC bias is being provided to the coils during transmit and receive per GEMS System Service Manual.	Correct MRI system DC bias problem.

3-6 Coil Element Information

Mode	Rcv 1	Rcv 2	Rcv 3	Rcv 4	Rcv 5	Rcv 6	Rcv 7	Rcv 8	PIN 1	PIN 2	PIN 3	PIN 4	PIN 5	PIN 6	PIN 7	PIN 8
8CTLTOP	ON	ON	ON	ON	ON	ON	ON	ON	-5	-5	-5	-5	+7	+7		
8CTL MID	ON	ON	ON	ON	ON	ON	ON	ON	+7	-5	-5	-5	-5	+7		
8CTLBOT	ON	ON	ON	ON	ON	ON	ON	ON	+7	+7	-5	-5	-5	-5		
8CTL123	ON	ON	ON	ON	ON	ON			-5	-5	-5	+7	+7	+7		
8CTL234			ON	ON	ON	ON	ON	ON	+7	-5	-5	-5	+7	+7		
8CTL345	ON	ON			ON	ON	ON	ON	+7	+7	-5	-5	-5	+7		
8CTL456	ON	ON	ON	ON			ON	ON	+7	+7	+7	-5	-5	-5		
8CTL12	ON	ON	ON	ON					-5	-5	+7	+7	+7	+7		
8CTL23			ON	ON	ON	ON			+7	-5	-5	+7	+7	+7		
8CTL34					ON	ON	ON	ON	+7	+7	-5	-5	+7	+7		
8CTL45	ON	ON					ON	ON	+7	+7	+7	-5	-5	+7		
8CTL56	ON	ON	ON	ON					+7	+7	+7	+7	-5	-5		

SECTION 4 – MAINTENANCE

4-1 Coil Care

Detach the coil connector from the scanner before attempting to clean the coil. Do not touch the connectors with bare fingers. Never press a sharp object against the surface of the connector. Do not reattach the connector after cleaning the coil until the coil has dried completely. Electrical shock may result if the coil is attached to the system during cleaning or when it is wet.

Do not spray or pour cleaning solution directly on coil. Do not submerge coil in solution. The coil contains sensitive electronics components that could be damaged by the solution. The coil cannot be sterilized and should be cleaned only according to the procedure outlined in this section.

The following solutions are recommended for the coil and pad surfaces: (1) a ten percent bleach solution (some discoloration may occur), (2) one ounce commercial dishwashing liquid mixed with one gallon of water or (3) warm water. Apply cleaning solution to a soft cotton cloth and proceed to clean. To prevent soiling of the coil, the user should place a cotton sheet over the coil before positioning the patient. If the coil is soiled, clean the coil as described above.

4-2 Special Care Requirements

Prior to returning a coil for service, use a ten percent bleach solution (as described above) to eliminate risk of exposure to potentially infectious materials.

SECTION 5 – REPLACEMENT

Simple removals that are clearly obvious are not described here.

Unless otherwise noted, the steps for re-assembly are simply the reverse order of the steps described for disassembly.

5-1 External Cable Replacement

Note: The wrist strap **MUST** be worn all times when handling Coils.

Step 1 – Remove coil cover.

Step 2 – Remove the two screws that hold down the cable assembly.

Step 3 – Disconnect SMB plugs and the DC Plug on the cable end from the interface boards.

Step 4 – Remove the old cable assembly.

Step 5 – Fix a new cable assembly on the CTL coil and reattach the cable lock with two screws. Connect all eight SMB plugs and the DC Plug to the corresponding receptacles on the interface board (see *Figure 11*).

Step 6 – Close the coil cover.


Figure 10: External cable assembly.


Figure 11: Fixing the cable assembly on the interface board.

SECTION 6 – RENEWAL PARTS

6-1 Field Replaceable Units

FIELD REPLACEABLE UNITS – TABLE 6-1

Part Name	GE Part #	USAI Part #
Coil with Coil ID Cable	2377426-2	U1-165140
Cable Assembly	2377426-7	U1-155131
Adapter Block	2225545-9	U1-155137
Anterior Bridge	2418267	U1-155132
Phantom Kit	2377426-3	2422887
Table Locator Kit	2415553	N/A

6-2 Other Replaceable Accessories

OTHER REPLACEABLE ACCESSORIES LIST – TABLE 6-2

Part Name	GE Part #	USAI Part #
Patient Comfort Pad Kit	2377426-6	U1-155133


Figure 12: Parts of the coil.

SECTION 7 – APPENDIX**7-1 SNR Data Sheet**

Use the table provided below to record the calculated signal to noise ratio (SNR) data obtained from the Functional Checks section (see section 3).

Date	Configuration	Element/ Roi	Signal Mean	Noise Std. Dev.	SNR	Limit/LSL	Pass/Fail
	8CTL12	1					
	8CTL12	2					
	8CTL34	3					
	8CTL34	4					
	8CTL56	5					
	8CTL56	6					

Comments

Exam No. _____

Site Name: _____

Date: _____

Magnet Strength: _____

Date	Configuration	Element/ Roi	Signal Mean	Noise Std. Dev.	SNR	Limit/LSL	Pass/Fail
	8CTL12	1					
	8CTL12	2					
	8CTL34	3					
	8CTL34	4					
	8CTL56	5					
	8CTL56	6					

Comments

Exam No. _____

Site Name: _____

Date: _____

Magnet Strength: _____

GE Healthcare Coils
USA INSTRUMENTS

1515 Danner Drive
Aurora, Ohio 44202
USA

www.gehealthcare.com

Imagination at work