

Troubleshooting the patient table

Table arm board

Symptom	Failure point	Possible cause	Troubleshooting	Resolution
The armboard does not latch.	Armboard clutches	The cable is too tight to release the clutch.	Remove the handle cover and loosen the adjustment screw until the clutch can be locked. If loosening the screw resolves the failure, there is no need to replace the armboard. If the screw does not resolve the condition, replace the armboard.	Loosen the adjustment screw or replace the armboard.
	Handle mechanism	The handle cannot return to the home position.	N/A	Replace the armboard.
		The arms, crimps, or cables are broken or worn out.	Remove the handle cover to examine the arms and cables for breaks or wear. If confirmed, replace the armboard.	
The armboard will not release with the handle.	Release cable	There is a gap between the arm and the baseplate pin.	Remove the handle cover to examine the gap between the arm and the baseplate pin; if there is a gap, turn the adjustment screw to reduce the gap to 0 mm. If the armboard releases after the adjustment, there is no need to replace it. Otherwise, replace the armboard.	Tighten the adjustment screw or replace the armboard.
	Armboard clutches	The armboard clutch shoulder screw is loose.	Remove the armboard and examine whether the clutch shoulder screw (large slotted brass screw) is too loose. If it is, replace the armboard.	Replace the armboard.

Table lowering

Symptom	Failure point	Possible cause	Troubleshooting	Resolution
The table does not lower with dock or table down pedal.	Safety lever	The safety lever is lowered.	Make sure the safety lever is not lowered. With the table mechanically and electrically docked, lower the table using a down pedal.	Raise the safety lever.
	Cradle position	The cradle is not in the home position.	Make sure the safety lever is not lowered. With the table mechanically and electrically docked, lower the table using a down pedal.	Press the home button on the magnet or pull the egress handle to move the cradle to the home position.
	Cradle interlock	The interlocks are holding the down link bar above the cylinder down lever.	Make sure the bracket is not stuck and holding the link bar up. Interference from cables can prevent the bracket from releasing.	Press the home button on the magnet or pull the egress handle to move the cradle to the home position. Adjust the flipper interlock cable per the flipper interlock adjustment procedure.
	Down pedal link bracket	The down link bar is resting on top of the cylinder down lever.	Make sure the down pedal linkage inside the dock is not broken and is installed properly.	Adjust the dock push rod to create a 1/16 to 1/8 inch gap between the down link bar and the cylinder down lever. Make sure the routing of the high pressure line does not allow the link bar to contact it.
			Make sure the down link bar is fully retracted by the down link retract spring when the down pedal is not depressed.	Replace the down link retract spring.
The table lowers slower with the dock down pedal than the table down pedal.	Push rod	The push rod is not adjusted correctly.	Examine the push rod position per the flipper interlock adjustment procedure.	Adjust the push rod until the gap between the down link bar and the cylinder down lever is 1/16 to 1/8 inch.

Table raising

Symptom	Failure point	Possible cause	Troubleshooting	Resolution
The table does not raise with the reciprocating pump.	Cylinder down lever	The cylinder down lever is stuck in the open position.	Examine the dock push rod setting to make sure it is not pushing on the cylinder down lever when both down pedals are not depressed.	Adjust the push rod until the gap between the down link bar and the cylinder down lever is 1/16 to 1/8 inch.
	Reciprocating pump linkages	The reciprocating pump linkages are worn too much.	Examine the reciprocating pump linkages for wear and make sure all the pins are installed.	Replace any worn linkage(s).
	Hydraulic system	There is too much air in the hydraulic system.	Examine the hydraulic system for too much air.	If REQUIRED, bleed air out of the hydraulic system.
		There is not enough hydraulic oil in the hydraulic system.	Compare the oil fill level to the oil fill line on the hydraulic reservoir with the table in the down position.	Add oil as required to the oil fill line.

Symptom	Failure point	Possible cause	Troubleshooting	Resolution
The table does not raise with the rotary pump (the motor operates but the table does not raise).	Cylinder down lever	The cylinder down lever is stuck in the open position.	Examine the cylinder down lever to make sure it returns home with neither down pedal being depressed.	Examine the cylinder down lever for component interference and release the cylinder down lever from the open position.
			Examine the dock push rod setting to make sure it is not pushing on the cylinder down lever when both down pedals are not depressed.	Adjust the push rod until the gap between the down link bar and the cylinder down lever is 1/16 to 1/8 inch.
	Rotary pump coupler	The dock motor coupler and rotary pump coupler are misaligned.	Make sure the switch inside of the dock is operating and installed correctly.	- Align the rotary pump. See Aligning the rotary pump. - Do the dock and table alignment procedure.
		The dock motor coupler and/or the rotary pump coupler are worn.	Visually examine each coupler for signs of significant wear.	Replace either coupler as required.
	Hydraulic system	There is not enough hydraulic oil in the hydraulic system.	Compare the oil fill level to the oil fill line on the hydraulic reservoir with the table in the down position.	Add oil as required to the oil fill line.

Long drive motor failures

Symptom	Failure point	Possible cause	Troubleshooting	Resolution
The cradle does not move.	Cable connection	The table is not electrically docked.	N/A	Electrically dock the table.
		Cable(s) are disconnected.	Examine the following cables and make sure that they are all connected: - Motor power cable (if the motor does not spin at all) - Motor encoder cable (if the motor does not spin at all) - Electromagnetic clutch cable (if the motor will spin but the cradle does not move)	Connect any cables that are disconnected.
		There is debris in a cable connection.	Examine all connections of the following cables to make sure there is no debris or other contamination present: - Cables at the I/O panel on top of the ISC cabinet - Dock connectors	Remove any debris or contamination from cable connectors.
	Rigibelt	The screws have loosened or backed out of the rigibelt mounting block.	Examine the upper and lower rigibelt drive assembly screws to see if they are loose or have backed out.	Replace the original screws with the locking patch M5screws, or adjust the torque as necessary. See Examining and replacing the dockable table rigibelt block screws .

Symptom	Failure point	Possible cause	Troubleshooting	Resolution
The cradle only moves approximately 100 mm.	Secondary latch	The secondary latch adjustment is not within specifications.	N/A	Adjust the secondary latch.
		The secondary latch has failed.	N/A	Replace the secondary latch.
The cradle only moves approximately 2 mm.	Primary latch	The primary latch adjustment is not within specifications.	N/A	Adjust the primary latch.
		The primary latch has failed.	N/A	Replace the primary latch.
		The tuning file is incorrect.	N/A	Confirm the software revision.
The cradle is stuck at approximately the 2000-2380 mm position.	Table and bridge alignment	The cradle is making contact with the tabletop or bridge.	N/A	Examine and re-align the tabletop and bridge.
		The rigibelt connection block is hitting the tabletop or bridge cover plate.	N/A	
The long drive motor is spinning but the cradle does not move and only the lower half of the electromagnetic clutch is rotating.	Electromagnetic clutch power	The electromagnetic clutch has no 22.5V power when the table is electrically docked.	Examine the following components: • electromagnetic clutch cable • DCPS	Restore power to the electromagnetic clutch.
	Electromagnetic clutch	The electromagnetic clutch has 22.5V power but the clutch is not engaged and the cradle does not move when the table is electrically docked.	N/A	Replace the electromagnetic clutch.
		The electromagnetic clutch is engaged but slipping is observed and the cradle does not move. Alignment is not within specification.	N/A	
The lower half of the electromagnetic clutch does not rotate.	Drive belt	The motor spins but the drive belt is broken.	N/A	Replace the drive belt.
		The motor spins but the drive belt is loose.	N/A	Tighten the belt to 93.5-106 N with the sonic tool or 3 mm per the pen tool.
	Motor	If the motor does not spin at all and all cable connections are clean, the motor has failed.	N/A	Replace the motor.
The cradle only moves approximately 20 mm when a cradle movement button is pushed.	String encoder	The encoder string is broken.	N/A	Replace the encoder.
		The string encoder has failed.		
	Encoder string mounting	The screw that holds the encoder string to the interface block is broken.	Manually move the cradle fully into the bore. Examine for the encoder string between the table and the bridge.	Replace the screw and interface block and attach the encoder string to the interface block.

Symptom	Failure point	Possible cause	Troubleshooting	Resolution
The cradle cannot be pushed or pulled in or out of the bore manually when the table is only mechanically docked and not electrically docked.	Rigibelt	The rigibelt drive has failed.	N/A	Replace the rigibelt drive.

Long drive noise

Symptom	Failure point	Possible cause	Troubleshooting	Resolution
Cradle movement makes noise.	Table/bridge	The table or bridge is misaligned.	N/A	Examine and re-align the tabletop and bridge.
	Object interference	There is an unexpected object on the tabletop or bridge surface.	N/A	Remove any objects from the tabletop and/or bridge surface.
	Cradle wheel	The cradle wheel(s) are worn out or otherwise failed.	N/A	Replace any worn or failed wheel(s).
	Long drive belt	The drive belt tension is too high.	N/A	Adjust the drive belt tension or replace the drive belt if required.
	Rigibelt	The rigibelt is jumping teeth.	N/A	Replace the rigibelt drive.
	Gearbox	The gearbox failed.	N/A	Replace the motor and gearbox assembly.
	Primary latch	The primary latch is making contact with the bottom of the cradle.	Examine the primary latch.	Adjust the primary latch. Replace the latch if required.
	Secondary latch	The secondary latch is making contact with the bottom of the cradle.	Examine the secondary latch.	Adjust the secondary latch. Replace the latch if required.
	Cradle flap	The cradle flap is scrubbing with the tabletop or the tabletop rear cover.	N/A	Replace the cradle flap.
	String roller	The encoder string roller has friction with the surfaces that hold the roller shaft.	N/A	Add grease between the encoder string roller and the surfaces that hold the roller shaft.

Table docking

Symptom	Failure point	Possible cause	Troubleshooting	Resolution
The table does not mechanically dock.	Dock hook	The dock hook tension is too high.	See if the dock hook attempts to engage but cannot rotate over center.	Do the dock hook tension adjustment.
		The dock hook has low to no tension.	Examine the dock hook to make sure it is not broken.	Replace the dock hook if broken.
			See if the dock hook attempts to engage but rotates back to the undocked position after failing to engage.	Do the dock hook tension adjustment.
	Dock link screw	The dock link screw is set incorrectly.	Make sure the front of the dock hook weldment barely touches the bumper tube when the table is undocked.	Adjust the dock link screw to within specifications.
		The bushing is worn at the screw contact point.	N/A	Replace the worn bushing.
	Rotary pump oil catch pan interference	The table is too low relative to the dock, causing the rotary pump oil catch pan to interfere with the dock nose cover.	Inspect for interference. Examine the alignment of the dock bosses and the table front casting bosses with the table docked.	Do the table alignment procedure to adjust caster height if required.
	Mechanical dock linkage	The mechanical dock linkage is jammed due to pressing the mechanical and electrical dock pedals simultaneously.	See if the dock hook weldment returns to barely touching the bumper tube when the table is undocked.	Release the jammed linkage by pulling the mechanical linkage forward and/or by pushing the electrical sled backward. Use care to not touch the electrical dock connector pins.
	Enclosure interference	There is interference between the lower base enclosure and dock enclosure.	Make sure there is no interference.	Adjust the mounting of an enclosure to avoid interference.
The table does not electrically dock and the electrical dock pedal requires little to no force to press and the electrical sled does not move forward.	Electrical dock linkage	The electrical dock linkage is broken.	Examine the electrical dock linkage hook. It is located under the front of the table near the rotary pump.	Replace the broken linkage.
		The electrical dock linkage is stuck.	Examine the electrical linkage to see if it is below the bell crank pin. The linkage is located under the front of the table.	Press the electrical dock pedal down, then press the mechanical dock pedal down. While holding the mechanical dock pedal down, release the electrical dock pedal.
The table does not electrically dock and the electrical dock pedal requires a high amount of force as the electrical sled starts to move forward.	Table alignment	The table electrical sled pins are hitting the dock and are not aligned to the holes on the docking connector.	Examine the alignment of the dock bosses to the table front casting bosses with the table docked.	Do the table alignment adjustment procedure.
	Connector return spring	The rocker clamps on either side of the electrical sled slipped early and hit the dock.	Make sure nothing is catching or holding the electrical sled from sliding forward and resolve any interference.	Rotate the rocker clamps back to their retracted position. Replace any broken connector return springs above the sled.

Symptom	Failure point	Possible cause	Troubleshooting	Resolution
The table does not undock and there is no pedal resistance.	Undock linkage components	With the table mechanically and electrically undocked, the undock pedal fails to activate the undock components.	Make sure the white delrin undock hook block is not broken or missing.	If required, replace the delrin undock hook block.
			Make sure the flipper interlock is set correctly, allowing the undock linkage to function when the cradle is at the home position.	Adjust the flipper interlock cable per the flipper interlock adjustment procedure.
The table electrically undocks but does not mechanically undock.	Undock linkage components	With the table mechanically and electrically docked, the undock pedal fails to activate the mechanical undock components.	Make sure the white delrin undock hook block is not broken or missing.	If required, replace the delrin undock hook block.
	Cradle position	The cradle is not in the home position.	Make sure the cradle is in the home position.	Press the home button on the magnet or pull the egress handle to move the cradle to the home position.
	Cradle interlock	The interlocks are holding the electrical dock link bar.	Make sure the flipper is not stuck and holding the link bar up. Interference from cables can prevent the flipper from releasing.	Resolve any interference so the flipper can move freely.

Port cover

Symptom	Failure point	Possible cause	Troubleshooting	Resolution
The P1/P2/P4 port cover will not open or close.	Object interference	There is an object in the gap between the port cover and the cradle cover.	N/A	Clean the gap between the port cover and the cradle cover.
	Cable	Cables interfere with the brass bar of the port cover.	N/A	Re-route the interfering cables.
	Cradle front/rear cover	The cradle front/rear cover is broken.	N/A	Replace the cradle front/rear cover.
	Port cover	The port cover is broken.	N/A	Replace the port cover.
	Screws	The screws on the port are loose.	N/A	Tighten the port screws.
	Cradle hardtop	The hardtop cover interferes with the port cover.	Make sure the screws on the cradle that limit the movement of the hardtop are installed.	Replace any missing limit screws.

Touch-and-Go (TnG)

Symptom	Failure point	Possible cause	Troubleshooting	Resolution
For any Touch-and-Go (TnG) failures, refer to Troubleshooting Touch-and-Go (TnG) .				

Cradle motion failures

Symptom	Possible error code/message	Failure point	Possible cause	Troubleshooting	Resolution
For more information, see Cradle motion hardware and theory of operation .					

Symptom	Possible error code/message	Failure point	Possible cause	Troubleshooting	Resolution
The cradle does not move.	N/A	Cable connection	Cable(s) are disconnected.	Make sure the dock motor runs. The dock motor and long drive motor share the PHPS 60V supply. Reseat cable runs M3311 and M3315 by temporarily bypassing filters and pigtail cables, and connecting directly to the PHPS.	Reseat cable runs M3358, M3369, and M3303.
		M3303 cable	M3303 may be damaged when pulled through the overhead cable trough, especially during installation.	Examine whether M3303 is damaged. See Cradle motion hardware and theory of operation .	Replace M3303. See Interconnect cables .
The cradle exhibits uncontrolled acceleration when the Cradle In or Cradle Out buttons are pressed.	N/A	Cable connection	Cable(s) are disconnected.	Reseat cable runs M3369, M3303, and PHPS J39.	Reseat cables.
		Cable connection (J3)	Damaged pin.	Examine J3 on the right side of the magnet for damaged pins.	Replace M3303. See Interconnect cables .
		M3303 cable	M3303 may be damaged when pulled through the overhead cable trough, especially during installation.	Temporarily bypass the J39 pigtail cable inside the Penetration (PEN) cabinet by connecting M3303 directly to PHPS J39. Examine whether M3303 is damaged. See Cradle motion hardware and theory of operation .	If necessary, replace M3303. See Interconnect cables .
The cradle visibly shakes during cradle travel into or out of the magnet bore.	N/A	Drive belt	Improper installation.	Examine whether the drive belt is installed correctly.	Reinstall the drive belt. See Installation and Upgrade Manuals .
		Cable connection	Damaged pin.	Examine all connections in the cable path from Scan Room Interface (SRI) J21 to PHPS2 J39 for damaged pins. Try routing cable Run M3311 directly to PHPS2 J36 and temporarily bypass the filter at PHPS J36.	Replace cable(s) that have damaged pin(s).
The cradle moves a short distance in one direction and stops and the position displayed on the In-Room Display (IRD) is not shown or is not changing.	N/A	Cable connection	Cable(s) are disconnected.	This symptom usually means that position encoder signals A or B are missing. Review the messages in the log. Reseat position encoder cable at Table Interface Box (TIB) J2. Examine and reseat cable run M3351 at TIB J5 and SRI-4 J2. Complete the actions above for the symptom #id_2024764_ph_w1p_4f5_3dc .	Reseat cables.
		Encoder	The encoder is broken.	After checking the cable connections, use the SRI real-time data diagnostic to make sure the encoder is sending signals. See Doing the Scan Room Interface (SRI) real-time data diagnostic .	Replace the position encoder.

Symptom	Possible error code/message	Failure point	Possible cause	Troubleshooting	Resolution
The cradle does not always respond when the Cradle In or Advance to Scan buttons are pressed, but the position displayed on the In-Room Display (IRD) updates relative to the LPCA location.	2268821: Cradle didn't start moving when requested. 2268879: Cradle didn't start moving when requested.	Cable connection	Cable(s) are disconnected.	This symptom suggests a problem with the PHPS2 and possibly the connection between the SRI J21 and PHPS2 J39. Complete the actions above for the symptom #id_2024764_ph_j2w_qf5_3dc .	Reseat cables.
		Cable connection (J3)	Damaged pin.		Replace M3303. See Interconnect cables .
		M3303 cable	M3303 may be damaged when pulled through the overhead cable trough, especially during installation.		If necessary, replace M3303. See Interconnect cables .
		PHPS2	The PHPS2 has failed.	See if the above action solved the problem. If not, then replace the PHPS2 and make sure that the SRI motion test is successful. See Doing the Scan Room Interface (SRI) motion test .	Replace the PHPS2.
		SRI	The SRI has failed.	See if the above action(s) solved the problem. If not, then replace the SRI.	Replace the SRI.
The log messages indicate a problem with the cradle velocity.	N/A	Cable connection	Cable(s) are disconnected.	There may be an issue with A, B, and/or Z signals from the long drive motor to the SRI-4. Examine and reseat cable run M3351 at TIB J5 and SRI-4 J2. Examine whether the SRI real-time data diagnostic is successful. See Doing the Scan Room Interface (SRI) real-time data diagnostic . Complete the actions above for the symptom #id_2024764_ph_ipj_wg5_3dc .	Reseat cables.
The cradle moves freely using the Cradle In and Cradle Out buttons on the left and right enclosure panels, but the magnet enclosure display does not update.	N/A	Cable connection	Cable(s) are disconnected.	Examine whether the magnet enclosure display test is successful. See Doing a magnet enclosure display test . A successful test suggests the system is not detecting the Home switch signal. See if the Home switch LED changes statuses when the cradle is moved to the Home position. Review the status of the Home switch with the SRI real-time data diagnostic. See Doing the Scan Room Interface (SRI) real-time data diagnostic . End-of-travel (EOT) and cradle Egress switches may also be affected. If the switch signals are not detected, temporarily bypass the filter (remove from path) at TIB J1. See Cradle motion hardware and theory of operation .	Reseat cables at TIB J1, TIB J5, and SRI J2.
The cradle responds to the Cradle In and Cradle Out buttons on one left or right enclosure panel, but not the other.	N/A	Cable connection	Cable failure.	Swap the left and right enclosure panels. If the issue persists on the same side, but not the other, then replace the Y-cable in run M3337. The left and right front enclosure panels share motion signal lines to SRI-4 J4.	Replace the Y-cable run M3337. Refer to part number 5170095 in Interconnect cables .

Symptom	Possible error code/message	Failure point	Possible cause	Troubleshooting	Resolution
The landmark advances to within 1 mm of the target after pressing Advance to Scan and then stops.	N/A	Cable connection	Cable(s) are disconnected.	Press the Advance to Scan button again to advance to 0 mm. The symptom may be due to the position encoder signals intermittently reaching SRI-4. Reseat cable run M3351 at TIB J5 and SRI-4 J2; usually reseating this cable resolves the issue. If not, then reseat the cable from the position encoder at TIB J2.	Reseat cables.
The system halts during an exam, when the patient is not moving, and the message log indicates the cradle moved off the scan plane.	2261552: The table moved off of scan plane while scanning.	Cable position	Cable is too close to the magnet bore or RF drive cables.	Make sure the patient is not moving and update the route of cable run M3351.	Move cable run M3351 farther away from the magnet bore and the RF drive cables. The RF drive cables are located between the dual drive transmit/receive switch (DTRSW) and the body coil.
The log messages indicate a problem with Home or EOT.	N/A	Brass plunger	Brass plungers inside the LPCA sometimes get stuck, blocking the optical switch light beam.	The optical switch LED indicator is orange when the light beam is unblocked.	Reset the brass plungers inside the LPCA.
The log message indicates the cradle is attached while the table is undocked and/or lowered.	2268941: The SRI detected that the cradle is attached while the table is undocked and/or lowered.	Plastic spring	Over time, the spring remains slightly compressed, blocking the light beam path of the LPCA cradle detect switch.	N/A	Replace the spring. See Replacing the LPCA cradle sensing components .
Under heavy load, the cradle travels a short distance and reports an error, whereas under light or no load, the cradle moves normally.	2268819: Velocity tolerance exceeded.	Long drive assembly belt	The belt does not follow the path shown in the <i>Installation Manual</i> . See Installation and Upgrade Manuals .	Examine whether the SRI motion test fails. See Doing the Scan Room Interface (SRI) motion test . See if the cradle moves freely under heavy load when detached from the LPCA. Use a flashlight and view the underside of the long drive assembly and inspect the belt routing path through the long drive assembly. Make sure the path matches that shown in the <i>Installation Manual</i> . See Installation and Upgrade Manuals .	Update the belt routing path as described in the <i>Installation Manual</i> . See Installation and Upgrade Manuals .
		Assembly bolts	The long assembly bolts were incorrectly installed near the rollers. (The short bolts should have been installed, instead.)	Use a flashlight and view the underside of the long drive assembly and illuminate the space around the belt rollers. If the assembly bolts are too long, they will contact the rollers, thereby impeding roller rotation.	See if the shorter assembly bolts have been installed in a different location. Swap the long bolts with the shorter bolts.

Symptom	Possible error code/message	Failure point	Possible cause	Troubleshooting	Resolution
During scanning, the cradle is intermittently driven out of the magnet bore, sometimes to the Home position, unprompted, at high speed.	2261552: The table moved off of scan plane while scanning.	Egress switch	The screw securing the wire on the Egress switch is loose.	See if the <code>sriEvents.xml</code> file contains "egress handle" and "egress move" entries for each event. Examine whether the wire on the Egress switch is loose or damaged.	Apply blue Loctite to the screw threads and tighten the screw on the wire.

5690015-2EN | Revision 7 | en-US | [gehealthcare.com](https://www.gehealthcare.com)

©2025 GE HealthCare Proprietary

GE is a trademark of General Electric Company used under trademark license.