SERVICE MANUAL

Invivo Corporation
3545 SW 47th Avenue
Gainesville, FL 32608

Document TR0228S

Revision 12

GE 1.5T 8 Channel Neurovascular Array

GE Catalog Part Number: M3087JB

Invivo Part Number: 800228

© Invivo Corporation, 2002

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation “**damage in shipment**” written on **all** copies of the freight or express bill **before** delivery is accepted or “signed for” by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately** upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14** day period.

File a report with

- Call 1–800–548–3366 and use option 8.
- Fill out a report on <http://egems.med.ge.com/edq/home.jsp>
- Contact your local service coordinator for more information on this process.

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TABLE OF CONTENTS

Section 1 – Introduction	8
1-1 Product Identification and Shipping List	8
1-2 Compatibility	8
1-3 Related Documentation	8
1-4 Environmental Requirements	9
1-5 Theory of Operation	9
Section 2 – Setup and Calibration	12
2-1 Coil Installation	12
2-1-1 Special Install Notes	12
2-1-2 Installing the Coil	12
2-2 Installation Functional Checks	13
2-3 Customer Notification	13
2-4 Periodic Quality Assurance Check	13
Section 3 – Functional Checks	14
3-1 Scanner Verification	14
3-2 Coil Imaging Performance Verification	14
3-2-1 QA Phantom and Coil Setup with Legacy Phantom	14
3-2-2 QA Phantom and Coil Setup with Unified Phantom	16
3-2-3 SNR Analysis	18
3-2-4 Plot Results	21
3-2-5 Example Images	27
3-3 External Cable Check	31
3-3-1 RF Signal Channels	31
3-3-2 PIN Diode Control Channels	32
3-4 PIN Diodes Check	33
3-5 Mechanical Hardware Check	33
Section 4 – Maintenance	34
4-1 Coil Care	34
4-2 Special Care Requirements	34
Section 5 – Replacement	35
5-1 Disassembly of Coil	35
5-2 External Cable Replacement	35
5-3 Mechanical Hardware Replacement	36
5-3-1 Superior Latch Assembly FRU Replacement	36
5-3-2 Right / Left Latch Assembly FRU Replacement	36
Section 6 – Renewal Parts	38
6-1 Field Replaceable Units	38
6-2 Other Replaceable Accessories	38
Section 7 – Appendix	39
7-1 Schematic	39
7-2 Coil Configuration	39

SECTION 1 – INTRODUCTION

1-1 Product Identification and Shipping List

This is a service manual for the GE SIGNA 1.5T 8 Channel Neurovascular Array.

NOTE: Coils do not ship with phantoms. Phantoms come in a unified phantom set with the MR system.

SHIPPING LIST – TABLE 1-1

Description	GE Part #	INVIVO #	Qty
1.5T 8 Channel Neurovascular Array	2317115-2	103045	1
Operator Manual	2317115-4	500078	1
Service Manual	2317115-5	TR0228S	1
NV Array Positioner	5344806	N/a	1
Wedge Pad	2317115-8	101407	2
HNC Pad	2317115-6	102955	1

1-2 Compatibility

This coil is compatible with the following hardware configurations


- Signa Horizon and LX Horizon 1.5T with Excite Technology (8-channel)
- Excite II (11.x), Excite III (12.x), HDx/HDxt/HDi that are upgrades.

1-3 Related Documentation

- Signa 11.x Service Methods CD, 2333500-1
- 8 Channel Neurovascular Array Operator's Manual, 2317115-4


The neck section consist of 5 elements, arranged as detailed in Figures 2 through 5. Table 1 describes the chest and neck element names.

Figure 1-5-2

Front Elements

Figure 1-5-3

Back Elements

Figure 1-5-4

Right Side Element

Figure 1-5-5


Left Side Element

RC	Right Chest
LC	Left Chest
A/PN	Anterior/Posterior Neck
CC	Center Chest
L/RN	Left/Right Neck



There are two modes of operation that utilize the two coil sections. Mode one is the head only mode where each of the 8 head elements is directed to the output. Mode two utilizes a switchable 8-to-3 combiner at the output of the head section and the five elements in the neck section for a total of 8 signals at the output. See the block diagram on the following page for more detail.

1.5T 8 Channel Neurovascular Array Block Diagram

MODE 1

MODE 2

SECTION 2 – SETUP AND CALIBRATION

2-1 Coil Installation

2-1-1 Special Install Notes

1. The mirror assembly is installed by sliding [in the superior direction] the grooved section of the mirror assembly into the grooved section atop the anterior coil. To remove, press the release button on the mirror assembly while sliding the mirror assembly in the inferior direction.
2. The upper mirror is designed to “break away” should it come into contact with the bore. The upper mirror can be reset by rotating it back into position until the latches are secured.
3. Remove the protective sheet of plastic from the surface of each mirror.
4. The 8 Channel NV Coil requires 10.0 M4 or greater software.
5. In order to use the 8 Channel NV coil, sites with 10.0 M4 or greater software require the installation of the 8 Channel NV Array QA Tool (2376441 rev. SW002). See direction 2382643 for procedure, as found in the shipping contents of the coil and the Service Methods CD.

2-1-2 Installing the Coil

2-1-2-1 Configuration File Installation

In order to use the 8 Channel NV coil, config files on the system must be updated.

1. Using the config file manager, remove all currently installed configs for the 8 Channel NV array by MRI Devices (Invivo).
2. Using the config file manager, install the following coil names:

8NVHEAD_A	8NVARRAY_A
8NVHEAD ONE RCVR	8NVANGIO_A
8NVHEAD NECK_A	8NVANGIO LAT
8NVHEAD NECK	8NVANGIO MED
8NVHService	8NVAService

2-1-2-2 Customer Protocol Modification

Update the system protocol files using CoilCorrect.sh. See section 9-2 of “Host Software Install Release 10.X”, direction OW2SCA8G, as found on the Service Methods CD. Translate the config names as shown in Table 2-1-2-2.

CONFIG NAMES TRANSLATION CHART – TABLE 2-1-2-2

10.0 M4 Name	New Name
8NVHEAD	8NVHEAD_A
8NVHEAD ONE RCVR	8NVHEAD ONE RCVR
8NVHEAD LONG	8NVHEAD NECK_A
8NVHEAD LONG 3	8NVHEAD NECK
8NVARRAY	8NVARRAY_A
8NVARRAY 7	8NVANGIO_A
8NVANGIO	8NVANGIO LAT
8NVANGIO 3	8NVANGIO MED

2-1-2-3 System Reboot and Save Info

1. Reboot the system to allow all coil configuration changes to take place.
2. Perform a “save info” to preserve the new settings in case they need to be reloaded in the future.

2-2 Installation Functional Checks

1. Perform system level Signal to Noise Check. Refer to Service Methods CD; System Level Procedures; Functional Checks; Signal to Noise Check.
2. Perform Section 3 – Coil Imaging Performance Verification.

2-3 Customer Notification

1. Leave the customer information kit, 2385392-100 and customer letter 2393321 at the site or deliver to the customer.
2. Send “8NVArraySupport@med.ge.com” confirmation that the coil has been installed. Provide the system ID, customer contact name, and phone number in the e-mail.

2-4 Periodic Quality Assurance Check

On a periodic basis, such as during installation and planned maintenance, perform the quality assurance checks as outlined below to ensure the coil is operating properly.

1. Check external cable for cracks or cuts.
2. Check plastic screws for tightness, damage, or missing.
3. Perform Section 3-2 – Coil Imaging Performance Verification and record data values in Data Sheet (Table 7-1).
4. If SNR is trending downwards, perform a system level signal to noise check. Also check the interconnect pins for damage or abuse.

SECTION 3 – FUNCTIONAL CHECKS

3-1 Scanner Verification

Perform system level Signal to Noise Check. Refer to Service Methods CD; System Level Procedures; Functional Checks; Signal to Noise Check.

3-2 Coil Imaging Performance Verification

3-2-1 Tools Required

TOOLS REQUIRED USING LEGACY PHANTOM – TABLE 3-2-1

Description	GE Part #	Qty
Legacy 8 Ch NV Phantom Kit	2317115-3	1

TOOLS REQUIRED USING UNIFIED PHANTOM – TABLE 3-2-2

Description	GE Part #	Qty
17cm TLT Head Spherical Phantom	46-265826G6	1
10cm Spherical Phantom	46-317586G1	1
Large Cylindrical unified phantom	5342679	1
Phantom positioner, NV	5344806	1

3-2-2 QA Phantom and Coil Setup with Legacy Phantom

1. End the previous exam if this has not already been done.
2. Attach the base of the 8 channel neurovascular array to the cradle in the same position on the cradle as the standard head coil as shown in Figure 3-2-2-1 and remove the top.
3. Note the serial number of the coil.
4. Place phantom positioner on coil as shown in Figure 3-2-2-2.
5. Place phantom on positioner as shown in Figure 3-2-2-3.
6. Attach top to coil as shown in Figure 3-2-2-4.
7. Establish an axial landmark on the line on top of the coil. Advance to scan.

Figure 3-2-2-1

Remove Coil Top

Figure 3-2-2-2

Place Phantom Positioner on coil

Figure 3-2-2-3

Place Phantom on Positioner

Figure 3-2-2-4

Place Top on Coil

Figure 3-2-2-5

Correct Bottle Placement

NOTE: For correct bottle placement, insert the bottom of each phantom bottle fully, until it touches the curved inside surface of the main phantom. Refer to Figure 3-2-2-5.

3-2-3 QA Phantom and Coil Setup with Unified Phantom

Visual	Description
 Figure 3-2-3-1	1. Place the base of the NV array coil on the patient table and latch it to the table as shown in Figure 3-2-3-1.
 Figure 3-2-3-2	2. Place the phantom positioner in the coil as shown in the Figure 3-2-3-2.
 Figure 3-2-3-3	3. Place the phantoms in the positioner as shown in the Figure 3-2-3-3.
	4. Carefully align the connector pins on the base of the Coil with the top part of the coil and latch both the coil sections into place as in Figure 3-2-3-4. (The scanner will not operate if the coil sections are not latched)

Figure 3-2-3-4	correctly).
----------------	-------------

3-2-4 MCQA / SNR Test

Continue for 14x Software Versions or Greater

1. On the Common Service Desktop, select Image Quality, and then the MCQA tool.
2. Make the proper coil selection and press the "Start Scan" button.
3. Further information on how to run the MCQA tool can be found on 14x/15x service CD or on web at: Troubleshooting/Image Quality/multicoil Quality Assurance.

Continue for 11x and 12x Software Versions

1. From the service browser, select Image Quality -> Coil SNR -> [Click here to start this tool](#)
2. Select the coil type you are testing, enter the serial number and select Start Scan.

Note: Additional documentation can be accessed by selecting File -> Help -> 8NV

Note: The SNR tool requires that **8NVAservice** and **8NVHservice** service configuration files be installed.


The image quality check uses two consecutive acquisitions for the same slice. The signal is measured at the center of the phantom with a ROI of about 60% of the phantom size. The noise is measured on the subtracted image of the two consecutive signal images using the same ROI. SNR is calculated by dividing the Signal mean by the noise standard deviation. All measurements are performed on composite images. Interim images are acquired for reference only.

Proceed to TABLE 3-2-6-2 SNR LIMITS

3-2-5 SNR Analysis 11.x and 12.x only

1. In the command window enter **snr_for_8nvarray_QA** and press enter. A screen like that shown in Figure 3-2-5-1 should appear.


Figure 3-2-5-1 SNR Analysis Screen

The image shows a software interface for SNR analysis. At the top, there are two input fields: "Exam number" with the value "50000" and "Coil serial number" with the value "none". Below these are two buttons: "Calculate SNR" and "Plot Data". The main area of the screen is titled "SNR Results" and contains a large, empty rectangular box for displaying the results.

2. In the Exam number window, enter the exam number as shown in Figures 3-2-5-2 and 3-2-5-3 (The Delete key must be used instead of Backspace).
3. In the Coil serial number window, enter serial number and bay number as shown in Figures 3-2-5-4 and 3-2-5-5 (Optional).


Figure 3-2-5-2

The screenshot shows a software interface with a label "Exam number" in blue text. To its right is a text input field containing the number "50000" in red. Below the input field is a cyan button with the text "Calculate SNR" in orange. At the bottom of the window is a horizontal bar.

Highlight Digits for Exam Number

Figure 3-2-5-3

The screenshot shows a software interface with a label "Exam number" in blue text. To its right is a text input field containing the number "51640" in red. Below the input field is a cyan button with the text "Calculate SNR" in orange. At the bottom of the window is a horizontal bar.

Type In Exam Number

Figure 3-2-5-4

The screenshot shows a software interface with a label "Coil serial number" in blue text. To its right is a text input field containing the word "none" in red. Below the input field is a cyan button with the text "Plot Data" in orange. At the bottom left of the window, the letters "ts" are visible.

Highlight Coil Serial Number

Figure 3-2-5-5


The screenshot shows a software interface with a label "Coil serial number" in blue text. To its right is a text input field containing the text "bay4a U12345" in red. Below the input field is a cyan button with the text "Plot Data" in orange.

Enter Coil Serial Number and Bay Number

- Click the Calculate SNR button as shown in Figure 3-2-5-6.

Figure 3-2-5-6

Exam number 51640

Calculate SNR

Click Calculate SNR Button


- Results will appear as shown in Figure 3-2-5-7.

Figure 3-2-5-7

SNR Results			
8 Channel Neurovascular Array Quality Assurance Test			
serial number bay4a U12345			
8NVHEAD mode *****			
phantom center =(127,129)			
Individual Element NEMA SNR from composite images			
ROI	signal	noise	snr
1	758.9	4.73	226.8
2	871.0	4.82	255.5
3	889.4	4.91	256.1
4	867.1	4.88	251.1
5	857.1	5.12	236.8
6	876.3	5.13	241.6
7	900.9	4.87	261.5
8	802.5	4.66	243.7
Composite SNR Using NEMA Method			
ROI	signal	noise	snr
C	825.7	5.03	222.0
8NVARRAY head area *****			
phantom center =(127,129)			
Individual Element NEMA SNR from composite images			
ROI	signal	noise	snr
1	873.9	4.54	272.3
2	965.1	4.35	313.5
3	988.0	4.02	347.9
4	1022.0	4.33	333.5
5	1069.4	4.15	364.3
6	1119.4	4.53	349.8
7	1083.4	4.78	320.2
8	938.8	4.93	269.2
Composite SNR Using NEMA Method			
ROI	signal	noise	snr
C	942.5	4.37	291.4
8NVARRAY neck area *****			
phantom center =(124,129)			
Composite SNR Using NEMA Method			
ROI	signal	noise	snr
C	722.5	5.28	184.9
8NVARRAY chest mode *****			
phantom center =(122,134)			
Individual Element NEMA SNR from composite images			
ROI	signal	noise	snr
1	237.0	4.70	71.3
2	211.0	4.35	68.7
3	283.6	4.84	82.9

SNR Results

6. Table 3-2-6-1 shows the test results that correspond to each test.

Note: The results are also stored in a file in the /usr/g/service/data directory. The filename will be similar to 8nvarray_51897_c7ga0522.coilqa where 51897 is the exam number.

7. Record the coil serial number and the associated filename in Table 3-2-5-3 for future reference.

Note: The analysis tool can be used to process multiple runs by entering a new exam number and coil serial number and pressing the Calculate SNR button again. It is not necessary to exit the analysis tool after each run.

8. Proceed to the next section for instructions to plot the data.

3-2-6 Plot Results 11x and 12x only

Note: Perform steps 1 through 3 only if the SNR tool is not already running.

Note: Table 3-2-6-1 explains the output of the tool. Table 3-2-6-2 provides SNR acceptance limits. Table 3-2-6-3 provides a data sheet for recording file names for future reference.

1. Open a command window on the tools desktop.
2. Navigate to the /usr/local/bin/CoilsTools directory.

3. In the command window enter **snr_for_8nvarray_QA** and press enter. A screen like that shown in Figure 3-2-6-1 should appear.
4. Click on the plot data button as shown in Figure 3-2-6-1.
5. Click on the “List 8NVARARRAY data files” button as shown in Figure 3-2-6-2.
6. Click on the name of the file to be plotted as shown in Figure 3-2-6-3.
7. Click on the plot data button as shown in Figure 3-2-6-4 to process the file.
8. Example of a nominal coil is shown in Figure 3-2-6-5.

Figure 3-2-6-1


Coil serial number

bay4a U12345

Plot Data

Start Plot Tool

Figure 3-2-6-2

8NVARARRAY Coil Data

Data path /usr/g/service/data

List 8NVARARRAY data files

Plot data

8nvarray_51612_d32g5855.coilQA
8nvarray_51640_d3cg4858.coilQA
8nvarray_51895_d47h5234.coilQA
8nvarray_51948_d4984436.coilQA
8nvarray_51951_d49a1517.coilQA
8nvarray_51957_d49k4043.coilQA
8nvarray_51958_d49k5845.coilQA
8nvarray_51959_d49l1206.coilQA
8nvarray_51960_d49l3301.coilQA

List Files

Figure 3-2-6-3

8NVARRAY Coil Data

Data path


8nvarray_51612_d32g5855.coilQA
8nvarray_51640_d3cg4858.coilQA
8nvarray_51895_d47h5234.coilQA
8nvarray_51948_d4984436.coilQA
8nvarray_51951_d49a1517.coilQA
8nvarray_51957_d49k4043.coilQA
8nvarray_51958_d49k5845.coilQA
8nvarray_51959_d49l1206.coilQA
8nvarray_51960_d49l3301.coilQA

Select File

Figure 3-2-6-4

Process File

Figure 3-2-6-5

Results Compared to Upper and Lower Limits

Note: Red lines show upper and lower acceptance limits. Blue line shows test results.

TABLE 3-2-6-1 EXAMPLE OUTPUT OF TOOL

8 Channel Neurovascular Array Quality Assurance Test			
serial number bay4a U12345			
8NVHEAD mode *****			
phantom center =(127,129)			
Individual Element NEMA SNR from composite images			
ROI	signal	noise	snr
1	758.9	4.73	226.8
2	871.0	4.82	255.5
3	889.4	4.91	256.1
4	867.1	4.88	251.1
5	857.1	5.12	236.8
6	876.3	5.13	241.6
7	900.9	4.87	261.5
8	802.5	4.66	243.7
Composite SNR Using NEMA Method			
ROI	signal	noise	snr
C	825.7	5.03	222.0
8NVARRAY head area *****			
phantom center =(127,129)			
Individual Element NEMA SNR from composite images			
ROI	signal	noise	snr
1	873.9	4.54	272.3
2	965.1	4.35	313.5
3	988.0	4.02	347.9
4	1022.0	4.33	333.5
5	1069.4	4.15	364.3
6	1119.4	4.53	349.8
7	1083.4	4.78	320.2
8	938.8	4.93	269.2
Composite SNR Using NEMA Method			
ROI	signal	noise	snr
C	942.5	4.37	291.4
8NVARRAY neck area *****			
phantom center =(124,129)			
Composite SNR Using NEMA Method			
ROI	signal	noise	snr
C	722.5	5.28	184.9
8NVARRAY chest mode *****			
phantom center =(122,134)			
Individual Element NEMA SNR from composite images			
ROI	signal	noise	snr
1	237.0	4.70	71.3
2	211.0	4.35	68.7
3	283.6	4.84	82.9

Using ROIs shown in Figure 3-2-2-8 signal is measured in each of the two composite images of the head area and averaged. The same ROI for each element is used to measure noise in the filtered difference image.

Composite image results

Using ROIs shown in Figure 3-2-2-8 signal is measured in each of the two composite images of the head area and averaged. The same ROI for each element is used to measure noise in the filtered difference image.

Composite image results

Composite image results

Composite image results measured on a circular ROI on each of 3 bottles in chest area.

TABLE 3-2-6-2 SNR LIMITS

			USL	Nominal	LSL
8NVHEAD mode	head	elem1	268	229	189
		elem2	301	257	213
		elem3	305	260	215
		elem4	309	264	218
		elem5	307	262	217
		elem6	313	267	221
		elem7	311	265	219
		elem8	275	235	194
		comp	270	231	191
8NVARRAY mode	head	elem1	327	279	231
		elem2	388	331	274
		elem3	407	348	288
		elem4	426	364	301
		elem5	429	366	303
		elem6	423	361	299
		elem7	394	336	278
		elem8	334	285	236
		comp	352	300	248
	neck		230	196	162
	chest 1		88	75	61
	chest 2		81	69	56
	chest 3		94	80	66

TABLE 3-2-6-3 DATA SHEET

Coil Serial Number	Filename

3-2-7 Example Images

Figure 3-2-7-1
Receiver 1
8NVHEAD Mode Localizer Image



Figure 3-2-7-2 8NVHEAD Mode Signal Images



Figure 3-2-7-5 8NVARRAY Mode Signal Images Slice 2 (neck)
Receiver 1 **Receiver 2** **Receiver 3**



Figure 3-2-7-6 8NVARRAY Mode Signal Images Slice 3 (head)
Receiver 1 Receiver 2 Receiver 3

3-3 External Cable Check


Each MCP signal connection has a DC blocking capacitor in the cable connector. Therefore, continuity of the RF signal conductors cannot be determined in the usual way, by measuring from pin to pin using a multimeter. The external cable check must be performed in two parts, section 3-3-1 for the RF signal conductors and 3-3-2 for the PIN diode control conductors.

Use Figure 3-3 for pin identification while performing the procedures in sections 3-3 and 3-4.

System and Coil Connectors Pin Identification Figure 3-3
Bendix System Connector

SubD Coil Cable Connector

3-3-1 RF Signal Channels

The following procedure will rule out a short to GND for each of the 16 signal pins detailed in Table 3-3-1.

Note: This check does not determine continuity for each of the MCP signal connections.

Check the coil cable using the following procedure. Referring to Section 5-2, External Cable Replacement, disconnect the SubD coil cable connector. Using a digital multimeter, measure each of the 16 signal pins designated in Table 3-3-1 with respect to GND. Verify that the signal pins are not shorted to GND.

EXTERNAL CABLE EXPECTED READINGS – TABLE 3-3-1

Signal Channel	SubD Coil Cable Connector	GND
MCP1	A1	2, 5, 7, 10, 12, 14, 17, 22, 24 and shield of coaxial contacts
MCP2	1	
MCP3	16	
MCP4	6	
MCP5	9	
MCP6	23	
MCP7	13	
MCP8	A2	
Signal Channel	Bendix System Connector	GND
MCP1	B1	A12, A14, A16, A18, A20, A22, A24 and A26
MCP2	A2	
MCP3	B3	
MCP4	A4	
MCP5	B34	
MCP6	A35	
MCP7	B36	
MCP8	A37	

If one or more of the pins are shorted to GND, replace the coil cable.

3-3-2 PIN Diode Control Channels

With the coil cable disconnected, use a multimeter to determine continuity between the designated pins for each control channel in Table 3-3-2. Also verify that the control pins are not shorted to GND.

EXTERNAL CABLE EXPECTED READINGS – TABLE 3-3-2

Control Channel	SubD Coil Cable Connector	Bendix System Connector
MC 1	3	B12
MC 2	4	B14
MC 3	19	B16
MC 4	20	B18
MC 5	8	B20
MC 6	21	B22
MC 7	11	B24
MC 8	25	B26
+15 V	15	A25, B25
GND	2, 5, 7, 10, 12, 14, 17, 22, 24 and shield of coaxial contacts	A12, A14, A16, A18, A20, A22, A24 and A26

If one or more of the readings are > 3 ohms, replace the coil cable. If all readings are < 3 ohms, proceed to the **Diode Check**.

3-4 PIN Diodes Check

With the coil cable connected to the coil, use a multimeter to measure the diode drop across the coil connector pins detailed in Table 3-4. Use Figure 3-3 for pin identification while performing this procedure.

PIN DIODE EXPECTED READINGS – TABLE 3-4

Bendix System Connector		Reading
Positive Lead	Negative Lead	
B12	GND [A12]	1.6 ± 0.2
GND [A12]	B12	Open
B14	GND [A14]	1.6 ± 0.2
GND [A14]	B14	Open
B16	GND [A16]	1.6 ± 0.2
GND [A16]	B16	Open
B18	GND [A18]	1.6 ± 0.2
GND [A18]	B18	Open
B20	GND [A20]	1.6 ± 0.2
GND [A20]	B20	Open
B22	GND [A22]	1.5 ± 0.2
GND [A22]	B22	Open
B24	GND [A24]	1.5 ± 0.2
GND [A24]	B24	Open
B26	GND [A26]	1.5 ± 0.2
GND [A26]	B26	Open

If one or more of the readings indicate a short, replace the coil. If all of the above readings are correct, redo 3-2 *Coil Imaging Performance Verification*. If the problem still exist, replace the 1.5T 8 Channel Neurovascular Array.

3-5 Mechanical Hardware Check

Not Applicable

SECTION 4 – MAINTENANCE

4-1 Coil Care

Detach coil connector from scanner before attempting to clean. Do not reattach after cleaning until coil has dried completely. Having the coil attached to the system during cleaning or when it is wet may result in electrical shock.

Do not spray or pour cleaning solution directly on coil. Do not submerge coil in solution. The coil contains sensitive electronics components that could be damaged by the solution.

The GE SIGNA 1.5T 8 Channel Neurovascular Array, Baseplate and Head Pads may be cleaned by wiping with a cloth dampened with a solution of 30% isopropyl alcohol and 70% tap water or a 10% Bleach solution.

4-2 Special Care Requirements

None

SECTION 5 – REPLACEMENT

Simple removals that are clearly obvious are not described here.

Unless otherwise noted, the steps for re-assembly are simply the reverse order of the steps described for disassembly.

5-1 Disassembly of Coil


In order to perform the following procedures:

- 5-3-1 Superior Latch Assembly FRU Replacement
- 5-3-2 Right / Left Latch Assembly FRU Replacement

it will be necessary to remove the anterior coil cover from the anterior coil section.

1. Loosen, but do not remove, the 17 screws that fasten the anterior coil cover to the anterior coil section. The 17 fastener locations are detailed in Figure 5-1.
2. Carefully lift the anterior coil section away from the anterior coil cover.
3. Place the anterior coil section aside, being careful to keep the screws stored in their respective holes.
4. During FRU replacement and cover reinstallation, assure that the electronics in the anterior coil section are not disturbed.

Figure 5-1 Anterior Coil Fastener Locations

5-2 External Cable Replacement

Remove the 4 coil connector mounting screws located at the superior end of the coil. Remove the coil connector by pulling away from the coil in the superior direction.

Cables bearing this Coil ID label are compatible with software versions 10.X and newer. With software versions 11.X and newer, the system will generate an error when an incompatible coil is prescribed.

5-3 Mechanical Hardware Replacement

5-3-1 Superior Latch Assembly FRU Replacement


In order to perform the Superior Latch Assembly FRU Replacement, it will be necessary to remove the anterior coil cover from the anterior coil section. Follow the procedure detailed in section 5-1 Disassembly of Coil, then continue below.

1. Remove the 4 screws that fasten the superior latch assembly to the anterior coil cover. The 4 fastener locations are detailed in Figure 5-3-1-1.
2. Replace the superior latch using the Superior Latch Assembly FRU.

Figure 5-3-1-1 Superior Latch Assembly Fastener Locations

Figure 5-3-1-2 Superior Latch Assembly


5-3-2 Right / Left Latch Assembly FRU Replacement

In order to perform the Right / Left Latch Assembly FRU Replacement, it will be necessary to remove the anterior coil cover from the anterior coil section. Follow the procedure detailed in section 5-1 Disassembly of Coil, then continue below.

1. Open the Superior Latch Assembly FRU package, and find the FRU Replacement Instructions, illustrated in Figure 5-3-1.

2. Follow the FRU Replacement Instructions while performing the Superior Latch Assembly replacement.

**Figure 5-3-1 Superior Latch Assembly
FRU Replacement Instructions**

SECTION 6 – RENEWAL PARTS

6-1 Field Replaceable Units

FIELD REPLACEABLE UNITS LIST – TABLE 6-1

Description	GE Part #	INVIVO #
1.5T 8 Channel Neurovascular Array	2317115-2	103045
Phantom Kit	2317115-3	103366
Cable	2317115-7	103025
Right / Left Latch Kit	2317115-10	105014
Superior Latch Kit	2317115-11	104975
17cm TLT Head Spherical Phantom	46-265826G6	N/a
10cm Spherical Phantom	46-317586G1	N/a
Large Cylindrical unified phantom	5342679	N/a
Phantom positioner, NV	5344806	N/a

Cables bearing this Coil ID label are compatible with software versions 10.X and newer. With software versions 11.X and newer, the system will generate an error when an incompatible coil is prescribed.

6-2 Other Replaceable Accessories

OTHER REPLACEABLE ACCESSORIES LIST – TABLE 6-2

Description	GE Part #	INVIVO #
Wedge Pad	2317115-8	101407
HNC Pad	2317115-6	102955

SECTION 7 – APPENDIX

7-1 Schematic

Refer to Block Diagram under Section 1-5 Theory of Operation.

7-2 Coil Configuration

Check the system coil configuration file for the 8NVHEAD and 8NVARRAY coil configurations. If 8NVHEAD and 8NVARRAY have not been installed, select 8NVHEAD and 8NVARRAY from the list of available coils, and install.

Revision History

Rev	Date	Author	Primary Reason for Change
A	10-24-02	L. Hyler	Initial Release
B	03-22-03	L. Hyler	Update QA Procedure – New Phantom and Positioner
C	03-24-03	L. Hyler	Refine Tables 3-3-1 and 3-3-2 Update Figures 3-2-1-3 and 3-2-1-4
D	03-24-03	L. Hyler	Change directory names in Sections 3-2-1 and 3-2-3
E	03-26-03	L. Hyler	Add mirror instructions to 2-1-1 Special Install Notes
F	03-26-03	L. Hyler	Update Service Methods CD in Section 1-3
G	03-31-03	L. Hyler	In 3-2-2, Change snr_8nvarray_QA to snr_for_8nvarray_QA Update Section 6, Renewable Parts
01	04-01-03	L. Hyler	
02	04-14-03	L. Hyler	Edit Table 3-2-3-1 Figure Reference Replace Figures 3-2-3-2 and 3-2-3-3 Edit 3-2-3 Plot Results, Item 5, to match Figure 3-2-3-2 Button
03	05-06-03	L. Hyler	Change 8NVHEAD to 8NVARRAY in Figures 3-2-4-4, 3-2-4-5 and 3-2-4-6
04	07-09-03 09-11-03	L. Hyler	Add Figure 3-2-1-5 and Correct Bottle Placement NOTE Update TABLE 3-2-3-2 SNR LIMITS 3-2-2.1 correct command to snr_for_8nvarray_QA Update 2-1 Coil Installation [8chNV install_rev1.doc]
05	10-22-03	L. Hyler	Corrected Operator's Manual Part Number in Section 1-3 Updated Loop Designations in Figures 1-5-2 through 1-5-5 Revised TABLE 3-2-3-2 SNR LIMITS Redesigned section 3-2-1 for both version 10.X and 11.X s/w
06	11-19-03	L. Hyler	Add Coil ID Information in Sections 5-2 and 6-1 Update Photos of Phantom Positioner
07	09-02-04	L. Hyler	
08	07-06-05	D. Visscher	Rebrand MRI Devices with Invivo (DCN 1378)
9	6/19/08	V. Thacker	Customer Change Request CCN0021 (DCN3994)
10	3/28/08	R. Hamilton	Change text and images to illustrate new superior latch and mirror. (DCN 3754)
11	6/24/09	R. M. Scott	Add Coil Configuration (DCN 5847)
12	10/29/09	M. Vaughn	Update per Change Request CCN0027 (Unified phantom set)