

SIEMENS

MAGNETOM

MR

Troubleshooting Guide

Patient Handling Error Messages

For all Tim Systems
with Software Version syngo MR B17 (VB17A)

© Siemens, 2009

- For internal use only - All documents may only be used by authorized personnel for rendering services on Siemens Healthcare Products. Any document in electronic form may be printed once. Copy and distribution of electronic documents and hardcopies is prohibited. Offenders will be liable for damages. All other rights are reserved.

Print No.: MR00-000.840.15.03.02
Replaces: MR00-000.840.15.02.02

English
Doc. Gen. Date: 03.12

Copyright / Version / Disclaimer

Copyright

"© Siemens, 2009" refers to the copyright of a Siemens entity such as Siemens Aktiengesellschaft - Germany, Siemens Shenzhen Magnetic Resonance Ltd. - China, Siemens Shanghai Medical Equipment Ltd. - China, Siemens Medical Solutions USA Inc. - USA and/or Siemens Healthcare Diagnostics Inc. - USA.

Document Version

Siemens reserves the right to change its products and services at any time.

In addition, manuals are subject to change without notice. The hardcopy documents correspond to the version at the time of system delivery and/or printout. Versions to hardcopy documentation are not automatically distributed. Please contact your local Siemens office to order current version or refer to our website <http://www.healthcare.siemens.com>.

Disclaimer

Siemens provides this documentation "as is" without the assumption of any liability under any theory of law.

The service of equipment described herein is to be performed by qualified personnel who are employed by Siemens or one of its affiliates or who are otherwise authorized by Siemens or one of its affiliates to provide such services.

Assemblers and other persons who are not employed by or otherwise directly affiliated with or authorized by Siemens or one of its affiliates are not entitled to use this documentation without prior written authority.

History - Changes to previous version

Date	Author	Version	Changes
03.09	CS SD MR	1.0	Initial version
05.11	CS PS MR	2.0	Changes as described in Charm ID: MR_00378254
03.12	CS SD MR	3.0	Copyright and Content added.

Content of Patient Handling Error Messages

Component	Source Name	Total numbers
Patient Handling	PXX.mc	238

Patient Handling

(Message clas S/I/W/E)	Id	Source	Message Text	Explanation	Action	Information
Error	22	PXX.mc	Patient Table error: MDSD Software-Version incompatible Please save System Log-Files immediately by selecting System/Control/Tools. Do not close any message window, before the function has been finished successfully. Please call Siemens Service	Software version compatibility test between PTAB and MDSD Incompatibility error. Hostsoftware is not updated (no patch). Under normal circumstances, the error will not occur	- Update Hostsoftware Please call Siemens Service	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"MDSD SW-Vers incompatible"<<<
Error	28	PXX.mc	Patient Table error: Internal Error Please save System Log-Files immediately by selecting System/Control/Tools. Do not close any message window, before the function has been finished successfully. Please call Siemens Service.	Under normal circumstances, the error will not occur	Please call Siemens Service	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Internal Error PT005"<<<
Warning	29	PXX.mc	Patient Table warning: EEPROM Error CONT A4140 (Defaults will be used) Please save System Log-Files immediately by selecting System/Control/Tools. Do not close any message window, before the function has been finished successfully. Please call Siemens Service.	Error during access to parameter memory chip on Controller A4140. Possible Causes: - hardware failure - checksum failure - parameters where erased by command PATB-Software will use default parameters (i.e. length, light marker distance...)	Only if error occurs frequently! Exchange of: - A4140 (CONT) ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"EEPROM Error CONT A4140"<<<
Warning	30	PXX.mc	Patient Table warning: EEPROM Error PWRD A4110 (Defaults will be used) Please save System Log-Files immediately by selecting System/Control/Tools. Do not close any message window, before the function has been finished successfully. Please call Siemens Service.	Error during access to parameter memory chip on PWRD A4110. Possible Causes: - hardware failure - checksum failure PATB-Software will use default parameters	Only if error occurs frequently! Exchange of: - A4110 (PWRD) ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"EEPROM Error PWRD A4110"<<<
Warning	32	PXX.mc	Patient Table warning: EEPROM ERROR SUPF A4120 Defaults will be used Please save System Log-Files immediately by selecting System/Control/Tools. Do not close any message window, before the function has been finished successfully. Please call Siemens Service.	Error during access to parameter memory chip on SUPF A4120. Possible Causes: - hardware failure - checksum failure PATB-Software will use default parameters	Only if error occurs frequently! Exchange of: - A4120 (CONT) ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"EEPROM ERROR SUPF A4120"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	37	PXX.mc	Patient Table error: Communication Host fault Reboot Scanner, switch off/on Scanner. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	The communication to the host is also monitored by the CAN-Nodeguarding from the table.		In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Communication Host fault"<<<
Error	38	PXX.mc	Patient Table error: Communication MDSD fault To reset error on Patient Table, press Patient Table IN and OUT buttons in succession! If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	A serial interface ("private") exists between PTAB and MDSD. Frequent communication error reported.	Error occurs frequently? Check: - W4115 and all connectors -H3 on A4110 lights up? (Fuse F17) on A4110 -Fuse F8 on A4110 (PWRD) -Fuse F10 on A4110 (PWRD) -H24 on A4110 lights up? -Fuse F14 on A4110	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Communication MDSD fault"<<<
Error	39	PXX.mc	Patient Table error: IIC-Bus external fault Reboot Scanner, switch off/on Scanner. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Continuous failure on the external IIC-Bus. Connected to the bus are: - Operating panel PAN.. (A4161..A4163 or 4561..63 or A5163..64) - Displays - Switch distribution SWT (A4115) - Air Temp. Meas. Unit ATMU (A4170)	Error localizing (see circuit diagram overview): - Disconnect each bus subscriber separately (and switch off ptab separately), Disconnect front operating panel, connect dummy plug (D-Sub 15 male, jumper Pin 3 to Pin 6) to A4110 (PWRD)-X15, replace rear operating panel with jumper A4110-J11 - switch on ptab, reboot, check for failure	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"IIC-Bus external fault"<<<
Error	40	PXX.mc	Patient Table error: IIC-Bus internal fault Reboot Scanner, switch off/on Scanner. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Continuous failure on the internal IIC-Bus. The following components are communicating on this bus: - Controller CONT(A4140..A4140) - Power&Drive PWRD (A4110) - Supplied functions SUPF (A4120)	Error localizing (see circuit diagram overview): - disconnect SUPF (W2122), reboot, check failure, (perform counter-check) if necessary replace, - exchange CONT, reboot, check failure, (perform counter-check) if necessary replace, - exchange PWRD, reboot, check failure, (perform counter-check) if necessary replace. Ptab must be switched off if any pcb will be exchanged!! ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"IIC-Bus internal fault"<<<
Error	43	PXX.mc	Patient Table error: IIC-Bus : Fault at Switch Terminal Reboot Scanner, switch off/on Scanner. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Check the position of the switch S1.1 (RTT, NO RTT) on the PWRD A4110, the cable W4141 and the SWT A4115.		In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Fault at SWT-IIC-Bus"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	53	PXX.mc	Patient Table error: Error searching reference driving direction into scanner To reset error on Patient Table, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	During initializing S10 is in use for referencing the z-position. If the ptab moves from home position (home position switch activated) into the magnet bore and S10 does not respond after a fixed distance, this error will show up.	Control of S10: - Show status of S10 with service menu ptab, move table by hand approx. 5 cm around the home position, S10 must react by "a click" and a status change must show up on the display. - Check free movement of S10 in every high, control free horizontal movement of the table top, - For reference drive, check sufficient gap between magnet cover and ptab - Check power transmission for horizontal axis If this error occurs together with horizontal reference and S10 respectively, please check the adjustment of the table in the following way: 1. Visually examine the covers along the seam to the table to ensure they are LEVEL. 2. Monitor vertical movement(s) at the narrowest position between the table and covers. Remain at this position. 3. Power off the table, manually move the tabletop to the outer mechanical end stop. 4. Power on the table and monitor the subsequent automatic tabletop movement for referencing: During the first reference phase, the tabletop moves toward the magnet or cover. Check: The table should NOT touch the cover. Shortly before reaching the cover, the table should move away from the magnet. 5. If this is not the case: The reference position has to be changed by shifting the cam for S10 (under the tabletop). This enables the tabletop to move without obstruction and extend out slightly. NOTE: If the reference position is changed for removable tabletop systems, the adjustment at the trolley has to be corrected as well. 6. Additional information for those particularly interested: Reference movement: Starting point: The table is at the mechanical end stop. After Power-On, automatic reference search movement begins: The table slowly moves in the direction of the magnet (toward the cover) until S10 is deactivated, then slowly moves outward until S10 is activated, then outward somewhat more quickly (securely on the switch), and then a few tenths of a mm inward, which is barely noticeable at the tabletop, but can be seen at the drive output shaft (compensation for slack).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Ref.search Error inwards"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	54	PXX.mc	Patient Table error: Error searching Reference driving direction to home position To reset error on Patient Table, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	During initializing S10 is in use for referencing the z-position. If the ptab is not in the upper and outside position (S10 not active), the ptab will move (in a limited way). During moving to the outside position the ptab will "search" S10. An error will be reported when S10 is not activated during the "search" ride.	Control of S10: - Show status of S10 with service menu ptab, move table by hand approx. 5 cm around the home position, S10 must react by "a click" and a status change must show up on the display. - Check sufficient distance (10mm), from the switch point of S10 to the mechanical stop, for table top movement. - Control free horizontal movement of the table top - For reference drive, check sufficient gap between magnet cover and ptab - Check proper activation of S10 from the cam of the table top - Check power transmission for horizontal axis If this error occurs together with horizontal reference and S10 respectively, please check the adjustment of the table in the following way: 1. Visually examine the covers along the seam to the table to ensure they are LEVEL. 2. Monitor vertical movement(s) at the narrowest position between the table and covers. Remain at this position. 3. Power off the table, manually move the tabletop to the outer mechanical end stop. 4. Power on the table and monitor the subsequent automatic tabletop movement for referencing: During the first reference phase, the tabletop moves toward the magnet or cover. Check: The table should NOT touch the cover. Shortly before reaching the cover, the table should move away from the magnet. 5. If this is not the case: The reference position has to be changed by shifting the cam for S10 (under the tabletop). This enables the tabletop to move without obstruction and extend out slightly. NOTE: If the reference position is changed for removable tabletop systems, the adjustment at the trolley has to be corrected as well. 6. Additional information for those particularly interested: Reference movement: Starting point: The table is at the mechanical end stop. After Power-On, automatic reference search movement begins: The table slowly moves in the direction of the magnet(toward the cover) until S10 is deactivated, then slowly moves outward until S10 is activated, then outward somewhat more quickly (securely on the switch), and then a few tenths of a mm inward, which is barely noticeable at the tabletop, but can be seen at the drive output shaft (compensation for slack).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Ref.search Err. outwards"<<<
Error	59	PXX.mc	Patient Table error: Switch off MDSD failed Reboot Scanner! If error occurs frequently, STOP Button is not safe anymore! If the problem still persists, please call Siemens Service	The power circuit of the MDSD can be switched ON or OFF (shut down) by SW. If the disconnect is not sufficient this error will occur. If the error occurs frequently a safe shut down by pressing the STOP-button is not guaranteed!	Only if error occurs frequently! Successive exchange (in this order) of: - A4110 (PWRD) - A4140 (CONT) - A4320 (MDSD) ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch off MDSD failed"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	60	PXX.mc	Patient Table error: Switch Off PSY failed! Reboot Scanner! If error occurs frequently, STOP Button is not safe anymore!	The power circuit of the lifting motion can be switched ON or OFF (shut down) by Software. If the shut down is not sufficient this error will occur. If the error occurs frequently a safe shut down by pressing the STOP-button is not guaranteed!	Only if error occurs frequently! Successive exchange (in this order) of: - A4110 (PWRD) - A4140 (CONT) First exchange A4110(PWRD) ONLY. If this does not solve the Problem, then exchange the A4140 (CONT) too! ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Off PSY failed!"<<<
Error	61	PXX.mc	Patient Table warning: STOP button pressed by user To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	The operator has pressed one of the STOP-buttons of the system. If this has not happened: - One of the control panels is creating an error - A line fault in one of STOP-(opening contact) circuits is present STOP-circuits: * Intercom (on A4110-X20) * Operating panel front (on A4110-X15) * Operating panel rear (on A4110-X16 or Jumper J11) * Switch distribution SWT A4115/16 (on A4110-X10) * Optional surgery table ORT (on A4110-X14 or Jumper J12)	Check STOP-buttons (press and release!): - on each operating panel, before reacting the STOP button must be pressed a min of 0,5 mm Attention: Error detection (during this test a safe movement of the ptab is not guaranteed)!! For testing (jumper out buttons or circuits) set jumpers as follow: A4110-X15 Pin 3 with Pin 6 A4110-X16 Pin 3 with Pin 6 or J11 A4110-X20 Pin 17 with Pin 19 A4110-X10 Pin 24 with Pin 25 A4110-X14 Pin 3 with Pin 6 or J12	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"STOP-But. pressed by User"<<<
Error	62	PXX.mc	Patient Table warning: STOP activated automatically To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Due to another error the SW of the ptab has released the STOP-function	- Solve origin error	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"STOP activated by autom."<<<
Warning	63	PXX.mc	Patient Table warning: Table Carrier Collision Lift up table using UP control and remove hindrance! If the problem still persists, please call Siemens Service	The eject of the ptab carrier from the spindle nut is monitored by switch S30. Normally S30 is actuated and electrical closed. If the carrier is eject from nut S30 is not actuated and electrical open. The goal is to soften any collision damage during table down movement.	- Drive table up by pressing the UP button, if not possible use crank - Remove any hindrance below the ptab	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Table Carrier Collision"<<<
Warning	64	PXX.mc	Patient Table warning: Table Plate Collision Lift up table using UP control and remove hindrance! If trolley is in use lower carefully! If the problem still persists, please call Siemens Service	The eject of the ptab carrier is monitored by the switch S31 and S32 (series connection) Actuated S31/S32 is open The goal is to soften any collision damage during table down movement.	- Drive table up by pressing the UP button, if not possible use crank - Remove any hindrance below the ptab	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Table Plate Collision"<<<

(Message clas	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Warning	66	PXX.mc	Patient Table warning: Lifting gear getting hot Use spare, only if needed ! No immediate action necessary. Please inform Siemens Service at next visit.	For the table down movement a brake with a self locking mechanism on the spindle is used. During table down movement, especially with heavy loads, the mechanism gets warm, the temperature is monitored by sensor S23. Before the function is limited this warning is created	- Use table UP- DOWN movement only if absolutely necessary!! - Consider cool down phase	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Lifting gear getting hot"<<<
Warning	67	PXX.mc	Patient Table warning: Lifting gear over temp. Wait for cooling down! If the problem still persists, please call Siemens Service	For the table down movement a brake with a self locking mechanism on the spindle is used. During table down movement, especially with heavy load, the mechanism gets warm, the temperature is monitored by sensor S23. This error occurs if the brake gets too hot (overtemp.). Only table UP movement is possible! Ready for MR-measurement	- Consider cool down phase	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Lifting gear over temp."<<<
Warning	68	PXX.mc	Patient Table warning: Sensor S23 is shorted. No immediate action is necessary! Please inform Siemens service at next visit!	Temperature sensor S23 (overtemp. spindle) is shorted (bypassed) detected by A4110 (PWRD).	- Check sensor S23, if necessary replace it (cable W4131 and S23 is one unit)	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Sensor S23 is shorted."<<<
Warning	69	PXX.mc	Patient Table warning: Sensor S23 is opened. No immediate action is necessary! Please inform Siemens service at next visit!	The temperature sensor S23 (overtemp. spindle) is not detected anymore. Probably a broken cable, or the faulty sensor is detected by A4110 (PWRD).	- Check sensor S23, if necessary replace it (cable W4131 and S23 is one unit)	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Sensor S23 is opened."<<<
Error	70	PXX.mc	Patient Table error: Vertical Controller Watchdog Reboot Scanner! If the problem still persists, please call Siemens Service	Under normal circumstances, the error will not occur	- Contact HSC in Erlangen - If error occurs frequently exchange A4140 (CONT) ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik)	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Vertical Controller Wdog"<<<
Error	71	PXX.mc	Patient Table error: No antivalence vertical encoder Reboot Scanner! If the problem still persists, please call Siemens Service	The signals of the position encoder are monitored against antivalence. An error will be generated if no antivalence exist.	Only if error occurs frequently! Successive exchange (in this order) of: - Vertical drive A4311 - A4140 (CONT) First exchange A4311(drive) ONLY. IF this does not solve the Problem, then exchange the A4140 (CONT) too! ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"No antivalence Vert.Enc."<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	72	PXX.mc	Patient Table error: Counting error 1 vertical encoder To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	The electrical impulse counting of the UP/DOWN movement is monitored against plausibility by a special hardware. An error will be generated if no plausibility exists. This error will be generated during vertical movement.	Only if error occurs frequently! Successive exchange (in this order) of: - A4311 (Vertical Drive) - A4140 (CONT) ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Count.Err.1 vert. enc."<<<
Error	73	PXX.mc	Patient Table error: Over voltage vertical powerstage To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Overvoltage is detected on the vertical final stage (A4110-PWRD) Possible causes: - Supply voltage too high - Too heavy load for slow down of the vertical movement, malfunction of spindle brake (energetic recovery) - faulty detection from A4110	If Error 81 (Y_SPINDLE_BRAKE_BREAKDOWN) occurs too, the most probably cause is a failure in spindle brake. If error (73) occurs frequently: - Check line voltage (230V) directly at table supply - Check spindle brake: (no power on Ptab, disconnect line) No load on the table, unbend belt of vertical drive, turn pulley of vertical drive by hand, the table must move downwards. A mechanical resistance is expected.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"OverVoltage Vert.PowStg"<<<
Error	74	PXX.mc	Patient Table error: Under voltage vertical power stage To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Undervoltage is detected on the vertical final stage (A4110-PWRD) Possible causes: - too low supply voltage - faulty hook up transistor on A4110 (PWRD) - faulty detection on A4110 If error 73 also occurs, consecutive error 74 is normal and without significance.	Check: - H22 on A4110 (on PWRD) lights up? - Fuse F12 on A4110 - Fuse F8 on A4110 - Supply voltage (230V) - Exchange A4110 (PWRD) if error occurs frequently ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Under VoltageVert.PoStg."<<<
Error	75	PXX.mc	Patient Table error: Short circuit phase A at vertical motor detected To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	The electrical impulse counting of the UP/DOWN movement is monitored against plausibility. The motor windings of the vertical drive are monitored against short circuit If faulty windings are detected, an error will be generated Possibility: Wrong detection on A4110	Only if error occurs frequently! Successive exchange (in this order) of: - Vertical drive A4311 - A4110 (PWRD) First exchange A4311(drive) ONLY. If this does not solve the Problem, then exchange the A4110 (PWRD) too! ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"PhA at Vert.Mot. shorted"<<<
Error	76	PXX.mc	Patient Table error: Short circuit phase B at vertical motor detected. To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	The motor windings of the vertical drive are monitored against short circuit. If faulty windings are detected, an error will be generated Possibility: Wrong detection on A4110	Only if error occurs frequently! Successive exchange (in this order) of: - Vertical drive A4311 - A4110 (PWRD) First exchange A4311(drive) ONLY. If this does not solve the Problem, then exchange the A4110 (PWRD) too! ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"PhB at Vert.Mot. shorted"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	77	PXX.mc	Patient Table error: Vertical position loop error overflow To reset table stop, press Patient Table IN and OUT buttons in succession! Check vertical movement for obstacles! If the problem still persists, please call Siemens Service	During the vertical movement the axis control is monitored against the target value. A too big difference between target and actual value (in the position) will lead to an error. Possible causes: - Tight movement of the mechanics - Wrong values from the position encoder cause a change in the motor control	- Check switching noise of motor brake If failure occurs with heavy loads during UP movement: - Check smooth movement of the mechanics (too heavy load, no gripping effect), "cracking" noise during vertical movement indicates twisted guideways. If failure occurs with heavy loads during DOWN movement: - Check brake torque of spindle brake (no load on the table, ptab mains supply OFF), release tooth belt (refer to REP manual) turn pulley, left/right against a small resistance for UP/ DOWN movement. The mechanical torque to move the table UP or DOWN must be about equal. - Check the release of the motor brake: With released tooth belt (reconnect ptab mains supply), move the table UP by pressing the UP button, the brake will release by a "click" noise and the motor axle will turn with out squeak noise. Error occurs without load and UP movement: - Exchange vertical drive A4311 ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Vert. Pos.Loop.Overflow"<<<
Error	78	PXX.mc	Patient Table error: Vertical over speed Error To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	spare (place holder) for future use		In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Vertical over speed Err."<<<
Error	79	PXX.mc	Patient Table error: Vertical commutation find error Reduce weight on table! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	During the first switch ON or loadware restart the internal motor initialization (commutation) of the vertical drive failed. Possible causes: - Tight movement of the mechanics - Position encoder failed - Error on vertical final stage at A4110 - Error on drive A4311	Error occurs during vertical movement with load: - Check smooth movement of the mechanics (too heavy load, gripping effect), "cracking" noise during vertical movement indicates twisted guideways. - Check brake torque of spindle brake (no load on the table, PTAB mains supply OFF), release tooth belt (refer to REP manual) turn pulley, left/right against a small resistance for UP/ DOWN movement. The mechanical torque to move the table UP or DOWN must be equal. Error occurs without load and during vertical movement: - Exchange vertical drive A4311 - if no effect exchange A4110 ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Verti. Commu. Find Error"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	80	PXX.mc	Patient Table error: Over temperature vertical drive To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	During vertical movement the motor temperature will be calculated (forecasted). Switch off by high motor temperature Possible causes: - To many moving cycles (20 times UP/DOWN) with heavy load - Tight/jammed mechanics - Tight spindle brake - Unreleased motor brake - Error on vertical final stage	Only if vertical movement is not in heavy use: - Check smooth movement of the mechanics (too heavy load, gripping effect), "cracking" noise during vertical movement indicates twisted guideways. - Check brake torque of spindle brake (no load on the table, ptab mains supply OFF), release tooth belt (refer to REP manual) turn pulley, left/right against a small resistance for UP/ DOWN movement. The mechanical torque to move the table UP or DOWN must be about equal. - Check the release of the motor brake: With released tooth belt (reconnect ptab mains supply), move the table UP by pressing the UP button, the brake will release by a "click" and the motor axle will turn with out squeaking noise. Only if error "Vertical Position Loop Error overflow" occurs frequently and no warning "Lifting Gear Over Temperature" is present: Exchange A4110 (PWRD) ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Over Temp. Vert. Drive"<<<
Error	81	PXX.mc	Patient Table error: Breakdown Spindle Brake Reboot Scanner! If the problem still persists, please call Siemens Service and remove patient from the patient table	The brake effect of the spindle brake is monitored during DOWN movement. An error occurs if the brake effect is to low. If error occurs frequently: Attention ptab might crash down! Remove patient from the ptab (system)!	Check brake torque of spindle brake (no load on the table, ptab mains supply OFF), release tooth belt very carefully (beware of table carrier moving down! Refer to REP manual) turn pulley, left/right against a small resistance for UP/ DOWN movement. The mechanical torque to move the table UP or DOWN must be about equal.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Breakdown Spindle Brake"<<<
Error	82	PXX.mc	Patient Table error: Breakdown of vertical drive brake Reboot Scanner! If the problem still persists, please call Siemens Service and remove patient from the patient table	During switch on and restart loadware the effect of the motor brake is monitored. An error occurs if the effect of the motor brake is not sufficient. Attention: The lift axis (UP/DOWN) is no more double safe against crash down! If error occurs frequently: Attention! PTAB might crash down! Remove patient from the table!	Check braking effect of motor brake: - Release tooth belt (refer to REP manual) turn pulley by hand is almost impossible, heavy torque (brake torque >3Nm). To turn the motor shaft manually a torque of 3Nm is necessary. If brake effect is too low: - Exchange vertical drive A4311 ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Vert. Drive Brake Fail"<<<
Informational	83	PXX.mc	Patient Table information: Vertical Drive Retraction Failure1 (Information) Check Weight on Table ! If the problem still persists, please call Siemens Service	During switch on and restart loadware the effect of the motor brake is monitored. Among other things the vertical Axis is moved up a very small distance (motor brake is opened) and than the same distance down. An error occurs if this downward movement is not possible.	Check braking effect of spindle brake: - Release tooth belt (refer to REP manual) turn pulley at spindle (downwards) by hand is possible with a max. brake torque of 10 Nm, typically 4-6 Nm, depending on load) If brake effect is too high: - make sure, that during power on and restart loadware weight on table is less than specified - Check the release of the motor brake: a "click" noise during test at power on and restart loadware must be present - Exchange vertical drive only, if there is no "click" noise; otherwise - Exchange spindle	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Vert.Dr.Retraction Fail1"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	84	PXX.mc	Patient Table error: Vertical Load Overweight To reset table stop, press Patient Table IN and OUT buttons in succession! Check weight on table If the problem still persists, please call Siemens Service	The vertical axis could not move UP during the first switch on or restart loadware. Possible causes: - Too heavy load - Tight/jammed mechanics - Motor brake cannot be released	If no overload: - Check smooth movement of the mechanics (too heavy load, no gripping effect), "cracking" noise during vertical movement indicates twisted guideways. - Check the release of the motor brake: A "click" noise during test at power on and restart loadware must be present	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Vertical Load Overweight"<<<
Error	85	PXX.mc	Patient Table error: Breakdown Vertical Belt Reboot Scanner! If problem still persists please call Siemens Service and remove patient from the patient table.	Broken toothbelt of the vertical drive Attention: The lift axle (UP/DOWN) is no more doubly safe against crash down! If error occurs frequently: Attention ptab might crash down! Remove patient from the ptab (system)!	Check toothbelt of vertical drive	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Breakdown Vertical Belt"<<<
Error	87	PXX.mc	Patient Table error: Counting error 3 vertical encoder To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	The electrical impulse counting of the UP/DOWN movement is monitored against plausibility by software. An error will be generated if no plausibility exist. This error will generated during movement of vertical Axis.	Only if error occurs frequently! Successive exchange (in this order) of: - A4311 (Vertical Drive) - A4140 (CONT) ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Count.Err.3 vert. enc."<<<
Informational	88	PXX.mc	Patient Table information: Vertical Drive Retraction Failure2 (Information) Check Weight on Table !	During the switch on and Restart Loadware the effect of the motor brake is monitored. An error occurs if this downward movement is not possible. Among other things the vertical Axis is moved up a very small distance (motor brake is opened) and than the same distance down. An error occurs if this downward movement is not possible.	Check braking effect of Spindle brake: - Release tooth belt (refer to REP manual) turn pulley at spindle (downwards) by hand is possible with a max. brake torque of 10 Nm, typically 4-6 Nm, depends on load) If brake effect is too high: - make sure, that during power on and restart loadware weight on table is less than specified - Check the release of the motor brake: A "click" noise during test at power on and restart loadware must be present - Exchange spindle brake	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Vert.Dr.Retraction Fail2"<<<
Error	89	PXX.mc	Patient Table error: Internal Error: Watchdog MDSD Reboot Scanner! If the problem still persists, please call Siemens Service	Internal error MDSD, should never occur	If error occurs frequently: - Exchange A4320 (MDSD) Contact HSC in Erlangen ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Internal Error: WD MDSD"<<<
Error	90	PXX.mc	Patient Table error: Horizontal Controller Watchdog Reboot Scanner! If the problem still persists, please call Siemens Service	Internal error MDSD, should never occur	If error occurs frequently: - Exchange A4320 (MDSD) Contact HSC in Erlangen ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Horizontal Cont. Wdog"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	91	PXX.mc	Patient Table error: No antivalence horizontal encoder Reboot Scanner! If the problem still persists, please call Siemens Service	The signals of the MDSD position encoder are monitored against antivalence. An error will be generated if no antivalence exist.	Only if error occurs frequently! - Exchange A4320 (MDSD) ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"No antivalenceHoriz.Enc."<<<
Error	92	PXX.mc	Patient Table error: Counting error horizontal encoder Reboot Scanner! If the problem still persists, please call Siemens Service	The electrical impulses counted from the position encoder of the horizontal drive (MDSD) is monitored against plausibility. An error will be generated if no plausibility exist.	Only if error occurs frequently! - Exchange A4320 (MDSD) ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Count. Err. horiz. enc."<<<
Error	93	PXX.mc	Patient Table error: Over temperature horizontal power stage To reset table stop, press Patient Table IN and OUT buttons in succession! Check horizontal movement for obstacles! If the problem still persists, please call Siemens Service	Overttemperature of the horizontal final stage in the MDSD detected. Possible cause: - Tight movement of the horizontal axis, many sequences of movement with heavy loads. - Misadjustment of the horizontal table alignment. For the IN/OUT movement a too high permanent force is necessary (misadjustment table driving uphill) - Failure in MDSD	- Check smooth movement of horizontal axis (IN/OUT): Move table to magnet position and press STOP button, move table by hand IN and OUT of the magnet, moving force must be <180 N - Check adjustment to make sure the table is level and not driving uphill	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Over Temp. Horiz.PowStg."<<<
Error	94	PXX.mc	Patient Table error: Under voltage MDSD auxiliary supply To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	24 Volt auxiliary supply of the MDSD measured too low (Mains on PWRD: +24V_DC_MDSD) Possible causes: - Short or overload on one of the 24V mains - Supply voltage too low (230V) - Contact problem of W4115 - Failure on the 24V mains of A4110 (PWRD) - MDSD failure On A4110 self resetting fuses F17 and F8 are introduced Self resetting is active after a proper cool down, switch OFF load, power OFF, phase. On A4110 self resetting fuses F17 and F8 are introduced Self resetting is active after a proper cool down, switch OFF load, power OFF, phase.	Check with and without MDSD connected (disconnect on A4110-X13): - LED A4110-H3 lights? - Check voltage: A4110-X4.1 to X4.2, target: 18,6V AC - Check voltage: A4110-X13.3 to X13.5, target: 18..27V DC	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Under Volt. MDSD Aux.Sup."<<<
Error	95	PXX.mc	Patient Table error: Over Voltage Horizontal Powerstage To reset table stop, press Patient Table IN and OUT buttons in succession! Check weight on table ! If the problem still persists, please call Siemens Service	Overvoltage on the horizontal final stage (in MDSD) detected. Possible causes: - Too high supply voltage - Too heavy load on the horizontal axis - Too heavy load for slow down of the horizontal movement (energetic recovery) - Detection failure in MDSD	If error occurs frequently: - Check mains supply (230V)	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"OverVoltage Horiz.PowStg"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	96	PXX.mc	Patient Table error: Under voltage horizontal power stage To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Undervoltage on the horizontal final stage (in MDSD) detected. Possible causes: - Contac problem on cable W4110 (PWRD) - Failed transistor on A4110 (PWRD) - Too low supply voltage - Detection failure in MDSD	Check: - A4110-H24 lights up ? (on PWRD) - Fuse F14 on A4110 - Fuse F10 on A4110 - Fuse F12 on A4110 If error occurs frequently: - Check W4115 - Supply voltage (230V) If error is permanent: Exchange A4110 (PWRD) Exchange A4320 (MDSD) ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Under Volt.Horiz.PowStg."<<<
Error	97	PXX.mc	Patient Table error: Short circuit phase A at horizontal motor detected To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	The motor windings of the horizontal drive in the MDSD are monitored against short circuit. If a fault is detected, an error will be generated.	If error occurs frequently: Exchange A4320 (MDSD) ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"PhA at HorizMot. shorted"<<<
Error	98	PXX.mc	Patient Table error: Short circuit phase B at horizontal motor detected To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	The motor windings of the horizontal drive in the MDSD are monitored against short circuit. If a fault is detected, an error will be generated.	If error occurs frequently: Exchange A4320 (MDSD) ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"PhB at HorizMot. shorted"<<<
Error	99	PXX.mc	Patient Table error: Horizontal position loop error overflow To reset table stop, press Patient Table IN and OUT buttons in succession! Check horizontal movement for obstacles! If the problem still persists, please call Siemens Service	During the horizontal movement the axis control is monitored against the target value. A too big difference between target and actual value (in the position) will lead to an error. Possible causes: - Horizontal movement not free, - Too heavy load ! - Tight/jammed mechanics	- Check for obstacles - Check for free movement of the mechanics (too heavy load, no gripping effect, any noise) - Check smooth movement of the mechanic (to heavy load, no gripping effect, any noise) - Check smooth movement of horizontal axis (IN/OUT): Move table to magnet position and press STOP button, move table by hand IN and OUT of the magnet, moving force must be <180 N - Check "moving in height", table high enough?	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Horiz. Pos.Loop.Overflow"<<<
Error	101	PXX.mc	Patient Table error: Horizontal commutation find Error To reset table stop, press Patient Table IN and OUT buttons in succession! Check horizontal movement for obstacles! If the problem still persists, please call Siemens Service	During the first switch ON the internal motor initialization (commutation) of the horizontal drive failed Possible causes: - Tight movement of the mechanics - Position encoder failed - Error on horizontal final stage	- Check horizontal movement for obstacles! - Check smooth movement of the mechanics (too heavy load, no gripping effect, any noise) Error occurs without load and during horizontal movement: - Exchange horizontal drive A4320 ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenaer Antriebstechnik).	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Horiz. Commu. Find Error"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	102	PXX.mc	Patient Table error: Nodeguarding MDSD failed Reboot Scanner! If the problem still persists, please call Siemens Service	Error occurs; If MDSD is connected to the Nodeguarding (via CAN-Bus to the host) and an error is detected.	Check CAN connection between MDSD and Host Contact HSC in Erlangen.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Nodeguarding MDSD failed"<<<
Error	103	PXX.mc	Patient Table error: Over temperature horizontal Motor To reset table stop, press Patient Table IN and OUT buttons in succession! Check horizontal movement for obstacles! If the problem still persists, please call Siemens Service	During horizontal movement the motor temperature will be calculated (forecasted) Switch off by high motor temperature Possible causes: - Too many moving cycles (20 times UP/DOWN) with heavy load without cool down phase - Tight/jammed mechanics - Misadjustment of the horizontal table alignment. For the IN/OUT movement a too high permanent force is necessary (misadjustment table driving uphill)	Only if horizontal movement was not in heavy use: - Check smooth movement of horizontal axis (IN/OUT): Move table to magnet position and press STOP button, move table by hand IN and OUT of the magnet, moving force must be <180 N - Check adjustment from the "moving in" height of the ptab - Check horizontal movement for obstacles!	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Over Temp. Horiz. Motor"<<<
Error	104	PXX.mc	Patient Table error: Unexpected reference pulse driving into scanner Reboot Scanner! If the problem still persists, please call Siemens Service	After the initial referencing the horizontal axis, the values of the reference switch S10 and the encoder from the MDSD, will be checked against plausibility. An error will be recorded if S10 is active (PTAB is not in the outer position)	- Check reference switch S10 for * actuation by foreign material i.e. wrong cable wiring * blocking * switch function (i.e. service mode ptab) - Check tooth belt/horizontal force transmission If ptab is on power, table top can be moved only couple of mm by hand - Drive table top 5 cm before the outer end position (z-direction) Block table movement to the inner position, try to move table (appr. 5 sec) to the inner position. No movement on the display visible. One of the "Horizontal Position Loop Error overflow" or "Over Temp. Horizontal Motor" error must be reported. If not, check the mech. power transmission	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Unexp.RefPuls into Scnr"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	104 continued				If this error occurs together with horizontal reference and S10 respectively, please check the adjustment of the table in the following way: 1. Visually examine the covers along the seam to the table to ensure they are LEVEL. 2. Monitor vertical movement(s) at the narrowest position between the table and covers. Remain at this position. 3. Power off the table, manually move the tabletop to the outer mechanical end stop. 4. Power on the table and monitor the subsequent automatic tabletop movement for referencing: During the first reference phase, the tabletop moves toward the magnet or cover. Check: The table should NOT touch the cover. Shortly before reaching the cover, the table should move away from the magnet. 5. If this is not the case: The reference position has to be changed by shifting the cam for S10 (under the tabletop). This enables the tabletop to move without obstruction and extend out slightly. NOTE: If the reference position is changed for removable tabletop systems, the adjustment at the trolley has to be corrected as well. 6. Additional information for those particularly interested: Reference movement: Starting point: The table is at the mechanical end stop. After Power-On, automatic reference search movement begins: The table slowly moves in the direction of the magnet (toward the cover) until S10 is deactivated, then slowly moves outward until S10 is activated, then outward somewhat more quickly (securely on the switch), and then a few tenths of a mm inward, which is barely noticeable at the tabletop, but can be seen at the drive output shaft (compensation for slack).	
Error	105	PXX.mc	Patient Table error: Missing reference pulse at home position Reboot Scanner! If the problem still persists, please call Siemens Service	After the initial referencing the horizontal axis, the values of the reference switch S10 and the encoder from the MDSD, will be checked against plausibility. An error will reported if S10 is not active in home position		In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Missing RefPuls at Home"<<<
Error	106	PXX.mc	Patient Table error: Reference at MDSD lost Reboot Scanner! If the problem still persists, please call Siemens Service	After the initial referencing the horizontal axis, the reference has to be obtained. If the reference has been lost of any causes this error will reported. Possible causes: - MDSD-Loadware was reset directly by host (via CAN Node ID4). - Other Failures in MDSD, reported by separate error messages.	- Look for other reported MDSD-Errors.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Reference at MDSD lost"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Informational	110	PXX.mc	Patient Table information: Counting error 2 vertical encoder Error was reset automatically	The electrical impulse counting of the UP/DOWN movement is monitored against plausibility by a special hardware. An error will be generated if no plausibility exist. This information will be generated if a problem occurs WITHOUT(!) use of vertical axis. Probably cause is signal disruption in succession of bad RF-Conditions.	Only if error occurs frequently! Successive exchange (in this order) of: - A4311 (Vertical Drive) - A4140 (CONT) ATTENTION: Any additional hints, regarding failures or errors descriptions of the exchanged pcb's, will be highly appreciated by JAT (Jenae)	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Count.Err.2 vert. enc."<<<
Error	111	PXX.mc	Patient Table error: Switch Test 01 failed Error during table top exchange Check left trolley hook! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combination of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: S48 is actuated (left hook unlocked) S46 is not actuated (left hook locked) Possible cause: - APC failed - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4152, W4142 minor possibility: - A4115/A4116 (SWT) failed	Check functions of S46 and S48 by using the service menu! Check additionally: - APC - W4152 - W4142	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 01 failed"<<<
Error	112	PXX.mc	Patient Table error: Switch Test 02 failed Error during table top exchange Check right trolley hook! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combination of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: S49 is actuated (right hook completely unlocked) S47 is not actuated (right hook completely locked) Possible cause: - APC failed - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4152, W4142 minor possibility: - A4115/A4116 (SWT) failed	Check functions of S47 and S49 by using the service menu! Check additionally: - APC - W4152 - W4142	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 02 failed"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	113	PXX.mc	Patient Table error: Switch Test 03 failed Error during table top exchange Table top does not rest on carrier plate, or hooks faulty! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTTsensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: Table top and carrier plate are already/still disconnected (S45), but the hooks are not locked (Trolley in close-up range of the table top) Possible cause: - APC failed - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4144, W4152, W4142 minor possibility: - A4115/A4116 (SWT) failed	Check functions of S45, S46, S47, S60 in APC by using the service menu! Check additionally: - APC - W4144 - W4152 - W4142	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 03 failed"<<<
Error	114	PXX.mc	Patient Table error: Switch Test 04 failed Error during table top exchange Check bracket of the patient comfort plugs! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: Plugs of patient comfort are connected (S61), but table top not in close-up range of carrier plate (S44) Possible cause: - bracket of patient comfort plugs is push up - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4144, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check functions of S44, S61 by using the service menu! Check additionally: - W4144 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 04 failed"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	115	PXX.mc	Patient Table error: Switch Test 05 failed Error during table top exchange Check table top and carrier plate for jammed items! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Se	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: S45 is actuated (same as table top on carrier) and S44 not actuated (same as table top not in close-up range of the carrier plate) (trolley not connected) Possible cause: - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4144, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check functions of S44 and S45 by using the service menu! Check additionally: - W4144 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 05 failed"<<<
Error	116	PXX.mc	Patient Table error: Switch Test 06 failed Error during table top exchange Check table top and carrier plate for jammed items! Check automatic plug control (APC)! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem st	Implausible switch combinations of the RTT-sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: APC connected (S60) and S44 is not ctuated (same as table top not in close-up range of the carrier plate) Possible cause: - APC failed - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4144, W4142, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check functions of S44 and S60 (in APC) by using the service menu! Check additionally: - APC - W4144 - W4142 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 06 failed"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	117	PXX.mc	Patient Table error: Switch Test 07 failed Error during table top exchange Check table top and carrier for inclination! Check table top and carrier plate for jammed items! Check automatic plug control (APC)! To reset table stop, press Patient Table IN and	Implausible switch combination of RTT sensors. All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: S45 is actuated (table top on carrier plate) and APC (S60) is not connected (trolley is docked in) Possible cause: - APC failed - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4144, W4142, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check functions of S45 and S60 (in APC) by using the service menu! Check additionally: - APC - W4144 - W4142 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 07 failed"<<<
Error	118	PXX.mc	Patient Table error: Switch Test 08 failed Error during table top exchange Check table top and carrier plate for jammed items! Check automatic plug control (APC)! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - Table top not completely released from trolley (S40), not completely released table top (S60, S46..S49), not correctly rested on carrier plate (S45). - Carrier plate too high (S20, S21). - Trolley too high, in slow motion area Possible cause: - APC failed - micro switch/sensor failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4144, W4152, W4142, W4141, W4143 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S45, S21, S60, S46..S49 by using the service menu. Check additionally: - APC - W4144 - W4152 - W4142 - W4141 - W4143	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 08 failed"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	119	PXX.mc	Patient Table error: Switch Test 09 failed Error during table top exchange Check table top and carrier plate for jammed items! Check automatic plug control (APC)! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: table top not in close-up range of the trolley (S41 actuated) but at least one-sided detected (S43 not actuated OR S42 actuated) Possible cause: - APC failed - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4152, W4142, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S41, S43, S42 by using the service menu. Check additionally: - APC - W4152 - W4142 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 09 failed"<<<
Error	120	PXX.mc	Patient Table error: Switch Test 10 failed Error during table top exchange Check table top and carrier plate for jammed items! Check automatic plug control (APC)! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - table top not in close-up range of the trolley (S41 actuated) - but both-sided detected (S43 not actuated AND S42 actuated) Possible cause: - APC failed - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4152, W4142, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check function of S41, S42, S43 by using the service menu. Check additionally: - APC - W4152 - W4142 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 10 failed"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	121	PXX.mc	Patient Table error: Switch Test 11 failed Error during table top exchange Check Trolley-tong for items, sheets! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - trolley docked (S40 actuated) - trolley not detected (S33 actuated) S33 (=high) wie freie Sicht aber auch S40 (=high) wie Trolley da Possible cause: - cable failed: W4149, W4150, W4141 minor possibility: - A4115/A4116 (SWT) failed - A4140 failed	Check the function of S40, S33 by using the service menu. Check additionally: - LED A4115-H33 follows sensor S40 . - W4149 - W4150 - W4141	Sensor S33 has only power (H33) if trolley is not docked (S40). Both sensors must not be permanently high at the same time. In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 11 failed"<<<
Error	122	PXX.mc	Patient Table error: Switch Test 12 failed Error during table top exchange Check if trolley is docked correctly! Check left hook! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - trolley is undocked (S40=low) - trolley is detected as obstacle by optical sensor (S33=low) - left hook is still locked (S46 not actuated, APC connected) Possible cause: - APC failed - micro switch/sensor failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4144, W4152, W4142, W4149, W4150, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S40, S33, S46, S60 (in APC) by using the service menu. Check additionally: - APC - W4144 - W4152 - W4142 - W4149 - W4150 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 12 failed"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	123	PXX.mc	Patient Table error: Switch Test 13 failed Error during table top exchange Check if trolley is docked correctly! Check right hook! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - trolley is undocked (S40=low) - trolley is detected as obstacle by optical sensor (S33=low) - right hook is still locked (S47 not actuated, APC connected) Possible cause: - APC failed - micro switch/sensor failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4144, W4152, W4142, W4149, W4150, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S40, S33, S47, S60 (in APC) by using the service menu. Check additionally: - APC - W4144 - W4152 - W4142 - W4149 - W4150 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 13 failed"<<<
Error	124	PXX.mc	Patient Table error: Switch Test 14 failed Error during table top exchange Check if trolley is docked correctly! Check the hooks! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - trolley is undocked (S40=low) - trolley is detected as obstacle by optical sensor (S33=low) - both hooks are still locked (S46 and S47 not actuated, APC connected) Possible cause: - APC failed - micro switch/sensor failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4144, W4152, W4142, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S40, S33, S46, S47, S60 by using the service menu. Check additionally: - APC - W4144 - W4152 - W4142 - W4149 - W4150 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 14 failed"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	125	PXX.mc	Patient Table error: Switch Test 15 failed Error during table top exchange Check if there are contaminations in the trolley docking tong! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - trolley is docked (S40 actuated) - lifting column in lowest position (S22 actuated) Possible cause: S50 ist defect, therefore the table didn't stop during lowering. - micro switch/sensor failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4150, W4143, W4142, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S40, S22, S50 by using the service menu. If Errors 63 or 64 also occurs, probably S50 is defect. Check additionally: - W4150 - W4143 - W4142 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 15 failed"<<<
Error	126	PXX.mc	Patient Table error: Switch Test 16 failed Error during table top exchange Table top does not rest correctly on trolley or carrier plate. Check if items are jammed Check if trolley is leveled If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - APC connected (S60), - table top is not detected (S45) - The two plates shall be disconnected, but the table top is only one-sided connected with the trolley (S42 or S43) Possible cause: - APC failed - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4144, W4152, W4142, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S45, S42, S43, APC by using the service menu. Check additionally: - APC - W4144 - W4152 - W4142 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 16 failed"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	127	PXX.mc	Patient Table error: Switch Test 17 failed Error during table top exchange Table top does not rest correctly on trolley or carrier plate. Check if items are jammed Check if trolley is leveled If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - APC connected (S60), - carrier plate is not detected (S45). - The two plates shall be disconnected, but the table top is not connected with the trolley (S42 and S43) Possible cause: - APC failed - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4144, W4152, W4142, W4141 minor possibility: - A4115/A4116 (SWT) failed*	Check the function of S45, S42, S43, APC by using the service menu. Check additionally: - APC - W4144 - W4152 - W4142 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 17 failed"<<<
Error	128	PXX.mc	Patient Table error: Switch Test 18 failed Error during table top exchange Table top does not rest correctly on trolley or carrier plate. Check if items are jammed! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - APC connected (S60) - the carrier plate is not detected (S45). - The two plates shall be disconnected, but the table top is not in the close-up range of the trolley (s41) Possible cause: - APC failed - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4144, W4152, W4142, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S45, S41, APC by using the service menu. Check additionally: - APC - W4144 - W4152 - W4142 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 18 failed"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	129	PXX.mc	Patient Table error: Switch Test 19 failed Error during table top exchange The hooks are locked but trolley is not detected. Check if items are jammed! Check if trolley is levelled! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - Trolley is not in close-up range (S41 actuated) - both hooks are locked (S46 and S47 actuated) Possible cause: - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - APC failed - cable failed: W4144, W4152, W4142, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S41, S46, S47 by using the service menu. Check additionally: - APC - W4144 - W4152 - W4142 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 19 failed"<<<
Error	130	PXX.mc	Patient Table error: Switch Test 20 failed Error during table top exchange Check if items are jammed! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - too low for table top exchange (S50) - S44 is actuated permanently and contact is open respectively - S50 is actuated permanently and contact is open respectively Possible cause: - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4144, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S50 and S44 by using the service menu. Check additionally: - W4144 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 20 failed"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	131	PXX.mc	Patient Table error: Switch Test 21 failed Error during table top exchange The carrier plate is too low for table top exchange! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: The carrier plate is too low for table top exchange (S22). Possible cause: - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - APC failed - cable failed: W4143, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S22 by using the service menu. Check additionally: - W4143 - W4141 maybe also: - APC - W4144 - W4152 - W4142 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 21 failed"<<<
Error	132	PXX.mc	Patient Table error: Switch Test 22 failed Error during table top exchange Check bracket of the patient comfort plugs! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - patient comfort plugs are connected (S61) - table top has been disconnected from the carrier plate (APC, S60) Possible cause: - bracket of the patient comfort plugs is open - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4144, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S44, S60 by using the service menu. Check additionally: - W4144 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 22 failed"<<<
Error	133	PXX.mc	Patient Table error: Switch Test 23 failed Error during table top exchange The automatic plug control (APC) at the carrier plate is defect or disturbed. Check if containments are in the APC To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - APC connected (S60) - table top rests on carrier plate (S45) - S62 (further switch on APC) not actuated. - S62 in APC disturbed	Check the function of S62 (in APC) by using the service menu. Check additionally: - APC - W4142	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 23 failed"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	134	PXX.mc	Patient Table error: Switch Test 24 failed Error at vertical position detection! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of vertical axis All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - S20 is actuated (stop up position) - S21 is not actuated (not in slow motion mode) Possible cause: - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4143, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S20, S21 by using the service menu. Check additionally: - W4143 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 24 failed"<<<
Error	135	PXX.mc	Patient Table error: Switch Test 25 failed Error at vertical position detection! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of vertical axis All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - S20 is actuated (stop up position) - S22 is actuated (lowest position) Possible cause: - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4143, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S20, S22 by using the service menu. Check additionally: - W4143 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 25 failed"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	136	PXX.mc	Patient Table error: Switch Test 26 failed Error at vertical position detection! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of vertical axis All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: S21 is actuated (in slow motion mode) S22 is actuated (lowest position) Possible cause: - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4143, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S21, S22 by using the service menu. Check additionally: - W4143 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 26 failed"<<<
Error	137	PXX.mc	Patient Table error: Switch Test 27 failed Error at vertical position detection! To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	Implausible switch combinations of vertical axis All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: - S20 and S21 are actuated (stop up position) - S22 is actuated (lowest position) Possible cause: - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4143, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S20, S21, S22 by using the service menu. Check additionally: - W4143 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 27 failed"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	138	PXX.mc	Patient Table error: Switch Test 28 failed Error at vertical position detection and/or at PTAB exchange system! Check that no items are jammed at the foot of the plates. To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	implausible switch combinations of vertical axis/RTT sensors All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: S21 is actuated (in slow motion mode) S50 is actuated (trolley is detected) Possible cause: - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4143, W4144, W4141 minor possibility: - A4115/A4116 (SWT) failed	Check the function of S50, S21 by using the service menu. service menu. Check additionally: - APC - W4143 - W4144 - W4141	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 28 failed"<<<
Error	139	PXX.mc	Patient Table error: Switch Test 29 failed Error at vertical position detection To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	unerlaubter Zustandsübergang (Hubachse) All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible: i.e. it is not possible that the PTAB is in the upper and the lower position Possible cause: - micro switch failed - micro switch blocked - cam failed, bended - cam not proper adjusted - cable failed: W4143, W4141, W4150 minor possibility: - A4115/A4116 (SWT) failed	- Check proper functions and installation of micro switches - Check function of S20, S21, S22, S50 - Check proper functions of cams - Check cables (on older systems) - Check the signals with the service menu: release switch and check against display	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 29 failed"<<<
Error	140	PXX.mc	Patient Table error: Switch Test 30 failed Error at vertical position detection and/or at reference switch off the horizontal movement To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible switch combination between S10 and lift axle; The PTAB is close to the home position (lift area) and S10 reports, "ptab not in outer position"	- Check reference switch S10 Call service menu on the display and look for S10 signal, try to move table (appr. 5cm) around the home position, S10 must switch ("click") and the change of signal must be reported in the display - Check for enough distance after	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 30 failed"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Warning	141	PXX.mc	Patient Table warning: Emergency decoupling was activated by User Close clutch again! To reset table stop, press Patient Table IN and OUT buttons in succession! Drive to home position. If the problem still persists, please call Siemens Service	Decoupling mechanism between MDSD and table top is activated. Emergency recovery of the patient	- close clutch again, - to reset, press Patient Table IN and OUT buttons in succession, - required reference drive (i.e. press on operating panel home position)	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Emcy.DecouplingActivated"<<<
Warning	142	PXX.mc	Patient Table warning: Switch Test 32 failed Invalid patient table top type (insert) detected, system not configured for insert. Please replace table top with correct type. If the problem still persits, consider patient evacuation, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Invalid Patient Table Top Type (insert)	Possible actions: - Replace table top with correct type or - Configure system for this table top type.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Invalid Pat.Table Type"<<<
Error	143	PXX.mc	Patient Table error: Switch Test 33 failed Error at horizontal position detection on PET-Insert table slider To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persits, consider patient evacuation, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	All switch combinations are checked for plausibility, wrong switch combinations are reported as an error Implausible switch combination between S63 and S64 on PET Insert Table; "Switches report PTAB Insert Slider in both outer and inner position"	- Check S63 and S64 - Check cabling	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Test 33 failed"<<<
Error	144	PXX.mc	Patient Table error: Inconsistent Host Command Reboot Scanner! If the problem stilll persists, please call Siemens Service	Inconsistent Host Command	1. Switch off/on PTAB electrically 2. Reboot Scanner 3. If error occurs again, create Savelog and inform HSC for further investigation	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Inconsistent Host Cmd"<<<
Informational	146	PXX.mc	Patient Table information: Low battery in trolley shuttle Replace with new battery! If the problem still persists, please call Siemens Service	The charge state of the batteries from the trolley remote control is monitored. Low charge state reported.	Exchange batteries from the trolley remote control	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Low Battery in Trolley"<<<
Warning	150	PXX.mc	Patient Table warning: Over Temperature SUPF Tunnel light was reduced automatically .	Heat sink on A4120 too hot Brightness of tunnel light was reduced automatically	Check: - ambience temperature (air conditioning) - cover aeration o.k.? - cover louvers not blocked ?	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Over Temperatur SUPF"<<<
Warning	151	PXX.mc	Patient Table warning: Heatsink at PWRD (down) getting hot. Tunnel light was reduced automatically.	Temperature warning heatsink A4110 (PWRD). Brightness of tunnel light was reduced automatically	Check: - ambience temperature (air conditioning) - cover aeration o.k.? - cover louvers not blocked ?	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"PWRD (down) getting hot"<<<
Error	152	PXX.mc	Patient Table error: Over Temperature PWRD (down) Wait for cooling down! After: reset table stop! If the problem still persists, please call Siemens Service	Heatsink on A4110 to hot. Automatic table stop activated	Check: - ambience temperature (air conditioning) - cover aeration o.k.? - cover louvers not covered ?	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Over Temp. PWRD (down)"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Warning	153	PXX.mc	Patient Table warning: Heatsink at PWRD (up) getting hot. Tunnel Light was reduced automatically.	Temperature warning heatsink A4110 (PWRD). Tunnel light was reduced automatically	Check: - ambience temperature (air conditioning) - cover aeration o.k.? - cover louvers not blocked ?	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"PWRD (up) getting hot"<<<
Error	154	PXX.mc	Patient Table error: Over Temperature PWRD (up) Wait for cooling down! After: reset table stop! If the problem still persists, please call Siemens Service	Heatsink on A4110 to hot. Automatic table stop activated	Check: - ambience temperature (air conditioning) - cover aeration o.k.? - cover louvers not blocked ?	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Over Temp. PWRD (up)"<<<
Warning	155	PXX.mc	Patient Table warning: Temperature sensor at SUPF is shorted. No immediate action necessary! Please inform Siemens Service at next visit!	Reported short circuit on temperature sensor from heatsink A4120 (SUPF)		In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Temp.Sens. SUPF shorted."<<<
Warning	156	PXX.mc	Patient Table warning: Temperature Sensor at SUPF is open. No immediate action necessary. Please inform Siemens service at next visit!	Disconnect of the temperature sensor from the heatsink A4120 (SUPF) detected.		In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Temp. Sens. SUPF opened"<<<
Warning	157	PXX.mc	Patient Table warning: Temperature sensor at PWRD_0 is shorted. No immediate action necessary. Please inform Siemens service at next visit!	Short circuit on the temperature sensor from the heatsink A4110 (PWRD) detected.		In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Temp.Sens. PWRD_0 short"<<<
Warning	158	PXX.mc	Patient Table warning: Temperature sensor at PWRD_0 is open. No immediate action necessary. Please inform Siemens service at next visit!	Disconnect of the temperature sensor from the heatsink A4110 (PWRD) detected.		In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Temp.Sens. PWRD_0 open"<<<
Warning	159	PXX.mc	Patient Table warning: Temperature Sensor at PWRD_1 is shorted. No immediate action necessary. Please inform Siemens service at next visit!	Short circuit on the temperature sensor from the heatsink A4110 (PWRD) detected.		In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Temp.Sens. PWRD_1 short"<<<
Warning	160	PXX.mc	Patient Table warning: Temperature sensor at PWRD_1 is open. No immediate action necessary. Please inform Siemens service at next visit!	Disconnect of the temperature sensor from the heatsink A4110 (PWRD) detected.		In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Temp.Sens. PWRD_1 open"<<<
Error	163	PXX.mc	Patient Table error: Gradient coil temperature sensor "PRI_X" is shorted. Reboot Scanner! If the problem still persists, please call Siemens Service	Short circuit from the corresponding temperature sensor of the gradient coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"G.Temp.Sens PRI_Xshorted"<<<
Error	164	PXX.mc	Patient Table error: Gradient coil temperature sensor "PRI_Y" is shorted. Reboot Scanner! If the problem still persists, please call Siemens Service	Short circuit from the corresponding temperature sensor of the gradient coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"G.Temp.Sens PRI_Yshorted"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	165	PXX.mc	Patient Table error: Gradient coil temperature sensor "PRI_Z" is shorted. Reboot Scanner! If the problem still persists, please call Siemens Service	Short circuit from the corresponding temperature sensor of the gradient coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"G.Temp.Sens PRI_Zshorted"<<<
Error	166	PXX.mc	Patient Table error: Gradient coil temperature sensor "SEC_X" is shorted. Reboot Scanner! If the problem still persists, please call Siemens Service	Short circuit from the corresponding temperature sensor of the gradient coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"G.Temp.Sens SEC_Xshorted"<<<
Error	167	PXX.mc	Patient Table error: Gradient coil temperature sensor "SEC_Y" is shorted. Reboot Scanner! If the problem still persists, please call Siemens Service	Short circuit from the corresponding temperature sensor of the gradient coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"G.Temp.Sens SEC_Yshorted"<<<
Error	168	PXX.mc	Patient Table error: Gradient coil temperature sensor "SEC_Z" is shorted. Reboot Scanner! If the problem still persists, please call Siemens Service	Short circuit from the corresponding temperature sensor of the gradient coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"G.Temp.Sens SEC_Zshorted"<<<
Error	169	PXX.mc	Patient Table error: Gradient coil temperature sensor "SOCKET_1" is shorted. Reboot Scanner! If the problem still persists, please call Siemens Service	Short circuit from the corresponding temperature sensor of the gradient coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"BC.Temp. Sens 1 shorted"<<<
Error	170	PXX.mc	Patient Table error: Gradient coil temperature sensor "SOCKET_2" is shorted. Reboot Scanner! If the problem still persists, please call Siemens Service	Short circuit from the corresponding temperature sensor of the gradient coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"BC.Temp. Sens 2 shorted"<<<
Error	171	PXX.mc	Patient Table error: Gradient coil temperature sensor "PRI_X" is open. Reboot Scanner! If the problem still persists, please call Siemens Service	Broken wire from the corresponding temperature sensor of the gradient coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"G.Temp.Sens PRI_X open"<<<
Error	172	PXX.mc	Patient Table error: Gradient coil temperature sensor "PRI_Y" is open. Reboot Scanner! If the problem still persists, please call Siemens Service	Broken wire from the corresponding temperature sensor of the gradient coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"G.Temp.Sens PRI_Y open"<<<
Error	173	PXX.mc	Patient Table error: Gradient coil temperature sensor "PRI_Z" is open. Reboot Scanner! If the problem still persists, please call Siemens Service	Broken wire from the corresponding temperature sensor of the gradient coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"G.Temp.Sens PRI_Z opened"<<<

(Message clas	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	174	PXX.mc	Patient Table error: Gradient coil temperature sensor "SEC_X" is open. Reboot Scanner! If the problem still persists, please call Siemens Service	Broken wire from the corresponding temperature sensor of the gradient coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"G.Temp.Sens SEC_X opened"<<<
Error	175	PXX.mc	Patient Table error: Gradient coil temperature sensor "SEC_Y" is open. Reboot Scanner! If the problem still persists, please call Siemens Service	Broken wire from the corresponding temperature sensor of the gradient coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"G.Temp.Sens SEC_Y opened"<<<
Error	176	PXX.mc	Patient Table error: Gradient coil temperature sensor "SEC_Z" is open. Reboot Scanner! If the problem still persists, please call Siemens Service	Broken wire from the corresponding temperature sensor of the gradient coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"G.Temp.Sens SEC_Z opened"<<<
Error	177	PXX.mc	Patient Table error: Gradient coil temperature sensor "SOCKET_1" is open. Reboot Scanner! If the problem still persists, please call Siemens Service	Broken wire from the corresponding temperature sensor of the body coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"BC.Temp. Sens 1 opened"<<<
Error	178	PXX.mc	Patient Table error: Gradient coil temperature sensor "SOCKET_2" is open. Reboot Scanner! If the problem still persists, please call Siemens Service	Broken wire from the corresponding temperature sensor of the body coil detected.	Install the spare sensors according to document "Replacement of Parts"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"BC.Temp. Sens 2 opened"<<<
Error	180	PXX.mc	Patient Table error: Self test of air temperature measurement unit failed Wait 15 minutes before starting test again. If the problem still persists, please call Siemens Service		If error occurs often change A4170 (ATMU) air temperature measurement unit	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Self test ATMU failed"<<<
Error	181	PXX.mc	Patient Table error: Self test of air temperature measurement not possible. Either unit or air is too warm. Wait 15 minutes before starting test again. If the problem still persists, please call Siemens Service	ATMU is tested by heating the sensor. The sensor has to cool down before test can be started again. If cabin air is warmer than 34°C test will not work.	Wait 15min. Start test once again!	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Too hot for selftest ATMU"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	183	PXX.mc	Patient Table error: Voltage 24V_MDSD is not ready (auxiliary voltage for horizontal drive) Reboot Scanner Switch off scanner, wait 30 sec, switch on scanner. If the problem still persists, please call Siemens Service	Breakdown of 24 Volt auxiliary supply +24V_DC_MDSD from A4110 detected (24 Volt auxiliary supply from MDSD) Possible causes: - Short circuit of mains supply - W4115 failed - Failed power supply circuit of 24 Volt auxiliary supply of A4110 (PWRD) - MDSD failed On A4110 self resetting fuses F17 and F8 are introduced Self resetting is active after a proper cool down, switch OFF load, power OFF, phase.	Check with and without MDSD connected (disconnect on A4110-X13): - LED A4110-H3 lights? - Check voltage: A4110-X4.1 to X4.2, target: 18,6V AC - Check voltage: A4110-X13.3 to X13.5, target: 18..27V DC	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Voltage 24V_MDSD not ready"<<<
Warning	184	PXX.mc	Patient Table warning: Voltage 24V_OR_T is not ready (auxiliary voltage for optional surgery table) Reboot Scanner Switch off scanner, wait 30 sec, switch on scanner. If the problem still persists, please call Siemens Service	Breakdown of 24 Volt auxiliary supply +24V_DC_OR_T from A4110 detected (24 Volt auxiliary supply for optional surgery table) Possible cause: - Short circuit of mains supply - W4115 failed - Failed power supply circuit of 24 Volt auxiliary supply of A4110 (PWRD) If surgery table is connected to X14: - cable failed, from the surgery table - Surgery table error Ignore error if no surgery table is connected! On A4110 self resetting fuses F15 and F8 are introduced Self resetting is active after a proper cool down, switch OFF load, power OFF, phase.	Check with and without surgery table connected (disconnect A4110-X14) - LED A4110-H1 lights? - Check voltage: A4110-X4.1 to X4.2, target: 18,6V AC - Check voltage: A4110-X13.3 to X13.5, target: 18..27V DC	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Voltage 24V_OR_T not rdy"<<<
Error	185	PXX.mc	Patient Table error: Voltage 24V_SW_T is not ready (supply for micro switches) Reboot Scanner Switch off scanner, wait 30 sec, switch on scanner. If the problem still persists, please call Siemens Service	Breakdown of 24 Volt auxiliary supply +24V_DC_SW_T from A4110 detected (24 Volt auxiliary supply for switch distribution SW_T, micro switches, except S33 und S40) Possible causes: - Short circuit from mains in one of the switches - W4141 failed, all cables to switches - A4115/16 (SW_T) failed - Failed power supply circuit of 24 Volt auxiliary supply of A4110 (PWRD) On A4110 self resetting fuses F16 and F8 are introduced Self resetting is active after a proper cool down, switch OFF load, power OFF, phase.	Check with and without connected load: - Check LEDs on A4110-H2 and A4115/4116-H32 disconnect plugs one at a time and check LED'S. With power ON ptab: - A4115/4116-X31 (RTT only, patient trolley) - A4115/4116-X35 (RTT only) - A4115/4116-X33 - A4115/4116-X37 - A4115/4116-X38 - A4115/4116-X39 - A4115/4116-X40 Power OFF ptab now and disconnect: - A4115/4116-X30 - A4110-X10	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Voltage 24V_SW_T not rdy"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	185 continued	PXX.mc			Power ON ptab and check LED's. Only for ptab with exchangeable table top is present: Change DIP switch on A4110 to "NO_RTT" (don't forget to switch back after test) Power ON ptab again When LED's glow brighter during disconnecting any of the connectors, there is a high possibility of a short circuit in one of these disconnected lines. Additional test's: - Check voltage: A4110-X4.1 to X4.2, target: ca. 18,6V AC - Check voltage: A4110-X10.10 to X10.11, target 18..27V DC	
Error	186	PXX.mc	Patient Table error: Voltage OUT_DC_Switches is not ready (supply voltage for optional sensors - for removable tabletop only) Reboot Scanner Switch off scanner, wait 30 sec, switch on scanner. If the problem still persists, please call Siemens Service	Breakdown of 24 Volt auxiliary supply OUT_DC_SWITICHES from A4110 detected (24 Volt auxiliary supply for optical sensor S33 only RTT removable table top) Possible causes: - Short circuit from mains in sensor S33 - W4141, W4149 failed - A4115/16 (SWT) failed - Failed power supply circuit of 24 Volt auxiliary supply of A4110 (PWRD) - MDSD failed On A4110 self resetting fuses F15 and F8 are introduced Self resetting is active after a proper cool down, switch OFF load, power OFF, phase.	Check with and without connected load: - Check LEDs on A4115 and 4116-H33 disconnect plugs from SWT in order. With power ON ptab: - A4115/4116-X41 (RTT only, removable table top) When LED's glows brighter during disconnecting any of these connectors, there is a high possibility of a short circuit in one of the disconnected lines. Additional test's: - Check voltage: A4110-X4.1 to X4.2, target: ca. 18,6V AC - Check voltage: A4110-X10.10 to X10.11, target 18..27V DC Power OFF ptab now and disconnect:	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"OUT_DC_Switches not rdy"<<<
Error	187	PXX.mc	Patient Table error: Voltage 14V_SWT not ready. Reboot Scanner. Switch off scanner, wait 30 sec, switch on scanner. If the problem still persists, please call Siemens Service	Breakdown of 14 Volt auxiliary supply +14V_DC_SWT from A4110 detected (14V-auxiliary supply for switch distribution, A4115 SWT generates internal 5 V) Possible causes: - Short circuit from mains(14V) - W4141 failed - A4115/16 (SWT) failed - Failed power supply circuit of 14 Volt auxiliary supply of A4110 (PWRD) On A4110 self resetting fuses F13 and F18 are introduced Self resetting is active after a proper cool down, switch OFF load, power OFF, phase.	Check: - LEDs A4110-H18 and A4115/4116-H30, -H31 lights? - W4141 - Voltage: A4110-X7.1 to X7.2, target: ca. 11V AC - Voltage: A4110-X10.24 to X10.22, target: 10..16V DC	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Voltage 14V_SWT not rdy"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	188	PXX.mc	Patient Table error: Voltage 14V_PAN_F not ready (supply voltage for front panels and display). Reboot Scanner. Switch off scanner, wait 30 sec, switch on scanner. If the problem still persists, please call Siemens Service	Breakdown of 14 Volt auxiliary supply +14V_DC_PAN_F from A4110 detected (14V-auxiliary supply for front display and panel) Possible cause: - Mains short circuit in one of these components - W4161, W4162, W4164 (W4561, W4562) failed - A4115/16 (SWT) failed - Failed power supply circuit of 14 Volt auxiliary supply of A4110 (PWRD) On A4110 self resetting fuses F11 and F18 are introduced Self resetting is active after a proper cool down, switch OFF load, power OFF, phase.	Check: - LEDs A4110-H17 lights? - Cable - Voltage: A4110-X71.1 to X7.2, target: ca. 11V AC - Voltage: A4110-X15.4 to X15.7, target: 10..16V DC It is possible to disconnect operator panel and display for error containment. The STOP circuit must be closed via a dummy plug (D-Sub 15 male, jumper pin3 to pin6) to drive the ptab (i.e.via remote) Only after switch OFF reconfigure	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Volt. 14V_PAN_F not rdy"<<<
Error	189	PXX.mc	Patient Table error: Voltage 14V_PAN_R not ready (supply voltage for rear panel and display) Reboot Scanner Switch off scanner, wait 30 sec, switch on scanner. If the problem still persists, please call Siemens Service	Breakdown of 14 Volt auxiliary supply +14V_DC_PAN_R from A4110 detected (14V-auxiliary supply for rear display and panel) Possible causes: - Mains short circuit in one of these components - W4163, W4165 (W4563) failed - A4115/16 (SWT) failed - Failed power supply circuit of 14 Volt auxiliary supply of A4110 (PWRD) On A4110 self resetting fuses F9 and F18 are introduced Self resetting is active after a proper cool down, switch OFF load, power OFF, phase.	Check: - LEDs A4110-H19 lights? - Cable - Voltage: A4110-X71.1 to X7.2, target: ca. 11V AC - Voltage: A4110-X16.4 to X16.7, target: 10..16V DC It is possible to disconnect operator panel and display for error containment. The STOP circuit must be closed via a dummy plug (D-Sub 15 male, jumper pin3 to pin6) to drive the ptab (i.e.via remote) Only after switch OFF reconfigure	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Volt. 14V_PAN_R not rdy"<<<
Warning	191	PXX.mc	Patient Table warning: Marking current situation as notable by user. This is a hint for debugging.	Marking current situation as notable by user. This is a hint for debugging and will lead to a log file entry.		In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Notable Mark by user"<<<
Error	192	PXX.mc	Patient Table error: Cannot reach Center! The required center position is outside of the valid range. Change patient positioning.	The area of moving the ptab is limited. The required center position is outside the valid range. This error will reported if the operator has marked a position outside of the moving area. Ptab will not start to move.	Change patient positioning.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Cannot reach Center!"<<<
Error	193	PXX.mc	Patient Table error: Cannot reach position! The required remote position is outside the valid range. Change patient positioning.	The area of moving the ptab is limited. The required remote position is outside the valid range. This error will reported if the host has requested movement to a position outside of the moving area. PTAB will not start to move.	Make sure, that host remote commands are working correct!	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Cannot reach position!"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	194	PXX.mc	Patient Table error: Button pressed at boot up Check panels in cabin To reset table stop, press Patient Table IN and OUT buttons in succession! If the problem still persists, please call Siemens Service	During power up and restart loadware must not be pressed any button for positioning table. Otherwise this error will be reported.	- Check panels for jammed moving buttons (UP, DOWN, SPEED, CENTER, HOME) To reset table stop, press Patient Table IN and OUT buttons in succession	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Button pressed at bootUp"<<<
Error	201	PXX.mc	OR Table error: Internal software error 1. Reset OR Table by pressing OR Table Stop button for approx. 10 seconds. 2. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Internal SW error, should not occur	1. Reset OR Table by pressing OR Table Stop button for approx. 10 seconds 2. Error report as detailed as possible about the source for error 3. Create Savelog and inform HSC for further investigation	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Internal Error"<<<
Error	202	PXX.mc	OR Table error: EEPROM malfunction 1. Reset OR Table by pressing OR Table Stop button for approx. 10 seconds. 2. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Contents of parameter-EEPROMS defective/not readable	1. Reset OR Table by pressing OR Table Stop button for approx. 10 seconds. 2. If error recurs, adjust OR Table again and store parameter 3. If error recurs despite save and restart, replace assembly B2100 4. Create Savelog and inform HSC for further investigation	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"EEPROM Error"<<<
Error	203	PXX.mc	OR Table error: Internal communication error 1. Reset OR Table by pressing OR Table Stop button for approx. 10 seconds. 2. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	IIC Communication does not function	1. Reset OR Table by pressing OR Table Stop button for approx. 10 seconds. 2. If error occurs continuously, check assemblies in the following order and replace if necessary: Cable W33, main operating console, assembly B2100 3. Create Savelog and inform HSC for further investigation	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"IIC-Bus internal fault"<<<
Error	205	PXX.mc	OR Table error: Host communication error 1. Reset OR Table by pressing OR Table Stop button for approx. 10 seconds. 2. If the problem still persists, reboot Scanner! 3. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Communication to the host is also monitored by the CAN-Node guarding from the table. Host connection to OR Table is relevant.		In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Communication Host fault"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	206	PXX.mc	OR Table error: Patient Table communication error 1. Reset OR Table by pressing OR Table Stop button for approx. 10 seconds. 2. If the problem still persists, reboot Scanner! 3. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	A serial interface ("private") exists between ORT and PTAB. Frequent communication error reported.	1. Reset OR Table by pressing OR Table Stop button for approx. 10 seconds. 2. Test PTAB (switched on, LW context..?) 3. If error occurs continuously, check assemblies and replace if necessary: Cable W40, PTAB (A4100), assembly B2100 4. Create Savelog and inform HSC for further investigation	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Communication PTAB fault"<<<
Error	208	PXX.mc	OR Table error: Internal D-switch error This error will also happen when the lifting column of the OR-table was manually pushed out of his locking position. Please perform the following actions: 1. Please turn the OR Table manually to the next locking position. 2. Press OR Table UP and DOWN buttons in succession to confirm error message. 3. If the problem still persists, the patient can not be transferred to the MR-system table anymore, consider patient evacuation (see description in the user manual). 4. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Invalid status of D-switch assembly (SD1...SD4) in OR or magnet position.	1. Check the switch adjustments (SD1...SD4) and cable W702) 2. Check the swivel joint assembly incl. "MV drive"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Invalid D-Position"<<<
Error	209	PXX.mc	OR Table error: Internal position K/T error 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. If the problem still persists, please use the "emergency mode" (see description in the user manual). Please check that the table will not collide with anything. 3. If the problem still persists, the patient can not be transferred to the MR-system table anymore, consider patient evacuation (see description in the user manual). 4. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Invalid position of axes K and T during swivelling or in the magnet position (encoder, switch SK, ST). Current switch position is shown in info field of the main operating console next to the error number.	1. Adjustment zero position/transfer position, check SK, ST, cable W130, W160 2. Check position and tolerance axes K and T	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Invalid Nullposition"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	210	PXX.mc	OR Table error: SU/SW switch error 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. Release and lift OR table top, remove possible obstacles and try to re-engage the OR table top again. 3. If the problem still persists, the patient can not be transferred to the MR-system table anymore, consider patient evacuation (see description in the user manual). 4. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Invalid status of locking switches SU1..SU4, SW1..SW4 (current switch position is shown in the info field of the main operating console next to the error number)	1. Check/adjust hook, claw assembly U and W 2. Check SU1-SW4, cable W50, W51, W60, W61 3. Check tabletop and "carrier frame C"	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Invalid U/W-Position"<<<
Error	211	PXX.mc	OR Table error: Invalid encoder signals 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. If the problem still persists, please use the "emergency mode". Please check that the table will not collide with anything. 3. If the problem still persists, the patient can not be transferred to the MR-system table anymore, consider patient evacuation (see description in the user manual). 4. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Anti-valence of encoder signals H, K, or T is not guaranteed. Axis is shown in the info field of the main operating console next to the error number.	1. Check/replace encoder H, K, or T 2. Check/replace cable W90, W120, W150 3. Check/replace control board B2100	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Encoder Antivalence Err."<<<
Error	212	PXX.mc	OR Table error: Battery malfunction 1. Press OR Table UP and DOWN buttons in succession to confirm error message. Note: In case of mains power failure, emergency backup power for the OR Table can not be guaranteed! 2. Please call Siemens Service.	Switch-off electronic, accu (self switch-off) does not function	Replace control board B2100	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Accu Switchoff Error"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	213	PXX.mc	OR Table error: Internal voltage supply malfunction 1 1. Shut down MR system and wait 1 minute before restart 2. If the problem still persists, consider patient evacuation (see description in the user manual). 3. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Switch-off intermediate circuit (internal voltage supply) UVW does not function.	Replace control board B2100	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switchoff-Error UVW"<<<
Error	214	PXX.mc	OR Table error: Internal voltage supply malfunction 2 1. Shut down MR system and wait 1 minute before restart. 2. If the problem still persists, consider patient evacuation (see description in the user manual). 3. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Switch-on intermediate circuit (internal voltage supply) UVW does not function. Possible subsequent error of other causes	1. Test safety lock (switches SK, ST, SD1..4, stop circuit) 2. Replace control board B2100	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switchon-Error UVW"<<<
Error	215	PXX.mc	OR Table error: Internal voltage supply malfunction 3 1. Shut down MR system and wait 1 minute before restart. 2. If the problem still persists, consider patient evacuation (see description in the user manual). 3. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Switch-off intermediate circuit (internal voltage supply) HKT does not function. Possible subsequent error of other causes.	Replace control board B2100	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switchoff-Error HKT"<<
Error	216	PXX.mc	OR Table error: Internal voltage supply malfunction 4 1. Shut down MR system and wait 1 minute before restart. 2. If the problem still persists, consider patient evacuation (see description in the user manual). 3. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Switch-on intermediate circuit (internal voltage supply) HKT does not function. Possible subsequent error of other causes.	1. Test safety lock (switches SD1..4, stop circuit) 2. Replace control board B2100	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switchon-Error HKT"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	217	PXX.mc	OR Table error: Auxiliary drive timeout error 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. If the problem still persists, consider patient evacuation (see description in the user manual). 3. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Mechanical end-stop auxiliary drive UVW was not detected within the expected time frame.	1. Test/adjustment auxiliary drive group UVW and its end stop. 2. Replace control board B2100.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"No Limit UVW"<<<
Error	218	PXX.mc	OR Table error: Overtemperature electronic 1. Please wait approximately 10 minutes for cooling down. 2. Press OR Table UP and DOWN buttons in succession to confirm error message. 3. If the problem still persists, consider patient evacuation (see description in the user manual). 4. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Overtemperature of control board B2100, temperature in the info display is shown in °C	1. Ask the customer how frequently the ORT has been moved and investigate if one of the motor drives were overloaded 2. Check temperature sensors 3. Replace control board B2100	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"MB Overtemperature"<<<
Error	219	PXX.mc	OR Table error: SU/SW switch error 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. Release and lift OR table top, remove possible obstacles and try to engage the OR table top again. 3. If the problem still persists, consider patient evacuation (see description in the user manual). 4. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Switch errors SU1, SU2, SW1, SW2 after locking/unlocking the hooks. Switch states are not as expected.	1. Check locking mechanism for obstacles. 2. Check / adjust hook, claw assembly U and W. 3. Check SU1-SW4, cable W50, W51, W60, W61. 4. Check auxiliary drives U, W.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Unexpected Movement HKT"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	220	PXX.mc	OR Table error: Unexpected OR Table movement 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. If the problem still persists, consider patient evacuation (see description in the user manual). 3. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Main axes HKT move without motorized drive. Can be a consequential error of emergency stop during movement.	1. Check main drives, their mountings and bearings. 2. Check encoders, their mountings and bearings.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Error U/W"<<<
Error	221	PXX.mc	OR Table error: OR Table position out of range 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. If the problem still persists, please use the "emergency mode" (see description in the user manual). Please check that the table will not collide with anything. 3. If the problem still persists, consider patient evacuation (see description in the user manual). 4. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Main axes (HKT) out of their permitted travel range.	1. Check main drives, their mounting and bearings. 2. Check encoders, their mounting, bearing and adjustment values.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Out of Range HKT-Pos."<<<
Error	222	PXX.mc	OR Table error: Watchdog fault 1. Reset OR Table by pressing OR Table Stop button for approx. 10 seconds. 2. If the problem still persists, consider patient evacuation (see description in the user manual). 3. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Error in control board B2100.	Replace control board B2100.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Watchdog Fault"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	223	PXX.mc	OR Table error: Motor drive power failure 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. Shut down MR system and wait 2 minutes before restart. 3. If the problem still persists, consider patient evacuation (see description in the user manual). 4. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Drive cannot acquire sufficient power.	1. Check drives H, K, T and their cables W80, W110, W140. 2. Check that the lift drives uniformly UP and DOWN with a load of approx. 150 kg. If not, replace control board B2100 3. Check self-locking of drives H, K, T.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Motor Open HKT"<<<
Error	224	PXX.mc	OR Table error: Short circuit in motor drive detected 1. Check for collisions or obstacles. 2. Press OR Table UP and DOWN buttons in succession to confirm error message. 3. Shut down MR system and wait 2 minutes before restart. 4. If the problem still persists, consider patient evacuation (see description in the user manual). 5. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Short-circuit detection of main drives, H, K or T.	Check drives H, K, T and their cables W80, W110, W140	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Motor Short HKT"<<<
Error	225	PXX.mc	OR Table error: Control board hardware error 1. Reset OR Table by pressing OR Table Stop button for approx. 10 seconds. 2. Shut down MR system and wait 2 minutes before restart. 3. If the problem still persists, consider patient evacuation (see description in the user manual). 4. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Internal error of control for main drives, H, K, or T.	Replace control board B2100.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Powerstage Fault"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	226	PXX.mc	OR Table error: Switch error of main axes 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. If the problem still persists, please use the "emergency mode" (see description in the user manual). Please check that the table will not collide with anything. 3. If the problem still persists, consider patient evacuation (see description in the user manual). 4. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Plausibility monitoring of switches SH, SK, ST for encoders H, K, T.	1. Adjustment/mounting of switches SH, SK, ST and learning of their positions. 2. Check the batches of drive 1. 3. Check the toothed belt tension of K.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Switch Error H/K/T"<<<
Error	227	PXX.mc	OR Table error: Detection of movement direction failed 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. If the problem still persists, consider patient evacuation (see description in the user manual). 3. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Detection of direction failed when moving the axes by the main drives.	1. Check polarity of electrical connections motor H, K, T (W80, W110, W140). 2. Replace control board B2100.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Speed Error H/K/T"<<<
Error	228	PXX.mc	OR Table error: UW Switches Motion Error 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. Press table down button to release locking mechanism, check for obstacles and try to engage the OR-table top again. 3. If the problem still persists, consider patient evacuation (see description in the user manual). 4. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Invalid status of switch SU/SW during vertical drive.	1. Check switch unit SU and SW, if necessary adjust 2. Check position of OR Table in OR-table carrier, correct if necessary 3. Check locking hooks and the mechanic	In case of PTAB stand alone mode the following text will appear at the patient table display: >>"OR-Table-Error 228"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	229	PXX.mc	OR Table error: Table Clamp Malfunction 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. Press table down button to release locking mechanism, check for obstacles and try to engage the OR-table top again. 3. If the problem still persists, consider patient evacuation (see description in the user manual). 4. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Invalid status of the locking switch S46....S49 of the OR-patient board.	Check switches S46...S49 and the actuating levers, adjust if necessary Check APC contact	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"OR-Table-Error 229"<<<
Warning	241	PXX.mc	OR Table warning: STOP button pressed by user 1. Release OR Table stop buttons. 2. Press OR Table UP and DOWN buttons in succession to confirm error message. 3. If the problem still persists, consider patient evacuation (see description in the user manual). 4. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	The operator has pressed one of the STOP-buttons of the system. If this has not happened, a circuit break in one of STOP-(opening contact) circuits is present.	1. Check the two OR Table stop buttons and their connections W170, W190. 2. Check PTAB stop buttons and their connections.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"STOP-But.pressed by User"<<<
Warning	242	PXX.mc	OR Table warning: OR Table stopped by Patient Table STOP button 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. Check if Patient Table is working properly. 3. If the problem still persists, consider patient evacuation (see description in the user manual). 4. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Table stop triggered through PTAB.	1. Reset table stop of PTAB. 2. Check PTAB stop buttons and their connections.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"STOP activated by auton."<<<
Warning	243	PXX.mc	OR Table warning: Table stop by OR Table 1. Release stop buttons of OR Table. 2. Press OR Table UP and DOWN buttons in succession to confirm error message. 3. If the problem still persists, please call Siemens Service.	Table stop by ORT (e.g. follow-up error).	1. Reset table stop ORT. 2. Check the two OR Table stop buttons and their connections W170, W190. 3. Check PTAB stop buttons and their connections.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"STOP activated by PTAB "<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Warning	244	PXX.mc	OR Table warning: Mains power failure 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. Continue in battery mode (emergency backup power). Be aware of limitations in the battery mode: a. The table top can not be transferred to the MR system. b. Once the mains power is available, the normal operation is returned. 3. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Detection of power failure (e.g. activate emergency shut-down for room)	1. Reset table stop and continue in accu mode. 2. Restore power supply. 3. If PTAB is supplied with line voltage, check cable W10, fuses F13 (ORT), F2 (PTAB). 4. Replace control board B2100.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Mains Power Fail"<<<
Warning	245	PXX.mc	OR-Table warning: Fixed Panel Fault 1. Please use the cable remote control Warning: The cable remote control does not work in the battery mode (emergency backup power). 2. Press OR Table UP and DOWN buttons in succession to confirm error message. 3. If the problem still persists, consider patient evacuation (see description in the user manual). 4. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Keyboard error on main operating console.	1. Check cable W33. 2. Replace operating console. 3. Replace control board B2100.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Fixed Panel Fault"<<<
Warning	246	PXX.mc	OR Table warning: Remote control failed 1. Please remove the remote control and use the main control of OR Table instead. 2. Press OR Table UP and DOWN buttons in succession to confirm error message. 3. If the problem still persists, consider patient evacuation (see description in the user manual). 4. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Keyboard error, remote control.	1. Check cables W30, W31, W35. 2. Replace remote control. 3. Replace control board B2100.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Mobile Panel Fault"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	247	PXX.mc	OR-Table error: PTAB has no reference 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. If the problem still persists, please use the "emergency mode" (see description in the user manual). Please check that the table will not collide with anything. 3. If the problem still persists, consider patient evacuation (see description in the user manual). 4. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	No reference position detected by PTAB.	1. Perform reference movement with PTAB. 2. Move PTAB table top into the magnet until OR Table transfer position.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"PTAB has no reference"<<<
Error	248	PXX.mc	OR Table error: Patient Table not ready 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. If the problem still persists, please use the "emergency mode" (see description in the user manual). Please check that the table will not collide with anything. 3. If the problem still persists, consider patient evacuation (see description in the user manual). 4. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	PTAB not ready for horizontal movement.	Check all prerequisites for horizontal movement.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"PTAB isn't ready to move"<<<
Error	249	PXX.mc	OR Table error: Patient Table not ready 1. Press OR Table UP and DOWN buttons in succession to confirm error message. 2. If the problem still persists, please use the "emergency mode" (see description in the user manual). Please check that the table will not collide with anything. 3. If the problem still persists, consider patient evacuation (see description in the user manual). 4. Please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	PTAB loadware is not running.	1. Reboot Scanner. 2. Create Savelog and inform HSC for further investigation.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"Accumulator missing"<<<

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Warning	250	PXX.mc	OR Table warning: Battery fault 1. Press OR Table UP and DOWN buttons in succession to confirm error message. Note: Battery charge can not be guaranteed 2. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Accu (for emergency operation) is not available.	1. Check/replace accu. 2. Check accu connection (e.g. polarity, fuse F, check and replace if necessary). 3. If fuse F1 is repeatedly destroyed, replace control board B2100.	In case of PTAB stand alone mode the following text will appear at the patient table display: >>>"PTAB is not running"<<<
Warning	2000	PXX.mc	Patient table warning: There is no front display attached. Please call Siemens Service.			None.
Error	2001	PXX.mc	Patient table error: Supplied functions not installed. Please call Siemens Service.			None.
Error	2002	PXX.mc	Patient table error: Power driver not installed. Please call Siemens Service.			None.
Error	2003	PXX.mc	Patient table error: Bus controller not installed. Please call Siemens Service.			None.
Error	2004	PXX.mc	Patient table error: Voltage error (14V DC PAN R). Please call Siemens Service.			None.
Error	2005	PXX.mc	Patient table error: Voltage error (14V DC PAN F). Please call Siemens Service.			None.
Error	2006	PXX.mc	Patient table error: Voltage error (14V DC SWT). Please call Siemens Service.			None.
Informational	2007	PXX.mc	Patient table information: Voltage error (OUT_DC SWITCHES). Please call Siemens Service.		Check trolley connection.	None.
Error	2008	PXX.mc	Patient table error: Voltage error (24V DC SWT). Please call Siemens Service.	check the position of the switch S1.1 (RTT, NO RTT) on the PWRD A4110, the cable W4141 and the SWT A4115.		None.
Error	2009	PXX.mc	Patient table error: Voltage error (24V DC ORT). Please call Siemens Service.			None.
Error	2010	PXX.mc	Patient table error: Voltage error (24V DC MDSD). Please call Siemens Service.			None.
Warning	2011	PXX.mc	Patient table warning: There is no front panel on right side attached. Please call Siemens Service.			None.
Warning	2012	PXX.mc	Patient table warning: There is no front panel on left side attached. Please call Siemens Service.			None.

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	2013	PXX.mc	Patient table error: The bore temperature sensor is not installed. Please call Siemens Service.	None.	Ref. general message text.	None.
Warning	2014	PXX.mc	Patient table warning: There is no rear panel attached. Please call Siemens Service.			None.
Warning	2015	PXX.mc	Patient table warning: There is no front control device on right side attached. Please call Siemens Service.			None.
Warning	2016	PXX.mc	Patient table warning: There is no front control device on left side attached. Please call Siemens Service.			None.
Warning	2017	PXX.mc	Patient table warning: There is no rear control device attached. Please call Siemens Service.			None.
Warning	2018	PXX.mc	Patient table warning: There is no rear display attached. Please call Siemens Service.			None.
Informational	2020	PXX.mc	Patient Table Information: Table Stop occurred			
Informational	2021	PXX.mc	Patient Table Information: Table Stop released			
Error	3000	PXX.mc	Patient table error: While scanning, the Patient table left the "sleep" mode. This could cause noise in the images. Retry. If the problem still persists, please call Siemens Service	Patient table CAN woke up while scanning.	* Start Test Tools and check the status of all CAN units. * Evaluate previous error messages. * Start sequence again. * If not successful, reboot Scanner to solve a possible software problem. * If not successful, reset hardware by switching off/on Scanner.	None.
Error	3001	PXX.mc	Patient table error: No communication with unit Patient table. Retry, reboot Scanner, switch off/on Scanner. If the problem still persists, please call Siemens Service	Communication with Patient table not possible. Patient table communication timeout occurred.	* Start Test Tools and run CAN tests. * Evaluate previous error messages. * Reboot Scanner to solve a possible software hang-up. * If not successful, reset hardware by switching off/on Scanner.	None.
Error	3002	PXX.mc	Patient table error: Patient table Device fault. Retry, reboot Scanner, switch off/on Scanner. If the problem still persists, please call Siemens Service	Patient table device fault	* Start Test Tools and run CAN tests. * Evaluate previous error messages. * Reboot Scanner to solve a possible software hang-up. * If not successful, reset hardware by switching off/on Scanner.	None.
Warning	3003	PXX.mc	Patient table warning: Patient table not operational. Retry, reboot Scanner, switch off/on Scanner. If the problem still persists, please call Siemens Service	Patient table not operational.	* Start Test Tools and run CAN tests. * Evaluate previous error messages. * Reboot Scanner to solve a possible software hang-up. * If not successful, reset hardware by switching off/on Scanner.	None.
Error	3004	PXX.mc	Patient table error: The patient table was reset unexpectedly. Retry, reboot Scanner, switch off/on Scanner. If the problem still persists, please call Siemens Service	Logon message received due to a reset of patient table CAN.	* Start Test Tools and check the status of all CAN units. * Search previous error messages, particularly for temperature and power faults. * Reboot Scanner and start sequence again.	None.

(Message clas	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Warning	4000	PXX.mc	"Move during Scan Device" warning: While scanning, the MDSD left the "sleep" mode. This could cause noise in the images. Retry. If the problem still persists, please call Siemens Service	MDSD CAN woke up while scanning.	* Start Test Tools and check the status of all CAN units. * Evaluate previous error messages. * Start sequence again. * If not successful, reboot Scanner to solve a possible software problem. * If not successful, reset hardware by switching off/on Scanner.	None.
Error	4001	PXX.mc	"Move during Scan Device" error: No communication with unit MDSD. Retry, reboot Scanner, switch off/on Scanner. If the problem still persists, please call Siemens Service	Communication with MDSD not possible. MDSD communication timeout occurred.	* Start Test Tools and run CAN tests. * Evaluate previous error messages. * Reboot Scanner to solve a possible software hang-up. * If not successful, reset hardware by switching off/on Scanner.	None.
Error	4002	PXX.mc	"Move during Scan Device" error: MDSD Device fault. Retry, reboot Scanner, switch off/on Scanner. If the problem still persists, please call Siemens Service	MDSD device fault	* Start Test Tools and run CAN tests. * Evaluate previous error messages. * Reboot Scanner to solve a possible software hang-up. * If not successful, reset hardware by switching off/on Scanner.	None.
Warning	4003	PXX.mc	"Move during Scan Device" warning: MDSD not operational. Retry, reboot Scanner, switch off/on Scanner. If the problem still persists, please call Siemens Service	MDSD not operational.	* Start Test Tools and run CAN tests. * Evaluate previous error messages. * Reboot Scanner to solve a possible software hang-up. * If not successful, reset hardware by switching off/on Scanner.	None.
Error	4004	PXX.mc	"Move during Scan Device" error: The MDSD was reset unexpectedly. Retry, reboot Scanner, switch off/on Scanner.	Logon message received due to a reset of MDSD CAN.	* Start Test Tools and check the status of all CAN units. * Search previous error messages, particularly for temperature and power faults. * Reboot Scanner and start sequence again.	None.
Warning	5000	PXX.mc	OR Table warning: Bus communication error Reboot Scanner. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Not used until ORT@CAN bus!!! ORT woke up CAN bus via PTAB while scanning.	1. Start TestTools and check the status of all CAN units. 2. Evaluate previous error messages. 3. Start sequence again. 4. If not successful, reboot Scanner to activate a possible software. 5. If not successful, reset HW by switching OFF/ON Scanner. 6. Start TestTools and check the status of all CAN units. 7. Evaluate previous error messages. 8. Start sequence again. 9. If not successful, create Savelog and inform HSC for further investigation.	None.
Error	5001	PXX.mc	OR Table error: Communication error Reset OR Table by pressing OR Table Stop button for approx. 10 seconds. Reboot Scanner. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Not used until ORT@CAN bus!!! Communication with ORT not possible. ORT communication timeout occurred.	1. Start TestTools and run CAN tests. 2. Evaluate previous error messages. 3. Reboot Scanner to solve a possible software hang-up 4. If not successful, reset hardware by switching OFF/ON Scanner. 5. Create Savelog and inform HSC for further investigation	None.

(Message clas	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	5002	PXX.mc	OR Table error: OR Table device fault Reset OR Table by pressing OR Table Stop button for approx. 10 seconds. Reboot Scanner. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Not used until ORT@CAN bus!!! ORT device fault.	1. Start Test Tools and run CAN tests. 2. Evaluate previous error messages 3. Reboot Scanner to solve a possible software hang-up. 4. If not successful, reset hardware by switching OFF/ON Scanner. 5. Create Savelog and inform HSC for further investigation.	None.
Informational	5003	PXX.mc	OR Table information: This message is not used until ORT@CAN bus!!! OR Table device fault Reset OR Table by pressing OR Table Stop button for approx. 10 seconds. Reboot Scanner. If the problem still persists, please save System Log-Files by selecting "System/Control/Tools" and "Save System Log-Files". Please call Siemens Service.	Not used until ORT@CAN bus!!! ORT not operational.	1. Start Test Tools and run CAN tests. 2. Evaluate previous error messages. 3. Reboot Scanner to solve a possible software hang-up. 4. If not successful, reset hardware by switching OFF/ON Scanner.	None.
Error	5004	PXX.mc	OR Table error: This message is not used until ORT@CAN bus!!!	Not used until ORT@CAN bus!!! Logon message received due to a reset of ORT CAN.	1. Start Test Tools and check the status of all CAN units. 2. Search previous error messages, particularly for temperature and power faults. 3. Reboot Scanner and start sequence again.	None.
Warning	5100	PXX.mc	Patient Table warning: Gradient coil temperature sensor "PRI_X" is hot. Reboot Scanner! If the problem still persists, please call Siemens Service			
Error	5101	PXX.mc	Patient Table error: Gradient coil temperature sensor "PRI_X" is too hot. Wait for cooling down for approximately 10 minutes! Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "PRI_X" is over the error threshold during scanner boot.	1.a) Wait (approx. 10 minutes) for cooling down of the gradient coil and try again. Reduce the thermal load of the coil by appropriate measurement conditions. 1.b) Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts" 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	
Warning	5102	PXX.mc	Patient Table warning: Gradient coil temperature sensor "PRI_Y" is hot. Reboot Scanner! If the problem still persists, please call Siemens Service			

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	5103	PXX.mc	Patient Table error: Gradient coil temperature sensor "PRI_Y" is too hot. Wait for cooling down for approximately 10 minutes! Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "PRI_Y" is over the error threshold during scanner boot.	1.a) Wait (approx. 10 minutes) for cooling down of the gradient coil and try again. Reduce the thermal load of the coil by appropriate measurement conditions. 1.b) Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts" 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	
Warning	5104	PXX.mc	Patient Table warning: Gradient coil temperature sensor "PRI_Z" is hot. Reboot Scanner! If the problem still persists, please call Siemens Service			
Error	5105	PXX.mc	Patient Table error: Gradient coil temperature sensor "PRI_Z" is too hot. Wait for cooling down for approximately 10 minutes! Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "PRI_Z" is over the error threshold during scanner boot.	1.a) Wait (approx. 10 minutes) for cooling down of the gradient coil and try again. Reduce the thermal load of the coil by appropriate measurement conditions. 1.b) Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts" 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	
Warning	5106	PXX.mc	Patient Table warning: Gradient coil temperature sensor "SEC_X" is hot. Reboot Scanner! If the problem still persists, please call Siemens Service			
Error	5107	PXX.mc	Patient Table error: Gradient coil temperature sensor "SEC_X" is too hot. Wait for cooling down for approximately 10 minutes! Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "SEC_X" is over the error threshold during scanner boot.	1.a) Wait (approx. 10 minutes) for cooling down of the gradient coil and try again. Reduce the thermal load of the coil by appropriate measurement conditions. 1.b) Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts" 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	
Warning	5108	PXX.mc	Patient Table warning: Gradient coil temperature sensor "SEC_Y" is hot. Reboot Scanner! If the problem still persists, please call Siemens Service			

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	5109	PXX.mc	Patient Table error: Gradient coil temperature sensor "SEC_Y" is too hot. Wait for cooling down for approximately 10 minutes! Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "SEC_Y" is over the error threshold during scanner boot.	1.a) Wait (approx. 10 minutes) for cooling down of the gradient coil and try again. Reduce the thermal load of the coil by appropriate measurement conditions. 1.b) Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts" 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	
Warning	5110	PXX.mc	Patient Table warning: Gradient coil temperature sensor "SEC_Z" is hot. Reboot Scanner! If the problem still persists, please call Siemens Service			
Error	5111	PXX.mc	Patient Table error: Gradient coil temperature sensor "SEC_Z" is too hot. Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "SEC_Z" is over the error threshold during scanner boot.	1.a) Wait (approx. 10 minutes) for cooling down of the gradient coil and try again. Reduce the thermal load of the coil by appropriate measurement conditions. 1.b) Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts" 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	
Warning	5112	PXX.mc	Patient Table warning: Gradient coil temperature sensor "SOCKET_1" is hot. Reboot Scanner! If the problem still persists, please call Siemens Service			
Error	5113	PXX.mc	Patient Table error: Gradient coil temperature sensor "SOCKET_1" is too hot. Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "SOCKET_1" is over the error threshold during scanner boot.	1.a) Wait (approx. 10 minutes) for cooling down of the gradient coil and try again. Reduce the thermal load of the coil by appropriate measurement conditions. 1.b) Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts" 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	
Warning	5114	PXX.mc	Patient Table warning: Gradient coil temperature sensor "SOCKET_2" is hot. Reboot Scanner! If the problem still persists, please call Siemens Service			

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	5115	PXX.mc	Patient Table error: Gradient coil temperature sensor "SOCKET_2" is too hot. Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "SOCKET_2" is over the error threshold during scanner boot.	1.a) Wait (approx. 10 minutes) for cooling down of the gradient coil and try again. Reduce the thermal load of the coil by appropriate measurement conditions. 1.b) Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts" 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	
Warning	5200	PXX.mc	Patient Table warning: Gradient coil temperature sensor "PRI_X" is cold. Reboot Scanner! If the problem still persists, please call Siemens Service			
Error	5201	PXX.mc	Patient Table error: Gradient coil temperature sensor "PRI_X" is too cold. Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "PRI_X" is under the error threshold during scanner boot.	1. Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts". 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	
Warning	5202	PXX.mc	Patient Table warning: Gradient coil temperature sensor "PRI_Y" is cold. Reboot Scanner! If the problem still persists, please call Siemens Service			
Error	5203	PXX.mc	Patient Table error: Gradient coil temperature sensor "PRI_Y" is too cold. Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "PRI_Y" is under the error threshold during scanner boot.	1. Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts". 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	
Warning	5204	PXX.mc	Patient Table warning: Gradient coil temperature sensor "PRI_Z" is cold. Reboot Scanner! If the problem still persists, please call Siemens Service			

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	5205	PXX.mc	Patient Table error: Gradient coil temperature sensor "PRI_Z" is too cold. Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "PRI_Z" is under the error threshold during scanner boot.	1. Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts". 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	
Warning	5206	PXX.mc	Patient Table warning: Gradient coil temperature sensor "SEC_X" is cold. Reboot Scanner! If the problem still persists, please call Siemens Service			
Error	5207	PXX.mc	Patient Table error: Gradient coil temperature sensor "SEC_X" is too cold. Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "SEC_X" is under the error threshold during scanner boot.	1. Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts". 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	
Warning	5208	PXX.mc	Patient Table warning: Gradient coil temperature sensor "SEC_Y" is cold. Reboot Scanner! If the problem still persists, please call Siemens Service			
Error	5209	PXX.mc	Patient Table error: Gradient coil temperature sensor "SEC_Y" is too cold. Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "SEC_Y" is under the error threshold during scanner boot.	1. Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts". 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	
Warning	5210	PXX.mc	Patient Table warning: Gradient coil temperature sensor "SEC_Z" is cold. Reboot Scanner! If the problem still persists, please call Siemens Service			

(Message class)	Id	Source	Message Text	Explanation	Action	Information
S/I/W/E						
Error	5211	PXX.mc	Patient Table error: Gradient coil temperature sensor "SEC_Z" is too cold. Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "SEC_Z" is under the error threshold during scanner boot.	1. Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts". 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	
Warning	5212	PXX.mc	Patient Table warning: Gradient coil temperature sensor "SOCKET_1" is cold. Reboot Scanner! If the problem still persists, please call Siemens Service			
Error	5213	PXX.mc	Patient Table error: Gradient coil temperature sensor "SOCKET_1" is too cold. Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "SOCKET_1" is under the error threshold during scanner boot.	1. Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts". 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	
Warning	5214	PXX.mc	Patient Table warning: Gradient coil temperature sensor "SOCKET_2" is cold. Reboot Scanner! If the problem still persists, please call Siemens Service			
Error	5215	PXX.mc	Patient Table error: Gradient coil temperature sensor "SOCKET_2" is too cold. Reboot Scanner! If the problem still persists, please call Siemens Service	Temperature value from gradient coil sensor "SOCKET_2" is under the error threshold during scanner boot.	1. Check the cooling system. 2. Reboot Scanner, if the problem still persists continue. 3. Check the signal lines between the gradient coil NTC's, X2 and the patient table electronics, A4120 SUPF/X50 (refer to wiring diagrams). Measure and compare the NTC resistances. 4. Install the spare sensors according to document "Replacement of Parts". 5. Reboot Scanner 6. If error occurs again, create Savelog and inform HSC for further investigation	