

TABLE OF CONTENTS

TABLE OF CONTENTS 1

1- BORE VENTING 2

2- BORE LIGHTING 2

REVISION HISTORY 3

Description - Checks for bore venting and lighting.

1- BORE VENTING

Bore ventilation air is available for the comfort of the patient in the scanner. The body coil has a blower box, located in the RF screen room, that blows air into the magnet bore. Verify proper functioning of bore ventilation using the following procedure:

1. With bore venting on (FAN ON magnet enclosure button), check that air is freely flowing through all nozzles of the Patient Comfort Module. If not, check for obstructions and remove them with proper nonferrous tools.
2. Check for leaks by placing tape over nozzles. Feel for air leaks at the end of the bore. Repair air hoses as needed.

2- BORE LIGHTING

Bore lighting is provided by fiber optics that channel light from a single light source at the rear end ring of the magnet enclosure to the inside of frosted light tubes that are built into the bore enclosure. The following procedure verifies proper operation of the bore light system:

1. With power applied to bore light power supply, check that the bore light system is on. If not, check the power supply or the bulb.

For power supply information, refer to the procedure for Magnet Enclosure P.S. Module Calibration (+12 Vdc Bore Light Power Supply).

For bore light replacement, refer to the procedure for Bore; Halogen Light Replacement, Section 1- Bore Light.

REVISION HISTORY

REV	DATE	AUTHOR	PRIMARY REASONS FOR CHANGE
0	June 15, 1998	J. Saperstein	Initial conversion from Toolbook to Word formats.