GE Medical Systems

Technical Publications

**Direction 2209282
Revision 4**

Signa® Release 8.x Software Options Installation By MOD

Copyright© 1999, 2000 by General Electric Co.

Operating Documentation

WARNING

- THIS SERVICE MANUAL IS AVAILABLE IN ENGLISH ONLY.
- IF A CUSTOMER'S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER'S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR OR PATIENT FROM ELECTRIC SHOCK, MECHANICAL OR OTHER HAZARDS.

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1 OVERVIEW

This procedure describes how to install software options on the Signa 8.x system using the MOD drive. The Software Option package contains the appropriate software option MOD.

2 OPTION SOFTWARE INSTALLATION

1. Insert the MOD disk in the SGI computer MOD drive.
2. Select the **[Service Tools Desktop]**.
3. Select **[Install...]** from the Service Desktop Manager.
4. After Type Root Password is displayed, type **operator** as the password.

2-1 For Release 8.3

1. After the Guided Install window opens (30 sec. for Octane; 60 sec. for Indigo), select **Option Key tab**.
2. Select **MOD** box at Install Option Key portion of window, then select **Install**.
3. When installation process is complete, either load the next MOD, or when all option keys have been loaded, select **File ->Quit** to exit the Guided Install tool.

2-2 For Release 8.2.5

1. From the "Install" window, select **[Option Install]**
2. From the "Option Install" window, select **[Change Status]** next to the appropriate Option Name to be installed.
3. At the "Select Device" window, select the appropriate drive device, i.e., **[MOD]**.
4. When installation process is complete, either load the next MOD option disk and repeat steps 2 & 3, or in the "Install Keys" window, select **[Quit]**.
5. Wait about 20 seconds for various windows to clear, then from the "Option Install" window, select **[Quit]**.
6. At the "Status Changed" window, select **[OK]**.
7. At the "Install" window, select **[Quit Install]**.

3 REBOOT COMPUTER

1. Reboot the computer.
2. The installation is now complete.

4 SIGNA 8.3 OPTION CATALOG MOD CONTENTS

What follows is a list of available catalog option MOD(s) for the Signa 8.3 and later systems and a list of what software options are activated with each MOD.

EchoSpeed Option Key MOD 2209991

- Turns on the EchoSpeed Option.

M1090LD- Proactive Service

- Proactive Service - Activates Proactive Diagnostics and Dialout.

M1090LF - Phase Contrast Vascular

- Phase Contrast Vascular Imaging - CVI

M1090LG - Cine

- Cine - Images usually of moving anatomy scanned and arrange for playback rapidly in movie format.

M1090LJ - IVI

- Interactive Vascular Imaging - IVI

M1090LL - 3D Surface

- 3D Surface/Voxtool

M1090LM - Research

- Research - Allows for access to low level control of scan parameters.

M1090LN - Research PSD

- Research PSD - Advanced PSD development.

M1090LP - Echo Plus

- DW EPI - Diffusion Weighted Echo Planer Imaging

M1090LT - SmartPrep

- SmartPrep - Automated triggering for contrast enhanced MRA

M1090MW - MS/MA

- Multi-Slice/Multi-Angle Oblique

M1090NN - FuncTool

- Func Tool - Raw Data Analyses tool.

M1090PL - ConnectPro

- HIS RIS - Provides DICOM modality work list connectivity to hospital or radiology information system.

M1090KS - Lx UG Options Europe/Asia

- Spectroscopy/PROBE - Single voxel proton spectroscopy with PROBE-S and PROBE-P
- DW EPI - Diffusion Weighted Echo Planer Imaging
- Flair EPI
- SmartPrep - Automated triggering for contrast enhanced MRA
- SSFSE - Single Shot Fast Spin Echo

M1090CR - Cardiac All

- Bloodsupp - Blood Suppression
- Fastcine - Fast gated multi-phase images of the heart that show the heart as it contracts and expands during cardiac cycle in a movie format.
- Tagging - The placement of "tags" in an MR image is used to assess myocardial wall functions visually.
- Sgdperf - SGD enhancements which increases slew rates and duty cycles for higher resolution thin slice imaging.

M1033JA - iDRIVE PRO

- iDRIVE Pro Interactive imaging uses a fast GRE pule sequence and a new user interface to acquire continuous 2D images at a high frame rate and provide on-image navigation. Not available on systems with a single i860 AP.

M1033JB - SmartPrep2000

- Smart Prep - Automated triggering for contrast enhanced MRA
- SmartPrep2000 Upgrade - Provides SmartPrep automated table motion.

M1033JH - SmartPrep2000 Upgrade

- SmartPrep2000 Upgrade - Provides SmartPrep automated table motion.

M1033JD - Probe2000

- Spectroscopy/PROBE - Single voxel proton spectroscopy with PROBE-S and PROBE-P
- Probe2000 Upgrade - Multi voxel proton spectroscopy (PROBE - SI)
- Multi-Nuclear Spectroscopy - Provides access to additional spectroscopy PSD which may be used in proton mode.

M1033JJ - Probe2000 Upgrade

- Probe2000 Upgrade - Multi voxel proton spectroscopy (PROBE - SI)
- Multi-Nuclear Spectroscopy - Provides access to additional spectroscopy PSD which may be used in proton mode.

M1033JF - ScanTools2000

- Cine - Images usually of moving anatomy scanned and arrange for playback rapidly in movie format.
- Phase Contrast Vascular Imaging
- Echo Planar Imaging (EPI)
- 3D Surface/Voxtool
- Flair EPI
- iDRIVE - Basic Real Time Interactive Imaging package.

M1033JL - LX Tools (former MiniTools)

- Fast Spin Echo & FLAIR
- Fast Gradient Echo
- Time of Flight
- Interactive Vascular Imaging
- SPECIAL - Spectral / Spatial Fat Suppression. Chemical saturation using a inversion pulse (180° pulse). This key will always suppress just fat and can not be changed to suppress water.
- SSFSE - Single Shot Fast Spin Echo
- Three Plane Localizer
- E3DTOF
- FSE-XL
- Clariview - Image Filtering Package
- Sgdperf - SGD (PSD) enhancements which increases slew rates and duty cycles for higher resolution thin slice imaging.

M1033J - ScanTools2000 Upgrade

- Cine
- Phase Contrast Vascular Imaging
- Echo Planar Imaging
- 3D Surface/Voxtool
- Flair EPI
- iDRIVE Pro Interactive imaging uses a fast GRE pulse sequence and a new user interface to acquire continuous 2D images at a high frame rate and provide on-image navigation.

REVISION HISTORY

REV	DATE	AUTHOR	PRIMARY REASONS FOR CHANGE
3	Oct, 1999	Piepenburg	Interleaf Version
4	Feb 21, 2000	Hawthorne	Converted to Word. Added section 3 explaining what options are turned on by what MOD.