

TABLE OF CONTENTS

table of contents.....	1
1- Scan Preparation.....	2
2- System Restoration	3
revision history.....	4

Description - This is a basic test to check the functionality of the PIN Switch Driver Board. It is designed to use the error detecting software to determine miscalibration while scanning. There is no additional hardware required to run this test.

1- SCAN PREPARATION

1. Move switch on the ISE CERD (MR2A30A13) to RF OUT DISABLE or disconnect connector J3 RF IN at RF/Pen cabinet I/F panel (MR1A7).
2. Prepare to scan using protocol #10.1 shown in the Service Protocol list under Signa Help, or for alternate proprietary procedure, use Table 1.

TABLE 1
DD/TR BOARD SCAN PROTOCOL

Note

This is the alternate proprietary procedure available for GE use, and to sites with a valid Advanced Service Package Limited License.

1. Click on **[New Pt]**, and enter
Id: **geservice**
Name: **dd & tr**
Weight (Lb): **111**
Set Patient Protocols to **Service**.

2. Depending on your Release,

8.0: Select **[other]** to go to the protocol list on the system. Then select the protocol by name: **DD/TR Board**.

8.1: In the protocol field, type **o.10.1** (o=Other, 1=series number) to load the protocol.

3. Select **[Research Modes]**, then select **[Setup Params]**. Increase TG to **200**.

4. **[Save Series]**, then **[Prepare to Scan]**.

If you need help, see the User Interface Tutorial.

3. At Scan Protocol screen, press **[Manual Prescan]**.

4. Open message display and select **[View Log]**. This displays any error messages caused by the head coil scan.

5. Record any errors. Select **[Close]**.

6. Select **[Done]**. Select **[New Series]** and **dd & tr** service protocol. At Patient Position, select Coil and Body.

7. Select **[Save Series]** and **[Manual Prescan]**.

8. Open message display and select **[View Log]**. This displays any error messages caused by the body coil scan.
9. Record any errors. Select **[Close]**.
10. Select **[Done]**. Select **[New Series]** and **dd & tr** service protocol. At Patient Position, select Coil and a Receiver Only Surface Coil.
11. Select **[Save Series]** and **[Manual Prescan]**.
12. Open message display and select **[View Log]**. This will display any error messages caused by the Surface Coil scan.
13. Record any errors. Select **[Close]**.
14. Select **[Done]**. Select **[New Series]** and **dd & tr** service protocol. At Patient Position, setup for any Spectro scan. Press **[Save Series]**, and **[Manual Prescan]**.

Note

Performing a Spectro Scan - Perform a spectroscopy scan only if the Spectro option is installed. Change the NUC CV to 1.0 (hydrogen).

15. Open message display and select **[View Log]**. This will display any error messages caused by the spectro coil scan.
16. Record any errors. Select **[Close]**.
17. Troubleshoot any errors that may have occurred in the previous scans.

2- SYSTEM RESTORATION

Move switch on ISE CERD (MR2A30A13) to RF OUT NORMAL, or reinstall RF IN connector J3 at RF/Pen cabinet I/F panel.

REVISION HISTORY

REV	DATE	AUTHOR	PRIMARY REASONS FOR CHANGE
0	July 14, 1998	J. Saperstein	Initial conversion from Toolbook to Word.