

TABLE OF CONTENTS

TABLE OF CONTENTS	1
1- MANUFACTURER	2
2- MODEL, TYPE REFERENCE.....	2
3- STANDARDS	2
4- NOTES AND REQUIREMENTS	2
5- CLASSIFICATION OF EQUIPMENT.....	2
6- ENVIRONMENTAL CONDITIONS	3
7- LINE VOLTAGE.....	3
8- OUTPUT	3
9- COOLING CONDITIONS	3
10- IMPORTANT WARNING STATEMENTS	3
11- IMPORTANT SYMBOLS.....	4
REVISION HISTORY.....	6

Description - This section provides an overview of the Agency Standards Specifications/Technical Description.

1- MANUFACTURER

Erbtec Engineering, Inc.
2760 29th St.
Boulder, Colorado 80301 USA

2- MODEL, TYPE REFERENCE

RF Penetration II Cabinet:

- 2160200 1.5T Configuration
- 2160200-3 1.0T Configuration
- 2160200-2 1.5T Upgrade Configuration
- 2160200-4 1.0T Upgrade Configuration

System Support Module: 2160200-43

3- STANDARDS

The RF/PEN II cabinet has been tested for compliance with:

- US Standard UL-2601 covering Medical and Dental Equipment
- Canadian Standard C22.2 No. 601.1 covering Medical Electrical Equipment
- International Standard IEC 601-1 covering Medical Electrical Equipment

In addition to the above standards, the system support module has been tested for compliance with EMC Standard IEC 601-1-2 covering Electromagnetic Compatibility.

4- NOTES AND REQUIREMENTS

Compliance to EMC requirement IEC 601-1-2 is to be proven in end use application.

This product is only to be operated on the secondary of an IEC 601-1 compliant isolation transformer.

Electrical supply connection shall not differ from external ratings label.

Mains socket-outlet must properly mate with mains plug.

For availability of circuit diagrams, descriptions, and field replaceable parts, contact General Electric Medical Systems.

5- CLASSIFICATION OF EQUIPMENT

- Protection Class 1
- Stationary Equipment
- Not AP/APG Proof
- Continuous Operation
- To be operated from IEC601-1 approved supply
- Type B Equipment
- Ordinary Equipment

6- ENVIRONMENTAL CONDITIONS

See Table 1 - Operating Environment

7- LINE VOLTAGE

208Vac 3N PE $\pm$ TBD, 50/60Hz, 10.5kVa

8- OUTPUT

- 50 ohm impedance
- Head Mode: 2kW/800W pulsed (Max Linear)
- Body Mode: 20kW/8kW pulsed (Max Linear)

9- COOLING CONDITIONS

For continued safe operation, position cabinet to allow maximum air flow in and around it. Do not block air intake or exhaust ports.

10- IMPORTANT WARNING STATEMENTS

For safe operation of the equipment, and the safety of the service engineer, follow all warnings.



Identifies a condition or action for which a specific hazard is known to exist which *will cause severe personal injury, death, or substantial property damage* if the instructions are ignored.

Example:



PERSONAL INJURY HAZARD!! TURN OFF ALL POWER TO CABINET TO AVOID ELECTROCUTION.



Identifies a condition or action for which a specific hazard is known to exist which *may cause severe personal injury, death, or substantial property damage* if the instructions are ignored.

Example:



POSSIBLE FIRE HAZARD! FOR CONTINUED PROTECTION AGAINST FIRE, REPLACE FUSE ONLY WITH ONE OF SAME TYPE AND RATING.

10- IMPORTANT WARNING STATEMENTS (continued)



Identifies a condition or action for which a potential hazard may exist *which will or can cause minor personal injury or property damage* if the instructions are ignored.

Example:



Possible equipment damage. Do not block airflow.

11- IMPORTANT SYMBOLS



Attention, consult accompanying documents



Equipotentiality



Protective Earth (ground)



Alternating Current



Three-phase Alternating Current with neutral conductor



Type B Equipment



Dangerous Voltage



Off (only for a part of EQUIPMENT)

11- IMPORTANT SYMBOLS (continued)

On (only for a part of EQUIPMENT)

REVISION HISTORY

REV	DATE	AUTHOR	PRIMARY REASONS FOR CHANGE
0	June 15, 1998	Erbtec Engineering	Initial Toolbook conversion to Word. (KK)