

TABLE OF CONTENTS

TABLE OF CONTENTS	1
1- IMAGE QUALITY PROBLEM	2
2- FAST RECEIVER DIAGNOSTICS	2
REVISION HISTORY	3

Description - This material covers troubleshooting the Fast Receiver Module.

1- IMAGE QUALITY PROBLEM

1. Perform a normal scan.
2. Perform an EPI scan with a phantom, using the suspect coil.
3. During the scan, verify the activity of the LEDs on the front panel of the Fast Receiver Module:
 - UART LED: Flashes during scan.
 - LO LED: On all the time.
 - A/D CLK LED: On and off during scan.
 - DAT BAK LED: On and off during scan.
 - Bottom Row, 100 kHz, 200 kHz, 300 kHz, 400 kHz, and 500 kHz: Appropriate LED is lit at download of scan based on TSP time of scan.
4. Disconnect Digital Interface cable and 2.5 mHz RF In cable from Receiver 0, and do a normal scan. Reconnect Digital Interface cable and 2.5 mHz RF In cable to Receiver 0.
5. With an oscilloscope, examine the analog signals below:
 - SEQ (digital signal)
 - WN (digital signal)
 - I/Q MON (analog signal, 1 kohm)
 - I/Q OUT (analog signal, 50 ohm)

2- FAST RECEIVER DIAGNOSTICS

See the procedure Fast Receiver Module Diagnostics.

REVISION HISTORY

REV	DATE	AUTHOR	PRIMARY REASONS FOR CHANGE
0	June 3, 1998	J. Saperstein	Initial Conversion from Toolbook to Word.
1	Oct 13, 1999	M. Keber	Added correct proprietary heading to document.