

TABLE OF CONTENTS

TABLE OF CONTENTS	1
1- INTRODUCTION	2
1-1 Process Fork Error.....	2
1-2 Default Limit on the Total Number of Processes Per User	2
2- INCREASING THE DEFAULT LIMIT FROM 250 TO 400 PER USER	2
REVISION HISTORY	3

1- INTRODUCTION

These instructions apply to all MR Signa LX systems equipped with an interface to Focused Ultrasound (FUS) option and operating with 9.1 M4 and M4A software.

1-1 Process Fork Error

On 9.1 scanners that interface with the FUS option, while prescribing scans more than 110 series, a process fork error occurs with the message "resource temporarily unavailable". Subsequent operating system commands executed as **signa** user fails with a "fork failed" error message.

1-2 Default Limit on the Total Number of Processes Per User

The default limit for the maximum number of processes allowed per user on the IRIX operating system in the 9.1 scanners is configured to be 250 processes per user. If scans are prescribed to have more than 110 series, there will be more than 250 processes running as **signa/sdc** user on the system.

2- INCREASING THE DEFAULT LIMIT FROM 250 TO 400 PER USER

On a 9.1 scanner with the FUS option, login as **root** user and from a command window follow the steps below to change the default limit on the number of processes allowed per user:

1. To check the current set limit, type the command **systune maxup** from a command window as **root** user as shown in the example below:

```
{sdc@t4} [2] systune maxup  
maxup = 250 (0xfa)
```

2. If the current set limit is NOT set to 400 and shows 250 as displayed from the output of Step 1 above, proceed to change the limit to 400 with the command **systune maxup 400** and confirm the change, when prompted, by typing **y** as shown in the example below. (This step is not necessary if the output of Step 1 above shows 400):

```
{sdc@t4} [3] systune maxup 400  
maxup = 250 (0xfa)  
Do you really want to change maxup to 400 (0x190)? (y/n) y  
{sdc@t4} [4]
```

3. Restart the system for the changes made in Step 2 above to take effect.

Note:

This procedure needs to be re-executed following every Load from Cold software installation of the 9.1 scanners with FUS.

REVISION HISTORY

REV	DATE	AUTHOR	PRIMARY REASONS FOR CHANGE
1	5/28/2009	Maya Rai	Initial Version.