

TABLE OF CONTENTS

TABLE OF CONTENTS	1
1- OVERVIEW	2
2- PATIENT TRANSPORT (PT 1) - SHEET 2	2
2-1 Hydraulics	2
2-2 Pedal Functions.....	2
2-3 Cradle Interlock.....	2
2-4 Emergency Release Lever	3
REVISION HISTORY.....	4

Description - This section contains theory for the Patient Transport. Sheet numbers listed at the end of some subsection titles refer to the sheet number within the Patient Handling Block Diagrams.

1- OVERVIEW

Hardware covered in this section includes the Patient Transport for Signa systems.

2- PATIENT TRANSPORT (PT 1) - SHEET 2

The patient transport is a lightweight, portable, patient table that can be rigidly attached to the magnet enclosure by means of the dock assembly, described in the magnet enclosure theory. The dock assembly provides a secure clamping station to which the transport hook may be affixed during patient transfer. The transport is held in place by a single 500-lb force provided by an over-the-center cam mounted on the transport. The table is capable of being raised and lowered whether it is attached to the magnet enclosure, or at a remote location. During operation, three caster locks prevent the movement of the unit during patient transfer. One caster has a directional lock for straight line tracking.

2-1 Hydraulics

A single-acting hydraulic pump, located at the docking end of the table, provides hydraulic pressure to raise the table. The hydraulic pump receives its power, either from an electric motor located in the dock, or from a foot pedal located at the rear of table. The table is lowered by gravity when the down valve, located at the base of the hydraulic cylinder on the table, is opened. This valve is mechanically linked to a foot pedal on the table, and when table is docked, to a set of pedals located on either side of the dock assembly.

2-2 Pedal Functions

Table functions are controlled by four foot pedals located at the rear of the table:

First pedal, at left, activates the docking mechanism to lock the table to the magnet enclosure. Second pedal detaches the table from the magnet enclosure. Third pedal operates the hydraulic pump (by means of mechanical linkage connected to the foot pedal) to raise the table. Fourth pedal lowers the table by opening the down valve.

2-3 Cradle Interlock

The table is also equipped with a cradle interlock, which prevents the table from being undocked or lowered unless the cradle is in the proper position on the table top. Inadvertent movement of the cradle is also protected by three mechanical latches built into the transport.

2-4 Emergency Release Lever

There is an emergency release lever on the dock assembly that allows the operator to free the transport in the event of a mechanical or hydraulic failure. When pulled, the dock hook is pulled completely away from the transport hook, freeing the transport. A spring returns the dock hook to its normal position.

REVISION HISTORY

REV	DATE	AUTHOR	PRIMARY REASONS FOR CHANGE
0	AUG 10, 1998	W. Pasciak	Initial release converted to Word