

D-1 Probe-S SNR Protocol

1. Run any proton head localizer to ensure the phantom is centered (up/down/in/out/left/right).
2. Set up the scan prescription per Table D-1.
3. For *TwinSpeed*, add the GradMode field, WHOLE or ZOOM.

TABLE D-1
PROBE-S SNR PROTOCOL

<u>PATIENT REGISTER .</u>	[New Pt]	<u>USER CVs SCREEN</u>	
<u>PATIENT INFORMATION</u>		CV3 scan mode	[1.00]
Patient Id	geservice	CV4 total # of scans	[32.00]
Patient Name	probes snr	CV15 1=fast version	
Weight (Lb)	111 - IMPORTANT	(fewer disdaqs)	[0] (default)
	[Landmark]	NOTE:	
Landmark	[>] [Nasion]	CV15 does not appear if VLTC was performed.	
		<u>SCANNING RANGE</u>	
<u>PATIENT PROTOCOLS</u>	[Patient Position]	FOV	[24]
		Start	S10
<u>PATIENT POSITION</u>		End	R10
Patient Position	[>] [Supine]	Start	L10
Patient Entry	[>] [Head First]	End	I10
Coil	[...] [Head] [Accept]	Start	A10
<u>IMAGING PARAMETERS</u>		End	P10
Plane	[>] [Axial]	Table Delta	0.00 (default)
Mode	[>] [MRS]	<u>ACQUISITION TIMING</u>	
Pulse Seq	[...] [Probe-S]	Freq	[256]
	[Accept]	Phase	[160]
Imaging Options	[Ext. Dyn. Range]	NEX	[2.00]
Psd Name	no entry	Freq Dir	[>] [A/P]
Protocol	no entry	Auto Center Freq	[>] [Water]
<u>SCAN TIMING</u>			Autoshim selected
* of Echoes	1 (default)		Contrast not selected
TE	[30]	[Save Series]	
TR	[2000]	[Prepare to Scan]	
		[Research Operations] [Display CVs]	
		CV Name	CV Value
		tempC	(greenest value on MRS Phantom)
			[Accept]
		[Research Operations] [Download]	

D-2 Probe-P SNR Protocol

1. Run any proton head localizer to insure the phantom is centered (up/down/in/out/left/right).
2. Set up the scan prescription per Table D-2.
3. For *TwinSpeed*, add the GradMode field, WHOLE or ZOOM

TABLE D-2
PROBE-P SNR PROTOCOL

<u>PATIENT REGISTER</u>	[New Pt]	<u>USER CVs SCREEN</u>	
<u>PATIENT INFORMATION</u>		CV3 scan mode	[1.00]
Patient Id	geservice	CV4 total # of scans	[32.00]
Patient Name	probep snr	CV15 1=fast version	
Weight (Lb)	111 - IMPORTANT	(fewer disdaqs)	[0] (default)
	[Landmark]	NOTE:	
Landmark	[>] [Nasion]	CV15 does not appear if VLTC was performed.	
		<u>SCANNING RANGE</u>	
<u>PATIENT PROTOCOLS</u>	[Patient Position]	FOV	[24]
		Start	S10
<u>PATIENT POSITION</u>		End	R10
Patient Position	[>] [Supine]	Start	L10
Patient Entry	[>] [Head First]	End	I10
Coil	[...] [Head] [Accept]	Start	A10
<u>IMAGING PARAMETERS</u>		End	P10
Plane	[>] [Axial]	Table Delta 0.00 (default)	
Mode	[>] [MRS]	<u>ACQUISITION TIMING</u>	
Pulse Seq	[...] [Probe-P]	Freq	[256]
	[Accept]	Phase	[160]
Imaging Options	[Ext. Dyn. Range]	NEX	[2.00]
Psd Name	no entry	Freq Dir	[>] [A/P]
Protocol	no entry	Auto Center Freq	[>] [Water]
<u>SCAN TIMING</u>			Autoshim selected
* of Echoes	1 (default)		Contrast not selected
TE	[37]	[Save Series]	
TR	[2000]	[Prepare to Scan]	
		[Research Operations] [Display CVs]	
		CV Name CV Value	
		tempC (greenest value on MRS Phantom)	
		[Accept]	
		[Research Operations] [Download]	