
GEM Flex Coil, 1.5T and 3.0T Service Manual

Distributed by GE Healthcare

NeoCoil
Innovations to **sharpen** your image

NC021SCM Rev 17
2023-02
©2023 NeoCoil, LLC
All rights reserved.

Medical Device Regulation

This product conforms to the requirements of EU medical device regulation 2017/745 concerning medical devices when it bears the following CE Marking of Conformity; first year of CE Mark for 1.5T GEM Flex Coil is 2011 and first year of CE Mark for 3.0T GEM Flex Coil is 2012:

European Representative


EC REP EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

United Kingdom Representative

UK REP EMERGO CONSULTING (UK) LIMITED
c/o Cr360 – UL International
Compass House, Vision Park Histon
Cambridge
CB24 9BZ
United Kingdom

Swiss Representative

CH REP MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

CAUTION: The waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Contact an authorized representative of the manufacturer for information concerning the decommissioning of your equipment.

US Federal Law 21 CFR Part 801

CAUTION: Federal law restricts this device to sale by or on the order of a physician.

Contents

Magnetic Field Warning	v
Important Information	vi

Section 1: Introduction

Description	1
Compatibility	2
System Components	3
Accessories	5
Related Documentation	5
Damage in Transportation	5
Product Locator	6
Environmental Requirements	7
Coil Specifications	7
Accessory Specifications	7
Use of Symbols	8
Classification of ME Equipment	12
Coil Care	12
Visual Inspection	13

Section 2: Installation

Set-up Procedure	15
Personnel Requirements	15
Preliminary Requirements	15
Tools and Test Equipment	15
Coil and Phantom Setup	16
Coil Configuration Procedure	16
Multi-Coil Quality Assurance (MCQA) Test Procedure	18
Coil and Phantom Setup	18
QA for the GEM Flex Pediatric Interface	21

Section 3: Replacements

FRUs and Replacements	23
1.5T GEM Flex Coil Field Replaceable Units	23
3.0T GEM Flex Coil Field Replaceable Units	24
Consumable Components	24
Cable Replacement	25
Personnel Requirements	25
Preliminary Requirements	25
Tools and Test Equipment	25

Contents (continued)

Section 3: Replacements (continued)


Replacing Foam Pocket on Array	27
Personnel Requirements	27
Preliminary Requirements.....	27
Tools and Test Equipment.....	27
Antenna or Interface Replacement.....	34
Personnel Requirements	34
Preliminary Requirements.....	34
Tools and Test Equipment.....	34

Section 4: Troubleshooting

General Notes	35
Failed MCQA	35
Bias Fault Message.....	38
Coil Not Recognized	38
Artifacts.....	38

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ATTENTION!

Installation of this coil **MUST** be performed by a GE field engineer. Please contact your field engineer to arrange for installation.

Attempts to perform installation of this coil without proper system configuration can result in system errors.

1 Introduction

Description

The GEM Flex Coil is a unique system made up of high density 16 element antenna arrays that are formable and flexible. Made up of interchangeable parts, the GEM Flex Coil is a receive-only coil designed to provide high resolution, sub millimeter imaging. The versatility of the system enables many clinical applications for 8 and 16 channel systems. The array system is compatible with 2D, 3D, parallel and isotropic imaging, as well as coil signal intensity correction. The 1.5T and 3.0T GEM Flex Coils are compatible with GE 1.5 and 3.0 HD or DV Series MR scanners, respectively.

NOTE: The operator manual indicates that scanning 32 channel systems is possible, however GE has disabled this feature in the coil configuration files.

The **antenna arrays** provide coverage for 8 or 16 channel scanning, as well as allow for greater coverage of large anatomy. Each of the arrays features up to 16 elements to provide deep penetration SNR and improved S/I and R/L coverage.

The **interface** is the device into which an array is plugged. It can be used in conjunction with the immobilization device or alone.

The **GEM Flex Pediatric Interface** is installed in and can only be used with the NeoCoil Pediatric Stabilizer.

Compatibility

This system is intended for use on GE Magnetic Resonance Scanners with HD-Connectors (used on HD series scanners, 1.5T only) or P-Connectors (used on DV series scanners, 1.5T and 3.0T). See pictures below.

HD-Connector

P-Connector

The system compatibility of the GEM Flex Pediatric Interface is as shown in the table below:

System Type	Minimum Software
MR450w-GEM	DV25.0

System Components

Contents will vary depending on the configuration purchased by your site.


Description	Part	NeoCoil part #	GE part #
16 channel fixed interface for HD 1.5T MRI systems 	GEM Flex Interface 16ch Fixed, 1.5T HD-Connector	NC022203	5430000-17
16 channel fixed interface for 1.5T DV MRI systems 	GEM Flex Interface 16ch Fixed, 1.5T P-Connector	NC022204	5430000-5
16 channel fixed interface for HD 3.0T MRI systems 	GEM Flex Interface 16ch Fixed, 3.0T HD-Connector	NC032203	5430008-6
16 channel fixed interface for 3.0T DV MRI systems 	GEM Flex Interface 16ch Fixed, 3.0T P-Connector	NC032204	5430008-5
16 channel fixed interface for 1.5T GE 450w-GEM systems for use with the NeoCoil Pediatric Stabilizer 	GEM Flex Pediatric Interface, 16ch Fixed, 1.5T P-Connector	NC022210	5430000-25

System Components (continued)

Description	Part	NeoCoil part #	GE part #
16 channel large array for hip 	GEM Flex Coil 16-L Array, 1.5T Receiver Only GEM Flex Coil 16-L Array, 3.0T Receiver Only	NC030200 for 1.5T NC040200 for 3.0T	5430000-2 or 5430000-22 for 1.5T 5430008-2 or 5430008-12 for 3.0T
16 channel medium array for knee 	GEM Flex Coil 16-M Array, 1.5T Receiver Only GEM Flex Coil 16-M Array, 3.0T Receiver Only	NC023200 for 1.5T NC033200 for 3.0T	5430000-3 or 5430000-23 for 1.5T 5430008-3 or 5430008-13 for 3.0T
16 channel small array for elbow 	GEM Flex Coil 16-S Array, 1.5T Receiver Only GEM Flex Coil 16-S Array, 3.0T Receiver Only	NC024200 for 1.5T NC034200 for 3.0T	5430000-4 or 5430000-24 for 1.5T 5430008-4 or 5430008-14 for 3.0T

System Components (continued)

Accessories

Description	Part	NeoCoil part #	GE part #
Knee stabilization device for flat MRI tables 	GEM Flex Knee Stabilization for Flat MRI Tables	NC026406	5430000-7
Knee stabilization device for curved MRI tables 	GEM Flex Knee Stabilization for Curved MRI Tables	NC026407	5430000-8
General purpose strap for GEM Flex Coil arrays 	GEM Flex General Purpose Strap	NC026403	NA
General purpose stabilization device for curved and flat tables 	GEM Flex Stabilization Device for Curved and Flat Tables	NC026402	5430000-13
Cover for the GEM Flex Interface Module 	GEM Flex Interface Cover	NC022504	NA

Related Documentation

- 1.5T GEM Flex Coil Instruction for Use, NeoCoil Part # NC021IFU
- 3.0T GEM Flex Coil Instruction for Use, NeoCoil Part # NC031IFU
- Instructions for Use Accessory Immobilization Devices, NeoCoil Part # NC074IFU
- Flex Coil Accessory Immobilization Devices Service Manual, NeoCoil Part # NC074SCM

Damage in Transportation

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or signed for by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier immediately upon discovery or, in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within a 14-day period.

To file a report with GE:

- Call 1-800-548-3366 and use option 6.
- Fill out a report at: <http://edq.health.ge.com/edq>
- Contact your local Service Coordinator for more information on this process.

Product Locator

Each serialized component is supplied with a product locator card (PLC). The number of PLCs for any given site will depend upon the MCat that's been ordered. Make sure to submit the complete number of installation cards for each catalog.

The Global Install Base Database (also known as the Product Locator System) tracks shipment, transshipment, and field location of the serialized models. There are two methods for submitting Product Locator information.

Method 1: For **U.S. ONLY** at this time, the preferred method of submitting information is the FE Site Verification website at: http://gib.gehealthcare.com/gib/gib_noaccess.jsp

The FE Site Verification has two components available on the web from the main menu:

- Installation/deinstallation product locator model and serial numbers
- Ability to add/modify the ship to address information

Method 2: Once the "Shipping Card" is filled out and submitted when shipped (extra cards are supplied for transshipments between storage and distribution points), and the "Installation Card" and extra shipment cards are attached:

- Verify that serial and model number on each rating plate matches the installation card numbers before removing the installation card.
- Process just the installation card and discard any extra shipment cards and labels. (There may be one or more shipment cards and bar code labels with the installation card. These shipment cards are used to trace the transfer of serialized units between various inventory storage and distribution points until the product reaches its final installation destination.)

NOTE: There could be as many as six of product locator cards, depending on the configuration ordered by the site. Make sure each applicable card is completed and submitted for proper warranty coverage and so that site history is maintained.

Environmental Requirements

Store the coil and accessories in an air conditioned scan room or equipment room until the GE field engineer can perform installation.

	Transport/Storage Conditions	Operating Conditions
Humidity	0% to 100% Non-condensing	30% to 75% Non-condensing
Temperature	-30° C to 55° C -22° F to 131° F	15° C to 21° C 59° F to 69.8° F
Atmospheric Pressure	50 kPa - 106 kPa	

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To maintain the highest possible image quality, and to protect patients and clinicians from potential injury or illness, it is important to keep the GEM Flex Coil clean and well maintained.



See NC074SCM for additional accessory specifications.

Use of Symbols

Coil users should understand the following symbols used for caution and warning explanations:

Symbol	Description
	Class II equipment
	Type BF applied
	CAUTION: Indicates a hazard that may cause a minor injury or equipment damage.
	WARNING: Indicates a hazard that may cause serious injury.
	Medical Device
	Do not loop cable.
	Dispose of RF coil properly.
	Consult Instructions for Use
	Instructions for Use are available in an electronic format.

Use of Symbols (continued)

Symbol	Description
	Place foam cover over the interface prior to use.
	Temperature limit
	Humidity limitation
	Atmospheric pressure limitation
	Receive only RF coil
	Keep dry
	Keep away from sunlight
	Fragile, handle with care
	Catalogue number
	Serial number
	Manufacturer
	Date of manufacture

Use of Symbols (continued)

Symbol	Description
	Country of manufacture
	Unique device identifier
	Importer
	Distributor
	Authorized representative in the European Community
	CE Marking of Conformity
	ETL Recognized Component (Canada and U.S.)
	UK Responsible Person
	UK Conformity Assessed
	Swiss Authorized Representative

Classification of ME Equipment

Classification & Ratings

CLASS II equipment in accordance with the international safety standard IEC 60601-1.

Type BF APPLIED PART complying with the specified requirements of the international safety standard IEC 60601-1.

This equipment is acceptable for Continuous Operation, having ordinary protection against ingress of liquids or fluids (IPX0).

All components are considered non-sterile and are not for use in the presence of flammable anesthetics.

Coil Care

WARNING: Electric shock may occur if the coil is attached to the system during cleaning or when it is still wet. Detach antenna cables from the scanner and interface before attempting to clean. Do not touch connectors with bare fingers. Never press sharp objects against connector surface. Do not reattach connector after cleaning until the coil has dried completely.

CAUTION: Do not spray or pour cleaning solution directly onto the antenna array, interface, or stabilization device. Do not submerge any of the parts in the solution. Parts of the equipment contain sensitive electronic components that may be damaged by the solution. Under no circumstances should the GEM Flex Coil be placed into any type of sterilizer. All parts should be cleaned only according to the procedure described below.

All of the parts of the GEM Flex Coil must be thoroughly cleaned and stored using the following procedure:

- Wipe with a cloth that has been dampened in a solution of 10% bleach and 90% tap water, or 30% isopropyl alcohol and 70% tap water.
- Should any part of the coil need to be returned for service, wipe it down with a 10% bleach solution (as described above) to minimize risk of exposure to potentially infectious agents.
- Allow the surfaces to dry completely.
- Store the GEM Flex Coil in an air conditioned scan or equipment room that meets the temperature and humidity parameters listed under *Environmental Requirements* on page 8.

Visual Inspection

After cleaning, visually inspect the following to ensure that the device and accessories remain safe and effective:

- Surfaces of the device and accessories are dry
- Soil and contamination have been removed

THIS PAGE LEFT
INTENTIONALLY BLANK.

2 Installation

Set-up Procedure

Personnel Requirements

Required Persons	Preliminary Reqs	Procedure
1	See Table Below	1 hour

Preliminary Requirements

Tools and Test Equipment

Equipment	Quantity	Field Strength	Phantom	GE Part #
16-S Array	1	1.5T & 3.0T	Small Cylindrical Unified Phantom	5342680
16-M Array	1	1.5T	Large Cylindrical Unified Phantom	5342679
		3.0T	Large Cylindrical Unified Phantom, SiOil	5342679-2
16-L Array	2	1.5T	TL Unified Phantom	5343347
		3.0T	TL Unified Phantom, SiOil	5343347-2
GEM Flex Pediatric Interface	1	1.5T	TL Unified Phantom	5343347

Coil and Phantom Setup

NOTE: Record the serial number of the coil(s) being used.

NOTE: The coil system cable must be directed towards the magnet.

CAUTION: All RF coil related tests must be run on a well-calibrated system that has passed all system tests, particularly white pixel and correlated noise.

Coil Configuration Procedure

Make sure that the latest software is loaded on the system prior to installing the GEM Flex Coil(s). Reference Compatibility table for minimum software requirements. Refer to GE Service procedure "Installing Software Updates" for how to load a Service Pack.

NOTE: Each array must be configured separately.

1. Plug the interface and array into the system connector.
2. Check the Coil Configuration Manager (HD Scanner) or the Coil Database Explorer (DV Scanner) to ensure the following configurations are installed.

Antenna	Coil Mode	Description
GEM Flex Coil 16-S Array	GEM Flex Small	16 channels
GEM Flex Coil 16-M Array	GEM Flex Medium	16 channels
GEM Flex Coil 16-L Array	GEM Flex LG Full	16 channels
GEM Flex Coil 16-L Array	GEM Flex LG Ctr	Medial 10 channels

Coil Configuration Procedure (continued)

NOTE: On the 450w-GEM Scanners, the following Coil Modes are available with the GEM Flex Pediatric Interface:

Antenna	Coil Mode	Description
GEM Flex Coil 16-L Array	Head Only	Head Only
GEM Flex Coil 16-M Array	Spine Only	Spine Only
GEM Flex Coil 16-M Array	Up Spine 32cm	Upper Spine 32cm
GEM Flex Coil 16-M Array	Bot Spine 32cm	Bottom Spine 32cm
GEM Flex Coil 16-M Array	Up Spine 24cm	Upper Spine 24cm
GEM Flex Coil 16-M Array	Bot Spine 24cm	Bottom Spine 24cm
GEM Flex Coil 16-M Array	Low Spine 24cm	Lower Spine 24cm
GEM Flex Coil 16-L Array	Head and Up Spine 32cm	Head and Upper Spine 32cm
GEM Flex Coil 16-M Array		

- Proceed with MCQA test procedure.

Multi-Coil Quality Assurance (MCQA) Test Procedure

Required Persons	Preliminary Reqs	Procedure	Finalization
1	0 min	20min	0min

Coil and Phantom Setup

NOTE: Record the serial number of the coil(s) being used. operating time, and FE number.

NOTE: The coil system cable must be directed towards the magnet.

NOTE: Each array must be tested separately.

CAUTION: All RF coil related tests must be run on a well-calibrated system that has passed all system tests, particularly white pixel and correlated noise.

1. Remove any other surface coils or pads from the table.
2. Place the appropriate phantom on the table in the R/L and S/I directions. When using two brick phantoms, ensure that the phantoms are aligned with and touching each other.
3. Lay the 16-L array lengthwise on top of the phantoms, with the spine of the array (cable exit) on top of the intersection of the phantoms. Ensure cable exits to the right of the bore; orientation should be maintained as shown in Figure 1. If testing a 3.0T 16-L array, secure the array to the phantom using a strap. See Figure 1.

Position the 16-S and 16-M arrays so that the spine of the array (cable exit) is resting on the table. See Figure 1.

4. Plug the connector from the antenna into the interface, ensuring that the connection is secure.
 - A green LED light on the interface will illuminate when the array is connected.

Multi-Coil Quality Assurance (MCQA) Test Procedure (continued)

5. Route the system cable down the center of the table, taking care to avoid looping the cable.
6. Landmark on the center of the array. See arrow in Figure 1.


Figure 1 Landmarking the 16-L array (3.0T array shown) on the left and the 16-S or 16-M array on the right

7. From the MR service desktop window, Image quality tab, select Multi-Coil QA Tool.
8. Confirm that the serial number in the MCQA tool matches the equipment being tested.

Multi-Coil Quality Assurance (MCQA) Test Procedure (continued)

9. Check for PASS/FAIL in the last column in the Status section. See Figure 2.


Figure 2 Sample user interface for MCQA

IMPORTANT!: The sample screen above is for reference only. Values shown in the sample above may not match the values seen at the customer site.

When the test is completed the results display the integrated Signal-to-Noise Ratio (SNR) as a function of the receive channel number. The PASS/FAIL status will show PASS when all coil elements are functioning properly. If the MCQA indicates any problems go to the Troubleshooting section.

QA for the GEM Flex Pediatric Interface

The GEM Flex Pediatric Interface is installed within the NeoCoil Pediatric Stabilizer.

IMPORTANT!: Leave the Pediatric Interface in the Pediatric Stabilizer to perform QA testing.

1. Remove any other surface coils or pads from the table.
2. Place the NeoCoil Pediatric Stabilizer on the table oriented in a head-first position.
3. Place your hand in the gap, with your fingers under the cover, and gently pry cover loose. See Figure 3.

Figure 3 Cover on Pediatric Stabilizer

4. Look through the window in the base of the Stabilizer and make note of the serial number that is applied to the Pediatric Interface.
5. Replace the cover back on to the base of the Stabilizer.
6. Lay the 16-M array on the flat section of the Stabilizer.
7. Plug the connector from the antenna into the interface, ensuring that the connection is secure.
 - A green LED light on the interface will illuminate when the array is connected.
8. Route the system cable down the center of the table, taking care to avoid looping the cable, and plug into the connector at the foot end of the table.
9. Advance the table to the bore and landmark on the center of the array.

QA for the GEM Flex Pediatric Interface (continued)

10. Place the appropriate phantom on top of the 16-M array so that it is centered on the array.
See Figure 4.

Figure 4 QA Setup for Pediatric Interface

11. From the MR service desktop window, Image quality tab, select Multi-Coil QA Tool.
12. Enter the serial number for the GEM Flex Pediatric Interface in the MCQA tool.
13. Check for PASS/FAIL in the last column in the Status section. See Figure 5.


Figure 5 Sample user interface for MCQA

When the test is completed the results display the integrated Signal-to-Noise Ratio (SNR) as a function of the receive channel number. The PASS/FAIL status will show PASS when all coil elements are functioning properly.

IMPORTANT!: The sample screen above is for reference only. Values shown in the sample above may not match the values seen at your site.

3 Replacements

FRUs and Replacements

1.5T GEM Flex Coil Field Replaceable Units

Description	NeoCoil Part #	GE Part #	Replacement Time
GEM Flex 16ch System Cable, 1.5T HD-Connector	NC022006	5430000-21	45 min.
GEM Flex 16ch System Cable, 1.5T P-Connector	NC022007	5430000-15	45 min.
GEM Flex Coil 16-M Array, 1.5T Receive Only	NC023000	54300003 or 5430000-23	15 min.
GEM Flex Coil 16-S Array, 1.5T Receive Only	NC024000	54300004 or 5430000-24	15 min.
GEM Flex Coil 16-L Array, 1.5T Receive Only	NC030000	54300002 or 5430000-22	15 min.
GEM Flex Interface 16ch Fixed, 1.5T HD-Connector	NC022003	5430000-17	30 min.
GEM Flex Interface 16ch Fixed, 1.5T P-Connector	NC022004	5430000-5	30 min.
GEM Flex Pediatric Interface 16ch Fixed, 1.5T P-Connector	NC022210	5430000-25	15 min.
GEM Flex Knee Stabilization for Flat MRI Tables	NC026000	5430000-7	15 min.
GEM Flex Knee Stabilization for Curved MRI Tables	NC026001	5430000-8	15 min.

3.0T GEM Flex Coil Field Replaceable Units

Description	NeoCoil Part #	GE Part #	Replacement Time
GEM Flex 16ch System Cable, 3.0T P-Connector	NC032007	5430008-17	45 min.
GEM Flex Coil 16-M Array, 3.0T	NC033000	5430008-3 or 5430008-13	15 min.
GEM Flex Coil 16-S Array, 3.0T	NC034000	5430008-4 or 5430008-14	15 min.
GEM Flex Coil 16-L Array, 3.0T	NC040000	5430008-2 or 5430008-12	15 min.
GEM Flex Interface 16ch Fixed, 3.0T P-Connector	NC032004	5430008-5	30 min.
GEM Flex Knee Stabilization for Flat MRI Tables	NC026000	5430000-7	15 min.
GEM Flex Knee Stabilization for Curved MRI Tables	NC026001	5430000-8	15 min.

Consumable Components

Pads, straps, and covers are designed as customer consumable components. They are not intended to be covered by any GE service contract. If pads, and straps are replaced they should be billed to the customer.

Description	NeoCoil Part #	GE Part #	Replacement Time
GEM Flex 16-M Coil Covering Replacement Kit	NC023001	5430000-11	1 hour
GEM Flex 16-S Coil Covering Replacement Kit	NC024001	5430000-12	1 hour
GEM Flex 16-L Coil Covering Replacement Kit	NC030001	5430000-10	1 hour
GEM Flex General Purpose Strap and Interface Cover Kit	NC026006	5430000-9	15 min.
GEM Flex Stabilization Device for Curved and Flat Tables	NC026002	5430000-13	15 min.
GEM Flex Coil Balun Cover Kit	NC090000	5430008-16	5 min.

Cable Replacement

CAUTION: Wear an ESD strap to avoid damage to the array due to electrostatic discharge.

Personnel Requirements

Required Persons	Preliminary Reqs	Procedure
1	See Table Below	45 min.

Preliminary Requirements

Tools and Test Equipment

Item	Qty.	Part #
Flathead Screwdriver	1	NA
ESD KIT w/MONITOR (WRIST STRAP, 10FT) or equivalent	1	2220483

If cable fails the continuity check (refer to [page 35 in the Troubleshooting section](#)) or if bent or damaged pins are identified then the system cable should be replaced. To do so, execute the following steps:

1. Remove the anti-tamper labels.
2. Remove the screws on the cover of the interface, by puncturing the white dot over the screwhead.


Figure 6 Interface cover screw locations

Cable Replacement (continued)

- Loosen the two screws on the bottom of the receptacle bracket shown in Figure 7.

Loosen these two screws only.

Figure 7 Screw locations on receptacle bracket

- Remove screws from old cable and system cable connector (SCC) boards.


Figure 8 Screw locations on cable clamp and SCC boards

- Attach new cable by carefully routing the wires between the modules and screwing the SCC boards (outlined in Figure 8) in place.

IMPORTANT!: If a pin on any of the SCC boards is bent, return the cable for replacement.

- Re-attach the cable clamp of the system cable.
- Replace the cover.
- Tighten the screws that secure the receptacle bracket, otherwise connection problems are likely to occur.
- Apply new anti-tamper labels.
- Run MCQA procedure as described in Installation section.

NOTE: Do not write on or mark up the cable being returned.

Replacing Foam Pocket on Array

CAUTION: There is a high risk of possibility for damage to the feed boards or cable. Replacing the whole array is recommended.

CAUTION: Wear an ESD strap to avoid damage to the array due to electrostatic discharge.

CAUTION: Hands must be clean prior to service to avoid contamination of the sensitive components.

Personnel Requirements

Required Persons	Preliminary Reqs	Procedure
1	See Table Below	1 hour

Preliminary Requirements

Tools and Test Equipment

Item	Qty.	Part #
Flathead screwdriver	1	NA
Tweezers (optional)	1	NA
Wire cutters	1	NA
ESD Kit w/monitor (Wrist Strap, 10ft) or equivalent	1	2220483
Alcohol swab	1	46-183039P1
Foam pocket kit	1	See FRU table

Replacing Foam Pocket on Array (continued)

1. Remove and discard the screw cover decal from the bottom of the spine.
2. Remove and discard the two screw cover plates from the top of the spine. On the 3.0T 16-L array, remove the velcro patch and underlying label from the top of the spine.
3. Remove the screws from the spine - four on top and four on bottom for the 16-S, 16-M, and 1.5T 16-L arrays; and six on top and seven on the bottom for the 3.0T 16-L array - puncturing the white dots on top of the screwheads. See Figure 9 and Figure 10.


Figure 9 Top spine cover on 16-S, 16-M and 1.5T 16-L array (left) and 3.0T 16-L array (right)


Figure 10 Bottom spine cover on 16-S, 16-M and 1.5T 16-L array (left) and 3.0T 16-L array (right)

4. Remove spine cover from top of array. Figure 10
5. Remove the two screws in the cable clamp and take off cover. Figure 11


Figure 11 Screws in cable clamp

Replacing Foam Pocket on Array (continued)

6. Remove the middle spine piece, being careful of the wires. On the 3.0T 16-L array, remove the screw below the cable clamp prior to removed the middle spine piece. Figure 12

Figure 12 3.0T 16-L Screw Below Cable Clamp

7. On the 16-L array only, remove cable ties from the foam pocket leaving the ties wrapped around the wires.
8. Remove the two exposed feed boards on the 16-S and 16-M arrays, or the one feed board on the 1.5T 16-L coil, by unscrewing the feed board(s) (Figure 13) and carefully lifting the feed board(s) from the antenna array. The feed board on the 3.0T 16-L does not require removal.


Figure 13 Screw locations on 16-L feed board

9. Unscrew the screw(s)/washer(s) holding the lower spine piece to the antenna array.
10. Remove the lower spine piece. Figure 14

Figure 14 Lower spine piece removed

11. Remove foam pocket by gently inserting your hand into the pocket, slowly lifting the covers away from the array.

Replacing Foam Pocket on Array (continued)

12. Holding down the array with one hand inside the pocket, use your other hand to gently remove the pocket off of the array. Figure 15

Figure 15 Hand position while removing pocket

13. Replace left and right foam pockets.

NOTE: On the 16-S and 16-M arrays, the end of the array that shows the S/N goes into the pocket with the white flap, per Figure 16.

With the 16-L array, insert the longer end of the array into the longer foam pocket.


Figure 16 Pocket/Array alignment with call-out of S/N end of array

Replacing Foam Pocket on Array (continued)

14. Align lower spine piece ensuring bosses and posts are inserted into their respective holes in the array and the foam. Figure 17

Figure 17 Alignment of bosses and posts in lower spine

15. Install the screw(s)/washer(s) that holds the lower spine piece to the array. Figure 18

Figure 18 Replacing the lower spine to the array

16. Replace middle spine piece around cable and wires aligning on top of feed board(s). On the 3,0T 16-L array, install the one screw below the cable clamp area.
17. Align and insert the bosses (16-S and 16-M arrays only) into the holes in the top of the foam. Figure 19

Figure 19 Inserting bosses on top of foam pocket

18. On the 16-S, 16-M and 1.5T 16-L arrays, align and fully seat feed board(s) gently so as not to bend the delicate pins. Figure 3-15 left picture.

IMPORTANT!: If a pin on the feed board is bent return the entire array for replacement.

Replacing Foam Pocket on Array (continued)

19. Install the feed board(s) with screws (two each). Figure 20 right picture.

Pins on feed board

Figure 20 Installing the feed boards

20. On the 16-L only, replace the cable ties from where they were removed.
21. Place the middle spine piece in position ensuring posts are inserted into the holes in the foam; the bosses are aligned with the holes in the array and lower spine piece. Figure 21

IMPORTANT: Avoid crushing the wires.

Figure 21 Middle spine piece

22. Install the cable clamp using the two screws, making sure the coil clamp is in proper alignment. Be sure rating plate on the cable is up when the array is up. Figure 3-17


Figure 22 Cable clamp alignment, shown with spine piece already attached

Replacing Foam Pocket on Array (continued)

23. Gently squeeze the two spine pieces together while replacing the three screws that hold the lower spine piece to the middle spine piece on the 16-S, 16-M, and 1.5T 16-L arrays, Figure 3-18; or the seven screws on the 3.0T 16-L array.

Figure 23 Middle spine piece

24. Install the upper spine piece using four screws on top and one on the lower spine piece on the 16-S, 16-M, or 1.5T 16-L arrays; or the six screws on top and one on the lower spine piece for the 3.0T 16-L array.

25. Clean off residual adhesive from the lower and upper spine parts.

26. Apply the new screw cover decal (bottom) and the two new screw cover plates (top).

NOTE: Remove adhesive cover strip on screw cover plates prior to installation.

27. For the 3.0T 16-L array only, place the white label over the two exposed screws, within the recessed area, on the top spine piece. Then place the hook patch in the recessed area over the label. Figure 24


Figure 24 White label (outlined for better visibility) on 3.0T 16-L array

28. Run MCQA test.

Antenna or Interface Replacement

Personnel Requirements

Required Persons	Preliminary Reqs	Procedure
1	See Table Below	15 min.

Preliminary Requirements

Tools and Test Equipment

Item	Qty.	Part #
NA	NA	NA

Perform MCQA Procedure as described in the Installation section.

4 Troubleshooting

General Notes

When encountering a problem with the GEM Flex Coil try these general troubleshooting tips before conducting further investigation.

- The green system light on the interface is illuminated when the antenna array is plugged into it and the interface is plugged into the system.
- All connections have been re-seated.
- There are no bent pins on any of the connectors.

When troubleshooting the Pediatric System:

- First remove both antennae from the Pediatric Stabilizer and follow the procedures in this section.
- Evaluate the 1.5T GEM Flex Interface 16ch Fixed, P-Connector by following the procedure in this section.
- Finally, after eliminating the antenna(e) and GEM Flex Interface as the problem source, perform MCQA on the GEM Flex Pediatric Interface 16ch Fixed, 1.5T P-Connector.

Failed MCQA

When the MCQA analysis has reported a failure, the fault could be coming from either the antenna array, the interface or the system connector. To determine which is at fault, **run MCQA again with a different antenna** but with the same interface.


If the problem is resolved using a new antenna array, return the bad antenna array and order a new one.

If the problem remains, check the continuity of the interface cable with a multimeter using the applicable table and diagram, shown on the following pages. If the cable is causing the fault order a new one. You may need to disconnect each SCS board in order to access the soldered end of each wire within the interface. Follow steps for [“Cable Replacement” on page 25](#) in order to remove the cable for testing as well as replacing it afterwards.

If the problem still persists, and you have effectively ruled out the antenna and cable, return the interface for a replacement.

HD 16 Channel Cable

Inside of PA Module	System Connector	Diagram
DC WIRES	PIN	
BLACK	M-9	
BROWN	M-10	
BLUE	A-5	
GRAY	A-3	
CLEAR/BLACK	M-8	
CLEAR/RED	M-7	
BLACK	M-6	
WHITE	M-3	
RF WIRES	PIN	
WHITE/YELLOW	H-3	
CLEAR	D-4	
WHITE/BLACK	H-4	
GREEN	C-3	
WHITE	D-3	
WHITE/RED	H-1	
GRAY	D-2	
WHITE/ORANGE	H-2	
BROWN	C-1	
BLACK	D-1	
VIOLET	G-3	
ORANGE	G-1	
RED	C-2	
YELLOW	G-2	
WHITE/BROWN	G-4	
BLUE	C-4	

DV 16 Channel Cable

Inside of Interface	System Connector	Diagram
DC WIRES	PIN	
BLACK	1	
BROWN	2	
BLUE	9	
GRAY	7	
CLEAR/BLACK	5	
CLEAR/RED	6	
BLACK	19	
WHITE	18	
RF WIRES	PIN	
WHITE/YELLOW	R15	
CLEAR	R8	
WHITE/BLACK	R16	
WHITE	R7	
GREEN	R3	
WHITE/RED	R13	
GRAY	R6	
WHITE/ ORANGE	R14	
BROWN	R1	
BLACK	R5	
VIOLET	R11	
ORANGE	R9	
RED	R2	
YELLOW	R10	
WHITE/BROWN	R12	
BLUE	R4	

Bias Fault Message

If the message log contains a bias fault message, visually inspect the system cable and antenna for a bent pin.

If there is a bent pin in the system cable, return it for a new one. If, however, the bent pin is on the antenna cable return the antenna.

If there is no visibly bent pin, re-run the scan with a different GEM Flex Coil. If the bias fault is not seen with the second coil, return the coil that produced the error. If the bias fault happens again with the second GEM Flex Coil, return the system cable. *Do not write on or mark up the cable being returned.*

Coil Not Recognized

If the scanner does not recognize the coil after installation, reseal the connections and try the scan again. If this does not resolve the problem, visually inspect the system cable and antenna for a bent pin. Ensure that the scanner is at a compatible software level.

Artifacts

If there are artifacts when using the antenna-array, run MCQA and be certain that there is at least a three inch space between the coil and the bore wall.

THIS PAGE LEFT
INTENTIONALLY BLANK.

5 Contact Information

GE Healthcare Americas

USA 800-558-5102
Canada 800-668-0732

GE Healthcare Europe (33) 1-41-19-76-76

GE Healthcare Asia/Australia

China 86-21-62192228
Taiwan 886-2-2505-7900
Singapore 65-291-8528
Australia 61-2-9975-5501
Japan 81-120-48-2630
Korea 82-31-740-6119
India 91-80-845-2923

Website www.gehealthcare.com

NeoCoil, LLC
N27 W23910A Paul Rd.
Pewaukee, WI 53072 USA

Phone Number 1-888-674-2645

Website www.neocoil.com