GE Healthcare

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**WARNING
(EN)**

This service manual is available in English only.

- If a customer's service provider requires a language other than English, it is the customer's responsibility to provide translation services.
- Do not attempt to service the equipment unless this service manual has been consulted and is understood.
- Failure to heed this warning may result in injury to the service provider, operator or patient from electric shock, mechanical or other hazards.

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This page left intentionally blank.

Revision History

Rev	Date	Reason for Change
2	May, 2011	Updated for MR450w GEM.
1	July, 2010	Initial release.

This page left intentionally blank.

Table of Contents

Chapter 1 INTRODUCTION.....	15
1 1.5T and 3.0T 6-Channel Flex Coil Introduction.....	15
1.1 Shipping Contents.....	15
1.2 Damage in Transportation.....	15
1.3 Product Locator.....	16
1.4 Requirements and Specifications.....	16
1.5 Coil Description.....	17
1.6 Modes of Operation.....	18
Chapter 2 INSTALLATION.....	21
1 1.5T and 3.0T 6-Channel Flex Coil Setup for MCQA Test.....	21
1.3.1 Coil Configuration.....	21
1.3.2 Flex Coil Alignment.....	23
1.3.3 MCQA Test.....	26
Chapter 3 REPLACEMENTS.....	27
1 1.5T and 3.0T 6-Channel Flex Coil FRUs and Replacement.....	27
1.1 Introduction.....	27
1.2 Safety.....	27
1.3 Field Replaceable Units.....	27
1.4 Replacement.....	28
1.4.1 6-Channel Flex Phased Array Coil Half.....	28
1.4.2 Quick Disconnect (QD) Assembly.....	28
1.5 QD Assembly Pre-Amplifiers.....	29
1.6 Finalization.....	34
Chapter 4 TROUBLESHOOTING.....	35
1 1.5T and 3.0T 6-Channel Flex Coil Troubleshooting.....	35
1.1 Hardware Inspection.....	35
1.2 System Inspection.....	35
1.3 MCQA Analysis Report Overview.....	36

This page left intentionally blank.

Chapter 1 Introduction

1 1.5T and 3.0T 6-Channel Flex Coil Introduction

The 1.5T and 3.0T 6-Channel Flexible Phased Array Coil is compatible with the GE 1.5T MR450, MR450w, MR450w GEM and the GE 3.0T MR750, MR750w, and MR750w GEM whole body MRI systems. It is a receive only, multi-coil array designed to give the optimal Signal-to-Noise Ratio (SNR) and uniformity for high resolution examinations. The coil enables the imaging of multi-oblique slices of the head and other areas of the body. The multi-element design improves the coverage and resolution.

1.1 Shipping Contents

The table below includes the coil contents shipped for the specified catalog number.

Table 1-1: Shipping Contents by Catalog Number

Description for M0049SS	Part Number	Quantity
1.5T 6-Channel Flex Phased Array Coil with a Quick Disconnect Assembly	5363383	1
Operator Manual	5371823	1
Service Manual	5371824	1
Coil Configuration CD	5372128	1
Description for M0050SS	Part Number	Quantity
3.0T 6-Channel Flex Phased Array Coil with a Quick Disconnect Assembly	5363393	1
Operator Manual	5371823	1
Service Manual	5371824	1
Coil Configuration CD	5372128	1

1.2 Damage in Transportation

All packages should be closely examined at time of delivery. If damage is apparent, have the notation, "Damage in Shipment," written on all copies of the freight or express bill before delivery is accepted or signed for by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **must** be reported to the carrier immediately upon discovery or, in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14-day period.

To file a report:

- Call 1-800-548-3366 and select Option 6.
- Contact your local Service Coordinator for more information on this process.

1.3 Product Locator

NOTICE

The 6-Channel Flex Coil catalog contains three different, serialized components. Make sure to submit three Installation Cards for each catalog.

The Global Install Base Database (also known as the Product Locator System) tracks the shipment, trans-shipment, and field location of the serialized models. There are two methods for submitting Product Locator information.

Method 1: For U.S. ONLY at this time, the preferred method of submitting information is the FE Site Verification website at: <http://gein2.med.ge.com/gib>.

The FE Site Verification has three components available on the web site from the main menu:

- Installation/deinstallation product locator model and serial numbers
- Ability to add/modify the ship to address information
- Ability to obtain a copy of the FE training tutorial for using this web site. A downloadable copy is available on the Product Locator Support Central webpage: <http://supportcentral.ge.com/15563>.

Method 2: Once the Shipping Card is completed and submitted when shipped (extra cards are supplied for trans-shipments between storage and distribution points), and the Installation Card and extra shipment cards are attached:

- Verify that the serial and model number on each rating plate matches the installation card numbers before removing the installation card.
- Process the Installation Card and discard any extra shipment cards and labels. (There may be one or more shipment cards and bar code labels with the installation card. These shipment cards are used to trace the transfer of serialized units between various inventory storage and distribution points until the product reaches its final installation destination.)

NOTE: There could be as many as three product locator cards, one for each coil half and one for the Quick Disconnect (QD) assembly. Make sure each applicable card is completed and sent in for proper warranty coverage and that the site history is maintained.

1.4 Requirements and Specifications

Coil installation and configuration requires the services of a GE Service Representative. Store the coil in an air-conditioned Scan or Equipment Room.

Table 1-2: Environmental Requirements

Measurement Type	Transport and Storage	Operating Conditions
Barometric pressure	500 to 1060 hPa	700 to 1060 hPa
Humidity	10 to 100% noncondensing	10 to 90% noncondensing

Temperature	-40 to 65° C (-40 to 149° F)	15 to 21° C (59 to 69.8° F)
-------------	------------------------------	-----------------------------

Table 1-3: Coil Frequency Specifications

Field Strength	Coil Model	Frequency
1.5T	6-Channel Flex Phased Array Coil	63.86 MHz $\pm$ 50 kHz
3.0T	6-Channel Flex Phased Array Coil	127.73 MHz $\pm$ 50 kHz

1.5 Coil Description

The 1.5T and 3.0T 6-Channel Flexible Phased Array Coil (shown below) is a receive only, multi-coil array optimized for high resolution exams. The coil enables the user to image multi-oblique slices where the intended areas of interest to be included, but not limited to: head, knee, shoulder, foot, wrist, and other extremities.

Illustration 1-1: 6-Channel Flexible Phased Array Coil - Identical Halves

The 6-Channel Flexible Phased Array Coil consists of two flexible halves that enable large openings and flexibility. The flexible halves conform to the contours of the patient and accommodate the patient's anatomy. This maximizes signal integrity and clinical utility. Both coil halves are connected to the system through an intermediate adapter connector, called a Quick Disconnect (QD) assembly shown below.

NOTE: Coil setup flexibility allows the use of only one coil half if desired.

Illustration 1-2: Quick Disconnect Assembly

NOTE: The Quick Disconnect assembly connects directly to the Low Profile Carriage Assembly (LPCA) of the MRI system shown below.

Illustration 1-3: Low Profile Carriage Assembly

1.6 Modes of Operation

The 6-channel Flexible Phased Array Coil can be operated in multiple modes. Coil modes for neuro applications:

- **6ch Neuro Flex:** Both coil paddles are connected to the QD assembly, and all 6 channels are operational.

- **6ch Neuro Flex1:** Both coil paddles are connected to the QD assembly, but only the paddle connected to Port 1 on QD assembly and three channels are operational.
- **6ch Neuro Flex2:** Both coil paddles are connected to the QD assembly, but only the paddle connected to Port 2 on QD assembly and three channels are operational.
- **3ch Neuro Flex1:** Single paddle connected to Port 1 on QD assembly, and three channels are operational.
- **3ch Neuro Flex2:** Single paddle connected to Port 2 on QD assembly, and three channels are operational.

Coil modes for general purpose applications:

- **6ch GP Flex:** Both coil paddles are connected to the QD assembly, and all 6 channels are operational.
- **6ch GP Flex1:** Both coil paddles are connected to the QD assembly, but only the paddle connected to Port 1 on QD assembly and three channels are operational.
- **6ch GP Flex2:** Both coil paddles are connected to the QD assembly, but only the paddle connected to Port 2 on QD assembly and three channels are operational.
- **3ch GP Flex1:** Single paddle connected to Port 1 on QD assembly, and three channels are operational.
- **3ch GP Flex2:** Single paddle connected to Port 2 on QD assembly, and three channels are operational.

This page left intentionally blank.

Chapter 2 Installation

1 1.5T and 3.0T 6-Channel Flex Coil Setup for MCQA Test

1.1 Personnel Requirements

Personnel Requirements	Preliminary Reqs	Procedure	Finalization
1	Not Applicable	60 mins	Not Applicable

1.2 Preliminary Requirements

1.2.1 Tools and Test Equipment

Item	Qty	Effectivity	Part#	Manufacturer
Head TLT Sphere for 1.5T (green solution)	1	-	2135650	-
Head TLT Loader for 1.5T	1	-	2135652	-
Head TLT Sphere for 3.0T (pink solution)	1	-	2359877	-
Head TLT Loader for 3.0T	1	-	2360031	-
Curved Adapter Panel Assembly	1	GEM Systems	5395828	Sumi Motherson

1.2.2 Replacement Parts

Item	Qty	Effectivity	Part#	Manufacturer
Patient Security Strap (provided with system)	1	-	-	-


1.3 Procedure

1.3.1 Coil Configuration

Verify the coil configuration for the site. If the applicable 6-Channel Flex Coil is not configured for the site, select the coil from the configuration list:

1. From the Service Desktop, select **Configuration > Coil Database Explorer > Click Here to Start This Tool**.
2. On the *Coil Name* tab, select the pull-down menu for *Clinical Coil Name*.
3. Scroll through the list to verify the existence of the 10 coil modes for the 6-Channel Flex Coils (shown below).

Illustration 2-1: Clinical Coil Names

4. If the coil modes are not in the list, locate the Coil CD that shipped with the coil.
5. To install the 6-Channel Flex Coil configurations and protocols:
 - a. Insert the 6-Channel Flex Coil definitions disk into the default (typically upper) Host DVD drive.
 - b. Open a C-shell from the Service Desktop Manager, and type:

- **su**<Enter>
- <root password> <Enter>
- **mount /mnt/cdrom**<Enter>

NOTE: No action is required when this message displays: Mount : block device/dev/scd0 is write-protected, mounting read only.

- **mnt/cdrom/autorun**<Enter>


- c. When the *Coil Data Disk Installation* window displays, type **I** to install the CD contents.
- d. When the CD installation is complete, press <Enter> to exit the *Coil Data Disk Installation* window.
- e. Type **eject** into the original C-shell, and remove the CD from the Host DVD drive.
- f. Repeat [Step 1](#) through [Step 3](#) to verify the 6-Channel Flex Coil is in the list.

1.3.2 Flex Coil Alignment

The flex coil must be properly aligned prior to the MCQA test.

1. (For Systems with Flat GEM Table Only) Insert the Curved Panel Assembly into the magnet end of the table as shown.

Illustration 2-2: Curved Adapter Panel Assembly


2. Place one coil half onto the table top.

3. Place the phantom with loader onto the coil half.
4. Place the second coil half on top of the phantom loader.

NOTE: The flex coil cables should lay flat on the table top. Make sure these cables do **not** contain loops, such as those shown below.

5. Plug the flex coil Quick Disconnect (QD) into Port A on the Low Profile Carriage Assembly (LPCA).


Illustration 2-3: Flex Coil Connected to LPCA

6. Set the landmark using the alignment lights with the crosshairs marked on the anterior half of the coil.
 - a. At the magnet, press the [Alignment Light] button to turn on the laser alignment lights.

NOTE: Make sure the flex coil and phantom with loader are properly aligned (centered left to right) with the laser alignment lights.


Illustration 2-4: Setting Landmark on Flex Coil

- b. Move the cradle to align the coil landmarks with the laser alignment lights.
 - c. Press the [Landmark] button to save the alignment.
7. Install the patient security strap (provided with the system) to secure both coil halves and the phantom with loader to the table top as shown.

NOTE: Make sure the security strap is properly aligned (centered front to back) on the flex coil and phantom with loader.

Illustration 2-5: Flex Coil and Phantom with Loader

1.3.3 MCQA Test

To perform the Multi-Coil Quality Assurance (MCQA) test for a 1.5T or 3.0T system, locate the appropriate Service Methods DVD:

- Access the **Coil Vendor Manuals** in the Coils folder.
- Navigate to the Flex Coils table.
- Click on the link: *MCQA Tool*.

1.4 Finalization

No finalization steps.

Chapter 3 Replacements

1 1.5T and 3.0T 6-Channel Flex Coil FRUs and Replacement

1.1 Introduction

The two coil halves of the 6-Channel Flexible Phased Array Coil assembly are identical in design and manufacture. The coil halves are not matched pairs, and do not require tuning so individual coil halves can be ordered and replaced.

Each coil half plugs into an Intermediate Adaptor Connector, otherwise known as a Quick Disconnect (QD) assembly, which is then plugged into the Low Profile Carriage Assembly (LPCA) on the scanner. The QD assembly contains six imbedded pre-amps (one for each channel) that can be ordered separately.

Tools Required: 2 Flat-Blade Screwdrivers, 3/16 in. and 1/8 in.

1.2 Safety

NOTICE

To prevent damage to equipment components, observe the following Electro-Static Discharge (ESD) precautions:

- Work on a static-free mat.
- Wear an ESD wrist strap to ensure that any accumulated electro-static charge is discharged from your body to the ground.
- Create a common ground for the equipment to be worked on by connecting the static-free mat, ESD wrist strap and peripheral units to that piece of equipment.

For further information regarding ESD protection, refer to the Safety chapter in this publication.

1.3 Field Replaceable Units

Field Replaceable Units (FRUs) are listed below.

Table 3-1: FRU List

Description	Part Number
1.5T 6-Channel Flex Phased Array coil half (each coil half is identical; therefore, a single half is permitted for replacement)	5229623
1.5T DV 6-Channel Coil QD Assembly	5363382
3.0T 6-Channel Flex Phased Array coil half (Each coil half is identical so a single half may be replaced.)	5229646
3.0T DV 6-Channel Coil QD Assembly	5363382-3
6-Channel Flex Coil Accessories Spares Kit (Needed only when the coil is used on the Surgical Suite product; kit contains replacement parts to support the coil strain relief and mounting mechanism.)	5269901
1.5T Low Noise Pre-Amplifier (Channel location will be marked by the Field Engineer on the pre-amp with a permanent marker at the time of replacement.)	2409869

3.0T Low Noise Pre-Amplifier (Channel location will be marked by the Field Engineer on the pre-amp with a permanent marker at the time of replacement.)	2409871
Head TLT Sphere for 1.5T (Green Solution)	2135650
Head TLT Loader for 1.5T	2135652
Head TLT Sphere for 3.0T (Pink Solution)	2359877
Head TLT Loader for 3.0T	2360031

1.4 Replacement

The 6-Channel Flexible Phased Array Coil assembly consists of two coil halves and the QD assembly. There are no replaceable parts within each of the coil half assemblies.

1.4.1 6-Channel Flex Phased Array Coil Half

The two coil halves of the 6-Channel Flexible Phased Array Coil assembly are identical. Replace one coil half as needed.

1. Make sure the 6-Channel Flexible Phased Array Coil is moved out of the Magnet Room.
2. Disconnect the faulty 6-Channel Flexible Phased Array Coil connector from the QD assembly.
3. Remove the replacement 6-Channel Flexible Phased Array Coil from the shipping container.
4. Connect the new 6-Channel Flexible Phased Array Coil connector to the QD assembly.
5. Place the coil removed from QD assembly into the shipping container.
6. Set up the flex coil and perform the MCQA test as described in [Chapter 2, 1.5T and 3.0T DV 6-Channel Flex Coil Setup for MCQA Test](#).

1.4.2 Quick Disconnect (QD) Assembly

1. Make sure the 6-Channel Flexible Phased Array Coil is moved out of the Magnet Room.
2. Disconnect the 6-Channel Flexible Phased Array Coil first half connector from Port 1 of the QD assembly.

Illustration 3-1: QD Assembly - Port 1 and Port 2

3. Disconnect the 6-Channel Flexible Phased Array Coil second half connector from Port 2 of the QD assembly.
4. Remove the replacement QD assembly from the shipping container.
5. Connect the 6-Channel Flexible Phased Array Coil first half connector to Port 1 of the new QD assembly.
6. Connect the 6-Channel Flexible Phased Array Coil second half connector to Port 2 of the new QD assembly.
7. Place the QD assembly removed from the coil halves into the shipping container.
8. Set up the flex coil and perform the MCQA test as described in [Chapter 2, 1.5T and 3.0T 6-Channel Flex Coil Setup for MCQA Test](#).

1.5 QD Assembly Pre-Amplifiers

NOTICE


To prevent damage to equipment components, observe the following Electro-Static Discharge (ESD) precautions:

- Work on a static-free mat.
- Wear an ESD wrist strap to ensure that any accumulated electro-static charge is discharged from your body to the ground.
- Create a common ground for the equipment to be worked on by connecting the static-free mat, ESD wrist strap and peripheral units to that piece of equipment.

For further information regarding ESD protection, refer to the Safety chapter in this publication.


1. Move the QD assembly out of the Magnet Room.
2. Remove the four screws and bottom cover (shown below) from the QD assembly.

Illustration 3-2: QD Assembly - Four Screws and Bottom Cover


3. If necessary, use a 1/8 in. flat-blade screwdriver and gently pry the housing halves of the QD assembly apart as shown below.

Illustration 3-3: Separating Housing on QD Assembly

4. Remove the two screws and Lexan cover (shown below) from the QD assembly using a 3/16 in. flat-blade screwdriver.

Illustration 3-4: Quick Disconnect (QD) Assembly - Two Screws and Lexan Cover

**NOTICE**

Pay special attention while removing the pre-amp from the socket to avoid damage to the printed circuit, components, or pre-amp.


5. Remove the pre-amp from the QD assembly by prying under the pre-amp module on alternate sides of the module with a flat-blade screwdriver as shown.

Illustration 3-5: Removing Pre-Amp from QD Assembly


6. Refer to the illustration below to identify the socket and pin locations for the pre-amp module.

Illustration 3-6: Socket and Pin Locations on the Pre-Amp

7. Make a note of the channel (location) number on the pre-amp so the new pre-amp can be marked with the same number.

Illustration 3-7: Quick Disconnect - Preamp Locations for 1.5T and 3.0T

8. Remove the replacement pre-amp from the shipping container and static bag.
9. Install the new pre-amp into the QD assembly, and use a permanent marker to record the channel number on the new pre-amp.

Illustration 3-8: Quick Disconnect (QD) Assembly Pre-Amp - Marking Pre-Amp Number


10. Reinstall the Lexan cover and two screws on the QD assembly.
11. Place the pre-amp removed from the QD assembly into the static bag and into the shipping container.

12. Reinstall the bottom cover and four screws to the QD assembly.
13. Set up the flex coil and perform the MCQA test as described in [Chapter 2, 1.5T and 3.0T 6-Channel Flex Coil Setup for MCQA Test](#).

1.6 Finalization

Update the Global Installed Base (GIB) with the new serial number information of the product locator for each new assembly replaced.

Illustration 3-9: Quick Disconnect - ID Label Example

Chapter 4 Troubleshooting

1 1.5T and 3.0T 6-Channel Flex Coil Troubleshooting

The two coil halves of the 6-Channel Flexible Phased Array coil assembly are identical in design and manufacture. The foam sealed housing for the coil provides the required flexibility needed for neuro and general purpose imaging applications. This, however, restricts access to the internal components so troubleshooting the coil is performed at the coil half assembly level. The coil halves are not matched pairs, and do not require tuning so individual coil halves can be ordered and replaced.

To determine if the problem lies in the coil or the receive chain, try swapping the coil cables at the Quick Disconnect (QD) Assembly to see if the problem moves with the cable swap.

- If the problem moves, then the issue is with the coil half.
- If the problem does not move with the cable swap, use the MCRV Tool to determine where in the receive chain the problem lies.

Each coil half plugs into an Intermediate Adaptor Connector, otherwise known as a Quick Disconnect (QD) Assembly, which is then plugged into the Low Profile Carriage Assembly (LPCA) on the scanner. Imbedded pre-amps are provided in the QD Assembly.

1.1 Hardware Inspection

1. Inspect the foam housing of each coil half and look for cracks, discoloration or other damage.
2. Inspect the cables and look for signs of damage such as cuts in the insulation, kinks or twists, especially at the coil and connector entry points.
3. Inspect the cable connector and connector pins. Look for signs of damage to the pins such as corrosion, bent pins, or missing pins.
4. Make sure the field strength labeled on the coil is correct for the scanner it is used on.
5. Inspect the QD Assembly connectors for bent or missing pins.
6. Inspect the QD Assembly housing for cracks, discoloration or other damage.
7. Inspect the connectors on the LPCA and look for damaged or missing pins.
8. Make sure the coil cable connector is fully engaged into the QD Assembly connector.
9. Make sure the QD Assembly is fully engaged into the LPCA connector.

1.2 System Inspection

1. Verify the scanner's operating frequency is within the bandwidth of the coil. Refer to [Chapter 1, 1.5T and 3.0T 6-Channel Flex Coil Introduction](#) for the bandwidth specification for the coil. Corrections to system operating frequency may be necessary to improve coil performance.
2. Review the customer's images and look for shading or artifacts.

1.3 MCQA Analysis Report Overview

The following illustrations depict the type of report resulting from a faulty coil connection.

Illustration 4-1: Loss of Signal on Channel 5


Illustration 4-2: Loss of Signal on Channel 5


Illustration 4-3: Loss of Signal on Channel 3



Illustration 4-4: Loss of Signal on Channel 3

© 2010 General Electric Company.

GE Medical Systems, a General Electric Company, going to market as GE Healthcare.

www.gehealthcare.com