SERVICE MANUAL

Invivo Corporation
3545 SW 47th Avenue
Gainesville, FL 32608

Document TR0638S

Revision 1

3.0T GEM SMALL ANTERIOR ARRAY

GE Part Number: M7000GY

Invivo Part Number: 9896-031-12161

Manufactured By:

© Invivo Corporation, 2011

Distributed By:

GE Healthcare

Damage in Transportation

All packages should be closely examined at time of delivery. If damage is apparent, have notation “**damage in shipment**” written on **all** copies of the freight or express bill **before** delivery is accepted or “signed for” by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately** upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14** day period.

File a report with

- Call 1–800–548–3366 and use option 8.
- Fill out a report on <http://egems.med.ge.com/edq/home.jsp>
- Contact your local service coordinator for more information on this process.

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TABLE OF CONTENTS

SECTION 1 – INTRODUCTION.....	9
1.1 Product Identification and Shipping List	9
1.2 Shipping Contents	10
1.3 Compatibility	10
1.4 Related Documentation.....	10
1.5 Environmental Requirements.....	10
SECTION 2 – INSTALLATION PROCEDURE.....	11
2.1 Installation Finalization	12
SECTION 3 – SMALL ANTERIOR ARRAY MCQA TEST	11
3.1 Tools Required.....	12
3.2 MCQA Setup	12
3.3 MCQA Test.....	15
3.4 Periodic Quality Assurance Check	17
SECTION 4 – TROUBLESHOOTING	17
4.1 Mechanical Hardware Check	17
4.2 Electrical Hardware check	17
4.3 Coil ID Troubleshooting Tips.....	17
SECTION 5 – MAINTENANCE.....	18
5.1 Coil Care	18
SECTION 6 – RENEWAL PARTS.....	18
6.1 Field Replaceable Units	18
6.2 Customer Consumable Components.....	18
SECTION 7 – UPGARADING TO GEM SMALL ANTERIOR ARRAY	19
7.1 How to Upgrade.....	19

ATTENTION!

A GE field engineer **MUST** perform installation of this coil. Please contact your field engineer to arrange for installation.

Attempts to perform installation of this coil without proper system configuration will result in system errors.

SECTION 1 – INTRODUCTION

1.1 Product Identification and Shipping List

This is a service manual for the 3.0T GEM Small Anterior Array. This coil is identified by using the label below:


Illustration 1: Anterior Array Label

Illustration 2 shows the 3.0T GEM Small Anterior Array, Phantom Positioner, Anterior Pad, and Strap Set.


Illustration 2

1.2 Shipping Contents

SHIPPING LIST BOX 1– TABLE 1-2-1

Description	Invivo Part #	GE Part #	Qty
3.0T 32Ch Cardiac Coil - Anterior	4535-303-12291	5419534-2	1
Phantom Positioner, GEM Small Anterior Array (2 pieces)	4535-303-03571	5406404-11	1
Anterior Pad, GEM Small Anterior Array	4535-302-60711	5406404-10	1
Strap Set, GEM Small Anterior Array	4535-303-03591	5406404-12	
Operator Manual	4535-303-05151	N/A	1

1.3 Compatibility

This coil is compatible with the following hardware configurations:

Discovery MR750w Non-GEM*, and Discovery MR750w GEM 3.0T Systems

*In order to use the GEM Small Anterior Array on a MR750w Non-GEM system, the curved phantom positioner part 5366627-4 must be purchased and billed to the customer to run MCQA testing.

1.4 Related Documentation

- 5670000 Discovery MR750w and MR750w GEM 3.0T Service Methods DVD
- 1.5T & 3.0T GEM Small Anterior Array Operator's Manual, 4535-302-58991

1.5 Environmental Requirements

Storage

Store the coil in an air-conditioned scan room or equipment room. To store the coil, a storage space of 48.3 cm/19.0 in (width) x 16.5 cm/6.5 in (depth) x 101.6 cm/40.0 (length) is required.

Dimensions

Coil	48.3 cm x 16.5 cm x 101.6 cm	(19.0 in x 6.5 in x 18.0 in)
Phantom Positioner	50.6 cm x 50.6 cm x 17.7 cm	(20 in x 20 in x 7 in)

Weight

Coil	2.95 Kg	(6.5 lb)
Phantom Positioner	0.68 Kg	(1.5 lb)

SECTION 2: INSTALLATION PROCEDURE

There are no installation steps other than plugging in the coil and performing functional checks (MCQA test). Coil configurations will already be on the system software. The coil required configurations to use this coil are:

COIL CONFIGURATIONS– TABLE 2-1

System	Clinical Coil Name	Annotation Name	Description
MR750w-GEM	GE_32ChCardAntOnly	Small AA only 30	GEM Small AA
	Body 30 Small	Body 30 Small	GEM Chest/Body/Pelvis or GEM Whole Body
MR750w	GE_32Ch Cardiac Ant	32Ch Cardiac Ant	32Ch Cardiac Array by Invivo

These coil configurations can be confirmed present on the system in the “Coil Information” tab of Coil DB Explorer (more information on coil DB Explorer can be found on MR Service Methods under Coils -> Running Coil Database Explorer).

2.1 Installation Finalization

Complete Product Locator information for the coil. Failure to do this may result in failure of your site to receive future FMI's. Instructions are below.

- o Go to GIB: http://gib.gehealthcare.com/gib/gib_MainMenu.jsp
- o After the log in screen, click on [Global Installed Base entries]
- o Under “Serial Changes to Installed Base”, enter the system ID of the system on which this coil is being installed. Click [Submit]
- o Click on [New Entries]

SECTION 3: GEM SMALL ANTERIOR ARRAY MCQA TEST

Follow this procedure to perform the MCQA test. In order to run MCQA for this coil on MR750w Non-GEM, the curved phantom positioner part 5366627-4 must be purchased. For detail, please see the 32Ch Cardiac service manual for MR750w Non-GEM.

3.1 Tools Required


TOOLS REQUIRED - TABLE 3-1

Description	GE Part #	Qty
Unified TL Phantom Brick	5343347	1
GEM Small AA – Phantom Positioner (2 pcs)	5406404-11 (MR750w GEM flat table)	1

3.2 MCQA Setup

All images shown are on the GEM flat table.

TABLE 3-2 Steps for Coil and phantom setup

Visual	Description
	1. Remove any other surface coils (if present) from the cradle.
 Illustration 3  Illustration 4	2. Remove any pads present on the coil. The cable must point towards the magnet and plug into ports P1 or P2. Place lower half of phantom positioner on table. Ensure that the positioner is centered in the R/L direction. Insert coil and place top part of the positioner on the coil as shown in Illustration 3 and Illustration 4.


Illustration 5

3. Confirm the message on in-room monitor as shown in Illustration 5.

Use the IntelliTouch strip as shown in Illustration 6.

Illustration 6

Confirm that the message on the in-room monitor has disappeared as shown in Illustration 7.

Push here to use IntelliTouch strip

Illustration 7


Illustration 8


Illustration 9: MCQA SETUP WITH UNIFIED TL BRICK

4. Turn alignment button on.

Move table into bore until landmark lasers intersect with crosshairs on the coil as shown in Illustration 8.

Press 'Landmark' button

Place phantom on the coil. Ensure that the phantom is contained within the phantom positioner as shown in Illustration 9.

Press 'Advance to Scan' button.

3.3 MCQA Test

- 3.3.1 At the operator workstation, end the previous exam (if required).
- 3.3.2 From the Common Service Desktop, select the [Image Quality] tab, choose "Multi-Coil QA Tool", and then "Click here to start this tool" as shown in Illustration 10. Select 3.0T_DVw_32Ch_Cardiac_by_Invivo_Anterior_only.Unified.mctcfg as shown in Illustration 11.


Illustration 10: MCQA from the Service Desktop


Illustration 11

3.3.3 Further instructions on using the MCQA tool can be found on the Service Methods CD at: Troubleshooting -> System -> Multi-Coil Quality Assurance Tool.

For troubleshooting: The relationship between MCQA image order, coil Receive Pin and coil element are the same (IMG 1 = Rx Pin 1 = Element 1). Illustration 12 below shows the element layout.


Illustration 12: Element layout for GEM Small Anterior Array

3.4 Periodic Quality Assurance Check

QA checks should be performed during installation and during planned maintenance. Perform the quality assurance checks as outlined below to ensure the coil is operating properly.

3.4.1 Check external cable for cracks or cuts.

3.4.2 Perform Section 3: 3.0T GEM SMALL ANTERIOR ARRAY MCQA TEST

SECTION 4: TROUBLESHOOTING

Note: The Invivo Anterior Array is not to be opened for any reason!

4.1 Mechanical Hardware Check


- Verify that there are no loops in the cables.
- Verify that there are no metal or ferromagnetic objects close to the coil, patient or magnet (i.e., safety pin, hair pin).
- Verify that the coil is properly positioned.

4.2 Electrical Hardware check

Visually inspect all coaxial connectors and pins on the system plug. If there are any broken, deformed or recessed pins or coaxial connectors, replace the coil component.

Coil Connector Illustration 13

NOTE: The center pins of the System Connector *receive coaxial connectors* are extremely fragile [Connector Pin Out Illustration , R1-R8]. **Do not attempt to make any measurements at these center pins.**

For GND measurements, use any of the receive coaxial pins' outer shields (R1-R8).

RESISTANCE MEASUREMENTS – TABLE 4-2-1

System Connector		Reading
Positive Lead	Negative Lead	
1	2	< 8Ω
3	4	< 8Ω
5	GND	Open
6	GND	Open
20	GND	140 ±40Ω
21	GND	< 5Ω
22	GND	140 ±40Ω
23	GND	< 8Ω

If one or more of the readings are out of specification, replace the coil. If all of the above readings are correct, redo Section 3: GEM SMALL ANTERIOR ARRAY MCQA test. If the problem still exists, replace the 3.0T Anterior Array.

4.3 Coil ID Troubleshooting Tips

If the following error is displayed, there may be a problem with the Coil ID coil or system hardware:
The prescribed RF coil and detected RF coil do not match.

Troubleshooting instructions are available through the Service Methods DVD via the path:
Troubleshooting -> Diagnostics -> Hardware Location -> Magnet Room -> Coil ID Functional Diagnostic

The diagnostic will display the contents of the coil ID chip whenever the connector is inserted. Follow the troubleshooting instructions if the correct ID is not displayed. When the troubleshooting is complete, and if the correct ID is still not displayed, replace the coil.

Troubleshooting example:

Coil ID diagnostic started

Reading Coil ID's...

Port 1 (P1) Connected

Port 3 (P3) Connected

Port 3 (P3): 1 Coil Detected:

--> Coil Name: GEM PA

Coil is VALID

Coil ID: GEMsiggjNFhn3pJ6xvsAqJoBQvFU0qb

Serial Number: 5c000000117e4e0f

Field Strength: 30000

Body Transmit: ENABLED

Port 1 (P1): 1 Coil Detected:

--> Coil Name: GEM Small AA

Coil is VALID

Coil ID: INVlppAbXy7sgvQeY29ZizFcq18H2Cb7

Serial Number: 430000001189dd0f

Field Strength: 30000

Body Transmit: ENABLED

SECTION 5 – MAINTENANCE

5.1 Coil Care

WARNING!

Detach coil connector from scanner before attempting to clean. Do not reattach after cleaning until coil has dried completely. Having the coil attached to the system during cleaning or when it is wet may result in electrical shock.

CAUTION

Do not spray or pour cleaning solution directly on coil. Do not submerge coil in solution. The coil contains sensitive electronics components that could be damaged by the solution.

The 3.0T 32-Channel Cardiac Coil, posterior pads must be cleaned and stored using the following procedure:

- Wipe with a cloth that has been dampened in a solution of 10% bleach and 90% tap water, or 30% isopropyl alcohol and 70% tap water.
- **Do not pour any cleaning solution directly on the coil!**
- Let the coil housing and pads dry before use.
- **Under no circumstances should the coil be placed into any type of sterilizer.**

SECTION 6 – RENEWAL PARTS

6.1 Field Replaceable Units

FIELD REPLACEABLE UNITS LIST – TABLE 6-1

Description	GE Part #	Invivo #
3.0T 32Ch Cardiac Coil - Anterior	5419534-2	4535-303-14471
GEM Small Anterior Array Phantom Positioner (2 pieces)	5406404-11	4535-303-03601

6.2 Customer Consumable Components

CUSTOMER CONSUMABLE COMPONENTS – TABLE 6-2

Description	GE Part #	Invivo #
Anterior Pad, 32ch Cardiac coil anterior only	5406404-10	4535-303-03611
Strap Set, 32ch Cardiac coil anterior only	5406404-12	4535-303-03621

Note: Pads and straps are designed as customer consumable components. They are not intended to be covered by any GE service contract. If pads and straps are replaced they should be billed to the customer.

SECTION 7 – UPGRADING TO GEM SMALL ANTERIOR ARRAY

7.1 How to upgrade

It is possible to upgrade a 3.0T 32Ch Cardiac Coil for MR750w Non-GEM only to a 3.0T GEM Small Anterior Array. To do so, order the following parts and bill the customer accordingly.

(Note: 32Ch Cardiac for MR750 cannot be upgraded to 3.0T GEM Small Anterior Array)

Description	GE Part #	Invivo #
GEM Small Anterior Array Phantom Positioner (2 pieces)	5406404-11	4535-303-03601
Strap Set, 32ch Cardiac coil anterior only	5406404-12	4535-303-03621

Revision History

Rev	Date	Author	Primary Reason for Change	DCN
1	8-16-11	R. Law	Issued	10049