

MI&CT ESERVICE NOTE
SN07012409a

Class C information is Property of GE.

This short-term informational note is invalid 6 months after the issue date.

Date: January 23, 2007

Subject: Time Jump on Reboot

Effectivity: VCT and LS 5.X Systems running 06MW29.7
Discovery VCT Systems running 06MW34.6

Details/Issue:

New sites with the 06MW29.7 software have reported time changes upon reboot, and operators have noticed time changes on startup and reboot.

Cause:

The adjtime file is used to adjust for system time loss. If the adjtime file is corrupted, the system time could jump after each reboot. The adjtime file corruption typically occurs when the setdate command is run before using reconfig to change the time zone.

Resolution:

If you...

- Have already installed a new system running 06MW29.7 or 06MW34.6 and required a time zone change or a have noticed a time change on startup or reboot, perform the steps in Section A.
- Are installing a new system running 06MW29.7 or 06MW34.6 and will require a time zone change, perform the steps in Section B.
Do not run setdate command before changing time zone.
- Are upgrading from 06MW03.4 or 06MW03.5 to 06MW29.7, follow the 06MW29.7 LFC instructions to avoid this issue.

Section A – To eliminate large time changes on reboot...

1. Reconfig to set the appropriate timezone
2. Reboot the system
3. Stop applications from coming up
4. Become root
 {ctuser@hostname} **su -**
 Password **#bigguy**
5. Remove the adjtime file
 {root@hostname} **rm /etc/adjtime**
6. Recreate an empty adjtime file
 {root@hostname} **touch /etc/adjtime**
7. Setdate to apply the correct time
8. Reboot the system

Section B – To prevent large time changes on reboot, perform these steps during Console Install

1. Reconfig to set the appropriate timezone
2. Reboot the system
3. Stop applications from coming up

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.

Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

4. Setdate to apply the correct time
5. Reboot the system

Affected Service Publication: (New Installation Manuals sent with 06MW29.7 SW)
LightSpeed 7.X Installation Manual (Book 2 of 2) Direction 5116411-100, Revision 09
LightSpeed 5.X Installation Manual (Book 2 of 2) Direction 2349182-100, Revision 14
BrightSpeed Elite, Edge, Excel Installation Manual (Book 2 of 2) Direction 5141176-100, Revision 4

Note: Discovery VCT Installation Manual Direction 5177652-100, Revision 1 is not affected.

PQR/SPR #: 13093409 (parent, FCTge25800), 13093780, 13094877, 13095035, 13095086, 13095087, 13095362, and FCTge26030.

Claire Elliott

MI&CT Service Methods Designer

MI&CT Engineering Manager Signature _____