

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.
Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

**MI&CT ESERVICE NOTE
SN07051613a**

Class C information is Property of GE.

This short-term informational note is invalid 6 months after the issue date.

Date: 05/16/2007

Subject: IMS Position Error 260134802 at Table Home position due to Firmware Fault.

Effectivity: VCT 64, VCT32 and RT16

Details:

The following IMS error may happen when the table is moved back to home position with heavy patient loaded.

```
Fri Oct 6 15:43:23 2006
Host : Table Sub-System      Ermes # : 260134802
Exception Class : Secondary   Severity : Pri/Soft
File : gtcb_movechk.cxx      Line# : 584
Function : No System Function Reported
Scan Type : None/Unknown Scan: 2485/2/0
```

```
IMS position is out of tolerance.
current position(um) : 298245065
target position (um) : 0
start position (um) : 298245163
velocity (um/s): 100000
time from start (ms) : 300
```

The error will not be cleared by hardware reset. To Recover the Technician needs to move the IMS toward gantry manually by hand or service switch on GTCB.

This problem is caused by firmware fault on VCT M3 software (05MW31.6). It was already resolved by M4 (06MW03.4) or later software.
But if the table characterize was done when the system was 05MW31.6 software, this problem may happen.

Resolution:

Re-characterize the IMS. See Service manual for the IMS Characterization Procedure.
The newly re-characterized IMS values will resolve the IMS Position Errors at Table Home position.

Affected Service Publication: N/A

PQR/SPR #: 13088711, 13101567

Takashi Uehara
Lead Service Integrator, GEYMS MI&CT Engineering

MI&CT Engineering Manager Signature _____ Mark Nyman _____.