

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.
Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

**MI&CT ESERVICE NOTE
SN07061215a**

Class C information is Property of GE.

This short-term informational note is invalid 6 months after the issue date.

Date: June 13, 2007

Subject: Missing Files for CT Belt Tension Check and Adjustment Procedure

Effectivity: Lightspeed systems with VCT Gantry - VCT, VCT XT, VCT Select, VCT Pro32, Discovery VCT

Details: The current procedure (CT Belt Tension Check and Adjustment) requires two files that were misplaced and unavailable to the Field.

1. BeltTension.dat (Must be renamed to BeltTension.exe)
2. Belt_Tension_Log.xls

The procedure indicates that these two files can be found in the "Downloadable Files" Directory but does not state where that specific directory is located.

Resolution: The procedure should have pointed the Field to the MI&CT Service Engineering website.

Temporary solution:

Acquire the required support files for the existing procedure as follows:

1. Log into the GE Intranet anytime prior to performing the CT belt Tension Check and Adjustment procedure.
2. Point your browser to the MI&CT Service Engineering website at:
http://www.ct.med.ge.com/tech_pubs/fct/system/fct_home.htm Or (Cut & Paste this address into your browser.
3. Scroll down on the left side of the page and locate the "[Downloads](#)" link.
4. Select "BeltTension_Ver_1.zip" and save it to your Laptop Desktop.
5. Extract the zip file and the two missing files needed for the existing procedure to your desktop.

Long Term Solution:

A completely new procedure is being developed for use in adjustment for the VCT Gantry drive belt and also the GT table that uses frequency measurements exclusively. This procedure will be made available as soon as parts it depends on become available to the Field globally and will be announce through a future Service Note. It will also become available upon the next release of the Service Methods DVD.

Affected Service Publication:

5193575-800 - LS 7.X Limited Access Service Info (Adv)
5193575-200 - LS 7.X Limited Access Service Info (Adv)
5116424-800 - LightSpeed 7.X Service Information (Adv)
5116424-200 - LightSpeed 7.X Service Information (Adv)
5178786-800 - Discovery VCT PET-LS7.X Service Methods (Adv)
5178786-200 - Discovery VCT PET-LS7.X Service Methods (Adv)

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.

Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

PQR/SPR #: 13075397, 13078737, 13079316, 13080798, 13078510

Mark Tercha

MI&CT Engineering

MI&CT Service Engineering Manager Reviewer: _____ Mark Nyman _____