

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.

Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

**MI&CT ESERVICE NOTE
SN07091114a**

Class C information is Property of GE.

This short-term informational note is invalid 6 months after the issue date.

Date: September 11, 2007

Subject: Use of Adapter part when replacing Anode Inverter tank

Effectivity: 1.x through 5.x (Non-HP60) Lightspeed CT and PET systems

Details: Part number 2171866 supplied with HV Tank (Anode)



Illustration 1

The anode Tank replacement procedures do not explain what these adapter bushings are used for. These adapters are sent with the replacement tank to accommodate the tanks use for both Highspeed and Lightspeed systems. The Adapters are only used when attaching the newer style Anode tank to the older Highspeed gantry due to the smaller bolt size found on the Highspeed gantry.

Resolution:

The same anode tank is used on the Lightspeed systems that use larger tank mounting bolts (3/8 inch). When the tank is installed on the systems indicated below, use of these Adapters **is not required and they should be discarded**.

- o All 1.x, 2.x, 3.x, 4.x, 5.x CT systems
- o All PET systems with 3.x, 4.x and 5.x CT systems within them.

The affected Service procedures will be updated to include this information upon their next release.

Note: These adapters (Part number 2171866) are not orderable.

Affected Service Publications:

- **CT 1.x (QX/I)** Document 2282386-100 Rev 16, Lightspeed 1.3 System Service Manual
 - o Anode Transformer Tank (Replacement)
- **CT 2.x Plus** Document 2243315-100 Rev 22, Lightspeed 2.X System Service Manual
 - o High Voltage Tank (Anode) Replacement
- **CT 3.x Lightspeed** Document 5114043 Rev 6, LightSpeed 3.X Service Information (Advanced)
 - o High Voltage Tank (Anode) , HP Anode / Cathode Inverter replacement
- **CT Lightspeed 5.X**
- 5193891-800, Rev 2, LightSpeed 5.X Pro16 Limited Access Service Methods – Adv
- 5196839 - LS RT16 and Xtra Ltd Access Service Info (Adv), Rev 3
- 5174162 - LightSpeed RT16 and Xtra Service Methods, Rev 2
- 5113142 - LightSpeed 5.X RT Service Methods, Adv, Rev 5

PQR/SPR #: 13101336

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE Healthcare is prohibited.
Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

Douglas R. Hofstetter
MI&CT Service Engineering
MI&CT Engineering Manager Reviewer: _____