

**MI&CT ESERVICE NOTE
SN07091912a**

Class C information is Property of GE.

This short-term informational note is invalid 6 months after the issue date.

Date: September 19, 2007

Subject: Protocols used in DMPR option are missing – any 07xWxx.x software
The root cause is that the protocols are multiple step and not single step.

Effectivity: Any CT/PET 07xWxx.x software with the DMPR option.

Details:

- Customers upgrading to any 07xWxx.x software that currently have multiple step protocols will no longer see them available in the list of reformat protocols available for DMPR.
- Any protocol built in Reformat for the intended use in DMPR must be a single step protocol. If the Next button is selected in a reformat protocol without any further prescription, the protocol is still marked a multiple step and will not be shown for use in DMPR.

Resolution:

- New single step protocols for use in DMPR will have to be built. If the same protocol name is to be used, delete the multiple step protocol, exit Reformat and then go back into Reformat and build the protocol and name it.
- A review of user protocols will need to be performed to make sure the reformat protocols are linked to their protocols using DMPR.
- Protocols intended for DMPR must be built in Reformat or if Volume Viewer is installed in the GENERAL category.

Affected Service Publication:

Any sites using 07MxWxx.x software installed with the DMPR option.

- 07MW11.10
- 07MW18.4
- 07BW08.3
- 07BW08.4
- 07MW11.4C
- 07MWHL36.4
- 07MWDVCT36.4

PQR/SPR #:

- none - Field issues reported to Karen Procknow.

Resa Lambert

CT Operator Console Service Methods – Service Engineering

MI&CT Engineering Manager Reviewer: _____ Mark Nyman _____