

MI&CT eSERVICE NOTE

SN07102512a

Class C information is Property of GE.

This short-term informational note is invalid 6 months after the issue date.

Date: October 25, 2007

Subject: GOC5 JBOD FRU Updates

Effectivity: All VCT, Pro32 and DVCT Systems

Details:

Due to obsolesce of the current Hard Disk Drive (HDD) and changes required by RoHS compliance for the Power Supply, two of the exiting FRUs are being replaced with newer models. This Service Note will explain the details and special consideration required when replacing these parts.

JBOD HDD: 5114536-10

The present Hard Disk Drive used in the scan disk array has gone End of Life (EOL). The Supplier cannot acquire these drives and have switched to using a new model. Global Parts and Repair Services will be supplying the current HDD until stock runs out, at which point the new HDD will be supplied. After Verification and Validation testing by Console Subsystem Engineering, no special considerations have been identified concerning the introduction of the new HDD. This includes the mixing of the current and new drives in the scan disk array.

JBOD PSU: 5114536-3

The Supplier of the JBOD PSU has made changes to the current PSU to make this assembly RoHS compliant. (RoHS is the initiative to reduce or eliminated the use of lead and other toxic chemicals in the construction of electrical assemblies.) Global Parts and Repair Services will be supplying the current PSU until stock runs out, at which point the new PSU will be supplied. Due to design changes with the new RoHS compliant PSU, special consideration has been identified concerning the mixing of the old verses new PSUs. The JBOD assembly utilizes two PSU in a redundant configuration. The new PSU cannot be installed into a JBOD chassis along with the current PSU model. Differences in the two PSUs electrical properties will cause reliability issues if the two PSU models are mixed in the redundant configuration. Consequently if a PSU is ordered and the current PSU model is no longer available, two of the new model PSU will need to be ordered.

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved.
 Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare
 Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE
 Healthcare is prohibited.
 Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

Resolution:

HDD:

	<u>Current HDD</u>	<u>New HDD</u>
GEHC Part #	5114536-10	5114536-11
Part Name	Hard Drive w/Mounting Sled	Hard Drive w/Mounting Sled
Vendor	Seagate ST37345LC FW002	Seagate ST37345LC FW001

Replacement Considerations

- None. Follow existing service documentation for replacing drives in scan disk array (JBOD).
- FE Assist Comments will be updated with part change-up information.

JBOD PSU:

Current PSU



New PSU



GEHC Part #	5114536-3	5114536-5
Part Name	Power Supply Module	Power Supply Module
Vendor	Bell Micro - nStor NSR-PWR	Bell Micro - nStor NSR-PWRP-GE

Replacement Considerations

- If system has the current 5114536-3 PSUs and there are 5114536-3 available, order one (1) 5114536-3 PSU.
- If system has the current 5114536-3 PSUs and there are no 5114536-3 available, order two (2) 5114536-5 PSUs.
- If system has the new 5114636-5 PSUs, order one (1) 5114536-5 PSU.
- Follow existing service documentation for replacing PSU in scan disk array (JBOD).
- FE Assist Comments will be updated with part change-up information.

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved.
Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare
Systems. Disclosure to persons not employed by GE Healthcare or not in possession of a Advanced Service Agreement with GE
Healthcare is prohibited.
Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

Affected Service Publication:

LightSpeed 7X Service Manuals

- GOC5 Illustrated Parts List – 5135498-100 (CEM-371322)
- Parts Replacement, Console, JBOD – (CEM-360989)

PQR/SPR #:

- No Quality Reports
- ECR/ECO 2049162

Michael L. Bauer
Console Subsystem SMDE
MI&CT Service Engineering