

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare. Disclosure to persons not employed by GE Healthcare or not in possession of an Advanced Service Agreement with GE Healthcare is prohibited. Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

**MI&CT eSERVICE NOTE
SN08060516a**

Class C information. Property of GE Healthcare

This short-term informational note is invalid 6 months after the issue date.

Date: June 5, 2008

Subject: Starting Console Diagnostics on 06MWxx.x Software

Effectivity: Systems running 06MWxx.x software (Includes PET_VCT_HP.80)

Details:

The script to run Console Diagnostics is different between the 06MWxx.x and 07MWxx.x software.

Resolution:

Use **ConsoleDiag** on 06MWxx.x software (Includes PET_VCT_HP.80).

Use **ConsoleDiags** on 07MWxx.x software.

Affected Service Publication:

LS 7.X Limited Access Service Methods – Adv, Direction 5193575

LS 5.X Pro16 Limited Access Service Info – Adv, Direction 5193891

LS 5.X RT Limited Access Service Information – Adv, Direction 5193855

LS 5.X (HP60) Limited Access Service Methods – Adv, Direction 5224328

LS 4.X Service Information – Adv, Direction 5114045

LS 3.X Service Information – Adv, Direction 5114043

Discovery VCT PET – LS 7.X Service Methods – Adv, Direction 5178786

Discovery ST, STE, RX HP60 Service Methods – Adv, Direction 5141123

Tracking Number (PQR/SPR/ECR):

PQR 13117118

Claire Elliott
MI&CT Service Engineering