

GE Healthcare MI&CT Class C eService Notes are copyright (c) protected to General Electric Company Inc with all rights reserved. Class C eService Notes issued by GE Healthcare are proprietary to GE Healthcare. Disclosure to persons not employed by GE Healthcare or not in possession of an Advanced Service Agreement with GE Healthcare is prohibited. Duplication and distribution of any class of eService Notes, in any form by entities other than GE Healthcare is strictly prohibited.

**MI&CT eSERVICE NOTE  
SN08101614a**

**Class C information. Property of GE Healthcare**

**Date:** October 16, 2008

**Subject:** Disconnection of the CT Table Foot Pedal Cables on GT and HP Tables Used in a [Surgical Environment].

**Effectivity:** All CT systems defined in the list below using the GT and HP Tables installed for use in a [Surgical environment]:

- Global Tables (GT) Tables- CT HD 750 Discovery, Lightspeed VCT/VCT32, RT16, RT Xtra
- H-Power (HP) Tables - Lightspeed Pro16, Lightspeed 5.X, and RT 4.

**Details:** Inadvertent use of the foot pedal option may cause a GT Table or HP Table to move unexpectedly.



**GT Table Style**



**HP Table Style**

**Resolution:** All GT and HP Tables installed in a [Surgical Environment] must have their foot pedals permanently disabled. Vertical motion may still be accomplished through the use of the push buttons found on the Gantry Control Panels. To proceed with the disconnection:

1. Inform your Customer prior to disconnection.
2. Escalate the issue through a CSO (Use CEIP).
3. Obtain and follow Engineering Instruction 5328450INS to disconnect the pedals.

**Affected Service Publication:**

LS 7.X VCT Installation Manual 5116411-100, LS 5.X Installation Manual 2349182-100  
LS 5.X Plus Installation Manual 5117716-100, LS 5.X (China) Installation Manual 5117716-100  
LS 5.X RT16 and RT Xtra Installation Manual 5154590-100

**Tracking Number (Complaint # or SPR #):** 13185731

Adam Nathan  
MICT Service Engineering