GE Healthcare

Direction 24239893TPH 1.5T HNS Coil Installation

Class A

Direction 24239893TPH
Rev 3
© 2008, General Electric Company,
All Rights Reserved

Revision History

Date	Rev	Change Description
5-Sept-08	0	Updated to meet compliance
27-Oct-08	1	Initial Release
19-March-09	2	• Updated Table of Contents with the “Scope” section • Updated “Damage in Transportation” statement • Extended the scope of installation to HDx and Discovery MR450 & MR750 • Updated the “Coil Config” section with a note related to pop-up box specific to Discovery MR systems • Removed the HDx designator from Section 5 • Update Table 5.3.1 with phantom part number • Added Section 6 - HNS Coil Setup for Reverse Ramped System
9-Sept-09	3	• Updated Language Policy per DOC0620154 • Added information per Phantom Elimination Program

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ATTENTION!

Installation of this coil MUST be performed by a GE field engineer. Please contact your field engineer to arrange for installation.

Attempts to perform installation of this coil without proper system configuration will result in system errors.

Table of Contents

1- DAMAGE IN TRANSPORTATION	11
2- Scope	12
3- COIL HARDWARE CONFIGURATIONS	13
4- COIL CONFIG.....	16
4-1 PERSONNEL REQUIREMENTS.....	16
4-2 OVERVIEW.....	16
4-3 PRELIMINARY REQUIREMENTS	16
4-3-1 Tools and Test Equipment	16
4-3-2 Safety.....	16
4-3-3 Required Conditions	16
4-4 PROCEDURE.....	16
5- HNS Coil MCQA Setup.....	21
5-1 PERSONNEL REQUIREMENTS.....	21
5-2 OVERVIEW	21
5-3 PRELIMINARY REQUIREMENTS	21
5-3-1 Tools and Test Equipment Utilizing Legacy Phantoms	21
5-3-2 Tools and Test Equipment Utilizing Unified Phantom Kit	21
5-3-3 Safety.....	22
5-3-4 Required Conditions	22
5-4 Procedure	22
6- HNS Coil Setup for Reverse Ramped System.....	24
6-1 Personnel Requirements	24
6-2 Overview.....	24
6-3 Tools and Test Equipment	24
6-4 Procedure	24

1- DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent write "Damage In Shipment" on ALL copies of the freight or express bill BEFORE delivery is accepted or "signed for" by a General Electric representative or hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14-day period.

To file a report:

- Call 1-800-548-3366 and use option 8.
- Fill out a report on <http://egems.med.ge.com/edq/home.jsp>
- Contact your local service coordinator for more information on this process.

Rev. September 17, 2006

2- SCOPE

This coil installation service manual is applicable to the HDx/HDxt Upgrade System, HDx/HDxt Forward Production System, Discovery MR450, and Optima MR450w.

The HNS Coil has 5 rating plates so there are 4 sets of Product Locator cards included with the shipment.

The Global Install Base Database (also known as Product Locator System) tracks shipment, trans-shipment, and field location of the serialized models. There are now two methods for submitting Product Locator information.

1. At this time, for U.S. ONLY, the preferred method of submitting information is the *FE Site Verification Web Site* at:

<http://gib.gehealthcare.com>

The FE Site Verification consists of three components that are available on the web from the main menu. They are:

- _ Install/deinstall product locator model and serial numbers
- _ Add/modify ship to address information
- _ Update CARES FE data for primary/secondary FEs

To obtain a copy of the FE training tutorial for using this website, a downloadable copy is available at the following:

Product Locator Support Central Page – <http://supportcentral.ge.com/15563>

Note: One “Shipping Card” is filled out and submitted when shipped (extra cards are supplied for trans-shipments between storage and distribution points), and the “Installation Card” and extra shipment cards are attached.

2. Verify the serial and model number on each rating plate matches the installation card numbers before removing the installation cards.

Note: There may be one or more shipment cards and bar code labels with the installation card. These shipment cards are used to trace the transfer of serialized units between various inventory storage and distribution points until the product reaches its final installation destination. Process just the installation card and discard any extra shipment cards and labels.

3. Enter the information from the installation cards into the GIB database for this customer.

3- COIL HARDWARE CONFIGURATIONS

The HNS coil has 5 components (or called units). The components are listed below:

- Posterior Head Component (or called Posterior Head Neck Unit)
- Face Component (or called Anterior Head Neck Unit)
- TL Component (or called TL unit)
- Horseshoe Component (or called adaptor block)
- Chest Component (or called Neck Chest unit)

The HNS coil may be configured in 9 different ways using its modular design of its components. The configurations are listed below:

- Posterior Head and Face
- Posterior Head, Face and TL
- Posterior Head and Horseshoe
- Posterior Head, Horseshoe and TL
- Posterior Head, Face and Chest
- Posterior Head, Face, Chest and TL
- Posterior Head, Horseshoe and Chest
- Posterior Head, Horseshoe, Chest and TL
- TL only

Each configuration of the HNS coil must be installed for the system to recognize it. However, we just need to use 5 configurations to setup the coil to get all coil configurations installed. Connecting all the coil components together may not result in the coil getting completely configured. Illustrations 1 through 5 show the valid coil component configurations during installation.

Illustration 1 – Setup 1; Posterior Head Component + Face Component


Illustration 2 – Setup 2; Posterior Head Component + Face Component + TL Component


Illustration 3 – Setup 3; Posterior Head Component + Horseshoe Component + TL Component


Illustration 4 – Setup 4; Chest Component Added To Setup 3


Illustration 5 – Setup 5; Posterior Head Component + Face Component + TL Component + Neck Component

4- COIL CONFIG

4-1 PERSONNEL REQUIREMENTS

Required Persons	Preliminary Reqs	Procedure	Finalization
1	Timing Not Applicable	30 Min	Timing Not Applicable

4-2 OVERVIEW

The HNS Coil shall be installed with the following procedure. All 5 coil setups that listed in Section 3 shall be performed in order to get all configurations installed.

4-3 PRELIMINARY REQUIREMENTS

4-3-1 Tools and Test Equipment

Item	Qty	Effectivity	Part #	Manufacturer
None				

4-3-2 Safety

No special safety requirement

4-3-3 Required Conditions

None

4-4 PROCEDURE


1. Go to the scan desktop and select the [New Patient] button.
2. Refer to Section 3 and assemble the coil in one of the 5 valid hardware setup configurations and make sure the coil components have a snug fit with each other.
3. Plug the connector, from the Head Neck Unit (HNU) of the coil, into the appropriate Port connector on the LPCA Port.

Note: The appropriate coil and connector labels can be seen on the coil.

Note: For Upgrade Systems, the system will automatically detect the new coil and a pop-up message box will appear stating the coil has not yet been configured for this system. See Illustration 6.

Note: For HDx Forward Production Systems, Discovery MR450 Systems, and Optima MR450w Systems, if the coil is not configured on the system a pop-up box will appear. See Illustration 6. If the box does not appear, then proceed to step 5.

Illustration 6 – New Coil Detected Pop-up

4. Select **[Yes]** to configure the new coil. If the process is successful a pop-up message will be displayed as shown in Illustration 7. If the configuration process is not successful then call the OnLine Center for coil configuration support. **Do not attempt to configure the coil manually.**

Illustration 7 – Configuration Successful Pop-up

5. Prescribe a new exam and enter patient information.
6. Click on the coil selection button and verify that the HNS coil is present under currently connected on the scanner user interface. For example, if you had used Setup 1 you would see the clinical coil modes as shown in Illustration 8.

Illustration 8a - Verification that the head coil is connected and that it is recognized by the system. Illustration 8a applies to Discovery MR450 and Optima MR450w Systems.



Illustration 8b - Verification that the head coil is connected and that it is recognized by the system.

Note: In forward production systems, the presence of coil entries like those shown in Illustration 8 is a confirmation that the coil installed properly. If you do not see the HNS coil modes, **do not attempt to install the coil manually**, call the OnLine center for coil configuration support.

Note: After all 5 of the coil Setups have been configured, you will see all the coil entries in the Coil Names window as shown in illustration 9.

Illustration 9 - HNS coil modes present under currently connected

Coil Selection

Coil Type :

- ☒ Currently Connected
- ☐ Head/Neck
- ☐ Upper Extremity
- ☐ Chest/Torso/Pelvis
- ☐ Spine
- ☐ Lower Extremity/Pv Run-Off
- ☐ Other
- ☐ Body

Coil
Body Coil by GE
HNS by GE (PosteriorHead, Face, Ch
HNS by GE (PosteriorHead, Face, Ch
HNS by GE (TL)

Configurations	Extents
HNS CTLBOT	
HNS CTL45	
HNS CTL34	
HNS CTL345	
HNS CTL56	
HNS CTL456	

Accept Cancel

7. After successfully configuring the coil proceed to Section 5, HNS Coil MCQA Setup.

5- HNS COIL MCQA SETUP

5-1 PERSONNEL REQUIREMENTS

Required Persons	Preliminary Reqs	Procedure	Finalization
1	Timing Not Applicable	30 Min	Timing Not Applicable

5-2 OVERVIEW

To ensure modularity in design and in diagnostics the coils must be tested in the following 3 configurations:

1. The TL unit
2. Posterior Head Chest Horseshoe TL
3. Posterior Head Face

The MCQA tool can ONLY run on the above mentioned individual coil configurations and NOT all coils configurations at the same time. MCQA will not run and will report that the configuration is invalid if an attempt is made to scan with an unsupported combination of coils.

5-3 PRELIMINARY REQUIREMENTS

5-3-1 Tools and Test Equipment Utilizing Legacy Phantoms

Item	Qty	Effectivity	Part #	Manufacturer
Head Phantom	1	Required	2424360 (Part of phantom kit 2418078 or 2418236)	Dielectric
Neck Phantom	1	Required	2422888 (Part of phantom kit 2418078 or 2418236)	Dielectric
HNU Phantom Positioner	1	Required	24220112 (Part of phantom kit 2418078 or 2418236)	Dielectric
TL Phantom	2	Required	2424361 (Part of phantom kit 2418078 or 2424359)	Dielectric

5-3-2 Tools and Test Equipment Utilizing Unified Phantom Kit

Item	Qty	Effectivity	Part #	Manufacturer
TL Unified Phantom	2	Required	5343347	Dielectric
Large Cylindrical Unified Phantom	1	Required	5342679	Dielectric
Phantom Positioner, HNS	1	Required	5344844	Medical Comfort

5-3-3 Safety

No special safety requirement.

5-3-4 Required Conditions

The appropriate coil configuration file for the HNS coil must be installed to run this tool. Refer to Section 4 – Coil Config.

5-4 Procedure

1. Record the coil serial #s for the Posterior Head, Chest and TL components in Table 1

Table 1 – Coil Component Serial Numbers

Posterior Head Component  A white, U-shaped medical device with a blue cable attached to the side.	Serial # _____
Chest Component  A yellow, U-shaped medical device with a blue cable attached to the top.	Serial # _____
TL Component  A white, rectangular medical device with a blue cable attached to the side.	Serial # _____

2. Perform MCQA setup following the steps in Chapter 2 of 1.5T Head Neck Spine (HNS Service Methods, Direction 2423983TPH.
3. Perform the Surface coil MCQA Test. Click here to go the MCQA test. [Multi-Coil Quality Assurance Tool](#)

Test instructions can be found in the following paths on the Service Methods DVD:

- Signa HDxt: Troubleshooting > Image Quality > Multi-Coil Quality Assurance Tool
- MR450: Troubleshooting > System > Multi-Coil Quality Assurance Tool

6- HNS COIL SETUP FOR REVERSE RAMPED SYSTEM

6-1 Personnel Requirements

Required Persons	Preliminary Requirements	Procedure	Finalization
1	N/A	15 mins	N/A

6-2 Overview

Procedure to make the coil compatible with reverse ramped magnets for HNS coil is provided.

6-3 Tools and Test Equipment


Item	Qty	Effectivity	Part#	Manufacture
Standard Tool Kit-Flathead screwdriver, 10in-lbs torque driver, 15in-lbs torque driver	1	N/A	N/A	N/A

6-4 Procedure


NOTICE

Observe anti-static precautions: Use wrist strap and ESD workstation (ESD mat connected to ground and place to connect ESD wrist strap).

Illustration 10 - Posterior Unit Interface board cable swap (HNS Rev 3)

The Rev 3 HNS coil can be made compatible with reverse ramped magnets for any system. To make the coil compatible the following steps are to be performed;

1. Put an ESD strap on and ground.
2. Unscrew the posterior cover and remove from coil. Retain the screws.
3. Swap cables E1-E2. Disconnect cables labeled E1 and E2. Connect the cable labeled E2 to location E1 on Interface sub-assembly. Similarly, connect the cable labeled E1 to location E2 on Interface sub-assembly.
4. Swap cables E3-E4. Disconnect cables labeled E3 and E4. Connect the cable labeled E3 to location E4 on Interface sub-assembly. Similarly, connect the cable labeled E4 to location E3 on Interface sub-assembly.
5. Swap cables E5-E6. Disconnect cables labeled E5 and E6. Connect the cable labeled E5 to location E6 on Interface sub-assembly. Similarly, connect the cable labeled E6 to location E5 on Interface sub-assembly.

Illustration 6-2 Posterior Unit- Interface board setup post-cable swap (HNS Rev 2)

The Rev 2 HNS coil can be made compatible with reverse ramped magnets for any system except DVMR system. To make the coil compatible the following steps are to be performed;

1. Put an ESD strap on and ground.
2. Unscrew the posterior cover and remove from coil. Retain the screws.
3. Swap cables E1-E2. Disconnect cables labeled E1 and E2. Connect the cable labeled E2 to location E2 on Interface sub-assembly. Similarly, connect the cable labeled E1 to location E1 on Interface sub-assembly.
4. Swap cables E3-E4. Disconnect cables labeled E3 and E4. Connect the cable labeled E3 to location E3 on Interface sub-assembly. Similarly, connect the cable labeled E4 to location E4 on Interface sub-assembly.
5. Swap cables E5-E6. Disconnect cables labeled E5 and E6. Connect the cable labeled E5 to location E5 on Interface sub-assembly. Similarly, connect the cable labeled E6 to location E6 on Interface sub-assembly.

www.gehealthcare.com

Imagination at work