

ACTION	MRC	MRSC
Change on the second tabcard (Language) the settings like in the picture below English – US.  Fig. 80:	☐	☐
First press OK, next Apply. Hint: Pressing “OK” in the “Regional and Language Options” Tab Card starts an automatic reboot, wait two steps with this action.	☐	☐
Open the Windows Explorer (It could be necessary to switch to “Advanced User” to open an Explorer)	☐	☐
Delete the file C:\MedCom\config\Mri\MRSystem\MrVCheck.ref	☐	☐
Press the OK Button and accept all Messages, which will be shown with “YES” or “OK”. A Reboot follows!	☐	☐