

Prerequisites for Scheduled Performance Verification of Chillers For Signa *TwinSpeed*

Revision 02/14/02

The following list indicates items that must be completed before contacting Ellis and Watts to schedule the Factory Authorized Performance Verification of the pertinent *Twin Speed* Cooling Cabinet. When this step is completed, then the PM schedule will be established for the year. The LTL 10 (M3088TE) is started up and serviced by GE, therefore, performance verification by an Ellis and Watts Service Provider is not required.

1. Appropriate line power is connected to the Chiller and correct phasing established.
2. Chiller is permanently mounted as required by GE and / or state / local code requirements.
3. **For units with Water Cooled condensers;** (LC 18MT/WC, See Ellis and Watts Tech. Pub. 471)
 - Piping has been flushed out to remove any debris before being attached to Chiller
 - Water supply to Chiller condenser has an inline filter or strainer.
 - Full Gradient and Shield Cooler (Usually this happens about a week before turnover to customer.)
4. **For Indoor Air-Cooled Condenser units;** (LC 20MT, See Ellis and Watts Tech. Pub. 472)
(LC 10MT, See Ellis and Watts Tech. Pub. 473)
 - Ambient temperature and humidity are controlled to meet the equipment room environmental requirements.
 - Full Gradient and Shield Cooler (For Dual Loop) or Full Gradient (For Single Loop) loading is available.(Usually this happens about a week before turnover to customer.)
5. **For Outdoor Condenser unit;** (LC 20MT/O See Ellis and Watts Tech. Pub. 470)
(LC 10MT/O See Ellis and Watts Tech Pub 474.)
 - All coolant piping connections are complete and leak tested (**Do not use tap water**).

NOTE: It is recommend that distilled water or nitrogen gas be used to pressurize and leak test. Full Gradient and Shield Cooler (For Dual Loop) or Full Gradient (For Single Loop) loading is available.(Usually this happens about a week before turnover to customer.)

Performance Verification is usually completed just before the system is turned over to the Customer. Ellis and Watts service department needs a two week notification to schedule the Local Service Provider. Please provide the name and address of the Imaging Center, FE contact name with the best method of contact, requested date of service, and Cryo System ID #. Contact information is listed below.

Ellis and Watts International Inc. is determined to be as helpful as possible with the installation of the chiller equipment that we provide to support the General Electric Medical Systems equipment, but, Ellis and Watts and our Approved Service Providers are not contracted to perform **any** of the installation of this equipment. If you desire for E&W or any of our ASP's to assist with the actual installation of the chiller equipment, then, you will need to contact the Ellis and Watts Service Department at 513-752-9000 x749 or E-mail at service@elliswatts.com to discuss scheduling and payment for this work, or start-up scheduling. Information is free so call any time to discuss questions about installation or operation of our chiller equipment.

Ellis and Watts International Inc.
4400 Glenwillow Lake Lane
Batavia, Oh. 45103