

NOTE:

Before using this document please make sure that it is the latest revision.

This cover sheet controls the revision status of this entire document.

The revision of all sheets of this document will be at the same level and approvals shall be performed electronically.

Revision #	Name	Date	Description	PCN#
A	Resa Lambert	10Oct2003	Design Release	
B	Resa Lambert	14Nov2003	Added IG1 ping	
C	Resa Lambert	14Nov2003	Updates: 3-1, 4-1, 4-5, 4-10, 4-16	
0				
1				
2				

Title	Preliminary Training or Troubleshooting for GRE Console			Rev
	First Made For: GRE Program			C.
	GE Medical Systems Milwaukee, WI	000000TST	Sheet 1 of 39	


TABLE OF CONTENTS

1. SCOPE.....	3
2. REQUIREMENTS.....	4
3. GENERAL INFORMATION	5
3-1 GRE — Information Checklist.....	6
4. GRE Subsystem Testing.....	9
4-1 DARC Node — Cable Checkout.....	9
4-2 DARC Node — Ping and Remote Shell	11
4-3 DARC Node — Serial Over LAN Checkout.....	12
4-4 DARC Node — Telnet Connection Failed	13
4-5 DARC Node — Reset.....	14
4-6 DARC Node — ifconfig eth3 Plus.....	21
4-7 DARC Node — DIP Diagnostic Help	22
4-8 DARC Node — DIP Diagnostic All	23
4-9 DARC Node — DIP Diagnostic Individual Testing.....	24
4-10 IG Node — Cable Checkout.....	29
4-11 IG Node — Ping and Remote Shell	29
4-12 IG Node — SOL connection from DARC Node to IG1.....	30
4-13 IG Node — Remote Shell to IG1 and RAC Board Information.....	31
4-14 IG Node — RAC Board Diagnostics	32
4-15 IG Node — racbist Diagnostics.....	33
4-16 SDDA — Cable Checkout	35
4-17 SDDA — View Error Message LOG.....	35
4-18 DARC Node — SDDA Utilities.....	35
4-19 SDDA — testHighSpeed	36
4-20 SDDA — reconfigScanDisk.....	38

Title		Preliminary Training or Troubleshooting for GRE Console			Rev C.
		First Made For: GRE Program			
	GE Medical Systems Milwaukee, WI	000000TST			
		Sheet	2	of	39

1. SCOPE

This engineering preliminary training procedure is used to perform various troubleshooting techniques on the various subsystems of the GRE system.



Title

Preliminary Training or Troubleshooting for GRE Console

First Made For:

GRE Program

Rev
C.

GE Medical Systems
Milwaukee, WI

000000TST

Sheet 3 of 39


2. REQUIREMENTS

Verify and record specific system hardware configuration.

- Open a shell. Type: **cat /usr/g/config/INFO <Enter>**
Record information on the GRE Information Checklist provided.

Verify the proper revision / type of Software has been loaded onto the GRE Console successfully or the DARC is having problems with the software load.

Ethernet Connections to the rear of the GRE Operator Console are installed and good.



Title	Preliminary Training or Troubleshooting for GRE Console First Made For: GRE Program	Rev C.
	GE Medical Systems Milwaukee, WI	000000TST Sheet 4 of 39


3. GENERAL INFORMATION

GRE	Global Recon Engine
DARC	Data Acquisition and Recon Control Node
SDDA	Scan Data Disk Arra
IG1	Image Generator #1 Node — the GRE Console is always shipped with #1
IG2	Image Generator #2 Node — future (listed as fail during boot-up if unpopulated)
IG3	Image Generator #3 Node — future (listed as fail during boot-up if unpopulated)
RAC	Recon Accelerator Card located in the IG Node
DIP	DAS Interface Processor located in the DARC Node — specific to DAS type
SOL	Serial Over LAN
rsh	remote shell
halt	Console is taken fully down with the only option being power OFF.

Title	Preliminary Training or Troubleshooting for GRE Console			Rev
	First Made For: GRE Program			C.
	GE Medical Systems Milwaukee, WI	000000TST	Sheet 5 of	39

Title		Preliminary Training or Troubleshooting for GRE Console			Rev
		First Made For: GRE Program			C.
	GE Medical Systems Milwaukee, WI	0000000TST Sheet 6 of 39			

setenv PRINTER_HOSTNAME	_____
setenv PRINTER_IPADDRESS	_____
setenv PRINTER no	
setenv O2_PRESENT n	
setenv GENTRY_TYPE h2 / hp/ wb	_____
setenv LANGUAGE en_US	
setenv DOSE y	
setenv NETWORK yes	
setenv USENIS yes	
setenv CHGDARCSUBNET yes	
setenv DOMAINNAME	_____
setenv NIS_SERVER	_____
setenv DARC_SUBNET	_____
setenv NEXT_PATIENT_EXAM_NUMBER 1	
setenv MOBILESYSTEM no	_____
setenv STREAK y	
setenv INROOMSCAN off	
setenv AAA n	
setenv AUTOMA_INDEX_TABLE_TYPE type2	_____
setenv RECONHARDWARE GRE / No O2 Pegasus	_____
setenv DASTYPE SDAS or MDAS 4 / 8/ 16	_____
setenv DETECTORTYPE Warp3	
setenv REFRESH_RATE 76	
setenv yp s	
setenv ypmaster \${SUITEID}_OC0	
setenv HOSTNAME SBC0	
setenv ALIAS \${SUITEID}_SBC0	
setenv UNSYSID \$GATEWAY_HOSTNAME	
setenv SERVID \$GATEWAY_HOSTNAME	
setenv tubeType \$TUBETYPE	
setenv ACRNEMA xxxx	
setenv REGEN no	
setenv HRCTTITLE	_____
setenv HRPORTRNUM	_____
setenv HRAETITLE	_____
setenv HRIPADDR	_____
{ctuser@hosatname}	




Title		Preliminary Training or Troubleshooting for GRE Console		Rev
		First Made For: GRE Program		C.
	GE Medical Systems Milwaukee, WI	000000TST		
		Sheet	8	of 39

4. GRE Subsystem Testing

This engineering preliminary procedure is used to perform a series of tests or scripts to verify the DARC Node, IG Node and SDDA (plus some internal assemblies) are operating properly. Other tests may be available but were not outlined in this document. Most of the information seen in the Unix Shells has been copied from actual GRE Consoles and not input via the computer keyboard by a user.

4-1 DARC Node — Cable Checkout

- 1. Verify the rear power cable is connected, switch is turned on and fan is operating.
- 2. Verify front power LED is illuminated.



- 3. Verify the Fiber Optic cable is connected to the DIP Board RX or RECV port (shown above – left most port) and it is in good condition.
- 4. Verify the 9-pin Sub D Scan XRAY Abort Cable is connected between the DIP Board and the SDDA (towards the right).
- 5. Verify the RHARD Interlock cable is connected to the 9-pin Sub D and directly below to the SDDA.
- 6. Verify the SCSI cable is connected to the SDDA and not kinked or in poor condition.

7. Verify Ethernet cable is in good condition and connected as follows:
- IG1 Ethernet 1 (bottom-left port) — DARC Node Ethernet 0 (top-left port)




8. Verify Ethernet cable is in good condition and connected as follows:
- HP (right-card – top port) Ethernet C — DARC Node Ethernet 3 (bottom-left port)
 -

(DARC Node eth3 shown here – top connection)



(HP Host Computer Ethernet C shown here – see label)



Title	Preliminary Training or Troubleshooting for GRE Console First Made For: GRE Program	Rev C.
	GE Medical Systems Milwaukee, WI	000000TST Sheet 10 of 39

4-2 DARC Node — Ping and Remote Shell

This subsection will verify the DARC Node can be communicated too or pinged. After verifying the communication line exists between the HP Host PC and the DARC Node a remote shell will be invoked to determine if the DARC Node is up.

Open a Unix Shell and type the following:

- {ctuser@hostname} **su - <Enter>**
- Password: **#bigguy <Enter>**
- [root@hostname]# **ping darc <Enter>**
- *PING darc (172.16.0.2) from 172.16.0.1 : 56(84) bytes of data.*
- *64 bytes from darc (172.16.0.2): icmp_seq=1 ttl=64 time=0.157 ms*
- *64 bytes from darc (172.16.0.2): icmp_seq=2 ttl=64 time=0.209 ms*
- *64 bytes from darc (172.16.0.2): icmp_seq=3 ttl=64 time=0.175 ms*
- *64 bytes from darc (172.16.0.2): icmp_seq=4 ttl=64 time=0.134 ms*
- **Select: Control-C to stop the ping**
- *--- darc ping statistics ---*
- *4 packets transmitted, 4 received, 0% loss, time 3000ms*
- *rtt min/avg/max/mdev = 0.134/0.168/0.209/0.031 ms*

- [root@hostname]# **rsh darc <Enter>**
- *Last login: Thu Oct 9 11:17:35 on ttyS1*
- *You have new mail.*
- [root@darc ~]# **exit <Enter>**
- *logout*
- *rlogin: connection closed.*
- [root@hostname]#

Title		Preliminary Training or Troubleshooting for GRE Console			Rev C.
		First Made For: GRE Program			
	GE Medical Systems Milwaukee, WI	000000TST			
		Sheet	11	of	39

4-3 DARC Node — Serial Over LAN Checkout

Verification of the Serial Over LAN check is the almost always the initial test performed on the DARC after power and cable verification. It has been placed after ping and remote shell checkout because they are very quick tests. With SOL Checkout the telnet connection is made and the SOL connection from the HP Host Computer to the DARC Node is then established. This test will verify the DARC Node was booted and is up.

Open a Unix Shell and type the following:

```
{ctuser@hostname} su - <Enter>
```

```
Password: #bigguy <Enter>
```

```
[root@hostname]# telnet localhost 623 <Enter>
```

```
Trying 127.0.0.1...
```

```
Connected to localhost.
```

*** Indicates HP Host Apps were loaded ***

*** Indicates cliservice proxy server is running ***

```
Escape character is '^['.
```

```
Server: darc <Enter>
```

```
Username: <Enter>
```

```
Password: <Enter>
```

```
Login successful
```

*** Indicates the SOL connection is established ***

```
dpccli> console <Enter>
```

```
GEMS Linux 1.7.9 *** Console Operating System Software Revision ***
```

```
See you on the darc side of the moon.
```

*** DARC Node is up ***

```
darc login: root <Enter>
```

*** see note ***

```
Password: #bigguy <Enter>
```

```
Last login: Thu Oct 9 11:11:15 on ttyS1
```

```
You have new mail.
```

```
[root@darc ~]# ~.
```

(tilde period to terminate telnet)

```
dpccli> exit <Enter>
```

```
Connection closed by foreign host.
```

```
[root@hostname]#
```

NOTE

If the darc login prompt is not present and localhost.localdomain login prompt is displayed then try reloading the DARC Node OS and DARC Application software per the procedure. Problem could be the DARC Apps did not load well. Always be sure the DARC OS matches the HP Host OS already loaded.

Title	Preliminary Training or Troubleshooting for GRE Console	Rev C.
	GE Medical Systems Milwaukee, WI	First Made For: GRE Program 000000TST Sheet 12 of 39

4-4 DARC Node — Telnet Connection Failed

If the telnet connection fails this subsection will guide you to the next steps in identifying the root cause (s).

Open a Unix Shell and type the following:

```
{ctuser@hostname} su - <Enter>
```

```
Password: #bigguy <Enter>
```

```
[root@hostname]# telnet localhost 623 <Enter>
```

```
Trying 127.0.0.1...
```

```
Failed to connect to localhost.
```

```
[root@hostname]# /etc/rc.d/init.d/cliservice <Enter>
```

```
[root@hostname]# telnet localhost 623 <Enter>
```

```
Trying 127.0.0.1...
```

```
Connected to localhost.
```

```
Escape character is '^['.
```

```
Server:
```

Discussion:

If the telnet connection is successful then perform SOL Checkout again.

If the telnet connection fails:

- Open a Unix Shell and halt the system
- Power Off the GRE Operator Console
- Wait 30 seconds.
- Power On the GRE Operator Console
- Retry the telnet command

If it still fails

- load the SOL CD if available for current version of Software.
- check or replace eth3 connection (DARC Node to HP Host ethA).
- reload the software package
- replace the HP Host PC or DARC Node as suggested by the OLC.

Title	Preliminary Training or Troubleshooting for GRE Console	Rev C.
	GE Medical Systems Milwaukee, WI	First Made For: GRE Program 000000TST Sheet 13 of 39

4-5 DARC Node — Reset

The DARC Node reset –c test should only be performed if the dpccli> console command failed to return the *OS revision information* and *See you on the darc side of the moon* message. Whenever this test is performed the user must verify the information on the screen as it scrolls past looking for failures.

NOTE

During the DARC Node reboot or reset a [Failure] may appear which is acceptable. One type of acceptable failure would be 'Stopping NIS' which may fail because NIS was never running. Record all failures or errors and contact OLC for further assistance if necessary.

Open a Unix Shell and type the following:

- {ctuser@hostname} su - <Enter>
- Password: #bigguy <Enter>

[root@hostname]# telnet localhost 623 <Enter>

Trying 127.0.0.1...

Connected to localhost.

Escape character is '^['.

Server: darc <Enter>

Username: <Enter>

Password: <Enter>

Login successful *** Indicates the SOL connection is established ***

dpccli> reset -c <Enter>

The DARC Node will be reset. Below is a partial display screen displayed during the reset – this is loaded with **OS 1.7.9** version (each OS version will be different):

Linux Software loading

Linux version 2.4.18-24SGI_XFS_1.2.0smp

(root@bob.ctt.med.ge.com) (gcc version 2.96 20000731 (Red Hat Linux 7.3 2.96-110)) #1 SMP Thu Feb 20 20:52:49 CST 2003

BIOS-provided physical RAM map:

BIOS-e820: 0000000000000000 - 000000000009b400 (usable)

BIOS-e820: 000000000009b400 - 00000000000a0000 (reserved)

BIOS-e820: 00000000000e8000 - 0000000000100000 (reserved)

Title Preliminary Training or Troubleshooting for GRE Console		First Made For: GRE Program		Rev C.
	GE Medical Systems Milwaukee, WI	0000000TST	Sheet 14 of	39

Bluser: 0000000000000000 - 000000000009b400 (usable)

user: 000000000009b400 - 00000000000a0000 P was found.

ACPI: RSDP located at physical address c00ff9b0

RSD PTR v0 [INTEL]

__va_rangeIC (acpi_id[0x0000] id[0x0] enabled[1])

CPU 0 (0x0000) enabledProcessor #0 Pentium 4(tm) XEON(400)
global_irq_base[0x30])INT_SRC_OVR (bus[0] irq[0x0] global_irq[0x2] polarity[0x0] triggerprocessor
Specification v1.4 command: ~.

Virtual Wire compatibility mode.

OEM ID: INTEL Product ID:5539.46 BogoMIPS

Memory: 1930604k/1966080k available (1694k kernel code, 30848k reserved, 1144k data, 188k init,
1048576k highmem)

Dentry cache hash table entries: 262144 (order: 9, 2097152 bytes)

Inode cache hash table entries: 131072 (order: 8, 1048576 bytes)

Mount cache hash ta done.

Checking 'hlt' instruction... OK.

POSIX conformance testing by UNIFIX

mtrr: v1.40 (20processor 1/1 eip 2000

Initializing CPU#1

masked ExtINT on CPU#1

ESR value before enabling v: Intel(R) Xeon(TM) CPU 2.80GHz stepping 07

Booting processor 2/6 eip 2000

Initializing CPU#25573.87 BogoMIPS

CPU: L1 I cache: 0K, L1 D cache: 8K

CPU: L2 cache: 512K

Setting 10 in the phys_id_present_map

...changing IO-APIC physical APIC ID to 10 ... ok.

Title		Preliminary Training or Troubleshooting for GRE Console		Rev
		First Made For: GRE Program		C.
	GE Medical Systems Milwaukee, WI	000000TST		
		Sheet	15	of 39

```

..T..... CPU clock speed is 2790.1331 MHz.

..... host bus clock speed is 132.1737 MHz.

cpu: 0, c88,D:8,S:51906,C:259533>

CPU1<T0:259520,T1:155696,D:12,S:51906,C:259533>

checking TSC synchro migration_task 1 on cpu=1

PCI: PCI BIOS revision 2.10 entry at 0xfdb85, last bus=4

PCI: Using) -> 16

PCI->APIC IRQ transform: (B0,I29,P1) -> 19

PCI->APIC IRQ transform: (B0,I31,P0) -> 17isapnp: No Plug & Play device found

speakup: initialized device: /dev/synth, node (MAJOR 10, MINOR 25)

Linux NET4.0 for Linux 2.4

Based upon Swansea University Computer Society NET3.039

Initializing RT netlink socket

apm: BIOS not found.

Starting kswapd

allocated 64 pages and 64 bhs reserved for the highmem bounces

VFS: Disk quotas v2_6.5.1

SGI XFS 1.2.0 with ACLs Unix98 ptys configured

Serial driver version 5.05c (2001-07-08) with MANY_PORTS MULTIPOINT SHARE_IRQ
SERIAL_PCI ISAPNP enabled

ttyS0 at 0x03f8 (irq = 4) is a 16550A

ttyS1 at 0x02f8 (irq = 3) is a 16550A

Real Time Clock Driver v1.10e

oprofile: can't get RTC I/O Ports

block: 1024 slots per queue, batch=256

Uniform Multi-Platform E-IDE driver Revision: 6.31

```

Title	Preliminary Training or Troubleshooting for GRE Console First Made For: GRE Program	Rev C.
	GE Medical Systems Milwaukee, WI	000000TST Sheet 16 of 39

Title		Preliminary Training or Troubleshooting for GRE Console			Rev	
		First Made For:			GRE Program	
		0000000TST			C.	
	GE Medical Systems Milwaukee, WI	Sheet 17 of 39				

Setting clock (localtime): Thu Oct 9 11:07:11 CDT 2003 [OK]

Activating swap partitions: [OK]

Setting hostname darc: [OK]

Mounting USB filesystem: [OK]

Initializing USB controller (usb-uhci): [OK]

Your system appears to have shut down uncleanly

Press Y within 1 seconds to force file system integrity check...

Checking root filesystem

[/sbin/fsck.xfs (1) -- /] fsck.xfs -a /dev/hda2

[OK]

Remounting root filesystem in read-write mode: [OK]

Finding module dependencies: [OK]

Checking filesystems

Checking all file systems.

[/sbin/fsck.xfs (1) -- /boot] fsck.xfs -a /dev/hda1

[/sbin/fsck.xfs (1) -- /usr/g] fsck.xfs -a /dev/hda6

[OK]

Mounting local filesystems: [OK]

Enabling local filesystem quotas: [OK]

Enabling swap space: [OK]

INIT: Entering runlevel: 3

Entering non-interactive startup

Setting network parameters: [OK]

Bringing up loopback interface: [OK]

Bringing up interface eth0: [OK]

Title	Preliminary Training or Troubleshooting for GRE Console First Made For: GRE Program	Rev C.
	GE Medical Systems Milwaukee, WI	000000TST Sheet 18 of 39

```

Bringing up interface eth1: [ OK ]
Bringing up interface eth2: [ OK ]
Bringing up interface eth3: [ OK ]
Starting system logger: [ OK ]
Starting kernel logger: [ OK ]
Starting portmapper: [ OK ]
Starting NFS statd: [ OK ]
Loading system font: [ OK ]
Initializing random number generator: [ OK ]
Mounting NFS filesystems: [ OK ]
Mounting other filesystems: [ OK ]
Starting dip: dip: found dip type 0x3229 i=2
dip: memtop      =0x78000000
dip: high_memory =0xf8000000
dip: __pa(highmem)=0x38000000
dip: DMA bus address 0x78000000 ; high mem at 0x38000000
dip: DIP base addr 0xe0000000 ioremap addr 0xf8b42000 reg_base=0x8000000
dip: Detected MDAS 16-Slice DIP: 2326523 B
[ OK ]
Starting igpower1: [ OK ]
Starting igpower2: Starting ig2:[FAILED]
[FAILED]
Starting igpower3: Starting ig3:[FAILED]
[FAILED]
Starting xinetd: [ OK ]

```

Title	Preliminary Training or Troubleshooting for GRE Console First Made For: GRE Program	Rev C.
	GE Medical Systems Milwaukee, WI	000000TST Sheet 19 of 39

Starting NFS services: [OK]

Starting NFS quotas: [OK]

Starting NFS mountd: [OK]

Starting NFS daemon: [OK]

Starting bmcscrip: [OK]

Starting dhcpd: [OK]

Starting console mouse services: [OK]

Starting crond: [OK]

Starting anacron: [OK]

Starting atd: [OK]

Starting cliservice: [OK]

GEMS Linux 1.7.9 ***** Console Operating System Software Revision *****

See you on the darc side of the moon. ***** DARC Node is up *****

darc login: **root** <Enter>

Password: **#bigguy** <Enter>

Last login: Thu Oct 9 11:05:08 on ttyS1

You have new mail.

root@darc ~]# ~. **(tilde period to terminate telnet)**

dpccli> **exit** <Enter>

Connection closed by foreign host.

[root@hostname]#

Title	Preliminary Training or Troubleshooting for GRE Console	Rev C.
	GE Medical Systems Milwaukee, WI	First Made For: GRE Program 000000TST Sheet 20 of 39

4-6 DARC Node — ifconfig eth3 Plus

This subsection will confirm the address at eth2 and the DARC Node Subnet address. Specific sites may require the DARC Subnet to be changed due to a conflict with the 172 base address which is initially loaded onto the HP Host and DARC Node from the Load From Cold software installation process.

Open a Unix Shell and type the following:

- {ctuser@hostname} **ifconfig** <Enter>

eth2 inet addr: 169 254.0.2 Bcast: 169.254.0.255 Mast: 255.255.255.0

- {ctuser@hostname} **su -** <Enter>
- Password: #bigguy <Enter>
- [root@hostname]# **rsh darc** <Enter>
- *Last login: Thu Oct 9 11:17:35 on ttyS1*
- *You have new mail.*
- [root@darc ~]# **ifconfig eth3** <Enter>
- inet addr: 172 16.0.2 Bcast: 172.16.0.255 Mast: 255.255.255.0

NOTE

The eth3 inet addr and Bcast may be different for specific sites that have experienced a 172 conflict. These specific sites will have a different address. One example is shown below at the end.

- [root@darc ~]# **exit** <Enter>
- [root@hostname]#

DARC Subnet Address Example:

[root@darc ~]# **ifconfig eth3** <Enter>

inet addr: 169 254.0.2 Bcast: 169.254.0.255 Mast: 255.255.255.0

Title		Preliminary Training or Troubleshooting for GRE Console			Rev C.
		First Made For: GRE Program			
	GE Medical Systems Milwaukee, WI	000000TST			
		Sheet	21	of	39

4-7 DARC Node — DIP Diagnostic Help

The DIP Board is located in the DARC Node. There are two connections to the DIP Board. A Fiber Optic Cable connects to RX or RECV port (left-most connection as viewed from the rear of the GRE Console). A 9-pin SubD Scan XRAY Abort Cable is also present which connects to the SDDA.

DIP Diagnostics are run via the DARC Node. A help menu is available as shown below. DIP Diagnostics may be run individually or as a group of all. Additionally, a number can be entered for iterations (DipDiag -a -i 5 <Enter>).

Test # 7 and test #8 will fail without a loopback cable installed between RX/RECV and TX/TXMT.

Loopback Cable installation: The Loopback cable should be known good before use in troubleshooting. Disconnect the Fiber Optic cable coming from the Gantry at the rear Console bulkhead connector. Connect one end of the loopback cable to the the rear Console bulkhead connector. Connect the other end of the loopback cable to the TX port of the DARC Node DIP Board. If Tests 7 and 8 fail – bypass the internal Console cable by connecting the loopback cable between the RX and TX ports of the DIP Board if possible.

Open a Unix Shell and type the following:

```
{ctuser@hostname} rsh darc <Enter>
```

```
Last login: Tue Oct 7 14:58:34 from oc
```

```
[ctuser@darc ~]$ cd /usr/g/bin <Enter>
```

```
[ctuser@darc bin]$ ./DipDiag -h <Enter>
```

```
usage: DipDiag [options]
```

```
options:
```

- a, --all,
- x, --execute, 1,2,3,...
- n, --do-not-execute, 1,2,3,...
- i, --iterations, N
- q, --quiet,
- p, --stop-on-error,
- h, --help,

Test numbers are as follows

- 1 - X-ray Enable Relay
- 2 - Test FEC Bit
- 3 - Test EDDR Bit
- 4 - Test Magic Number Register
- 5 - Memory Test
- 6 - Interrupt Test
- 7 - External LoopBack Test
- 8 - DMA Test

Without options the test is run in interactive mode.

```
[ctuser@darc bin]$
```

Title		Preliminary Training or Troubleshooting for GRE Console		Rev
		First Made For: GRE Program		C.
	GE Medical Systems Milwaukee, WI	000000TST		
		Sheet	22	of 39

4-8 DARC Node — DIP Diagnostic All

The DIP Board is located in the DARC Node. There are two connections to the DIP Board. A Fiber Optic Cable connects to RX or RECV port (left-most connection as viewed from the rear of the GRE Console). A 9-pin SubD Scan XRAY Abort Cable is also present which connects to the SDDA.

DIP Diagnostics are run via the DARC Node. A help menu is available. DIP Diagnostics may be run individually or as a group of all. Additionally, a number can be entered for iterations (DipDiag -a -i 5 <Enter>).

Test # 7 and test #8 will fail without the use of a loopback cable.

```
{ctuser@hostname} rsh darc <Enter>
Last login: Tue Oct 7 14:58:34 from oc
[ctuser@darc ~]$ cd /usr/g/bin <Enter>
```

```
[ctuser@darc bin]$ ./DipDiag -a <Enter>
12:05:40 opening LINUX device /dev/dip
creating H16 dip interface
DIP PN : 2326523 B ***Specific DIP Board identified here ***
test 1 Passed X-ray Enable Relay
test 2 Passed Test FEC Bit
test 3 Passed Test EDDR Bit
test 4 Passed Test Magic Number Register
fill mem
write to dip took 4.714 sec
verifying memory
read from dip took 19.469 sec
test 5 Passed Memory Test
test 6 Passed Interrupt Test
error - failed memcmp of pattern
test 7 Failed External LoopBack Test
error - DIP detected an error 0x69
test 8 Failed DMA Test
1 loops 2 failures *** See NOTE ***
2 tests failed
```

NOTE

Test #7 and test #8 will fail without a loopback cable installed between RX and TX of the DIP Board.

Title Preliminary Training or Troubleshooting for GRE Console		First Made For: GRE Program		Rev C.
	GE Medical Systems Milwaukee, WI	000000TST	Sheet 23 of 39	

4-9 DARC Node — DIP Diagnostic Individual Testing

The DIP Board is located in the DARC Node. There are two connections to the DIP Board. A Fiber Optic Cable connects to RX or RECV port (left-most connection as viewed from the rear of the GRE Console). A 9-pin SubD Scan XRAY Abort Cable is also present which connects to the SDDA.

DIP Diagnostics are run via the DARC Node. A help menu is available. DIP Diagnostics may be run individually. Additionally, a number can be entered for iterations (DipDiag -6 -i 5 <Enter>).

Test # 7 and test #8 will fail without the use of a loopback cable.

In the test case shown below only tests 1 through 6 were selected because a loopback cable was not installed on the Console. If a loopback cable is available select 12345678.

Open a Unix Shell and type the following:

```
{ctuser@hostname} rsh darc <Enter>
```

```
Last login: Tue Oct 7 14:58:34 from oc
```

```
[ctuser@darc ~]$ cd /usr/g/bin <Enter>
```

```
[ctuser@darc bin]$ ./DipDiag <Enter>
```

```
+-----+
|                                     |
|                               WARNING |
|                               ===== |
|                                     |
| - THIS TOOL IS FOR MANUFACTURING PROCESSES ONLY. |
| - DO NOT USE THIS TOOL WHILE SCANNING.           |
| - APPLICATION PROCESSES SHOULD BE SHUTDOWN WHILE |
|   USING THIS TOOL.                               |
|                                     |
+-----+
```

Do you want to continue? (Y/N): y <Enter>

12:05:40 opening LINUX device /dev/dip

creating H16 dip interface

info - **'DIP Board Existence Test' was successful**

info - **DIP board revision is '2326523-1B'**

```
+-----+
| DIP Diagnostics Main Menu |
+-----+
```

1. X-ray Enable Relay

Title	Preliminary Training or Troubleshooting for GRE Console	Rev C.
	GE Medical Systems Milwaukee, WI	First Made For: GRE Program 000000TST Sheet 24 of 39

- 2. Test FEC Bit
- 3. Test EDDR Bit
- 4. Test Magic Number Register
- 5. Memory Test
- 6. Interrupt Test
- 7. External LoopBack Test
- 8. DMA Test
- 0. Exit

Please make your selection: **123456 <Enter>** ***** No loopback cable is installed *****
+info - **'X-ray Enable Test' was successful**

+-----+
| DIP Diagnostics Main Menu |
+-----+

- 1. X-ray Enable Relay
- 2. Test FEC Bit
- 3. Test EDDR Bit
- 4. Test Magic Number Register
- 5. Memory Test
- 6. Interrupt Test
- 7. External LoopBack Test
- 8. DMA Test
- 0. Exit

Please make your selection: info - **'FEC Bit Test' was successful**

+-----+
| DIP Diagnostics Main Menu |
+-----+

- 1. X-ray Enable Relay
- 2. Test FEC Bit
- 3. Test EDDR Bit
- 4. Test Magic Number Register
- 5. Memory Test
- 6. Interrupt Test

Title		Preliminary Training or Troubleshooting for GRE Console			Rev C.
		First Made For: GRE Program			
	GE Medical Systems Milwaukee, WI	000000TST			
		Sheet	25	of	39

- 7. External LoopBack Test
- 8. DMA Test
- 0. Exit

Please make your selection: info - **'EDDR Bit Test' was successful**

+-----+
| DIP Diagnostics Main Menu |
+-----+

- 1. X-ray Enable Relay
- 2. Test FEC Bit
- 3. Test EDDR Bit
- 4. Test Magic Number Register
- 5. Memory Test
- 6. Interrupt Test
- 7. External LoopBack Test
- 8. DMA Test
- 0. Exit

Please make your selection: info - **'Magic Number Register Test' was successful**

+-----+
| DIP Diagnostics Main Menu |
+-----+

- 1. X-ray Enable Relay
- 2. Test FEC Bit
- 3. Test EDDR Bit
- 4. Test Magic Number Register
- 5. Memory Test
- 6. Interrupt Test
- 7. External LoopBack Test
- 8. DMA Test
- 0. Exit

Please make your selection: fill mem
write to dip took 4.714 sec

Title		Preliminary Training or Troubleshooting for GRE Console		Rev
		First Made For: GRE Program		C.
	GE Medical Systems Milwaukee, WI	000000TST		
		Sheet	26	of 39

verifying memory
read from dip took 19.469 sec
info - **'Memory Test' was successful**

+-----+
| DIP Diagnostics Main Menu |
+-----+

- 1. X-ray Enable Relay
- 2. Test FEC Bit
- 3. Test EDDR Bit
- 4. Test Magic Number Register
- 5. Memory Test
- 6. Interrupt Test
- 7. External LoopBack Test
- 8. DMA Test
- 0. Exit

Please make your selection: info - **'Interrupt Test' was successful**

+-----+
| DIP Diagnostics Main Menu |
+-----+

- 1. X-ray Enable Relay
- 2. Test FEC Bit
- 3. Test EDDR Bit
- 4. Test Magic Number Register
- 5. Memory Test
- 6. Interrupt Test
- 7. External LoopBack Test
- 8. DMA Test
- 0. Exit

Please make your selection: **0 <Enter>**

Exiting DIP Diagnostics.

Title		Preliminary Training or Troubleshooting for GRE Console			Rev C.
		First Made For: GRE Program			
	GE Medical Systems Milwaukee, WI	000000TST			
		Sheet	27	of	39

[ctuser@darc bin]\$ **exit** <Enter>
logout
rlogin: connection closed.
{ctuser@hostname}

Title		Preliminary Training or Troubleshooting for GRE Console		Rev
		First Made For: GRE Program		C.
	GE Medical Systems Milwaukee, WI	000000TST		
		Sheet	28	of 39

4-10 IG Node — Cable Checkout

1. Verify the rear power cable is connected, switch is turned on and fan is operating.
 2. Verify front power LED is illuminated.
 3. Verify Ethernet cabling is not kinked and in good condition and connected as follows:
- IG1Ethernet 1 (bottom-left port) — DARC Node Ethernet 0 (top-left port)

4-11 IG Node — Ping and Remote Shell

This subsection will verify the IG Node can be communicated too or pinged via the DARC Node. After verifying the communication line exists between the IG Node and the DARC Node a remote shell will be invoked to determine if the IG Node is up. This can be performed as ctuser or as root.

Open a Unix Shell and type the following:

- {ctuser@hostname} **su - <Enter>**
- Password: **#bigguy <Enter>**
- [root@hostname]# **rsh darc <Enter>**
Last login: Thu Oct 9 11:17:35 from oc
You have new mail.
- [root@darc ~]# **ping ig1 <Enter>**
- *PING ig1 (10.0.1.2) from (10.0.1.1): 56(84) bytes of data.*
- *64 bytes from ig1 (10.0.1.2) icmp_seq=1 ttl=64 time=0.157 ms*
- *64 bytes from ig1 (10.0.1.2) icmp_seq=2 ttl=64 time=0.209 ms*
- *64 bytes from ig1 (10.0.1.2) icmp_seq=3 ttl=64 time=0.175 ms*
- *64 bytes from ig1 (10.0.1.2) icmp_seq=4 ttl=64 time=0.134 ms*
- **Select: Control-C to stop the ping**
- *--- ig1 ping statistics ---*
- *4 packets transmitted, 4 received, 0% loss, time 3000ms*
- *rtt min/avg/max/mdev = 0.134/0.168/0.209/0.031 ms*
- [root@darc ~]# **rsh ig1 <Enter>**
- Password: **#bigguy <Enter>**
- Last login: Thu Oct 16 13:04:55 from darc
- [root@ig1 ~]# **ping ig1 <Enter>**
- *PING ig1 (10.0.1.1) from (10.0.1.2): 56(84) bytes of data.*
- *64 bytes from ig1 (10.0.1.2) icmp_seq=1 ttl=64 time=0.157 ms*
- *64 bytes from ig1 (10.0.1.2) icmp_seq=2 ttl=64 time=0.209 ms*
- *64 bytes from ig1 (10.0.1.2) icmp_seq=3 ttl=64 time=0.175 ms*
- *64 bytes from ig1 (10.0.1.2) icmp_seq=4 ttl=64 time=0.134 ms*
- **Select: Control-C to stop the ping**

Title	Preliminary Training or Troubleshooting for GRE Console	Rev C.
	GE Medical Systems Milwaukee, WI	First Made For: GRE Program 000000TST Sheet 29 of 39

- *--- ig1 ping statistics ---*
- *4 packets transmitted, 4 received, 0% loss, time 3000ms*
- *rtt min/avg/max/mdev = 0.134/0.168/0.209/0.031 ms*
- *[root@ig1 ~]# exit <Enter>*
- *logout*
- *rlogin: connection closed.*
- *[root@darc ~]# exit <Enter>*
- [root@ig1 ~]# exit <Enter>*
- *logout*
- *rlogin: connection closed.*
- *[root@hostname]#*

4-12 IG Node — SOL connection from DARC Node to IG1

The DARC Node controls the IG Node. The DARC Node must be up.

Open a Unix Shell and type the following:

{ctuser@hostname} su - <Enter>

Password: #bigguy <Enter>

[root@hostname]# rsh darc<Enter>

[root@hostname]# telnet localhost 623 <Enter>

Trying 127.0.0.1...

Connected to localhost.

Escape character is '^'.

Server: ig1 <Enter>

Username: <Enter>

Password: <Enter>

Login successful ***** Indicates the SOL connection from DARC to IG1 is established *****

dpccli> exit <Enter>

Title	Preliminary Training or Troubleshooting for GRE Console	Rev C.
	GE Medical Systems Milwaukee, WI	First Made For: GRE Program 000000TST Sheet 30 of 39

4-13 IG Node — Remote Shell to IG1 and RAC Board Information

The RAC Board resides in the IG Node, however, the DARC Node controls the IG Node. The DARC Node must be up to view the flash revision (fan_para-6-12) and the device driver information.

Open a Unix Shell and type the following:

```
{ctuser@hostname} su - <Enter>
```

```
Password: #bigguy <Enter>
```

```
[root@hostname]# rsh darc<Enter>
```

```
Last login: Wed Oct 8 11:20:22 from oc
```

```
[ctuser@darc ~]$ ls /usr/local/gemsrac <Enter>
```

```
fan_para-6-12.hexout power_up_cfg-6-11.hexout
```

***** fan_para-6-12.hexout is the filename we will use to see the actual information *****

```
[ctuser@darc ~]$ rsh ig1 <Enter>
```

```
Last login: Wed Oct 8 11:20:51 from darc
```

For convenience - here is theFLASH help menu:

```
usage [options] flash_rac hexfile
```

```
options:
```

```
-d device    Specify the device name (default  
/dev/gemsrac)
```

```
-e           Erase FLASH before programming.
```

```
-n           Do not write to the FLASH.
```

```
-f           Force programming, do not check the image  
first.
```

```
-v           Verify the FLASH against the hex input file.  
and do nothing else.
```

```
-l n         Write, erase, and/or verify FPGA load n
```

```
-o filename  Save binary output to file.
```

```
-S n         Program only the first n sectors (for debug)
```

```
[ctuser@ig1 ~]$ flash_rac -v -l 1 fan_para-6-12.hexout <Enter>
```

```
Board ID: 2316826 C
```

```
done reading file 'fan_para-6-12.hexout'
```

```
1305558 bytes
```

```
Calculated checksum byte is 0x7e
```

```
image size is padded to 1310720
```

```
verify OK
```

```
[ctuser@ig1 ~]$ cat /proc/driver/gemsrac <Enter>
```

```
RAC Device Driver Version 1.2
```

```
compiled Oct 3 2003 05:14:04
```

```
Board Physical Identifier: 2316826 C
```

```
PCI Command Reg : 0x117
```

Title		Preliminary Training or Troubleshooting for GRE Console		Rev
		First Made For: GRE Program		C.
GE Medical Systems Milwaukee, WI		000000TST		
		Sheet	31 of	39

```

PCI Error Status: 0x0
PCI Rev      : 0xe          or (0x8 with the 0x0 below)
PCI Subsys   : 0x1          or (0x0 with the 0x8 above)
DMA bus address : 0x78000000
[ctuser@ig1 ~]$ exit <Enter>
logout
rlogin: connection closed.
[ctuser@darc ~]$ exit <Enter>
logout
rlogin: connection closed.
{ctuser@hostname}

```

4-14 IG Node — RAC Board Diagnostics

The RAC Board resides in the IG Node, however, the DARC Node controls the IG Node. The DARC Node must be up to view the flash revision (fan_para-6-12) and the device driver information.

NOTE

RAC diagnostic Test 3 SRAM B Incr Pat will always fail the first pass run and must pass subsequent passes.

Open a Unix Shell and type the following:

```

{ctuser@hostname} cd /usr/g/bin <Enter>
{ctuser@hostname} ls rac* <Enter>
racdiags*
{ctuser@hostname} racdiags <Enter>

```

9 (counts from 0 to 9)

```

0 SRAM A Incr Pat      PASSED 10 loops
1 SRAM A Walking Pat   PASSED 10 loops
2 SRAM A Address Pat   PASSED 10 loops
3 SRAM B Incr Pat      PASSED 10 loops    *** This will fail on the initial pass ONLY ***
4 SRAM B Walking Pat   PASSED 10 loops
5 SRAM B Address Pat   PASSED 10 loops
6 FLASH Region 0       PASSED 10 loops
7 FLASH Region 1       PASSED 10 loops
13 DMA host<-image mem PASSED 10 loops
14 DMA host->projection mem PASSED 10 loops

```

10 Loops: All tests passed

```
{ctuser@hostname}
```

Title	Preliminary Training or Troubleshooting for GRE Console First Made For: GRE Program	Rev C.
	GE Medical Systems Milwaukee, WI	000000TST Sheet 32 of 39

4-15 IG Node — racbist Diagnostics

The racbist diagnostics are available in command line form. These are the same as racdiags except that the racdiags have been pre-set to run a series of 2 tests initially and not report their fail status. The tests that may initially may fail are # 0 and # 3. For illustrative purposes we will show an actual test file output which shows the initial failure of test #3. If racbist diagnostics are to be utilized it is recommended that test #0 and test #3 are run initially and then further diagnostics are run.

Also, racdiags do not run the complete list of racbist tests.

Open a Unix Shell and type the following:

```
{ctuser@hostname} su - <Enter>
```

```
Password: #bigguy <Enter>
```

```
[root@hostname]# rsh darc <Enter>
```

```
Last login: Thu Oct 16 13:40:38 from oc
```

```
You have new mail.
```

```
[root@darc ~]# rsh ig1 <Enter>
```

```
Password: #bigguy <Enter>
```

```
Last login: Thu Oct 16 13:04:55 from darc
```

```
[root@ig1 ~]# racbist <Enter>
```

0 SRAM A Incr Pat	PASSED 1 loops
1 SRAM A Walking Pat	PASSED 1 loops
2 SRAM A Address Pat	PASSED 1 loops
3 SRAM B Incr Pat	FAILED 1 of 1
4 SRAM B Walking Pat	PASSED 1 loops
5 SRAM B Address Pat	PASSED 1 loops
6 FLASH Region 0	PASSED 1 loops
7 FLASH Region 1	PASSED 1 loops
8 FLASH Region 2	PASSED 1 loops
9 FLASH Region 3	PASSED 1 loops
10 FLASH Region 4	PASSED 1 loops
11 FLASH Region 5	PASSED 1 loops
12 DMA host->image mem	PASSED 1 loops
13 DMA host<-image mem	PASSED 1 loops
14 DMA host->projection mem	PASSED 1 loops
15 DMA host<-projection mem	PASSED 1 loops

Title Preliminary Training or Troubleshooting for GRE Console		First Made For: GRE Program		Rev C.
	GE Medical Systems Milwaukee, WI	000000TST	Sheet 33 of 39	

```
[root@ig1 ~]# racbist --ex 0,3 <Enter> (run this initially as these may fail)
0 SRAM A Incr Pat      PASSED 1 loops
3 SRAM B Incr Pat      PASSED 1 loops
1 Loops: All tests passed
```

```
[root@ig1 ~]# racbist --ex 0,1,2,3,4,5,6,7,13,14 --it 10 <Enter>
9
```

```
0 SRAM A Incr Pat      PASSED 10 loops
1 SRAM A Walking Pat   PASSED 10 loops
2 SRAM A Address Pat   PASSED 10 loops
3 SRAM B Incr Pat      PASSED 10 loops
4 SRAM B Walking Pat   PASSED 10 loops
5 SRAM B Address Pat   PASSED 10 loops
6 FLASH Region 0       PASSED 10 loops
7 FLASH Region 1       PASSED 10 loops
13 DMA host<-image mem PASSED 10 loops
14 DMA host->projection mem PASSED 10 loops
```

10 Loops: All tests passed

```
[root@ig1 ~]# exit <Enter>
```

logout

rlogin: connection closed.

```
[root@darc ~]# exit <Enter>
```

logout

rlogin: connection closed.

```
[root@hostname]#
```

Title Preliminary Training or Troubleshooting for GRE Console		First Made For: GRE Program		Rev C.
	GE Medical Systems Milwaukee, WI	000000TST	Sheet 34 of	39

4-16 SDDA — Cable Checkout

1. Verify the rear power cable is connected, switch is turned on and fan is operating.
2. Verify front power LED is illuminated.
3. Verify the 9-pin Sub D Scan Abort cable is connected between the DIP Board and the SDDA (towards the right). Verify DIP errors are present during reboot when application software is started).
4. Verify the Interlock cable is connected to the 9-pin Sub D and directly above to the DARC Node.
5. Verify the SCSI cable is connected to the DARC Node and not kinked or in poor condition,(verify no SCSI - RAID errors are present during reboot when application software is started)
6. Verify the two Audio connections to the HP Host Computer are good.

4-17 SDDA — View Error Message LOG

Visually scan the Error Message Log for Scan Data Disk failures. Record the Error message, time, and error number for reference during the troubleshooting process.

4-18 DARC Node — SDDA Utilities

The DARC Node controls the SDDA. The DARC Node must be up.

Open a Unix Shell and type the following:

```
{ctuser@hostname} rsh (space) darc <Enter>
```

```
Last login: Tue Oct 7 14:58:34 from oc
```

```
[ctuser@darc ~]$ cd (space) /usr/g/scripts <Enter>
```

```
{ctuser@darc scripts} ./SDDAUtills(space) -U <Enter>
```

could be a lot of already mounted info – look for the OK below:

```
[Logged by SDDAUtills] Mounted /dev/md0 at /raw_data. [OK]
```

```
{ctuser@darc scripts} cat(space) /proc/scsi/scsi <Enter>
```

Attached devices:

Host: scsi0 Channel: 00 Id: 00 Lun: 00

Vendor: MITSUMI Model: CD-ROM SR244W Rev: T01A

Type: CD-ROM ANSI SCSI revision: 02

Host: scsi1 Channel: 00 Id: 00 Lun: 00

Vendor: SEAGATE Model: ST336753LW Rev: 0005

Type: Direct-Access ANSI SCSI revision: 03

Host: scsi1 Channel: 00 Id: 01 Lun: 00

Vendor: SEAGATE Model: ST336753LW Rev: 0005

Type: Direct-Access ANSI SCSI revision: 03

```
{ctuser@darc scripts} exit <Enter>
```

Title		Preliminary Training or Troubleshooting for GRE Console			Rev C.
		First Made For: GRE Program			
	GE Medical Systems Milwaukee, WI	0000000TST			
		Sheet	35	of	39

4-19 SDDA — testHighSpeed

The DARC Node mounts the SDDA drives. The SDDA houses two Scan Data Disk Drives and the Intercom / Prescribed Tilt Board. This test verifies the SCSI cabling between the DARC Node and the SDDA. It writes data to the drive, then reads it back and verifies a checksum looking for dropped packets. It ensures the disk is able to transfer data at a rate equal to or greater than 75 MB/sec. This test shall be run with Applications down.

Visually verify the front panel SDDA LEDs for HD1 and HD2 blink during the test.

Open a Unix Shell and type the following:

```
{ctuser@hostname} rsh darc <Enter>
```

```
Last login: Tue Oct 7 15:12:51 from oc
```

```
[ctuser@darc ~]$ cd /usr/g/bin <Enter>
```

```
[ctuser@darc bin]$ ls test* <Enter>
```

```
testHighSpeed
```

```
[ctuser@darc bin]$ ./testHighSpeed <Enter>
```

```
Starting pass_number 1
```

```
saving 1024 MB to disk
```

```
1024 MB in 10.480 sec 97.710 MB/sec
```

```
Starting pass_number 2
```

```
saving 1024 MB to disk
```

```
1024 MB in 10.030 sec 102.094 MB/sec
```

```
Starting pass_number 3
```

```
saving 1024 MB to disk
```

```
1024 MB in 9.340 sec 109.636 MB/sec
```

```
Starting pass_number 4
```

```
saving 1024 MB to disk
```

```
1024 MB in 8.150 sec 125.644 MB/sec
```

```
Starting pass_number 5
```

```
saving 1024 MB to disk
```

```
1024 MB in 8.440 sec 121.327 MB/sec
```

Average speed: 111.282 Mb/sec. testHighSpeed PASSED.

Title	Preliminary Training or Troubleshooting for GRE Console	Rev C.
	GE Medical Systems Milwaukee, WI	First Made For: GRE Program 000000TST Sheet 36 of 39

[ctuser@darc bin]\$ **exit** <Enter>
logout
rlogin: connection closed.
{ctuser@hostname}

Title		Preliminary Training or Troubleshooting for GRE Console			Rev C.
		First Made For: GRE Program			
	GE Medical Systems Milwaukee, WI	000000TST			
		Sheet	37	of	39

4-20 SDDA — reconfigScanDisk

The DARC Node mounts the SDDA drives. The SDDA houses two Scan Data Disk Drives and the Intercom / Prescribed Tilt Board. This test wipes out all Scan Data and reconfigures the Scan Disk. This test shall be run with Applications down.

Open a Unix Shell and type the following:

```
{ctuser@hostname} su - <Enter>
```

```
Password: #bigguy <Enter>
```

```
[root@hostname]# rsh darc <Enter>
```

```
Last login: Tue Oct 7 14:42:40 from oc
```

```
You have new mail.
```

```
[root@darc ~]# cd /usr/g/scripts <Enter>
```

```
[root@darc scripts]# ls reconfig* <Enter>
```

```
reconfigScanDisk.log      reconfigScanDisk
```

```
[root@darc scripts]# ./reconfigScanDisk <Enter>
```

```
Appears that /raw_data is mounted. Trying to unmount..OK
```

```
Checking if SCSI driver [ aic79xx ] is loaded..
```

```
SCSI driver is already loaded. Good.
```

```
Checking if primary disk is IDE..
```

```
Detected IDE disk as primary. Good.
```

```
Detecting other SCSI disks on the system..
```

```
Detected 2 SCSI disk(s) : /dev/sda /dev/sdb
```

```
Creating primary partition on /dev/sda
```

```
Creating primary partition on /dev/sdb
```

```
Creating entries in /etc/raidtab.gre..OK
```

```
Initializing RAID device..OK
```

```
Examining if I created /dev/md0 correctly..OK
```

```
Creating XFS filesystem..OK
```

```
Creating mount point directory : /raw_data ..OK
```

```
Verifying what I did..
```

```
Mounting RAID disks at /raw_data..MountDevice
```

```
/dev/md0 on /raw_data type xfs (rw,nodev)
```

```
Verifying mount , Status = 1
```

```
[ should be 1 ]
```

```
[Logged by SDDAUtils] Mounted /dev/md0 at /raw_data. [OK]
```

```
OK
```

```
Creating a file /raw_data/test and writing to it..OK
```

```
Deleting file /raw_data/test..OK
```

```
Verification completed successfully
```

Title	Preliminary Training or Troubleshooting for GRE Console First Made For: GRE Program	Rev C.
	GE Medical Systems Milwaukee, WI	000000TST Sheet 38 of 39

Changing ownership of /raw_data to ctuser [owner] and ctuser [group]..OK

Applying permissions [777] to /raw_data..OK

All done. Transcript saved in reconfigScanDisk.log

[root@darc scripts]# **df <Enter>**

```
Filesystem      1k-blocks    Used Available Use% Mounted on
/dev/hda2        4188164 1238768 2949396 30% /
/dev/hda1         27296   11488   15808 43% /boot
none            967508   16412  951096  2% /dev/shm
/dev/hda6       30647188   40080 30607108  1% /usr/g
oc:/usr/g/queue/prospective
                14935616 2052336 12883280 14% /usr/g/queue/prospective
oc:/usr/g        14935616 2052336 12883280 14% /usr/g/oc
/dev/md0         71673044 2605820 69067224  4% /raw_data
```

NOTE * If /raw_data is not present then perform the ./startSDDA and df again*****

[root@darc scripts]# **./startSDDA <Enter>**

[root@darc scripts]# **df <Enter>**

```
Filesystem      1k-blocks    Used Available Use% Mounted on
/dev/hda2        4188164 1238768 2949396 30% /
/dev/hda1         27296   11488   15808 43% /boot
none            967508   16412  951096  2% /dev/shm
/dev/hda6       30647188   40080 30607108  1% /usr/g
oc:/usr/g/queue/prospective
                14935616 2052336 12883280 14% /usr/g/queue/prospective
oc:/usr/g        14935616 2052336 12883280 14% /usr/g/oc
/dev/md0         71673044 2605820 69067224  4% /raw_data
```

[root@darc scripts]# **exit <Enter>**

rlogin: connection closed.

[root@bayhostname]#

Title		Preliminary Training or Troubleshooting for GRE Console			Rev C.
		First Made For: GRE Program			
	GE Medical Systems Milwaukee, WI	000000TST			
		Sheet	39	of	39