

H-POWER INTERCONNECT CABLES

NOTE:

Before using this document please make sure that it is the latest revision using GEMS PCICS.

**This cover sheet controls the revision status of the entire document.
All sheets of this document are at the same revision level.**

Revision #	Name	Date	Description
A	B. D. Black	7-Nov-01	Initial Draft.
B	B. D. Black	6-Dec-01	Peer Review for Smart ID/IF BD
C	B.D. Black	2-Jan-02	Review of H-Power Interconnects from 12/19 – Mechanical Drawings; Fan Harness
D	B. Kohl	20-Mar-02	Atlantis build 4/02
E	B. Kohl	3-Jun-02	Alpha release of Harness Drawings
F	B. Kohl	4-Jun-02	Peer Review with William Brennan
G	B. Kohl	27-Aug-02	Post-Endeavor and Discovery Builds
H	D. Moller-Gunderson	31-Jan-03	Revision
J	D. Moller-Gunderson	1-Jul-03	Revision; Changes made in cable 2343336 (Inverter to IF Board)
1	A. Okuneye	25-May-04	Revision

Electronic Release Process - NOT Design Released

TITLE: H-POWER Interconnect Cables – Rotating		FIRST MADE FOR: H-Power Func <i>Proprietary to GE Medical Systems</i>	REV 1
	GE Medical Systems Milwaukee, Wi. USA	2308721GSP	Sheet 1 of 20

The source file for this document was created using Microsoft WORD 7.0a and is stored on the CT Products Engineering network.

State: COMPLETE - Production parts may not be defined by this document; it is an uneditable finished draft.
See the GEMS ePDM System to determine the status of this document.

H-POWER INTERCONNECT CABLES

TABLE OF CONTENTS

1	SCOPE.....	3
1.1	OVERVIEW	3
2	REFERENCE INFORMATION.....	3
3	WIRING & CABLE	4
3.1	GEMS SUPPLIED CABLES	4
3.1.1	RUN #300: Signal Cable – ORP Power Supply	4
3.1.2	RUN #301: Signal Cable – ORP to Alignment Lights.....	4
3.1.3	RUN#302: ORP to Remote sensor: Only Used For Engineering Purposes	5
3.1.4	RUN#303: RCIB2 – ORP to JEDI.....	5
3.1.5	RUN #304: ORP to S/R COMM.....	6
3.1.6	RUN #305: RCIB1 – ORP to DAS	7
3.1.7	RUN #306: Service Terminal: Only Used For Engineering Purposes	8
3.1.8	RUN #307: Ethernet Connection	8
3.1.9	RUN #308: H-UDI Port on ORP - REMOVED.....	9
3.1.10	RUN #309: Real-Time Signals from JEDI to DAS.....	9
3.1.11	RUN #310: JEDI & H4 Power Board Interface	10
3.1.12	RUN #311: JEDI Inverter to JEDI HV Tank.....	10
3.1.13	RUN #312: HVDC to JEDI Inverter.....	11
3.1.14	RUN #313: 140kVDC Connection to X-Ray Tube	11
3.1.15	RUN #314: Stator Power.....	11
3.1.16	RUN #315: DAS RCIB	12
3.1.17	RUN #316: CCB Power Supply – REMOVED – Now belongs to DAS.....	13
3.1.18	RUN #317: AC Supply for DAS	13
3.1.19	RUN #318: AC Supply for H4 Power IF Board.....	14
3.1.20	RUN #319: JEDI Inverter AC Power.....	14
3.1.21	RUN #320: Heat Exchanger PS AC Input.....	14
3.1.22	RUN #321: Oil Pump Power.....	15
3.1.23	RUN #322: Pressure Limit Switch	15
3.1.24	RUN #323: Fan #1.....	15
3.1.25	RUN #324: Fan #2.....	16
3.1.26	RUN #325: Temperature Limit Switch	16
3.1.27	RUN #326: Smart ID Board.....	17
3.1.28	RUN #327: Power to ORP Power Supply: REMOVED.....	17
3.1.29	RUN#329: HX Power Supply to Fan Harness: REMOVED	17
3.1.30	RUN#330: Jedi Inverter to Auxiliary Box CAN Cable	18
3.1.31	RUN#331: Jedi to HV Tank.....	19
3.1.32	RUN#332: Jedi Inverter to Auxiliary Box.....	20

Electronic Release Process - NOT Design Released

TITLE: H-POWER Interconnect Cables – Rotating		FIRST MADE FOR: H-Power Func Proprietary to GE Medical Systems	REV 1
	GE Medical Systems Milwaukee, Wi. USA	2308721GSP	Sheet 2 of 20

The source file for this document was created using Microsoft WORD 7.0a and is stored on the CT Products Engineering network.

State: COMPLETE - Production parts may not be defined by this document; it is an uneditable finished draft.
See the GEMS ePDM System to determine the status of this document.

H-POWER INTERCONNECT CABLES

1 SCOPE

The purpose of this document is to define the electrical signals on all the interconnecting cables used on the rotating gantry of the H-Power system. The interconnect documents and drawings are referenced with part number 2308721.

1.1 OVERVIEW

The Rotating gantry of the H-Power system is composed of eight major physical assemblies. These are:

- ORP Board
- DAS
- Slip Ring
- Collimator
- JEDI Inverter
- Auxiliary Box
- JEDI High Voltage Tank
- Heat Exchanger
- X-Ray Tube

These assemblies are interconnected with various power and signal cables. This document identifies the physical wires used in each cable and the electrical signals on those wires. **This document refers to the components used on a 16-slice system.**

2 REFERENCE INFORMATION

The following documents form a part of this specification to the extent applicable. All documents are under revision control and any changes require a Product Change Notice (PCN). Where listed, document revisions are for information only. The latest revision of each document in effect at time of purchase shall apply.

TBD ala. Raw material specs, etc.?

Electronic Release Process - NOT Design Released

TITLE: H-POWER Interconnect Cables – Rotating		FIRST MADE FOR: H-Power Func <i>Proprietary to GE Medical Systems</i>	REV 1
	GE Medical Systems <i>Milwaukee, Wi. USA</i>	2308721GSP	Sheet 3 of 20

The source file for this document was created using *Microsoft WORD 7.0a* and is stored on the CT Products Engineering network.

H-POWER INTERCONNECT CABLES

3 WIRING & CABLE

----{ }---- Spade Lugs
 ----> Or <---- Male Connection
 ----< Or >---- Female Connection
 ----8 8---- Ring Terminals
 ----X X---- Ethernet

3.1 GEMS Supplied Cables

3.1.1 RUN #300: SIGNAL CABLE – ORP POWER SUPPLY

24V DAS Power Supply

P/N 2358011 – MDAS

ORP Board

P/N 2394982 - GDAS

	<u>Signal Name</u>	<u>J1</u>
----{ + OUT }----	+24V DC	----< 1 <---- P24_IN
----{ - OUT }----	GND	----< 2 <---- P24_RTN

Material : Spade Lugs to 2-position Mate-N-Lock Plug (AMP p/n 350777-1) with 2 sockets (AMP p/n 350536-2), to two spade tongues (AMP p/n 52947-1) shielded terminal. 20 AWG wire, 1 set blk/red twisted pair.

3.1.2 RUN #301: SIGNAL CABLE – ORP TO ALIGNMENT LIGHTS

ORP Board

P/N 2333306

Alignment Lights

<u>J2</u>	<u>Signal Name</u>	
----> 1 >----	+5V L_ON_P1	----> >---- LFT COR-1 --->RGT COR-1
----> 6 >----	P2_RTN	----> G >---- LFT COR-2 --->RGT COR-2
----> 2 >----	P1_RTN	----> >---- SAG-1 --->BOT INT-1
----> 7 >----	+5V L_ON_P3	----> G >---- SAG-2 --->BOT INT-2
----> 5 >----	P1_RTN	----> >----
----> 3 >----	+5V L_ON_P2	----> G >---- INT LFT-1 --->INT RGT-1
----> 8 >----	P4_RTN	----> >---- INT LFT-2 --->INT RGT-2
----> 4 >----	P3_RTN	----> G >---- EXT-1
----> 9 >----	+5V L_ON_P4	----> >---- EXT-2

Material : 9-position sub-D plug (AMP p/n 205204-4) with sockets (AMP p/n 66506-3) to seven 4-position mini MATE-N-LOK (AMP p/n 172167-1), twisted pairs "FT-4" with sockets (AMP p/n 770988-3). 22AWG Each a set of twisted pairs.

Electronic Release Process - NOT Design Released

TITLE:
H-POWER Interconnect Cables – Rotating

FIRST MADE FOR: H-Power Func
Proprietary to GE Medical Systems

REV
1



GE Medical Systems
Milwaukee, Wi. USA

2308721GSP

Sheet 4 of 20

H-POWER INTERCONNECT CABLES

3.1.3 RUN#302: ORP TO REMOTE SENSOR: ONLY USED FOR ENGINEERING PURPOSES

3.1.4 RUN#303: RCIB2 – ORP TO JEDI

<u>ORP Board</u>			<u>P/N 2333307</u>			<u>JEDI Inverter</u>		
16 pin Straight		15 pin Sub-D				15 pin Sub-D	16 pin Straight	
		<u>J5</u>	<u>Signal Name</u>			<u>J3</u>		
1	CAN_0V	----< 1 <-----	LGND			-----> 1 >----	CAN_0V	1
2	CAN_12V	----< 9 <-----	CAN_12V			-----> 9 >----	CAN_12V	2
3	CAN1_P	----< 2 <-----	CAN2_P			-----> 2 >----	CAN1_P	3
4	CAN1_N	----< 10 <-----	CAN2_N			-----> 10 >----	CAN1_N	4
5	FAULT_P	----< 3 <-----	FAULT2_P			-----> 3 >----	FAULT_P	5
6	FAULT_N	----< 11 <-----	FAULT2_N			-----> 11 >----	FAULT_N	6
7	RESET2_P	----< 4 <-----	RESET2_P			-----> 4 >----	RESET2_P	7
8	RESET2_N	----< 12 <-----	RESET2_N			-----> 12 >----	RESET2_N	8
9	EXPCMD_P	----< 5 <-----	EXPCMD_P			-----> 5 >----	EXPCMD_P	9
10	EXPCMD_N	----< 13 <-----	EXPCMD_N			-----> 13 >----	EXPCMD_N	10
11	XRON_P	----< 6 <-----	XRON_P			-----> 6 >----	XRON_P	11
12	XRON_N	----< 14 <-----	XRON_N			-----> 14 >----	XRON_N	12
13	DAS_TRIG_H	----< 7 <-----	DAS_TRIG_H			-----> 7 >----	DAS_TRIG_	13
14	DAS_TRIG_L	----< 15 <-----	DAS_TRIG_L			-----> 15 >----	DAS_TRIG_L	14
15	PSFB2	----< 8 <-----	PSFB2			-----> 8 >----	PSFB2	15
16	GND						GND	16

Material : Connection to ORP is 15-position right angle over-molded sub-D plug (AMP p/n 205206-3) with pins (AMP p/n 66506-3). Connection to Inverter is 15-position sub-D plug (AMP p/n 205206-3) with pins (66506-3) . Wire is 7 pair +1, 22 AWG with overall SHIELD, (GE P/N 46-243509P1) "FT-4"

Electronic Release Process - NOT Design Released

TITLE:
H-POWER Interconnect Cables – Rotating

FIRST MADE FOR: H-Power Func
Proprietary to GE Medical Systems

REV
1



GE Medical Systems
Milwaukee, Wi. USA

2308721GSP

Sheet **5** of **20**

H-POWER INTERCONNECT CABLES

3.1.5 RUN #304: ORP TO S/R COMM.

ORP Board

P/N 2333308

Slipring

15 pin Sub-D

15 pin Sub-D

<u>J6</u>	<u>Signal Name</u>	<u>J2</u>
----> 1 >----	SER_OUTBOUND	----> 1 >----
----> 2 >----	EXP_ENBL	----> 2 >----
----> 3 >----	GND	----> 3 >----
----> 4 >----	GND	----> 4 >----
----> 5 >----	GND	----> 5 >----
----> 6 >----	SER_OUTBOUND	----> 6 >----
----> 7 >----	EXP_ENBL	----> 7 >----
----> 8 >----	GND	----> 8 >----
----> 9 >----	GND	----> 9 >----
----> 10 >----	GND	----> 10 >----
----> 11 >----	SER_REF	----> 11 >----
----> 12 >----	SER_REF	----> 12 >----
----> 13 >----	SER_INBOUND	----> 13 >----
----> 14 >----	SER_INBOUND	----> 14 >----
----> 15 >----	SHIELD_GND	----> 15 >----

Material : Both connectors are High Density (3-Row) 15-pin sub-D, connection to ORP is receptacle (AMP p/n 748565-1) with sockets (AMP p/n 748610-2) and connection to S/R is plug (AMP p/n 748364-1) with pins (AMP p/n 748333-2) Wire is Belden #9728 22 AWG with overall SHIELD, "FT-4"

Electronic Release Process - NOT Design Released

TITLE:

H-POWER Interconnect Cables – Rotating

FIRST MADE FOR: H-Power Func

Proprietary to GE Medical Systems

REV

1



GE Medical Systems
Milwaukee, Wi. USA

2308721GSP

Sheet **6** of **20**

The source file for this document was created using *Microsoft WORD 7.0a* and is stored on the CT Products Engineering network.

State: COMPLETE - Production parts may not be defined by this document; it is an uneditable finished draft.
See the GEMS ePDM System to determine the status of this document.

H-POWER INTERCONNECT CABLES

3.1.6 RUN #305: RCIB1 – ORP TO DAS

ORP Board

P/N 2333309

DAS

15 pin Sub-D

J7

Signal Name

J31

15 pin Sub-D

----->	1	>-----	CAN_RTN	----->	1	>-----	ISO_GND
----->	2	>-----	CAN1_P	----->	2	>-----	RCIB_CAN_BUS+
----->	3	>-----	FAULT1_P	----->	3	>-----	RCIB_FLT_BUS+
----->	4	>-----	RESET1_P	----->	4	>-----	RCIB_CAN_RST+
----->	5	>-----	EXPCMD_P	----->	5	>-----	EXP_CMD+
----->	6	>-----	XRON_P	----->	6	>-----	
----->	7	>-----	DAS_TRIG_H	----->	7	>-----	E_TRIG+
----->	8	>-----	PSFB1	----->	8	>-----	CAN_PSOK_OUT
----->	9	>-----	CAN_+12V	----->	9	>-----	+12V_CAN
----->	10	>-----	CAN1_N	----->	10	>-----	RCIB_CAN_BUS-
----->	11	>-----	FAULT1_N	----->	11	>-----	RCIB_FLT_BUS-
----->	12	>-----	RESET1_N	----->	12	>-----	RCIB_CAN_RST-
----->	13	>-----	EXPCMD_N	----->	13	>-----	EXP_CMD-
----->	14		N/C	----->	14		N/C
----->	15	>-----	DAS_TRIG_L	----->	15	>-----	E_TRIG-

Material : Connection to ORP is 15-position Sub-D receptacle (AMP p/n 205205-2) with sockets (AMP p/n 66504-3) and connection to DAS is 15-position Sub-D plug (AMP p/n 205206-3) with pins (66506-3). Wire is 7 pair +1, 22 AWG with overall SHIELD (GE p/n 46-243509P1), "FT-4"

*Pins 6, 8, and 14 do not match up with the signals from the ORP, but they are signals that will not ever be used, therefore this problem is noted but not acted upon.

Electronic Release Process - NOT Design Released

TITLE:

H-POWER Interconnect Cables – Rotating

FIRST MADE FOR: H-Power Func

Proprietary to GE Medical Systems

REV

1



GE Medical Systems
Milwaukee, Wi. USA

2308721GSP

Sheet 7 of 20

H-POWER INTERCONNECT CABLES

3.1.7 RUN #306: SERVICE TERMINAL: ONLY USED FOR ENGINEERING PURPOSES

3.1.8 RUN #307: ETHERNET CONNECTION

ORP Board

Ethernet Connection

9 pin Ethernet
Cable

9 pin Ethernet
Cable

J9

Signal Name

----X 1 X----
----X 2 X----
----X 3 X----
----X 4 X----
----X 5 X----
----X 6 X----
----X 7 X----
----X 8 X----

E_TX0B_H
E_TX0B_L
E_RX1B_H

E_RX1B_L

----X 1 X----
----X 2 X----
----X 3 X----
----X 4 X----
----X 5 X----
----X 6 X----
----X 7 X----
----X 8 X----

Material : 7 pair +1, 22 AWG with overall SHIELD, "FT-4"

Electronic Release Process - NOT Design Released

TITLE:
H-POWER Interconnect Cables – Rotating

FIRST MADE FOR: H-Power Func
Proprietary to GE Medical Systems

REV
1



GE Medical Systems
Milwaukee, Wi. USA

2308721GSP

Sheet **8** of **20**

The source file for this document was created using *Microsoft WORD 7.0a* and is stored on the CT Products Engineering network.

State: COMPLETE - Production parts may not be defined by this document; it is an uneditable finished draft.
See the GEMS ePDM System to determine the status of this document.

H-POWER INTERCONNECT CABLES

3.1.9 RUN #308: H-UDI PORT ON ORP - REMOVED

3.1.10 RUN #309: REAL-TIME SIGNALS FROM JEDI TO DAS

JEDI Inverter

P/N 2333310

DAS

9 pin Sub-D

9 pin Sub-D

<u>J2</u>	<u>Signal Name</u>	<u>J32</u>
-----< 1 <-----	KV_MEAS_P	-----< 3 <-----
2		
-----< 3 <-----	KV_MEAS_N	-----< 7 <-----
4		
5		
6		
-----< 7 <-----	MA_MEAS_P	-----< 4 <-----
8		
-----< 9 <-----	MA_MEAS_N	-----< 8 <-----

Material : Both connectors are sub-D, the connector to the Jedi is a 9-position plug (AMP p/n 205204-4) with pins (AMP p/n 66506-3), while the connector to the DAS is a 9-position receptacle (AMP p/n 205203-3) with sockets (66504-3). Wire is Alpha #2212C 7 pair +1, 22 AWG with overall SHIELD, "FT-4"

Electronic Release Process - NOT Design Released

TITLE:

H-POWER Interconnect Cables – Rotating

FIRST MADE FOR: H-Power Func

Proprietary to GE Medical Systems

REV

1



GE Medical Systems
Milwaukee, Wi. USA

2308721GSP

Sheet **9** of **20**

The source file for this document was created using Microsoft WORD 7.0a and is stored on the CT Products Engineering network.

State: COMPLETE - Production parts may not be defined by this document; it is an uneditable finished draft.
See the GEMS ePDM System to determine the status of this document.

H-POWER INTERCONNECT CABLES

3.1.11 RUN #310: JEDI & H4 POWER BOARD INTERFACE

<u>JEDI Inverter</u>		<u>P/N 2343336</u>	<u>H-Power IF Board</u>
16 pin Straight	15 pin Sub-D <u>J5</u>	<u>Signal Name</u>	15 pin Sub-D <u>J2</u> Male
	----> 3 >----	T70_P	----< 3 <----
	----> 4 >----	T70_N	----< 4 <----
	----> 5 >----	FAN_P	----< 5 <----
	----> 13 >----	FAN_N	----< 13 <----
	----> 7 >----	Spare	----< 7 <----
	----> 15 >----	Spare	----< 15 <----
	----> 11 >----	T40_P	----< 11 <----
	----> 12 >----	T40_N	----< 12 <----
	----> 8 >----	Spare	----< 8 <----

Material : Connector to the Inverter is 15-position sub-D right-angle overmold receptacle (AMP p/n 205205-2) with sockets (AMP p/n 66504-3). Connector to the Power IF board is 15-position sub-D receptacle (AMP 205205-2) with sockets (AMP p/n 66504-3). Cable is 7 pair +1 22 AWG (GE p/n 46-315414P1) with overall SHIELD, "FT-4"

3.1.12 RUN #311: JEDI INVERTER TO JEDI HV TANK

<u>JEDI Inverter</u>		<u>P/N 2283352</u>	<u>JEDI HV Tank</u>
Ring Terminals Twisted Pair	<u>U12, U5</u>	<u>Signal Name</u>	Ring Terminal Connection
Stud diameter 6	----8 1 8----	Inverter/Tank Cable (HV)	----8 1 8---- Primary connection
	----8 2 8----	Inverter/Tank Cable (HV)	----8 2 8---- Diameter 8

Material : GE Part Number 2283352ADW

Electronic Release Process - NOT Design Released

TITLE: H-POWER Interconnect Cables – Rotating		FIRST MADE FOR: H-Power Func <i>Proprietary to GE Medical Systems</i>	REV 1
	GE Medical Systems Milwaukee, Wi. USA	2308721GSP	Sheet 10 of 20

The source file for this document was created using Microsoft WORD 7.0a and is stored on the CT Products Engineering network.

H-POWER INTERCONNECT CABLES

3.1.13 RUN #312: HVDC TO JEDI INVERTER

<u>Slipring</u>			<u>P/N 2363942</u>	<u>JEDI Inverter</u>		
Ring Terminal Twisted Pair	<u>J1</u>		<u>Signal Name</u>	<u>U1, U2</u>	Ring terminal Twisted Pair	
	-----8	1 8-----	HVDC	-----8 2 8-----	U1 for minus DC	
	-----8	2 8-----	HVDC_RTN	-----8 1 8-----	Bus: Diameter 8	
					U2 for plus DC	
					Bus: Diameter 6	

Material : 7 pair +1, 22 AWG with overall SHIELD, "FT-4"

3.1.14 RUN #313: 140kVDC CONNECTION TO X-RAY TUBE

<u>JEDI HV Tank</u>			<u>P/N 2337317</u>	<u>X-Ray Tube</u>		
3-conductor Oil-filled & sealed	<u>J2 HV Tank receptacle</u>		<u>Signal Name</u>	<u>J7</u>	3-conductor Receptacle	
	-----{	1 }-----	Large	-----{	3 }-----	
	-----{	2 }-----	Common	-----{	2 }-----	
	-----{	3 }-----	Small	-----{	1 }-----	

Material : 3 Conductor, shielded, 160kV HV Cable, "FT-4"

3.1.15 RUN #314: STATOR POWER

<u>Auxiliary Box</u>			<u>P/N 2317421</u>	<u>X-Ray Tube</u>		
Phoenix 3-pin connector	<u>J2 Rotor Board</u>		<u>Signal Name</u>	<u>J2</u>	3-pin flag Internal connection	
	-----{	1 }-----	"u" – Power, Phase A	-----{	1 }-----	
	-----{	2 }-----	"v" – Power, Phase B	-----{	2 }-----	
	-----{	3 }-----	"w" – Power, Phase C	-----{	3 }-----	
			GND			

Material : Shielded, 3-conductor 16 AWG Phoenix GE PN 2206319-2, "FT-4"

Electronic Release Process - NOT Design Released

TITLE: H-POWER Interconnect Cables – Rotating		FIRST MADE FOR: H-Power Func <i>Proprietary to GE Medical Systems</i>	REV 1
	GE Medical Systems Milwaukee, Wi. USA	2308721GSP	Sheet 11 of 20

The source file for this document was created using Microsoft WORD 7.0a and is stored on the CT Products Engineering network.

H-POWER INTERCONNECT CABLES

3.1.16 RUN #315: DAS RCIB

<u>Collimator</u>		<u>P/N 2333314</u>		<u>DAS</u>
15 pin Sub-D				15 pin Sub-D
	<u>J1</u>	<u>Signal Name</u>	<u>J30</u>	
	-----< 1 <-----	CAN_RTN	-----> 1 >-----	ISO_GND
J2 on collimator	-----< 2 <-----	CAN1_P	-----> 2 >-----	RCIB_CAN_BUS+
is a terminator	-----< 3 <-----	FAULT1_P	-----> 3 >-----	RCIB_FLT_BUS+
Connection	-----< 4 <-----	RESET1_P	-----> 4 >-----	RCIB_CAN_RST+
	-----< 5 <-----	EXPCMD_P	-----> 5 >-----	EXP_CMD+
	-----< 6 <-----	XRON_P	-----> 6 >-----	N/C
	-----< 7 <-----	DAS_TRIG_H	-----> 7 >-----	E_TRIG+
	-----< 8 <-----	PSFB1	-----> 8 >-----	CAN_PSOK_OUT
	-----< 9 <-----	CAN_+12V	-----> 9 >-----	+12V_CAN
	-----< 10 <-----	CAN1_N	-----> 10 >-----	RCIB_CAN_BUS-
	-----< 11 <-----	FAULT1_N	-----> 11 >-----	RCIB_FLT_BUS-
	-----< 12 <-----	RESET1_N	-----> 12 >-----	RCIB_CAN_RST-
	-----< 13 <-----	EXPCMD_N	-----> 13 >-----	EXP_CMD-
	-----< 14 <-----	XRON_N	-----> 14 >-----	XRON_N
	-----< 15 <-----	DAS_TRIG_L	-----> 15 >-----	E_TRIG-

14 not used by collimator

Material : Both connectors are 15-Position sub-D plug (AMP p/n 205206-3) with pins (AMP p/n 66506-3). Cable is 7 pair +1, 22 AWG (GE p/n 46-243509P1) with overall SHIELD, "FT-4"

***Pins 6, 8, and 14 do not match up with the signals from the CCB, but they are signals that will not ever be used, therefore this problem is noted but not acted upon.**

Electronic Release Process - NOT Design Released

TITLE:
H-POWER Interconnect Cables – Rotating

FIRST MADE FOR: H-Power Func
Proprietary to GE Medical Systems

REV
1



GE Medical Systems
Milwaukee, Wi. USA

2308721GSP

Sheet **12** of **20**

H-POWER INTERCONNECT CABLES

3.1.17 RUN #316: CCB POWER SUPPLY – REMOVED – NOW BELONGS TO DAS

3.1.18 RUN #317: AC SUPPLY FOR DAS

Slipring

P/N 2260623-2 – MDAS

15V DC PWR SUPPLY
(Fusebox)

P/N 2392420 – GDAS

Ring Terminal
Twisted Pair

J4

Signal Name

Faston push-pin
Twisted Pair

-----8	1	8-----	120VAC	-----{	1	}-----
-----8	2	8-----	0VAC	-----{	2	}-----

Material : Connector to the slipring is a ring terminal (AMP p/n 2-34855), connector to the DAS Fusebox is a Faston Insulated Terminator (AMP 640907-2). S/R to fusebox is a 12 AWG brn (GE p/n 2261425-4)/blu (GE p/n 2261425-21) twisted pair. "FT-4"

Electronic Release Process - NOT Design Released

TITLE:

H-POWER Interconnect Cables – Rotating

FIRST MADE FOR: H-Power Func

Proprietary to GE Medical Systems

REV

1



GE Medical Systems
Milwaukee, Wi. USA

2308721GSP

Sheet **13** of **20**

The source file for this document was created using Microsoft WORD 7.0a and is stored on the CT Products Engineering network.

State: COMPLETE - Production parts may not be defined by this document; it is an uneditable finished draft.
See the GEMS ePDM System to determine the status of this document.

H-POWER INTERCONNECT CABLES

3.1.19 RUN #318: AC SUPPLY FOR H4 POWER IF BOARD

<u>Slipring</u>		<u>P/N 2333312</u>	<u>H-Power IF Board</u>
Ring Terminal			4-Position Mat-N-Lok
Twisted Pair	<u>J3</u>	<u>Signal Name</u>	<u>J1</u> Twisted Pair
-----8 1 8-----		120VAC	-----< 1 <-----
		0VAC	-----< 3 <-----
		120VAC	-----< 2 <-----
-----8 2 8-----		0VAC	-----< 4 <-----

Material : The connector going to the slipring is a ring terminal and the connection going to the H-Power IF board is 4-Position Mate-N-Lok plug (AMP p/n 350779-1). 14 AWG, brn (GE p/n 2261425-3)/blu (GE p/n 2261425-21) twisted pair. "FT-4"

3.1.20 RUN #319: JEDI INVERTER AC POWER

<u>Slip Ring</u>		<u>P/N 2333311</u>	<u>JEDI Inverter (Aux Box)</u>
Ring Terminal			2-Pin Phoenix Male
Twisted Pair	<u>J4</u>	<u>Signal Name</u>	<u>J6</u> Twisted Pair
-----8 1 8-----		120VAC	-----< 1 <----- Via a fuse
-----8 2 8-----		0VAC	-----< 2 <----- Via a clamped gnd

Material : Before from the slipring to the fuse, please see run 317. The connector from the fuse to the Inverter is a 2-position Phoenix cap (Phoenix p/n 1808816). Crimp (Phoenix 3190645). S/R to fuse and fuse to Aux box is 16 AWG twisted pair brn (GE p/n 2261425-2)/blu(GE p/n 2261425-20). Connectors to the S/R are ring terminals (AMP 2-31894-2). "FT-4"

3.1.21 RUN #320: HEAT EXCHANGER PS AC INPUT

<u>H-Power IF Board</u>		<u>P/N 2334120</u>	<u>HX Power Supply</u>
2-pin Mat-N-Lok			2-Pin Spade
Twisted Pair	<u>J5</u>	<u>Signal Name</u>	<u>T1</u> Tongue
-----> 1 >-----		120VAC	-----{ 1 Scre
-----> 2 >-----		0VAC	-----{ 2 Scre

Material : "FT-4"

Electronic Release Process - NOT Design Released

TITLE: H-POWER Interconnect Cables – Rotating		FIRST MADE FOR: H-Power Func <i>Proprietary to GE Medical Systems</i>	REV 1
	GE Medical Systems Milwaukee, Wi. USA	2308721GSP	Sheet 14 of 20

The source file for this document was created using Microsoft WORD 7.0a and is stored on the CT Products Engineering network.

H-POWER INTERCONNECT CABLES

3.1.22 RUN #321: OIL PUMP POWER

<u>H-Power IF Board</u>	<u>P/N 2334120</u>	<u>X-Ray Tube Oil Pump</u>
4-Pos Mat-N-Lock In Line CAP Twisted Pair <u>J6</u>	<u>Signal Name</u>	4-Pos Mat-N-Lock Twisted Pair <u>J3</u>
----> 1 >----	120VAC	----< 1 <----
----> 2 >----	120VAC	----< 2 <----
----> 3 >----	0VAC	----< 3 <----
----> 4 >----	0VAC	----< 4 <----

Material : "FT-4"

3.1.23 RUN #322: PRESSURE LIMIT SWITCH

<u>H-Power IF Board</u>	<u>P/N 2334120</u>	<u>Pressure Limit Switch</u>
2-pos Phoenix Twisted Pair <u>J8</u>	<u>Signal Name</u>	2-Pin Molex Twisted Pair
----> 1 >----	OVT_N/OVP_P	----{ 1 }----
----> 2 >----	OVTP_N	----{ 2 }----

Via a switch

Material : "FT-4"

3.1.24 RUN #323: FAN #1

<u>H-Power IF Board</u>	<u>P/N 2334120</u>	<u>Fan #1</u>
4-pos Phoenix Twisted Pair <u>J11</u>	<u>Signal Name</u>	4-Pin Mat-N-Lock Twisted Pair <u>J1</u>
----{ 1 }----	48VDC	----{ 1 }----
----{ 2 }----	THERM_RTN	----{ 2 }----
----{ 3 }----	SPD_CTRL	----{ 4 }----
----{ 4 }----	GND	----{ 4 }----
	NC	----{ 3 }----

Material : "FT-4"

Electronic Release Process - NOT Design Released

TITLE:
H-POWER Interconnect Cables – Rotating

FIRST MADE FOR: H-Power Func
Proprietary to GE Medical Systems

REV
1



GE Medical Systems
Milwaukee, Wi. USA

2308721GSP

Sheet **15** of **20**

H-POWER INTERCONNECT CABLES

3.1.25 RUN #324: FAN #2

H-Power IF Board

P/N 2334120

Fan #2

4-pos Phoenix
Twisted Pair

J11

Signal Name

J1

4-Pin Mat-N-Lock
Twisted Pair

-----{ 1 }-----
-----{ 2 }-----
-----{ 3 }-----
-----{ 4 }-----

48VDC
THERM_RTN
SPD_CTRL
GND
NC

-----{ 1 }-----
-----{ 2 }-----
-----{ 4 }-----
-----{ 3 }-----

Material : "FT-4"

3.1.26 RUN #325: TEMPERATURE LIMIT SWITCH

H-Power IF Board

P/N 2334120

Temp Limit Switch

3-pos Phoenix
Twisted Pair

J7

Signal Name

J4

2-Pos Mat-N-Lock
Twisted Pair

-----> 1 >-----
-----> 2 >-----
-----> 3 >-----

OVTP_P
OVT_N
NC

-----{ 1 }-----
-----{ 2 }-----

Via a switch

Material : "FT-4"

Electronic Release Process - NOT Design Released

TITLE:

H-POWER Interconnect Cables – Rotating

FIRST MADE FOR: H-Power Func
Proprietary to GE Medical Systems

REV
1



GE Medical Systems
Milwaukee, Wi. USA

2308721GSP

Sheet **16** of **20**

The source file for this document was created using Microsoft WORD 7.0a and is stored on the CT Products Engineering network.

State: COMPLETE - Production parts may not be defined by this document; it is an uneditable finished draft.
See the GEMS ePDM System to determine the status of this document.

H-POWER INTERCONNECT CABLES

3.1.27 RUN #326: JEDI INVERTER TO SMART ID BOARD

<u>Jedi Inverter</u>	<u>P/N 2343336</u>	<u>Smart ID</u>
15 pos Sub-D Connector		9-Position Sub-D Female Connector
<u>J5</u>	<u>Signal Name</u>	<u>J1</u>
----> 1 >----	CAN_+15V_RTN	----> 1 >----
----> 9 >----	CAN_+15V	----> 6 >----
----> 2 >----	CAN_P	----> 2 >----
----> 10 >----	CAN_N	----> 7 >----
----> 3 >----	T70_P	----> 3 >----
----> 4 >----	T70_N	----> 8 >----
----> 6 >----	Spare	----> 4 >----
----> 14 >----	Spare	----> 9 >----
----> 8 >----	Spare	----> 5 >----

Material : 7 pair +1, 22 AWG with overall SHIELD, "FT-4"

3.1.28 RUN #327: POWER TO ORP POWER SUPPLY: REMOVED

3.1.29 RUN#329: HX POWER SUPPLY TO FAN HARNESS: REMOVED

Electronic Release Process - NOT Design Released

TITLE: H-POWER Interconnect Cables – Rotating		FIRST MADE FOR: H-Power Func <i>Proprietary to GE Medical Systems</i>	REV 1
	GE Medical Systems Milwaukee, Wi. USA	2308721GSP	Sheet 17 of 20

The source file for this document was created using Microsoft WORD 7.0a and is stored on the CT Products Engineering network.

H-POWER INTERCONNECT CABLES

3.1.30 RUN#330: JEDI INVERTER TO AUXILIARY BOX CAN CABLE

Jedi Inverter

P/N 2328759

Jedi Auxiliary Box

<u>J13</u>	<u>Signal Name</u>	
----> 1 >----	GND	----> 1 >----
----> 2 >----	P15V	----> 2 >----
----> 3 >----	M15V	----> 3 >----
----> 4 >----	GND	----> 4 >----
----> 5 >----	PP11	----> 5 >----
----> 6 >----	PP5	----> 6 >----
----> 7 >----	GND	----> 7 >----
----> 8 >----	P15V	----> 8 >----
----> 9 >----	PP12	----> 9 >----
----> 10 >----	PP6	----> 10 >----
----> 11 >----	PP13	----> 11 >----
----> 12 >----	PP7	----> 12 >----
----> 13 >----	PP14	----> 13 >----
----> 14 >----	PP8	----> 14 >----
----> 15 >----	GND	----> 15 >----
----> 16 >----	P15V	----> 16 >----
----> 17 >----	FORCE_DISCHARGE	----> 17 >----
----> 18 >----	INHIB_DISCHARGE	----> 18 >----
----> 19 >----	VCC	----> 19 >----
----> 20 >----	DISCHARGING_DC_BUS	----> 20 >----
----> 21 >----	PP15	----> 21 >----
----> 22 >----	PP9	----> 22 >----
----> 23 >----	M15V	----> 23 >----
----> 24 >----	R15V	----> 24 >----
----> 25 >----	GND	----> 25 >----
----> 26 >----	PP10	----> 26 >----

Electronic Release Process - NOT Design Released

TITLE:
H-POWER Interconnect Cables – Rotating

FIRST MADE FOR: H-Power Func
Proprietary to GE Medical Systems

REV
1



GE Medical Systems
Milwaukee, Wi. USA

2308721GSP

Sheet **18** of **20**

The source file for this document was created using *Microsoft WORD 7.0a* and is stored on the CT Products Engineering network.

State: COMPLETE - Production parts may not be defined by this document; it is an uneditable finished draft.
See the GEMS ePDM System to determine the status of this document.

H-POWER INTERCONNECT CABLES

3.1.31 RUN#331: JEDI TO HV TANK

HV Tank

P/N 2328760

Jedi Inverter

<u>J12</u>	<u>Signal Name</u>	<u>J12</u>
----> 1 >----	GND	----> 1 >----
----> 2 >----	CONT_TANK_OUT	----> 2 >----
----> 3 >----	TANK_TEMP_P	----> 3 >----
----> 4 >----	TANK_TEMP_N	----> 4 >----
----> 5 >----	P15V_PPC	----> 5 >----
----> 6 >----	TANK_ID_0	----> 6 >----
----> 7 >----	TANK_ID_1	----> 7 >----
----> 8 >----	TANK_ID_2	----> 8 >----
----> 9 >----	TANK_ID_3	----> 9 >----
----> 10 >----	0V_FAN_1	----> 10 >----
----> 11 >----	24V_FAN_1	----> 11 >----
----> 12 >----	24V_FAN_1	----> 12 >----
----> 13 >----	0V_FAN_1	----> 13 >----
----> 14 >----	INJECT_50mA	----> 14 >----
----> 15 >----	CA_MA_P	----> 15 >----
----> 16 >----	CA_MA_N	----> 16 >----
----> 17 >----	AN_BLEEDER	----> 17 >----
----> 18 >----	CA_BLEEDER	----> 18 >----
----> 19 >----	CONT_TANK_IN	----> 19 >----
----> 20 >----	GND	----> 20 >----

Electronic Release Process - NOT Design Released

TITLE:

H-POWER Interconnect Cables – Rotating

FIRST MADE FOR: H-Power Func

Proprietary to GE Medical Systems

REV

1



GE Medical Systems
Milwaukee, Wi. USA

2308721GSP

Sheet **19** of **20**

The source file for this document was created using *Microsoft WORD 7.0a* and is stored on the CT Products Engineering network.

State: COMPLETE - Production parts may not be defined by this document; it is an uneditable finished draft.
See the GEMS ePDM System to determine the status of this document.

H-POWER INTERCONNECT CABLES

3.1.32 RUN#332: JEDI INVERTER TO AUXILIARY BOX

Jedi Inverter

P/N 2360922

Auxiliary Box

<u>J1</u>	<u>Signal Name</u>	<u>J1</u>
----> 1 >----	0V_15V_EXT	----> 5 >----
----> 2 >----	P_15V_EXT	----> 6 >----
----> 3 >----	0V_GATE	----> 7 >----
----> 4 >----	24V_GATE	----> 8 >----
----> 5 >----	0V_FAN_1	----> 9 >----
----> 6 >----	24V_FAN_1	----> 10 >----
----> 7 >----	0V_FAN_2	----> 11 >----
----> 8 >----	24V_FAN_2	----> 12 >----
----> 9 >----	0V_FAN_3	----> 13 >----
----> 10 >----	24V_FAN_3	----> 14 >----
----> 11 >----	0V_FAN_4	----> 15 >----
----> 12 >----	24V_FAN_4	----> 16 >----

Electronic Release Process - NOT Design Released

TITLE:

H-POWER Interconnect Cables – Rotating

FIRST MADE FOR: H-Power Func

Proprietary to GE Medical Systems

REV

1



GE Medical Systems
Milwaukee, Wi. USA

2308721GSP

Sheet **20** of **20**

The source file for this document was created using *Microsoft WORD 7.0a* and is stored on the CT Products Engineering network.

State: COMPLETE - Production parts may not be defined by this document; it is an uneditable finished draft.
See the GEMS ePDM System to determine the status of this document.