

# RESTRICTED SERVICE NOTE

Dec. 2, 1996

(GE EXCLUSIVE USE ONLY)

Page 1 of 6

FILE TITLE

Signa® Images/Artifacts

SN 65032

## PERIPHERAL SIGNAL ARTIFACTS


**Purpose:** To explain the source of some of the image artifacts that can occur on normally functioning systems. These artifacts are not due to defective hardware or software, but due to limitations in the MR physics and applications.

Peripheral signal artifacts can appear as either a bright spot or as a ribbon of signal smeared through the image. We sometimes refer to these kinds of artifacts as Startifacts (bright spot) or Annefacts (ribbon). Basically both of these artifacts have the same root cause. Signals are generated outside the desired field of view (FOV) and the receiver is able to detect them.

Startifact signals, which appear as a bright star close to the middle of the image, have their origins very far from isocenter. In that non-linear region the FID (free induction decay) signal coming off the RF 180 pulse or from a SAT pulse isn't crushed out and it aliases back into the image. (Aliasing is caused primarily by the fact that our gradient and magnetic coils have finite lengths, resulting in the gradient and  $B_0$  fields becoming increasingly non-linear away from iso-center.) A similar artifact also occurs in 3-D gradient PSDs for essentially the same reason.

Annefact appears in Fast Spin Echo (FSE) scans as a small bright ghost smeared through the image in the phase direction. The origin is also far from iso-center where the gradients are non-linear, and is combined with the  $B_0$  (magnetic) field's inhomogeneity to produce a spatially compressed (therefore very bright) spin echo signal which is aliased back into the image. Uncompensated eddy currents in this area cause phase errors in the compressed signal and smear it through the image.

**Example:** One of the most common complaints occurs during Sagittal (and Coronal) spine scans using a 4-channel phase array coil with Fast Spin Echo (FSE). The artifact occurs when phase encoding is in the S/I (Z-axis) direction and No Phase Wrap is turned on. (The reason for swapping phase and frequency (SPF), to phase encode in the S/I direction, is to reduce CSF flow related artifacts. No Phase Wrap (NP) is selected to prevent the anatomy, which extends outside the Field of View (FOV), from aliasing back into the image in the S/I direction.) The image below has 2 distinct artifacts, a Startifact (oval) and an Annefact (rectangle).



----- Image Parameters -----  
Image ..... 3 / 9                      Tr 3000.00 ms  
Echo ..... 1/1                         Te 102.00 ms  
**Plane            Sagittal**  
**SPF....Yes**  
**PSI) name:    fse**  
  
Slice Thickness. 4.00 mm   NEX 2.00  
Slice Spacing... 1.00 mm   **Coil CTLBOT**  
Scan Time.. 3: 18.10   ETL.. 8  
BandWidth.. 16.0 Khz  
  
Scan Matrix   256 x 256   **FOV 26.00 cm**  
Scan Options: FC ST **NP**  
SAT Pulses : SI

**Workaround** By selecting the receive coils that match the imaging FOV (i.e. LS45, LS56, CS12, etc.), you can lessen the likelihood that you will pickup the peripheral signals that are generated outside the field of view.

# RESTRICTED SERVICE NOTE

Dec. 2, 1996

(GE EXCLUSIVE USE ONLY)

Page 2 of 6

FILE TITLE

Signa® Images/Artifacts


SN 65032

Figures One and Two show simulated and ideal plots of frequency versus distance along the Z-axis for the magnetic and gradient fields. Figure Three combines both these fields, and is important in understanding what is happening when aliasing occurs.

In Figure Three there is a peak in the simulated plot sometimes referred to as the "turn around point". Note that some frequencies are repeated as one moves beyond the turn around point. This repetition in frequencies is the root cause of these artifacts - essentially there are multiple points in space with the same  $B_0$ .

**Note:** The turn around point will vary depending on magnet (i.e. S1, S2, S3, etc.) and gradient coil types. These plots shown here are for discussion purposes only and do not show the full complexity of the issue. Other factors that contribute to the issue are  $B_1$ ,  $T_2^*$  and eddy currents.

FIGURE ONE -  $B_0$  FREQUENCY VS. DISTANCE (Z-axis only)  
[Typical Scanning Range +/- 24cm]



# RESTRICTED SERVICE NOTE

Dec. 2, 1996

(GE EXCLUSIVE USE ONLY)

Page 3 of 6

FILE TITLE

Signa® Images/Artifacts

SN 65032

FIGURE TWO - GRADIENT FIELD FREQUENCY VS. DISTANCE (Z-axis only)





FIGURE THREE - B0 PLUS GRADIENT FIELD FREQUENCY VS. DISTANCE (Z-axis only)



# RESTRICTED SERVICE NOTE

Dec. 2, 1996

(GE EXCLUSIVE USE ONLY)

Page 4 of 6

FILE TITLE

Signa® Images/Artifacts

SN 65032

## OTHER EXAMPLES:

Annefact. 5.4 system. FSE optimization lightened the artifact, but didn't completely remove it.

**Workaround:** Use CS12 and do not swap phase and frequency.




```
|----- Image Parameters -----|
Image ..... 6 / 24      Tr      2962.96 ms
Echo ..... 2/2         Te       102.00ms
Plane  Sagittal       Ti       0.00 ms
SPF....Yes          Contig.. 1          Flip.... 90 deg.
Slice Thickness.    3.00 mm    NEX  2.00
Slice Spacing...  1.00 mm      Coil CTLTOP
ScanTime.. 6:25.29  ETL.. 8    RawNo..9218
BandWidth..32.0 Khz  Ave Sar..0.70  Pk Sar 1.39

Study Number.... 1029          Series 2
Scan Matrix   512 x 256        FOV  24.00 cm
Scan Options: ST PG NP
SAT Pulses :
PSD name:      fse
```

Startifact and Annefact with 3dgre scan. Signal alaising from signals received from coils 3&4.

**Workaround:** Use CS1 2 instead of CTLTOP



```
|----- Image Parameters -----|
Image ..... 28 / 60    Tr       42.00 ms
Echo ..... 1 / 0      Te       15.00 ms
Plane   Axial         Ti       0.00 ms
SPF.... No           Contig.. 1    Flip ..... 5 deg.
Slice Thickness.  2.00 mm  NEX  1.00
Slice Spacing...  0.00 mm  Coil CTLTOP
ScanTime.. 4:21.62  ETL.. 0    RawNo..48157
BandWidth.. 16.0 Khz  Ave Sar..0.47  Pk Sar..0.94
Scan Matrix   256 x 128      FOV  20.00 cm
Scan Options: FC ST
SAT Pulses :
a: Loc:50.0:Width:80.0
p:Loc:6493.6:Width:80.0
PSD name: 3dgrass
```

# RESTRICTED SERVICE NOTE

Dec. 2, 1996

(GE EXCLUSIVE USE ONLY)

Page 5 of 6

FILE TITLE

Signa® Images/Artifacts

SN 65032

Startifact with FSE and 4 coil array;


**Workaround:** Don't swap phase and frequency.



```
|----- Image Parameters -----|
Image ..... 15/28      Tr      3000.00 ms
Echo ..... 1/2        Te       16.00 ms
Plane Sagittal Ti      0.00 ms
SPF Yes    Contig..1    Flip    90 deg.
Acq 1/2    Fract Echo..2 Protocol: T SPINE
Slice Thickness.      3.00 mm    NEX    2.00
Slice Spacing ... 1.00 mm    Coil CTLMID
ScanTime .. 6:30.10    ETL.. 8    RawNo..35335
Bandwidth..32.0 Khz    Ave Sar..0.99    Pk Sar..1.98
Scan Matrix.... 512 x 256    FOV    38.00 cm
Scan Options: ST NP
SAT Pulses : SI
a:Loc:22.7:Width:75.0
PSD name:      fse
```

Large Startifact and Annefact with CTLTOP

**Workaround:** Uses CS12 instead of CTLTOP to eliminate signals from coils 3-4 alaising into image.



```
|----- Image Parameters -----|
Image ..... 7 / 20      Tr      450. 00 ms
Echo ..... 1 / 1        Te       15.00 ms
Plane      Unknown      Ti      0.00 ms
SPF....No    Contig.. 1    Flip.... 20 deg.
Acq.... 1/1    Fract Echo..1 Protocol:CERVICAL SPINE
Slice Thickness.  3.00 mm    NEX    3.00
Slice Spacing...  1.00 mm    Coil CTLTOP
ScanTime.. 8:45.70    ETL.. 0    RawNo..3078
BandWidth..16.0 Khz    Ave Sar..0.58    Pk Sar..1.16
Scan Matrix  256 x 192    FOV    20.00 cm
Scan Options: FC ST
SAT Pulses : s:Loc:36.3:Width:80.0
i:Loc:6480.3:Width:80.0    a:Loc:33.7:Width:80.0
PSD name:      2dfast
```

# RESTRICTED SERVICE NOTE

Dec. 2, 1996

(GE EXCLUSIVE USE ONLY)

Page 6 of 6


FILE TITLE

Signa® Images/Artifacts

SN 65032

Startifact with MEMP (2d SE) scan using CTLTOP

**Workaround:** Use CS12 to prevent alaising from signals picked up by coils 3&4.



```
|----- Image Parameters -----|
Image .....4/ 2      Tr    800.00 ms
Echo .....1/ 1      Te    18.00 ms
Plane .....Unknown   Ti     0.00 ms
SPF....No      Contig.. 1      Flip.... 90 deg.
Acq.... 1/0      Fract Echo..0      Protocol:
Slice Thickness. 5.00 mm      NEX 3.00
Slice Spacing... 1.00 mm      Coil CTLTOP
ScanTime.. 7:50.45   ETL.. 0      RawNo..63491
BandWidth.. 16.0 Khz   Ave Sar. .068   Pk Sar. 1.35
Scan Matrix   256 x 192      FOV 20.00 cm
Scan Options: ST
SAT Pulses :   i: Loc:6509.6:Width:0.0
PSD name: memp
```

Startifact with Fast Spin Echo in the Extremity Coil

**Workaround:** Don't swap phase and frequency.



```
|----- Image Parameters -----|
Image   14/32      Tr   3400.00 ms
Echo    2/2      Te  115.00 ms
Plane   Unknown   Ti    0.00 ms
SPF....Yes      Contig.. 1      Flip.... 90 deg.
Slice Thickness. 4.00 mm      NEX 2.00
Slice Spacing... 1.00 mm      Coil EXTREM
ScanTime.. 7:22.10   ETL.. 8      RawNo..2566
BandWidth.. 16.0 Khz
Scan Matrix  256 x 256   FOV 16.00 cm
Scan Options:      NP
SAT Pulses :      None
PSD name: fse
```

## Acknowledgements

I want to thank Glen Reynolds of Systems Engineering for his time and input.

John A. Johnson  
MR SERVICE ENGINEERING