GE HealthCare

ICE RS-DLP Adapter Installation Manual

9687257-1EN Rev. 3

General Service Documentation.

©2025-2026 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

Reproduction and/or distribution is prohibited.

This manual is a reference for the installation of optional ICE RS to DLP Adapter. The optional ICE RS to DLP Adapter is available in the following products: Vivid Explorer with software v211, Vivid Advanced with software v211.

Revision History

Table -1 Revision History

Revision	Date (YYYY-MM-DD)	Reason for change
Rev. 1	2025-08-26	Initial release.
Rev. 2	2025-10-15	Revision up
Rev. 3	2026-02-24	Update product name

Please verify that you are using the latest revision of this document. Information pertaining to this document is maintained on MyWorkshop. If you need to know the latest revision, contact your distributor, local GE HealthCare Sales Representative or in the USA call the GE HealthCare Ultrasound Clinical Answer Center at 1 800 682 5327 or 1 262 524 5698.

Language Disclaimers

WARNING

- (English)
(EN)
- This Service Manual is available in English only.
- If a customer's service provider requires a language other than English, it is the customer's responsibility to provide translation services.
 - Do not attempt to service the equipment unless this Service Manual has been consulted and is understood.
 - Failure to heed this Warning may result in injury to the service provider, operator or patient from electric shock, mechanical or other hazards.

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Legal Notes

The contents of this publication may not be copied or duplicated in any form, in whole or in part, without prior written permission of GE HealthCare.

GE HealthCare makes no representations or warranties with respect to the information herein. In addition, the information is subject to change without notice. Every precaution has been taken in the preparation of this document. Nevertheless, GE HealthCare assumes no responsibility for errors, omissions, or any damages, including special or consequential, resulting from the use of this information. GE HealthCare will issue updates to this information periodically, as needed. If there are any questions regarding the information contained in this manual, please contact your GE HealthCare Representative.

Damage in Transportation

Follow this procedure if damage is apparent:

1. Write "Damage in Shipment" on ALL copies of the freight or express bill BEFORE delivery is accepted or "signed for" by a GE HealthCare representative or hospital receiving agent.
2. Report the damage to the carrier.
 - Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier.
 - A transportation company will not pay a claim for damage if an inspection is not requested within this 14-day period.

Certified Electrical Contractor Statement - For USA Only

All electrical installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations, and testing shall be performed by qualified GE HealthCare personnel. In performing all electrical work on these products, GE HealthCare will use its own specially trained field engineers. All of GE HealthCare's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE HealthCare equipment shall only utilize qualified personnel (i.e. GE HealthCare field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

Omission and Errors

If there are any omissions, errors or suggestions for improving this documentation, contact the GE HealthCare Ultrasound Global Documentation Group with specific information listing the system type, manual title, part number or direction number, revision number, page number and suggestion details.

Mail the information to:

GE Medical Systems (China) Co., Ltd.

Service Documentation

No.19, Changjiang Road,

WuXi National Hi-Tech Dev.Zone,

214028 Jiangsu

P.R. China

TEL: +86 510 85225888

FAX: +86 510 85226688

GE Vingmed Ultrasound AS

Service Documentation

Strandpromenaden 45

3191 Horten

NORWAY

GE HealthCare employees should use Post-Market Quality Management (PQM) to report service documentation issues.

Table of Contents

Chapter 1 - Optional ICE RS-DLP Adapter Installation Instruction

Optional ICE RS-DLP Adapter Installation Instruction 22

 Introduction 22

 Tools 22

 Manpower 22

 Preparations 22

 Installation Procedure 22

Chapter 1

Optional ICE RS-DLP Adapter Installation Instruction

In this chapter

Optional ICE RS-DLP Adapter Installation Instruction.22

Optional ICE RS-DLP Adapter Installation Instruction

Introduction

This instruction describes the Optional ICE RS-DLP Adapter installation for ultrasound system.

Tools

None.

Manpower

5 minutes.

Preparations

None.

Installation Procedure

READ and UNDERSTAND these instructions thoroughly before proceeding with the installation. Perform each step in sequence and check it off when completes. If a problem occurs after completion, go back through the procedure and check for implementation errors before contacting your Ultrasound Regional Engineer.

CAUTION

Possible operational damage. Failure to strictly follow ESD (Electrostatic Discharge) precautions during this upgrade may cause constant or intermittent operational abnormalities. Strictly follow all precautions.

NOTE

Not all features described in this instruction may be available or cleared for sale in all markets. Please contact your local GE HealthCare Ultrasound representative to get the latest information.

Install the Optional ICE RS-DLP Adapter

The ICE RS-DLP Adapter is designed to provide a connection between the DLP connector to the RS connector, connect the Vivid Explorer and Vivid Advanced to ICE Catheters with RS connector.

1. Connect the ICE RS-DLP Adapter to the DLP port and turn the knob to lock the adapter.

NOTE

Be careful not to connect it in reverse.

Figure 1-1 Connect the ICE RS-DLP Adapter to the DLP port

2. Connect the ICE Catheters with RS connector to the Adapter and flip the connector locking lever to the right.

Figure 1-2 Connect the ICE Catheters with RS connector

**Functional Check**

After installation, please check if the ICE probes can work properly.

C

certified electrical contractor statement - for USA only xvii

D

damage in transportation xvi

F

Functional Check 23

I

Installation 22

L

language disclaimers iv

legal notes xv

M

Manpower 22

O

omission and errors xviii

P

Preparations 22

R

revision history iii

T

Tools 22

transportation
damage xvi