



 Applicable to ISO 13485 QMS

TECHNICAL BULLETIN

(This SN is classified as class A)

Technical Bulletin Unique Identifier	DOC2705276
Title	Service Note for FMI 67940 release
Subject	This Service Note announces the release of FMI 67940
Applicability (e.g., affected product, process)	SIGNA Pioneer SIGNA Premier SIGNA Architect SIGNA Artist Discovery MR750w 3.0T Discovery MR750 3.0T Optima MR450w 1.5T SIGNA Voyager SIGNA Creator/Explorer SIGNA MR355/MR360/MR380 Brivo MR355/Optima MR360 1.5T Signa HDxt
Date	11-July-2022
Originating Business (select one)	☐ GEHC Digital (Imaging and Care Area Solutions, Value Based Care) ☐ Clinical Care Solutions (Ultrasound, Lifecare Solutions) ☒ Imaging (Magnetic Resonance, Molecular Imaging and Computed Tomography, X-ray, Interventional, Surgery, Womens Health, Global Services) ☐ Other: _____
Originating Modality	MR
Details	The FMI67940 provides: • PX29.1_R02 SP01 to SIGNA Pioneer running on PX29.1_R02 • RX29.1_R02 SP01 to SIGNA Premier running on RX29.1_R02 • DV29.1_R02 SP01 to SIGNA Architect/ SIGNA Artist/ Discovery MR750w 3.0T/ Discovery MR750 3.0T/ Optima MR450w 1.5T running on DV29.1_R02 • VX29.1_R02 SP01 to SIGNA Voyager running on VX29.1_R02

* Printed copies are uncontrolled unless otherwise identified *

Before using this document, consult MyWorkshop for the latest revision.

GE Healthcare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.



 Applicable to ISO 13485 QMS

	• SV29.1_R02 SP01 to SIGNA Creator/Explorer running on SV29.1_R02 • HD29.1_R02 SP01 to 1.5T Signa HDxt running on HD29.1_R02 • PX26.4_R02 SP01 to SIGNA Voyager (China Mainland only) running on PX26.4_R02 • SV25.4_R02_SP01 to SIGNA Creator/Explorer/MR355/MR360 running on SV25.4_R02 • HD29.1_R02 with SP01 to 1.5T Signa HDxt running on HD28 • DV28.1_R07 to 1.5T SIGNA Artist running on DV28.1 • VX28.0_R07 to 1.5T SIGNA Voyager running on VX28.0 • SV25.3_R07 to 1.5T SIGNA Creator/Explorer/MR355/MR360/MR380 running on SV25.3 • PX26.3_R05 to 1.5T SIGNA Voyager (China Mainland only) running on PX26.3 • SV20.2_R02 to 1.5T Brivo MR355/Optima MR360 (China Mainland) running on SV20.2 • SV23.2_R02 to 1.5T Brivo MR355/Optima MR360 (China Mainland) running on SV23.2 The FMI67940 provides the software with quality fixes and improvements to sites.
Tracking Number	N/A
Affected Service Publications	N/A
Resolution	See attached
Technical Bulletin POC	Qin Yu (SSO: 503314108)



Applicable to ISO 13485 QMS

Effectivity

SIGNA Pioneer
 SIGNA Premier
 SIGNA Architect
 SIGNA Artist
 Discovery MR750w 3.0T
 Discovery MR750 3.0T
 Optima MR450w 1.5T
 SIGNA Voyager
 SIGNA Creator/Explorer
 SIGNA MR355/MR360/MR380
 Brivo MR355/Optima MR360
 1.5T Signa HDxt

Overview

This Service Note announces the release of FMI 67940 (OnRequest FMI).

Details

1. The release of FMI 67940

The FMI67940 provides:

- PX29.1_R02 SP01 to SIGNA Pioneer running on PX29.1_R02
- RX29.1_R02 SP01 to SIGNA Premier running on RX29.1_R02
- DV29.1_R02 SP01 to SIGNA Architect/ SIGNA Artist/ Discovery MR750w 3.0T/ Discovery MR750 3.0T/ Optima MR450w 1.5T running on DV29.1_R02
- VX29.1_R02 SP01 to SIGNA Voyager running on VX29.1_R02
- SV29.1_R02 SP01 to SIGNA Creator/Explorer running on SV29.1_R02
- HD29.1_R02 SP01 to 1.5T Signa HDxt running on HD29.1_R02
- PX26.4_R02 SP01 to SIGNA Voyager running on PX26.4_R02 (for China Mainland)
- SV25.4_R02 SP01 to SIGNA Creator/Explorer, SIGNA MR355/MR360 running on SV25.4_R02 (for China Mainland)



 Applicable to ISO 13485 QMS

- HD29.1_R02 with SP01 to 1.5T Signa HDxt running on HD28
- DV28.1_R07 to 1.5T SIGNA Artist running on DV28.1
- VX28.0_R07 to 1.5T SIGNA Voyager running on VX28.0
- SV25.3_R07 to 1.5T SIGNA Creator/Explorer/MR355/MR360/MR380 running on SV25.3
- PX26.3_R05 to 1.5T SIGNA Voyager running on PX26.3 (for China Mainland)
- SV20.2_R02 to 1.5T Brivo MR355/Optima MR360 running on SV20.2 (for China Mainland)
- SV23.2_R02 to 1.5T Brivo MR355 /Optima MR360 running on SV23.2 (for China Mainland)

The FMI67940 (OnRequest FMI) provides the software with quality fixes and improvements. Please refer to the “Fixed Issue List” in section4 for details.

2. FMI Kit Contents

FMI67940A Kit PX29.1 R02 SP01 for sites running PX29.1 R02, Part No: 5841744

Level	PN#	Description	Qty
1	5841744	FMI67940A PX29.1_R02 Service Pack SP01	
2	5904171	PX29.1_R02 Service Pack 1	1
2	5917420	FMI67940A-F MR29.1 R02 SP01 FMI Instruction	1

FMI67940B Kit RX29.1 R02 SP01 for sites running RX29.1 R02, Part No: 5841745

Level	PN#	Description	Qty
1	5841745	FMI67940B RX29.1_R02 Service Pack SP01	
2	5904170	RX29.1_R02 Service Pack 1	1
2	5917420	FMI67940A-F MR29.1 R02 SP01 FMI Instruction	1

FMI67940C Kit DV29.1 R02 SP01 for sites running DV29.1 R02, Part No: 5841746

Level	PN#	Description	Qty
1	5841746	FMI67940C DV29.1_R02 Service Pack SP01	
2	5904169	DV29.1_R02 Service Pack 1	1
2	5917420	FMI67940A-F MR29.1 R02 SP01 FMI Instruction	1

FMI67940D Kit VX29.1 R02 SP01 for sites running VX29.1 R02, Part No: 5841747

Level	PN#	Description	Qty
1	5841747	FMI67940D VX29.1_R02 Service Pack SP01	
2	5896559	VX29.1_R02 Service Pack 1	1
2	5917420	FMI67940A-F MR29.1 R02 SP01 FMI Instruction	1

GE Healthcare Technical Bulletin



GE Healthcare

Applicable to ISO 13485 QMS

FMI67940E Kit SV29.1 R02 SP01 for sites running SV29.1 R02, Part No: 5841748

Level	PN#	Description	Qty
1	5841748	FMI67940E SV29.1_R02 Service Pack SP01	
2	5896602	SV29.1_R02 Service Pack 1	1
2	5168656	Common Part:DVD-R	2
2	5917420	FMI67940A-F MR29.1 R02 SP01 FMI Instruction	1

FMI67940F Kit HD29.1 R02 SP01 for sites running HD29.1 R02, Part No: 5914658

Level	PN#	Description	Qty
1	5914658	FMI67940F HD29.1_R02 Service Pack SP01	
2	5896601	HD29.1_R02 Service Pack 1	1
2	5917420	FMI67940A-F MR29.1 R02 SP01 FMI Instruction	1

FMI67940G Kit PX26.4 R02 SP01 for sites running PX26.4 R02, Part No: 5914666

Level	PN#	Description	Qty
1	5914666	FMI67940G PX26.4_R02 Service Pack SP01	
2	5904557	PX26.4_R02 Service Pack 1	1
2	5917421	FMI67940G/H PX26.4_R02 SP01 & SV25.4_R02 SP01 FMI Instruction	1

FMI67940H Kit SV25.4 R02 SP01 for sites running SV25.4 R02, Part No: 5914764

Level	PN#	Description	Qty
1	5914764	FMI67940H SV25.4_R02 Service Pack SP01	
2	5903723	SV25.4_R02 Service Pack 1	1
2	5168656	Common Part:DVD-R	2
2	5917421	FMI67940G/H PX26.4_R02 SP01 & SV25.4_R02 SP01 FMI Instruction	1

FMI 67940J Kit HD29.1_R02 with SP01 for sites running HD28, Part No: 5914765

Level	PN#	Description	Qty
1	5914765	FMI67940J HD29.1_R02 for HD28 sites	
2	5881227	MR HD29.1_R02 Software Collector	1
3	5881262	HD29.1_R02 USB Installer with Label	1

* Printed copies are uncontrolled unless otherwise identified *

Before using this document, consult MyWorkshop for the latest revision.

GE Healthcare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.

Approved

GE Healthcare Technical Bulletin



GE Healthcare

Applicable to ISO 13485 QMS

3	5861614	Backup USB with Label	1
3	5824702	EULA License Information	1
3	5896601	HD29.1_R02_software ServicePack01	1
2	5165713	DVD Media for Linux OS - 50 pack	0.1
2	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1
2	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
2	6339528-299	HOPE OM	1
2	5860821-399	MR 29.1 Privacy and Security manual	1
2	5851463-199	29.0+ Cybersecurity Changes	1
2	5917422	FMI67940J HD29.1_R02 FMI Instruction	1

FMI67940K Kit DV28.1 R07 for sites running DV28.1, Part No: 5914766

Level	Part #	Description	Qty
1	5914766	FMI67940K DV28.1_R07	
2	5917687	Artist IPM DV28.1_R07 Software Collector	1
3	5917688	Artist IPM DV28.1_R07 USB Installer with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	46-302200P7	SERIAL NO MODEL NO DESCRIPTION and MANUFACTURED DATE 2.25 X 1.40 PLATE	1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
2	5917423	FMI67940K/L/M MR28 R07 FMI Instruction	1

FMI67940L Kit VX28.0 R07 for sites running VX28.0, Part No: 5914767

Level	Part #	Description	Qty
1	5914767	FMI67940L VX28.0_R07	
2	5904578	Voyager IPMAir VX28.0_R07 Software Collector	1
3	5917501	Voyager IPMAir VX28.0_R07 USB Installer with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
2	5917423	FMI67940K/L/M MR28 R07 FMI Instruction	1

FMI67940M Kit SV25.3 R07 for sites running SV25.3, Part No: 5914768

Level	Part #	Description	Qty
-------	--------	-------------	-----

* Printed copies are uncontrolled unless otherwise identified *

Before using this document, consult MyWorkshop for the latest revision.

GE Healthcare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.

GE Healthcare Technical Bulletin



GE Healthcare

Applicable to ISO 13485 QMS

1	5914768	FMI67940M SV25.3_R07	
2	5903727	SV25.3_R07 Software Collector	1
3	5917554	SV25.3_R07 USB Installer with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
2	5917423	FMI67940K/L/M MR28 R07 FMI Instruction	1

FMI 67940N Kit PX26.3_R05 for sites running PX26.3, Part No: 5914769

Level	PN#	Description	Qty
1	5914769	FMI67940N PX26.3_R05	
2	5903761	Voyager PX26.3_R05 Software Collector	1
3	5904585	Voyager PX26.3_R05 USB with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
2	5917424	FMI67940N PX26.3_R05 FMI Instruction	1

FMI67940P Kit SV20.2_R02 for sites running SV20.2 Part No: 5877355

Level	Part #	Description	Qty
1	5877355	FMI67940P SV20.2_R02	
2	5904589	SV20.2_R02 software Collector	1
3	5917680	SV20.2_R02 USB Installer with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
2	5920203	FMI67940P/Q SV20.2_R02 & SV23.2_R02 FMI Instruction	1

FMI67940Q Kit SV23.2_R02 for sites running SV23.2 Part No: 5877357

Level	Part #	Description	Qty
1	5877357	FMI67940Q SV23.2_R02	
2	5904590	SV23.2_R02 software Collector	1
3	5917681	SV23.2_R02 USB Installer with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1

* Printed copies are uncontrolled unless otherwise identified *

Before using this document, consult MyWorkshop for the latest revision.

GE Healthcare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.



 Applicable to ISO 13485 QMS

3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
2	5920203	FMI67940P/Q SV20.2_R02 & SV23.2_R02 FMI Instruction	1

3. Training

For sites upgrading from HD28 to HD29.1, 0.5 day of applications training is required.

4. Fixed Issues List

SPRs fixed in PX29.1_R02 SP01/RX29.1_R02 SP01/DV29.1_R02 SP01/VX29.1_R02 SP01/SV29.1_R02 SP01/HD29.1_R02 SP01/PX26.4_R02 SP01/SV25.4_R02 SP01

SPR#	Brief Description
HCSDM00685848	Scan DISCO, end exam before all images generated, the DT and time on image are incorrect
HCSDM00679478	System not doing imaging calibration when coil calibration happens.
HCSDM00679477	Intermittently cooling start error occurred during the scan and scan stopped
HCSDM00679475	Acqs Before Pause missing for 2D FSE application
HCSDM00679354	Distcorr, maxtensor, muse option keys are not compatible
HCSDM00679300	Change PV coil driver coil_type from 1 to 0 due to PV coil is Non MEMS coil
HCSDM00679280	Cube T2 flair shading on pediatric head
HCSDM00679279	DISCO FatSat generates garbage images
HCSDM00679175	Scanner was not scannable due to nbPowerMonType in MRconfig.cfg was set to 0
HCSDM00678136	Multi Viewer not responding and filming not responding
HCSDM00677870	3DFRFSE not part of scantools
HCSDM00675534	Grafidy's data processing code does not initialize the variables after malloc()
HCSDM00681500	Image cut and dark banding artifact observed for 32Ch Body LT Bone coil mode
HCSDM00681114	Enable power save mode for SSSPS-Lite and XFD-PS
HCSDM00680224	Dolphin card and PERC removal changes

GE Healthcare Technical Bulletin



GE Healthcare

Applicable to ISO 13485 QMS

HCSDM00680264	MXH930 code support and backwards compatibility for IXH610
HCSDM00679430	Remove Run Diagnostics and Stop Diagnostics from GI > ICN Configuration
HCSDM00682491	Fix issue with vre_dacq restart
HCSDM00681928	IXH610 system does not report correct error when link is down
HCSDM00682559	NFS fails to read dacq channels on vre_dacq initialization
HCSDM00681912	Host did not update ermes database properly in .PC service pack
HCSDM00682471	Coil changes in a ".TAR" file for Service Pack installation
HCSDM00685695	Image Header isn't launching
HCSDM00685839	Dolphin1 ICN memory corruption on bootup
HCSDM00679287	Fix Mobile AE Title cannot be set properly.
HCSDM00686277	SSFSE refocus flip angle enable on UI using typein ssfse_refocus
HCSDM00686899	Dolphin changes causing ICN crash

SPRs fixed in HD29.1_R02 and SP01

SPR#	Brief Description
HCSDM00657412	Display of po+H57:H167p-up message during the change in Patient position
HCSDM00666198	Image artifacts when Slice ASSET is used with NoSlabWrap Factor larger than 1
HCSDM00663091	incorrect error message is observed when user tries to add an unsupported satband
HCSDM00659647	ARES magnet - disable in R02 branch
HCSDM00654947	build infrastructure updates for MR29.1_R02
HCSDM00657759	DWI RPG ORIG images not generated completely
HCSDM00656387	RF Amp addition for Signa Artist & Signa Voyager
HCSDM00641443	[Point Release/BaliA] Need to add system tools support for ARES magnet
HCSDM00662970	Stream spec update to synchronize error messages from mainline & pipeline update to v3.0.1
HCSDM00657075	[MR29.1_R02] enable CSA16KW configuration for RAP tool
HCSDM00648225	Bridge Type i18n translation.
HCSDM00662997	MPRAGE Protocol edit issue with incorrect error message
HCSDM00661494	Refactoring GRAF3 commands for CAM to ICE upgrade

*** Printed copies are uncontrolled unless otherwise identified ***

Before using this document, consult MyWorkshop for the latest revision.

GE Healthcare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.



Applicable to ISO 13485 QMS

HCSDM00656168	The SV25.3 OM can't install on 29.1_R01 SW.
HCSDM00654871	[Atlantis/Qin R02] FSE + DynTG + FS pulsegen failed
HCSDM00662364	[eDelivery-PS] Support additional locale for Portuguese and Spanish
HCSDM00654555	[R02] No tuning CFH waveform appears on the second shimming Volume during breast imaging
HCSDM00654940	[Bali-A] AGP IRF driver code not returning correct q_phase values for pcal
HCSDM00655571	S/W tools support for new product addition : Voyager AIR
HCSDM00655800	[Motus Table on Premier - Impacted by Pioneer] Guide install & Config update for Premier with Motus because of Pioneer bridge
HCSDM00656160	[Double promote] Table configuration is marked as RED and present a different value from the saved value in SaveINFO media.
HCSDM00656163	[Double promote] Empty Dock install page in ICW
HCSDM00659568	[R02] Enable HOEC protocols in BRM config (HOPE, Explorer/Creator)
HCSDM00655333	Host Tools Support for ARES Magnet in Bali SW
HCSDM00654430	[Atlantis/Qin R02] AIRx Brain OML landmark is slightly different.
HCSDM00654431	[Atlantis/Qin R02] FGRE Cardiac TimeCourse DL Recon Artifact
HCSDM00655794	Pioneer and Voyager Motus ICW changes to support - Need to update the Patient Table Install check list similar to Premier MOT
HCSDM00662783	[MR29.1 R02] SPT result file does not contain details for HeadSNR on nonHCI system
HCSDM00654872	[Atlantis/Qin R02] Poor oZTEo uniformity for Spine anatomy on Artist
HCSDM00656196	Error in restoring Voyager channel config value from PX26.2 to MR29.1_R01
HCSDM00656256	[MR29.1_R02] Bali-A Security TR OS packages upgrades
HCSDM00654819	CSD to show RAP tool on right side window when clicked from Calibration tab RAP Tool
HCSDM00655666	[Atlantis R02] Remove FIESTA TC protocol from GE library for all 3T configurations
HCSDM00654818	EPT page in ICW Maintenance mode is missing a Check Status button
HCSDM00654821	Strong artifact of DWI CMB image due to wrong HOEC calibration data
HCSDM00655939	Adding 450w non-gem support
HCSDM00654556	[R02] DISCO FatSat generates garbage images
HCSDM00654831	HOS upgrade mode changes for ProDrive Shim Supply in Premier
HCSDM00654833	iLinq error page is not displayed when iLinq is configured and network not

*** Printed copies are uncontrolled unless otherwise identified ***

Before using this document, consult MyWorkshop for the latest revision.

GE Healthcare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.

GE Healthcare Technical Bulletin



GE Healthcare

Applicable to ISO 13485 QMS

	available.
HCSDM00656401	[Bali A] Inaccurate SAR Display limit when SAR derating
HCSDM00656396	Grad Spec changes in about MR
HCSDM00659536	[Atlantis/Qin] - TDEL is 0 for all wash-in phases in DISCO with SmartPrep
HCSDM00657095	CSE 3.1.21.1 RFV1 V1 Coreload Integration
HCSDM00657664	3D TOF MgT protocol grays out
HCSDM00657106	[double promote]Platform IRD TDR code implement
HCSDM00658020	ICW changes for Mulan & HD
HCSDM00656416	[29.1 R02] Update max. B1rms limits for Pioneer-Motus
HCSDM00661285	[ATLANTIS FRT] Scan crash (found from Kibana)
HCSDM00660744	[Atlantis/Qin R02] Ghosting observed on 0.5 NEX Coronal 2D FSE FS scans on Abdomen
HCSDM00657580	APx addition in scantools for Optima 450w nongem
HCSDM00659942	Increase CV12 max from 120 to 180 for CUBE PD for 1.5T same as 3T
HCSDM00656859	Shading Artifacts of PSMDE on 1.5T
HCSDM00659535	Image cut off at back for Lspine exam on 750w
HCSDM00659503	Support Routine generate_fliptable failure with SSFSE PSD
HCSDM00660355	Defensive fix option key not getting restored
HCSDM00657468	"Total Phase" is wiped out for disco
HCSDM00662718	IIS check error on Voyager Air config
HCSDM00662387	[Atlantis FRT] ifcc crash in Atlantis build (found from Kibana Dashboard)
HCSDM00660610	Demote CSA 16KW for R02 from Signa Artist & Signa Voyager
HCSDM00659502	Gradient Fault on Ax T2 FSE scan on c-spine scan with ART and high resolution
HCSDM00659723	90 flip with FGRE does not provide T2* contrast.-Protocol update
HCSDM00660235	Extend Lspine CUBE T2 mode flipmod(130 80 130 80) to pelvis
HCSDM00659306	[Undo] APx addition in scantools for Optima 450w nongem
HCSDM00657451	[eDelivery-PS] Handling RIaS logger file permission issues
HCSDM00659140	Enabling MIP & MPR Inline PP Tasks for MENSA Application
HCSDM00662644	[R02] SSFSE TRF with Multinex Pause mode get blank images
HCSDM00661084	cube stir under voltage on explorer

*** Printed copies are uncontrolled unless otherwise identified ***

Before using this document, consult MyWorkshop for the latest revision.

GE Healthcare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.

GE Healthcare Technical Bulletin



GE Healthcare

Applicable to ISO 13485 QMS

HCSDM00661338	Device disconnected from IPM cloud after software upgrade
HCSDM00659359	3DMERGE with Spatial Sat has gradient over-current fault
HCSDM00663310	[Qin II] Update parallel imaging file for adding 8ch Flex L/M/S coils to Voyager system
HCSDM00660979	[eDelivery-PS] - tool tip i18n translations
HCSDM00659719	DST Flickering Issue cause SCN crash
HCSDM00661041	[Atlantis/Qin R02] oZTEo MR Bone uniformity issue seen on Hip and Spine on 3T, and Head/Neck on 1.5T
HCSDM00660220	[Atlantis FRT]Blue Screen SCN crash
HCSDM00661100	Zero size pdf file created in USB during ProtocolExchange export for gep/adult wholebody Multistation protocol
HCSDM00661346	[MR29.1_R02] promote ge protocol for 450w non-GEM
HCSDM00658242	[Double promote] Month and Day of Date/Time tab is showing wrong values
HCSDM00662818	ICW Issue on CAM to ICE Upgrade
HCSDM00662917	[MR29.1 R02] On a staged system, IIS may not see Sysgain results as valid
HCSDM00657200	[eDelivery-PS] Logger and Notification Defect fixes and I18n translation strings implementation
HCSDM00662361	Option key migration issue fix
HCSDM00657203	[eDelivery]RSD_INFO.xml date format mismatch,and entries in decline group are not proper. Notification tooltip is not proper
HCSDM00657201	[eDelivery-PS] RIaS : i18n issue fix in splash screen
HCSDM00657205	[eDelivery-PS]RIaS: Exception handling incase of save info failure during full software download
HCSDM00657206	[eDelivery-PS]RIaS: Replace hostrestart with reboot in RIaS workflow
HCSDM00657202	[eDelivery-PS]build comparison was not proper while restoring rsd_info.xml
HCSDM00657204	[eDelivery-PS] Paused download is not resuming intermittently on reboot.
HCSDM00664998	[ATLANTIS FRT]scn crash(from Kibana)
HCSDM00663549	build infra : software version files promote
HCSDM00665042	ICW upgrade mode issue after eDelivery
HCSDM00664009	Product Configuration dat file update for Pioneer & Voyager

*** Printed copies are uncontrolled unless otherwise identified ***

Before using this document, consult MyWorkshop for the latest revision.

GE Healthcare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.

GE Healthcare Technical Bulletin



GE Healthcare

Applicable to ISO 13485 QMS

HCSDM00665205	System loses network connectivity when the mobile site host name is different from system host name.
HCSDM00664280	Scan Crash with BRAVO application when slices per slab entered as 2048
HCSDM00666220	NotifyEvent messages are not getting sent to the ExsiClient after setting NotifyEvent all=on
HCSDM00664368	Protocol update required for 450W Non Gem- Static Configuration
HCSDM00659574	VX29.1 R01 Service Pack configure file promote
HCSDM00660281	[Qin II] Promote the Main Line coil driver/coil DB to MR29.1_R01 branch
HCSDM00651603	[IO]Auto Parameters Promote for 1.5T-withdrawn-nothing changed
HCSDM00652362	Cabmon3 IO: Continuously monitor Platform ICC flow mode
HCSDM00659096	Correct graf3 command order
HCSDM00656194	Soften 33ch limitation to 36ch for Voyager/Pioneer
HCSDM00654835	feetFirstPDHSEnabled = "n0" shown in Premier MRconfig.cfg
HCSDM00659903	[R02] Change SPT order to put CNT first
HCSDM00656914	[Bali-a] Expand receiver TNS counters from 16 to 32 bits
HCSDM00656915	Mitigate 5V_A Spurs
HCSDM00659361	MAGiC T2 value accuracy enhancement with satCrusher optimization
HCSDM00659120	[R02] Enable HOECC for 1.5T systems
HCSDM00662782	Add additional build label information in the Release Notice email
HCSDM00660239	[MR29.1_R02] Bali-A Security TR - Disable Firefox login password saving feature
HCSDM00661037	[Atlantis FRT] 32ch Torso Coil acceleration limits is too low
HCSDM00657207	[eDelivery-PS] Uninstallation of service pack from GI and update rsd_info.xml

SPRs fixed in DV28.1_R07, VX28.0_R07, SV25.3_R07, PX26.3_R05, SV20.2_R02, SV23.2_R02

SPR#	Brief Description
------	-------------------

* Printed copies are uncontrolled unless otherwise identified *

Before using this document, consult MyWorkshop for the latest revision.

GE Healthcare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.



 Applicable to ISO 13485 QMS

HCSDM00680369	Scan DISCO, end exam before all images generated, the DT and time on image are incorrect
HCSDM00679356	Update nbPowerMonType from 0 to 1 when upgrading from SV25.2 to MR28
HCSDM00678609	PVA coil was not allowed to scan on non-MEMS config
HCSDM00676257	Cube T2 flair shading on pediatric head
HCSDM00690269	DISCO FatSat generates garbage images, turn off run_eval_from_get_flag for DISCO
HCSDM00679290	MDE shading image quality issue.
HCSDM00679291	No tuning waveform appears on the second shimming Volume during breast imaging
HCSDM00679678	Multi Viewer not responding and filming not responding
HCSDM00663418	Enable power save mode for SSSPS-Lite and XFD-PS on Mulan SKYLake
HCSDM00678650	Blue Screen
HCSDM00682331	Remove Run Diagnostics and Stop Diagnostics from GI > ICN Configuration
HCSDM00682383	Dolphin card and PERC removal changes merge to VX28 branch
HCSDM00682415	MXH930 code support and backwards compatibility for IXH610
HCSDM00686631	GEN7-DL-P ICN Support for VX28
HCSDM00686613	High Channel Scans failing in PRICE software with GEN7-DL
HCSDM00688153	AW Unagi release integration into VX28 branch
HCSDM00690019	ICN-ICE Datapath diag failed
HCSDM00691175	ICN1 Hard Disk Stress the tools for this test case is fail to use
HCSDM00688913	Enable ICN Firmware upgrade in Gen7 and Gen7-DL
HCSDM00691165	ICN Image Reconstruction iterations=1 When run this test fail
HCSDM00692065	Update fixed R1/R2 value in signal acquisition for 1.5T system with MEGASWITCH
HCSDM00691488	vreicetestcl doesn't work on PRICE system
HCSDM00679289	Voyager mobile couldn't Connect with PACs due to AE title out of sync
HCSDM00675387	Auto transfer by series cannot transfer images after moving the cradle
HCSDM00695825	Table configuration is showing red after restore info in MR380 system
HCSDM00682333	update R07 mapping to support OM installation

5. Installation

The installation procedure could be found in the FMI Instructions of each FMI kit.

File Revision

GE Healthcare Technical Bulletin



GE Healthcare

Applicable to ISO 13485 QMS

Rev	Page	Description of Content Changed	Reason for Change or Change Control Number	Author Name (SSO)
1	All	Initial Release	Initial Release	Qiao Zhanpeng (503211725)
2	All	Added content for FMI67940J, K, L, M, N, P and Q.	For FMI67940 phase 2 release.	QIN YU (503314108)

*** Printed copies are uncontrolled unless otherwise identified ***

Before using this document, consult MyWorkshop for the latest revision.

GE Healthcare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.