

GE HealthCare
Technical Bulletin



Technical Bulletin Unique Identifier	DOC2842958
Title	Service Note for FMI 67943 release
Subject	This Service Note announces the release of FMI 67943
Applicability (e.g., affected product, process)	1.5T SIGNA Voyager 1.5T SIGNA Creator/Explorer 1.5T SIGNA MR355/MR360 1.5T SIGNA HDxt 1.5T SIGNA Artist 1.5T Optima MR450w 3.0T SIGNA Architect 3.0T Discovery MR750w 3.0T Discovery MR750 3.0T SIGNA Pioneer 3.0T SIGNA Premier
Date	02-Jun-2023
Relevant Product Business (select one)	☒ Imaging ☐ Patient Care Solutions ☐ Ultrasound ☐ Other: _____
Originating Modality	MR

*** Printed copies are uncontrolled unless otherwise identified ***

Before using this document, consult MyWorkshop for the latest revision.

GE HealthCare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.

GE HealthCare

Technical Bulletin



Details	This Business Mandatory FMI67943 provides software MR29.1_R04 to IB sites: • VX29.1_R04 to SIGNA Voyager running on VX29.1 • SV29.1_R04 to SIGNA Creator/Explorer running on SV29.1 • HD29.1_R04 to SIGNA HDxt running on HD29.1 • DV29.1_R04 to Architect/MR750w/MR750/Artist/MR450w running on DV29.1_R01/R02 • PX29.1_R04 to SIGNA Pioneer running on PX29.1_R01/R02 • RX29.1_R04 to SIGNA Premier running on RX29.1 • PX26.4_R04 to SIGNA Voyager running on PX26.4 • SV25.4_R04 to SIGNA Creator/Explorer, SIGNA MR355/MR360 running on SV25.4 • DV28.4_R04 to SIGNA Architect running on DV28.4 • DV27.3_R04 to SIGNA Architect running on DV27.3 • DV26.4_R04 to Discovery MR750 & MR750W running on DV26.4 Details of the quality fixes please refer to the “Unique Issue List” in section 4.
Tracking Number	N/A
Affected Service Publications	N/A
Resolution	See attached
Technical Bulletin Owner	Qin Yu (SSO: 503343456)

(This SN is classified as class A)

*** Printed copies are uncontrolled unless otherwise identified ***

Before using this document, consult MyWorkshop for the latest revision.

GE HealthCare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.

GE HealthCare

Technical Bulletin



Effectivity:

1.5T SIGNA Voyager
 1.5T SIGNA Creator/Explorer
 1.5T SIGNA MR355/MR360
 1.5T SIGNA HDxt
 1.5T SIGNA Artist
 1.5T Optima MR450w
 3.0T SIGNA Architect
 3.0T Discovery MR750w
 3.0T Discovery MR750
 3.0T SIGNA Pioneer
 3.0T SIGNA Premier

Overview:

This Service Note announces the release of FMI 67943 (Mandatory FMI).

1. The release of FMI 67943:

This Business Mandatory FMI67943 provides software MR29.1_R04 to IB sites:

FMI 67943 ABC & FGH:

- VX29.1_R04 to SIGNA Voyager running on VX29.1
- SV29.1_R04 to SIGNA Creator/Explorer running on SV29.1
- HD29.1_R04 to SIGNA HDxt running on HD29.1
- DV29.1_R04 to Architect/MR750w/MR750/Artist/MR450w running on DV29.1_R01/R02
- PX29.1_R04 to SIGNA Pioneer running on PX29.1_R01/R02
- RX29.1_R04 to SIGNA Premier running on RX29.1

FMI 67943 DE & L:

- PX26.4_R04 to SIGNA Voyager running on PX26.4
- SV25.4_R04 to SIGNA Creator/Explorer, SIGNA MR355/MR360 running on SV25.4
- DV26.4_R04 to Discovery MR750 & MR750W running on DV26.4

FMI 67943 J:

- DV28.4_R04 to SIGNA Architect running on DV28.4
- DV27.3_R04 to SIGNA Architect running on DV27.3

Details of the quality fixes please refer to the “Unique Issue List” in **Section 4. Fixed Issue List**.

GE HealthCare

Technical Bulletin



2. FMI kit Content:

This FMI has 10 kits to address different software versions on multiple products.

FMI 67943A Kit VX29.1_R04 Part Number: 5914773

Level	Part Number	Description	Qty
1	5914773	FMI67943A VX29.1_R04	
2	5932876-1	MR VX29.1_R04 Software Collector	1
3	5932877-1	VX29.1_R04 USB Installer with Label	1
3	5824702	EULA License Information	1
3	5861614	Backup USB with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1
2	5939713	FMI67943 MR29.1 R04 FMI Instruction for Global	1

FMI 67943B Kit SV29.1_R04 Part Number: 5914774

Level	Part Number	Description	Qty
1	5914774	FMI67943B SV29.1_R04	
2	5932876-3	MR SV29.1_R04 Software Collector	1
3	5932877-3	SV29.1_R04 USB Installer with Label	1
3	5824702	EULA License Information	1
3	5861614	Backup USB with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1
2	5939713	FMI67943 MR29.1 R04 FMI Instruction for Global	1

FMI 67943C Kit HD29.1_R04 Part Number: 5914895

*** Printed copies are uncontrolled unless otherwise identified ***

Before using this document, consult MyWorkshop for the latest revision.

GE HealthCare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.

GE HealthCare

Technical Bulletin



Level	Part Number	Description	Qty
1	5914895	FMI67943C HD29.1_R04	
2	5932876	MR HD29.1_R04 Software Collector	1
3	5932877	HD29.1_R04 USB Installer with Label	1
3	5824702	EULA License Information	1
3	5861614	Backup USB with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1
2	5939713	FMI67943 MR29.1 R04 FMI Instruction for Global	1

FMI 67943D Kit PX26.4_R04 Part Number: 5914896

Level	Part Number	Description	Qty
1	5914896	FMI67943D PX26.4_R04	
2	5932876-2	MR PX26.4_R04 Software Collector	1
3	5932877-2	PX26.4_R04 USB Installer with Label	1
3	5824702	EULA License Information	1
3	5861614	Backup USB with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1
2	5939716	FMI67943 FMI Instruction for China	1

FMI 67943E Kit SV25.4_R04 Part Number: 5914897

Level	Part Number	Description	Qty
1	5914897	FMI67943E SV25.4_R04	
2	5932876-4	MR SV25.4_R04 Software Collector	1
3	5932877-4	SV25.4_R04 USB Installer with Label	1

*** Printed copies are uncontrolled unless otherwise identified ***

Before using this document, consult MyWorkshop for the latest revision.

GE HealthCare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.

GE HealthCare

Technical Bulletin



3	5824702	EULA License Information	1
3	5861614	Backup USB with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1
2	5939716	FMI67943 FMI Instruction for China	1

FMI 67943F Kit DV29.1_R04 Part Number: 5914898

Level	Part Number	Description	Qty
1	5914898	FMI67943F DV29.1_R04	
2	5939713	FMI67943 MR29.1 R04 FMI Instruction for Global	1
2	5932876-11	DV29.1 Software Collector	1
3	5824702	EULA License Information	1
3	5861614	Backup USB with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1
3	5932877-11	DV29.1 USB Installer with Label	1

FMI 67943G Kit PX29.1_R04 Part Number: 5914899

Level	Part Number	Description	Qty
1	5914899	FMI67943G PX29.1_R04	
2	5939713	FMI67943 MR29.1 R04 FMI Instruction for Global	1
2	5932876-7	PX29.1 Software Collector	1
3	5824702	EULA License Information	1
3	5861614	Backup USB with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1

*** Printed copies are uncontrolled unless otherwise identified ***

Before using this document, consult MyWorkshop for the latest revision.

GE HealthCare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.

Released

GE HealthCare

Technical Bulletin



3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1
3	5932877-7	PX29.1_R04 USB Installer with Label	1

FMI 67943H Kit RX29.1_R04 Part Number: 5914900

Level	Part Number	Description	Qty
1	5914900	FMI67943H RX29.1_R04	
2	5939713	FMI67943 MR29.1 R04 FMI Instruction for Global	1
2	5932876-5	RX29.1 Software Collector	1
3	5824702	EULA License Information	1
3	5861614	Backup USB with Label	1
3	5932877-5	RX29.1_R04 USB Installer with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1

FMI 67943J Kit DV28.4_R04 & DV27.3_R04 Part Number: 5914903

Level	Part Number	Description	Qty
1	5914903	FMI67943J DV28.4 R04 and DV27.3 R04	
2	5939714	FMI67943 FMI Instruction for China Architect only	1
2	5932876-13	DV28.4 Software Collector	1
3	5824702	EULA License Information	1
3	5861614	Backup USB with Label	1
3	5932877-13	DV28.4 USB Installer with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1
2	5932876-14	DV27.3 Software Collector	1
3	5824702	EULA License Information	1
3	5932877-14	DV27.3 USB Installer with Label	1
3	5861614	Backup USB with Label	1

*** Printed copies are uncontrolled unless otherwise identified ***

Before using this document, consult MyWorkshop for the latest revision.

GE HealthCare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.

Released

GE HealthCare

Technical Bulletin



3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1
FMI 67943L Kit DV26.4_R04 Part Number: 5914908			
Level	Part Number	Description	Qty
1	5914908	FMI67943L DV26.4 R04	
2	5939716	FMI67943 FMI Instruction for China	1
2	5932876-15	DV26.4 Software Collector	1
3	5824702	EULA License Information	1
3	5861614	Backup USB with Label	1
3	5932877-15	DV26.4_R04 USB Installer with Label	1
3	5165713	DVD Media for Linux OS - 50 pack	0.1
3	5661793	GE EQUIPMENT MAINTENANCE STICKER	1
3	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES, PRINTED BY DUPLITECH	1

3. Special Instructions:

- 1) Before on-site visit, please check the site **Product Software version** to make sure get the correct FMI kit.
- 2) FMI67943 ABC & FGH could be deployed via Remote Software Download (RSD) to sites as many as possible where RSD is permitted. Physical kits will be provided for sites where RSD is not performed, for example due to regulation limit or site environment limit.

3) For sites where ordered Option key Air Recon DL M7009SA, contact OLE to remove existing option key code on eLicense, then regenerate new option key code and download for site.

AIR Recon DL	M7009SA	airrcndl2d, ARDLdwepi
--------------	---------	--------------------------

NOTE: M7009SA has added new content in this upgrade, if the site does not remove the original M7009SA key code on the eLicense, it will not be able to import the new M7009SA key due to repeated naming, and OLE has the corresponding removal access.

The difference before and after the update on eLicense:

Original M7009SA

GE HealthCare

Technical Bulletin



M7009SA	AIR Recon DL M7009SA(MR)	N	airrcnd12d permanent 0 	Permanent
---------	--------------------------	---	--	-----------

New M7009SA

M7009SA	AIR Recon DL M7009SA(MR)	N	ARDLdwepl permanent 0  airrcnd12d permanent 0	Permanent
---------	--------------------------	---	---	-----------

- 4) If the customer needs to reverse the breast image, please let them contact Sales to set a zero-price order, after this install the follow option key.

Breast Image Orientation	M7030AA	BreastImgO
--------------------------	---------	------------

4. Fixed Issue List:

SPR # if applicable	Brief DescriptionD
HCSDM00726046	Merge PACS is intermittently having issues with Image Orientation (Patient) DICOM tag
HCSDM00712391	AutoVoice failed for Disco 1st phase with SmartPrep
HCSDM00712392	MAVRIC STIR minseqrfamp error with weight >160kg
HCSDM00712393	xTG failed issue for newborn baby
HCSDM00712400	Disable DISCO run_eval_from_get_flag for all 1.5T systems
HCSDM00712402	DW PROP DUO shading on Explorer
HCSDM00712501	Vibrant sp sat scan time fix
HCSDM00712561	Operator can't record customer voice
HCSDM00712563	Nothing is displayed in the Intensity Filter selection box after SSFSE setup
HCSDM00712564	setNavShimSupply crash when using NAV shim PS
HCSDM00712641	MP-RAGE PROMO image shift after the previous MP-RAGE PROMO scan aborted
HCSDM00712643	Overrange on MP-RAGE brain protocol
HCSDM00712644	PSD timeout when adding SSRF imaging option to Ax VIBRANT high spatial resolution GE protocol
HCSDM00712645	3dtof gradient trip due to missing ISI on ICE
HCSDM00712668	CUBE fineline reduction
HCSDM00712670	Slice ASSET for sagittal Vibrant has image stretch/shrink
HCSDM00712671	For FSE, Dynamic TG(CV38) & Blurring Cancellation (CV7) on at same time give noisy images

GE HealthCare

Technical Bulletin



SPR # if applicable	Brief DescriptionD
HCSDM00712673	3dtof doesn't show MT effect
HCSDM00712682	SSFSE abd T2 shading mitigation with PURE tuning for 3T
HCSDM00712688	Regular 3D SPGR/GRE show wrong slice resolution on the UI
HCSDM00712759	Update rdb_hdr_raw_pass_size in pfile header when CoilSwitch is ran in prescan
HCSDM00712798	3dradial ZTE scan failed with Resp. Trigger
HCSDM00712800	Prescan Cal PSD Encoding Center Shift
HCSDM00712966	Frequency Direction Saved in the Protocol is getting swapped in the scanning Session.
HCSDM00712967	Exsi Crash and Scan crash in SetCollectPutCvs
HCSDM00712968	ZTE protocol shows up advisory message for slice thickness during protocol load
HCSDM00713057	Intleave needs to be disabled in Spacing UI for MUSE
HCSDM00713058	C045 EPI protocol download failed
HCSDM00713060	cfx/y/zfs are not 3.21 when ART is ON with DW-EPI with XRMW
HCSDM00713061	epi2 timeout during HAST
HCSDM00713063	type-in PSD epialtoff and epiRTalt off fail to set pepolar =0 for hyperband scan
HCSDM00713064	dacq_adjust update for DW-EPI with Hyperband
HCSDM00713065	Severe A/P ghosting with Asset on GRE-EPI
HCSDM00713082	banding artifact worsened at site after 29.0 to 29.1 update. same issue in MR30.0
HCSDM00713350	IDEAL-IQ Echo Phase Correction Robustness Improvement
HCSDM00713353	HCSDM00: HAST:data acquisition error with A193 3D SAG VIBRANT Unilateral base ASPIR ASSET no turbo protocols
HCSDM00713444	MPRAGE PROMO download failed with minimum Recovery Time
HCSDM00713448	Dual Echo Artifact when cv6 is on
HCSDM00713608	Prone Position Axial Breast Image Rotate Feature
HCSDM00713706	fgre-xl 3plane not able to be saved
HCSDM00713725	Include pfov control in PUREMIX
HCSDM00713839	Update fixed R1/R2 value in signal acquisition for 1.5T system with MEGASWITCH
HCSDM00714191	subsyte update calibration Dir. in prescan for Invivo 8ch Wrist and 8ch Foot ankle coil, GE 8ch Foot ankle coil.
HCSDM00715219	Auto transfer by series cannot transfer images after moving the cradle
HCSDM00715220	Reformat sequence is automatically transferred without permission Issue
HCSDM00715638	Partial NEX over one should not be supported for BRAVO MP-RAGE

*** Printed copies are uncontrolled unless otherwise identified ***

Before using this document, consult MyWorkshop for the latest revision.

GE HealthCare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.

Released

GE HealthCare

Technical Bulletin



SPR # if applicable	Brief DescriptionD
HCSDM00717095	ghost artifacts in multi-slice Ax ungated Fiesta
HCSDM00717107	Acquisition failure in slice ASSET with R=3
HCSDM00717111	Flex long recon lag due to incorrect zFOV calculation
HCSDM00717787	Calibration coverage logic to determine if another cal is required is broken.
HCSDM00717901	Shoulder image dark line artifact from Atlantis pilot site
HCSDM00718008	Silenz MRA architect GE protocol has wrong HNU coil
HCSDM00719799	Integrate Option Key MR30. to MR29.1 build
HCSDM00720270	Series test page annotation for imaging calibration is incorrect for prescan options string.
HCSDM00720271	Incorrect series reused for imaging calibration.
HCSDM00720680	VX29.1_R02 (and RX29.1) Root Partition Full issue
HCSDM00724241	Coronal LAVA-Flex Robustness Improvement
HCSDM00722726	Breast Image Orientation based on Option Key
HCSDM00725856	Interleave still available with type-in for spacing on MUSE DWI.
HCSDM00728582	System cannot detect HEC via cable p/n 5729973 changed under SCR-0131764
HCSDM00726518	Performance degradation Issue with more number of slices in viewEdit/Setup Workflow.
HCSDM00713083	AIRx need extra option for Slice Orientation to meet with site's various slice order custom
HCSDM00713743	SSFSE with TRF protocol crash with calcLopeVarflip error
HCSDM00714190	Promote GE 8ch Foot Ankle and Invivo 8ch Foot Ankle coil to open ARC by customers' request

Revision	Version	Date	Author	Description
01	01	02-June-2023	Qin Yu	Initial Release

End of Document

*** Printed copies are uncontrolled unless otherwise identified ***

Before using this document, consult MyWorkshop for the latest revision.

GE HealthCare Confidential

General service documentation. Copyright General Electric company. GE proprietary. Reproduction and/or distribution is prohibited.

Released