Upgrading to Gen7/Gen7-DL Image Compute Node (ICN) on a Discovery MR750, MR750w, or Signa Architect System

Gen7/Gen7-DL Upgrade Manual

5829487
Revision 5
US English
© 2019-2022 General Electric Company
All rights reserved.

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Revision history

Record of reason for changes to this document.

Rev	Date	Description
5	January 2022	Routed in MyWorkshop with DOC2270583, revision 6. This revision adds support for the Dolphin 2 interface connection. The following updates were made: - Chapter 5, Section 5.2 Installing the ICN: - Added new step 1 - Added note to step 2 about changing between ICN styles - In Figure 5-2ICN cable connections (Gen 7 and Gen 7-DL, Dell R640), updated description in callout 1 - Added Figure 5-3 PCIe cable connection - Dolphin 2 card - Chapter 7, Section 7.1 Determining software installation: - Updated Table 7-2 Software upgrade steps to include an asterisk in the two SW 29 column headers for the Dolphin 2 footnote reference - Added footnote content for Dolphin 2 card and supported software
4	December 2020	- Updated title to include DV products. This document now applies to DV products as well as Architect. - Chapter 1 1.2.3 Updated harvest test information 1.2.4 Updated harvest list and part numbers 1.4 Updated product locator information 1.6 Added system and component environmental conditions - Chapter 5 5.2 Changed ISC to PGR - Chapter 7 7.1 Updated software upgrade steps table to include Rev 5 and earlier ICNs Revision 4 is referenced as DOC2270583, rev 5, in MyWorkshop.

Rev	Date	Description
3	June 2020	- Chapter 1 1.2.4 Added harvest ICN options and part numbers 1.3 Removed information not relevant to this manual 1.7 Removed non-magnetic tools - Chapter 5 - Added Replacing the Image Compute Node (ICN) network cables topic 5.2 Clarified installation steps and added Gen 7-DL relevant cabling information. - Chapter 7 - Changed chapter 7 to Determining software installation - Chapter 8 - Added Gen 7-DL and MR29 SW relevant information - Chapter 9 - Removed Logging in to the Image Compute Node (ICN) chapter Revision 3 is referenced as DOC2270583, rev 4, in MyWorkshop.
2	Nov 2019	Correct erroneous manual number on front page and footer throughout entire manual. Revision 2 is referenced as DOC2270583, rev 3, in MyWorkshop.
1	July 2019	Initial Release of part number 5829487. Revision 1 is referenced as DOC2270583, rev 2, in MyWorkshop.

Page intentionally left blank

Contents

Chapter 1 Getting started	13
1.1 General safety precautions	13
1.2 Material policy overview	13
1.2.1 Returned material.....	13
1.2.2 Disposal of material	13
1.2.3 Harvest test procedure.....	13
1.2.4 Harvest items	14
1.3 Product delivery instructions.....	15
1.4 Product locator.....	15
1.5 Damage in transportation.....	16
1.6 MR system shipping and storage requirements.....	17
1.7 Tools and materials	17
Chapter 2 Shutting down the host computer	19
2.1 Shutting down the host computer	19
Chapter 3 Applying LOTO.....	21
3.1 Applying LOTO - PGR PDU/Gradient Subsystem	21
Chapter 4 Removing the Image Compute Node (ICN)	25
4.1 Removing the ICN	25
Chapter 5 Installing the Image Compute Node (ICN)	27
5.1 Replacing the Image Compute Node (ICN) network cables	27
5.2 Installing the ICN	27
Chapter 6 Removing LOTO.....	29
6.1 Removing LOTO - PGR PDU/Gradient Subsystem.....	29
Chapter 7 Determining software installation	31
7.1 Determining software installation	31
Chapter 8 Configuring the Image Compute Node (ICN) - VRE	33
8.1 Configuring the ICN - VRE.....	33
Chapter 9 Doing a check scan.....	39
9.1 Doing a check scan.....	39

Page intentionally left blank

Chapter 1 Getting started

1.1 General safety precautions

Safety
Before working in any GE Healthcare MR suite or performing any GE Healthcare service procedure, you must: • Have read and understood all hazard conditions and safety requirements in the latest revision of the GE Healthcare <i>MR Service Safety Manual</i> (5452735). • Have successfully completed all relevant GE Healthcare Environmental Health and Safety (EHS) courses (or for non-GE employees, equivalent workplace training courses). • Comply with all site-specific training and workplace safety requirements. If you have any safety concerns at any time, do not begin work or immediately stop work and move to a safe location. Immediately contact your supervisor or site safety officer for instructions on how to proceed.

1.2 Material policy overview

1.2.1 Returned material

Based on GE Healthcare Policy 3.9, all requested material is to be promptly returned as directed, intact to either GoldSeal or Approved GE Recycling Center. Do not remove or swap any parts from the material to be returned.

1.2.2 Disposal of material

Regulations governing the disposal of this equipment or material vary from location to location and from country to country. To make sure you are compliant with all local, state, or national environment, health, and safety regulations, dispose of all items based on the local GEMS Policy and Procedures, and other applicable published guidance. Contact the GE Recycling Center for assistance.

1.2.3 Harvest test procedure

NOTICE

RISK OF DAMAGE

Equipment may be damaged prior to removal.

The equipment must be tested prior to the system shutdown/upgrade to allow appropriate harvest capabilities.

Before going on-site:

- Print the proper regional address labeling from *Harvest Return Addresses* (DOC2100599).
- Make sure you can download the appropriate *System Harvest Test Procedure* from the online documentation library using:
 - DOC1899393 for MR450, MR450W, MR750, MR750W, Artist, and Architect

- If you have not been contacted to test the system or cannot download the document, contact the local harvest team. For all other requests, contact the @HEALTH Global Harvest Team via email: GlobalHarvestTeam@ge.com.

After testing is complete, submit the test record for review. After these test records are reviewed, they will be processed for electronic signatures.

Package the equipment and follow your region's method of returning harvest equipment.

NOTE

Where possible, use packaging from the new upgraded component for the equivalent harvested parts.

1.2.4 Harvest items

Table 1-1 Return the following items for harvest after deinstallation:

Item	Part number
SUN 4100 (Gen2) ICN with VRF	5127452-2
or	or
SUN 4100 (Gen2) ICN with VRF and without Infiniband	5127452-6
or	or
SUN 4170 (Gen3) DV 16-Channel non-GEM	5337894
or	or
SUN 4170 (Gen3) DV 32-Channel non-GEM	5337894-2
or	or
SUN 4170M2 (Gen4) DV 16-Channel non-GEM	5911000
or	or
SUN 4170M2 (Gen4) DV 32-Channel 32 GB non-GEM	5911000-2
or	or
SUN 4170M2 (Gen4) DV 32-Channel 72 GB GEM	5911000-3
or	or
Dell R620 (Gen5) DV 16-Channel ICN with VRF	5921000
or	or
Dell R620 (Gen5) DV 32-Channel ICN with VRF	5921000-2
or	or
Dell R620 (Gen5) DV 32-Channel+ ICN with VRF	5921000-3
or	or
Dell R630 (Gen6) DV Value with Dolphin Card	5931000
or	or
Dell R630 (Gen6) Performance ICN with Dolphin 1 Interface Card	5931000-3
or	or
Dell R630 (Gen6) Advanced ICN with Dolphin 1 Interface Card	5931000-5
or	or

Table 1-1 Return the following items for harvest after deinstallation: (Table continued)

Item	Part number
Dell R640 (Gen7) Performance with Dolphin 1 Card or	5941000-3 or
Dell R640 (Gen7) Advanced ICN with Dolphin 1 Interface Card	5941000-5

1.3 Product delivery instructions

- The *Shipping Document* lists catalog numbers delivered. Review and confirm that order is delivered complete. Check impact on installation schedule if catalogs and/or packing boxes are missing and/or noted as shipped short.
- **NOTE**
For full systems and large components, two or more people are required to accept delivery of equipment due to weight and size.
- Labels that are attached to the outside of packing boxes summarize box contents.

Product Delivery Instructions (PDI) specify box contents, part numbers, and shipping procedure. The PDI is numbered according to the catalog number (for example, is for the Fixed Site Installation Kit). Lists of items included with each box are detailed by separate checklists, or a separate sheet that provides a summary of that box contents. Refer to PDI and packing lists for information specific to your shipment.

1.4 Product locator

During installation, the Install FE shall confirm that the model and serial numbers are part of the asset. Use either the gib web interface (http://gib.gehealthcare.com/gib/gib_entry.jsp) or the MyPLCReader app. This information will then be loaded directly to Siebel by the IB verification team. If components are missing, submit a Global Installation Quality (GIQ) request.

To install the MyPLCReader app, use the following procedure.

1. From the GE app icon, search for and Install the MyPLCReader app on your iPhone or iPad.

Figure 1-1 MyPLCReader app

2. Update the Configuration page.
 - a. The first time the app opens, the Configuration page will open.
 - b. Enter your SSO, Country and **if in US select the edit button to change the default email address** (by default, this app will send the email to the Europe IB Team Resource.EuropeanCentralAdminGIB@ge.com). The email to be used in US is usa.ib.verification@ge.com.
3. Start your first scan.
 - a. Go to the Menu, make the appropriate selection for the job you are performing: **New Install**, **Upgrade System**, or **Update System**.
 - b. Enter the applicable data, install end date, system ID, and RFS number.
 - c. Click **Scan Bar Codes** and start scanning each bar code on the PL Cards.
 - d. Click the **Email** button.
 - e. The email is sent to the default email address, which was set up earlier.

1.5 Damage in transportation

Examine all packages closely at delivery. If damage is apparent, do the following:

1. Make sure the notation *damage in shipment* is written on all copies of the freight or express bill before delivery is accepted or signed for by a General Electric representative or a hospital receiving agent.
2. Whether noted or concealed, you must report damage to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and hold the contents and containers for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14-day period.

1.6 MR system shipping and storage requirements

NOTE

Some content may not apply to your components.

1. The customer must provide a shipping dock or other staging area for unloading system components for full systems and large components (magnet, cabinets, patient table, and so on).
2. The customer is responsible for unloading and moving all system components from the delivery truck or van to final locations in the MR suite. This includes:
 - Door and other openings are sufficiently wide to allow passage
 - Sufficient space is provided for any required dollies or rigging
3. MR systems and components must be shipped and stored within the environmental requirements as listed in [Table 1-2 Transportation and storage environmental conditions on page 17](#).

Table 1-2 Transportation and storage environmental conditions

System equipment	Temperature range	Temperature change	Relative humidity	Humidity change
Electronic Cabinets and Equipment	-22 to 140°F (-30 to 60°C)	68°F/Hr (20°C/Hr)	10-80% non-condensing	30%/Hr
Surface Coil	-22 to 122°F (-30 to 50°C)	68°F/Hr (20°C/Hr)	10-95% non-condensing	30%/Hr

1.7 Tools and materials

Tools and test equipment			
Item	Quantity	Part number	Manufacturer
Brass Master Padlock with Brass Shackle	1	46-194427P320	-
Red LOTO Personal Lock Wrap	2	2393068	-
Red Warning LOTO Tag	1	46-194427P322	-
Digital Voltmeter (DVM)	1	46-194427P284	-


Page intentionally left blank

Chapter 2 Shutting down the host computer


2.1 Shutting down the host computer

1. To shut down the system, click the arrow on the **Tools High Level Access** tab.
2. Click **System Shutdown**.

Figure 2-1 System power off

3. Click **Yes** to confirm the system power off.
4. Select the reason for power off to confirm.

Page intentionally left blank

Chapter 3 Applying LOTO

3.1 Applying LOTO - PGR PDU/Gradient Subsystem

⚠ DANGER

ELECTROCUTION HAZARD!

HIGH VOLTAGE PRESENT.

USE PROPER LOCK OUT / TAG OUT PROCEDURES BEFORE servicing and/or maintaining.

Overview

The PGR cabinet contains the Power Distribution Unit (PDU) and eXtreme Gradient Amplifier (XGA). The XGA is part of the eXtreme Gradient Driver (XGD) subsystem. Prior to servicing and/or maintenance work on the PGR PDU or any element of the XGD, the steps in this Lock Out / Tag Out (LOTO) procedure must be completed by a LOTO-authorized GE Field Engineer (FE). Completing all the steps in this procedure ensures a safe environment and avoids equipment damage when working on these parts.

Table 3-1 Affected by LOTO

Name of Equipment	Number of Locks	Titles of Employees Authorized to perform LOTO	Titles of Affected Employees	How to Notify
- eXtreme Gradient Amplifier (XGA) - eXtreme Power Supply (XPS) - Power, Gradient, RF Cabinet Power Distribution Unit (PGR PDU) - Gradient Filter	1 per GE FE	GE FEs	Hospital personnel	Verbal, posted signs

NOTE

The following table describes the type, location, and magnitude of energy to be LOTO'd.

Energy source	Yes	No	Location of energy isolating means	Magnitude of energy
Electrical	x		PGR PDU circuit breaker panel, Main Disconnect Panel (MDP) circuit breaker panel	208 VAC, 420 VAC, 380-480 VAC
Pneumatic		x		
Hydraulic		x		
Gas/Water/Steam		x		
Chemical		x		
Mechanical Motion		x		
Gravity		x		
Springs		x		

Energy source	Yes	No	Location of energy isolating means	Magnitude of energy
Thermal		x		
Stored Energy	x		Time discharge of internal capacitors	208 VAC, 420 VAC
Air Under Pressure		x		
Oil Under Pressure		x		
Water Under Pressure	x		Heat Exchanger Cabinet (HEC), PGR Pump Circuit Breaker (PPMPCB) and RF amp quick disconnect	5.0 Bar (80 psi)
Gas Under Pressure		x		
Steam		x		
Other		x		

Types of equipment and/or methods selected to dissipate or isolate stored energy


Time method for dissipation of stored energy: 5 minutes.

Type(s) of equipment and/or method(s) used to ensure disconnection(s)

1. One personal lock and tag per FE at the site working on the de-energized equipment at the time of LOTO.
2. Voltmeter to read voltage for verification that energy has dissipated in the PDU.

Applying LOTO

1. To prepare for shutdown of equipment, notify affected personnel working in the area the LOTO is being performed.
2. At the Operator Workspace, shut down the host PC.
3. Locate and turn off the PGR PDU main breaker on the MDP that feeds the PGR cabinet.
Lock Out/Tag Out the PDU breaker on the MDP.

Figure 3-1 PDU breaker on new style MDP

4. Make sure all energy has been dissipated. Use a voltmeter to measure the line-to-line voltage between each INPUT phase voltage test point on the PDU.


The PDU main breaker must be in the On position to read accurate voltage on the test points.

NOTE

The test points are 100:1 ratio (480 VAC = 4.80 VAC, 208 VAC = 2.08 VAC). Therefore, a low voltage reading still indicates the presence of high voltage. **Your voltage reading for each of the three phases must be 0 VAC before proceeding.**

Input Voltage Range: 380 - 480 VAC = 3.8 - 4.8 VAC

Output Voltage Range: 208 VAC = 2.08 VAC

Figure 3-2 New style PDU test points

Page intentionally left blank

Chapter 4 Removing the Image Compute Node (ICN)

4.1 Removing the ICN

NOTICE


FAILURE TO PROPERLY POWER DOWN THE ICN CAN RESULT IN PREMATURE ICN FAILURE.

Never turn the ICN ON/OFF from the power button on the front. There is software running on ICN disk drives that can be corrupted if not shut down correctly.

Follow the procedure for removing power to the ICN chassis.

1. Open the PGR cabinet front door.
2. Loosen 2 captive screws under the screw cover and slide out the ICN.
3. Slide the ICN out of the cabinet until the ends of the rails touch the stoppers.


Figure 4-1 Slide the ICN out

4. Disconnect all the cables from the rear of the ICN, including the AC power plug(s). Note the locations for reinstallation.


5. Push the two release tabs down and move the ICN forward.

Figure 4-2 Slide rails and release tabs

6. Lift the ICN up to release it from the bracket.
There are four studs on each side of the ICN that fit into the rail.

Figure 4-3 Removing ICN from rails

Chapter 5 Installing the Image Compute Node (ICN)

5.1 Replacing the Image Compute Node (ICN) network cables

1. Remove the Ethernet cable (5488202-2) from port 2 of the Ethernet switch and discard.
2. Remove the Ethernet cable (5488202-3) from port 3 of the Ethernet switch and discard.
3. Install the Ethernet cable (5488202-8) into the Ethernet switch. The other end of this cable will connect to port 4 of the new ICN.

5.2 Installing the ICN

NOTICE

FAILURE TO PROPERLY POWER DOWN THE ICN CAN RESULT IN PREMATURE ICN FAILURE.

Never turn the ICN ON/OFF from the power button on the front. There is software running on ICN disk drives that can be corrupted if not shut down correctly.

Follow the procedure for removing power to the ICN chassis.


1. With the slide rails out, install the new ICN, lining up the studs on the ICN with the slots on the rails.
2. Connect all other cables that were previously removed to the rear of the new ICN.

NOTE


If changing between ICN styles (units with the Dolphin 1 interface with units with the Dolphin 2 interface) make sure the associated PICE cable is available. Some cabinets will contain both cables. Save the alternate cable with the cabinet for future use.

NOTE

Gen 7 and Gen 7-DL ICNs have the same cable connections.

Figure 5-2 ICN cable connections (Gen 7 and Gen 7-DL, Dell R640)

1	PCIe port (Dolphin 1 card), or PCIe P1 port (Dolphin 2 card)
2	iDRAC (management controller for ICN)
3	Raw data server (RDS) - optional, only for research sites
4	Main switch (port 4), only 1 network cable is needed for GEN 7 and Gen 7-DL ICN
5	Power

Figure 5-3 PCIe cable connection - Dolphin 2 card**NOTE**

Ethernet ports 1 & 3 are NOT used. Ethernet port 2 is used ONLY when the system is connected to a Raw Data Server (external to the MR System), it should be left free in all other cases. Not following this configuration may lead to VRE configuration failures and/or boot-up/TPS failures.

3. Push the two release tabs down and move the ICN all the way into the PGR cabinet until it locks in place.
4. Install 2 screws.
5. Restore all cabinet covers.

Chapter 6 Removing LOTO

6.1 Removing LOTO - PGR PDU/Gradient Subsystem

1. Alert necessary personnel that power is being restored.
2. Remove the lock and tag.
3. Switch on the MDP's PDU breaker.
4. Press the green **EMO RESET** button on the PDU.

Figure 6-1 EMO RESET button

5. Power up the host PC.

Page intentionally left blank

Chapter 7 Determining software installation

7.1 Determining software installation

Based on the Gen 7 Image Compute Node (ICN) compatibility, determine the next steps needed. Software Load From Cold (LFC) is required in some cases and other cases only require configuring the ICN. Refer to the below tables to determine the next actions.

Use the below task index in combination with [Table 7-2 Software upgrade steps on page 31](#) to determine the next steps for the upgrade.

Table 7-1 Task index

Task index	Task description	Document reference
1	Replace the ICN	4.1 Removing the ICN on page 25 and 5.2 Installing the ICN on page 27
2	Configure the ICN	8.1 Configuring the ICN - VRE on page 33
3	Upgrade software by performing LFC	Refer to <i>Installing host system software (quick start)</i> from the online documentation library.

The below table shows the next steps to be performed based on the task number per index.

Table 7-2 Software upgrade steps

Upgrade To ---->	Gen 7 ICN, SW 28	Gen 7 DL ICN, SW 28	Gen 7 ICN, SW 29*	Gen 7 DL ICN, SW 29*
Upgrade From				
Gen 5 and earlier ICN	Not Applicable	Not Applicable	1, 3	1, 3
Gen 6 ICN, SW 28	1, 2	Not Compatible	1, 3	1, 3
Gen 6 ICN, SW 29	Not Applicable	Not Applicable	1, 2	1, 2
Gen 7 ICN, SW 28	Not Applicable	Not Applicable	3	1, 3
Gen 7 ICN, SW 29	Not Applicable	Not Applicable	Not Applicable	1, 2

* Gen 7 ICN with the Dolphin 2 Interface requires a software level of DV29.1_R02_SP01 or DV29.0_R02_SP01 or later.

If the ICN is being upgraded to the Gen7-DL as part of the Air Recon DL option upgrade, make sure the airrcndl2d option key is installed and the system is rebooted for the key to initialize the Air Recon DL protocols. For more information refer to *Air Recon DL operation and FE troubleshooting theory* in the online documentation library.

Page intentionally left blank

Chapter 8 Configuring the Image Compute Node (ICN) - VRE

8.1 Configuring the ICN - VRE

Use this procedure if Load From Cold (LFC) with the required software has already been completed.

The following procedure describes manual ICN-VRE configuration. ICN compatibility is dependent on software revision. Refer to *PGR Cabinet* for compatibility. Volume Recon Engine (VRE) and Image Compute Node (ICN) are terms that are used interchangeably. VRE configuration installs the software onto the ICN. This procedure is needed after ICN replacement. This procedure may also be needed when troubleshooting issues with the ICN system. Refer to *Troubleshooting the Image Compute Node (ICN)* for more information.

1. **NOTE**

If you are replacing the ICN, power on the system and the host PC, and log in as root.

NOTE

(For DV29.0 or later) The root password is no longer **operator**. Passwords are managed according to *MR Cyber Security Manual* in the online documentation library.

Click **Guided Install** and then click **Start**.

2. Invoke Guided Install as root.


- a. When you power up the host computer and see the logon screen, click on **root login**.

Figure 8-1 Logon screen

- b. Type the user name: **root** and type the root password.

Figure 8-2 Welcome logon screen

NOTE

It is possible that the customer changed the default password. If you cannot log on , contact the customer for the correct password.

Logging on as the root user automatically shows the Linux desktop and launches the *Guided Install Starter* window.

- c. Click **Yes** to invoke the guided install.

NOTE

If **No** has already been selected, right-click on the Linux desktop and select **Start Guided Install**.


3. Select **ICN Configuration**, click **Configure ICN**, and click **OK**.
4. On the ICN configuration screen, click **View VRE Configuration file** to view the output status of the ICN.

NOTE

Information will only be shown if ICN was previously configured in the software. Note the serial number and check the output to verify the new serial number.

- Click **Configure ICN**.

Figure 8-3 Guided install – ICN configuration


- If prompted, insert the OS DVD, or USB as applicable, that came with the system software into the host PC and click **OK** to continue the ICN configuration (this step is only necessary if repository files are corrupted on the host computer).

ICN configuration progress and details appear in a separate window. At the end of ICN configuration, another dialog box is displayed with ICN configuration success or failure message.


7. If ICN configuration fails, review the dialog box and details in the separate window to determine next steps.

Figure 8-4 ICN progress details

8. When the **ICN Remote Console** window appears, click **OK**, and then **Continue**. This will install the Java security options.


Figure 8-5 ICN remote console

9. Once the *Java Security Warning* window appears select the following:
 - a. Select the **I accept the risk and want to run this application** checkbox.
 - b. Click on **Show Options** and select the **Do not show this again for apps from the publisher and the location above.** checkbox.

- a. Click **Run**.

Figure 8-6 Security Warning

The *Remote Console* window opens showing the progress of the installation.

10. At the end of VRE configuration, a dialog box displays with either a success or failure message. If VRE configuration fails, review the dialog box message and details in the separate window to determine next steps.

NOTE

On the ICN configuration screen, you can click **View VRE Configuration file** to view the output status of the ICN.

11. Click **OK** twice for the VRE configuration success and ICN configuration success messages.

NOTE

A reboot is required after complete ICN configuration.

12. Remove the OS DVD or USB SW installer.

Page intentionally left blank

Chapter 9 Doing a check scan

9.1 Doing a check scan

Make sure that you use the DQA phantom, so that the image can be analyzed for correct orientation.


NOTICE

Depending on site configuration and/or rules, you may want to temporarily disable Auto Archive or Auto Network.


1. Attach the head coil to the cradle and position the DQA phantom inside the coil.
2. Landmark on the center line of the phantom. The laser must be in the middle of the line. Press **Advance to Scan**.

Figure 9-1 Phantom orientation, split-top head coil

- Click the **New Worklist** icon.

Figure 9-2 New worklist icon

- Select a **Patient ID** workflow (other than **geservice**).

Figure 9-3 Entering a new work list

The screenshot shows a software interface for entering a new work list. At the top, there's a 'Patient Record' section with icons for adding, editing, and deleting records. Below this is a table with columns: Scheduled, Worklist, Name, Patient ID, Accession, SPS Description, Modality, Ref. Physician, Scheduled, Status, Local Protocol, Recomm., and Caution. The first row shows a patient scheduled for 01:58:31 PM, Local, with a Patient ID of 000000, Accession of 000000, SPS Description of MR, Modality of MR, Ref. Physician of 11-26-2007, Status of Started, and Local Protocol of 30.11 3-Plane 2D Localizer.

Below the table is a form with three main sections: Patient, Exam, and Other Information.

Patient Section:

- Name: Last Name (), First Name ()
- * Patient ID: A00000
- Scheduled Date: 11-26-2007
- Date Of Birth: MM-DD-YYYY ()
- Age: Year (), Month (), Weeks (), Days ()
- Sex: ()
- * Weight: Lbs (111.0), Kgs (50.35)

Exam Section:

- Ref. Physici...: Last Name (), First Name ()
- Accession: ()
- SPS Description: ()
- Exam Description: ()
- Procedure ID: ()
- Operator's Initial: ()
- Radiologist: ()
- Protocol: 30.11 3-Plane 2D Localizer
- Favorite Protocols: Select Favorite ()
- Show All Protocols... ()

Other Information Section:

- This is patient sensitive information, please adhere to patient confidentiality rules. Please verify the information using either the patient's medical record or directly with the patient to guarantee its accuracy.
- Allergies: NONE ()
- Pre-Med: NONE ()
- Pregnancy Status: No Entry ()
- History: ()

At the bottom right, there is a 'Reopen Exam' button.

- Enter for weight: **111**.
- Click **Show All Protocols...**
- When the **Protocol Selector** menu opens, click **Protocol Details**.

8. Select the Protocol Library: **GE**.

Figure 9-4 Setting protocols

9. Select **Adult**.
10. Click the **Template** tab.


11. Select **3 - Plane 2D Localizer** and click the arrow to move the selection to the box on the right.

Figure 9-5 Protocol selected

12. Click **Accept**.
13. Click **Save**.
14. Click **Start Exam**.
15. Enter the following image parameters.


Figure 9-6 Image parameters

Option	Description
Frequency field of view (FOV) (Freq. FOV)	30.0
Phase FOV	1.00
Slice Thickness	5.0

- 16. Click **Save RX**.
- 17. Click **Scan**.
- 18. Review the images to confirm that the orientation is correct.

Figure 9-7 Check scan image orientation

Axials:	
1	GE logo in upper right corner
2	A is displayed
Coronal & Sagittal:	
3	Comb in upper left corner
4	CS displayed in upper right corner

Page intentionally left blank