
Service Manual

for

16Ch Diagnostic Breast Coil - P-H16LE-030-01893

to be operated on

GE 3.0 T MR Systems

3.0T GEHC Part no.: 5772248

RAPID Biomedical GmbH

Manufacturer:

RAPID Biomedical GmbH
Kettelerstrasse 3-11
97222 Rimpf, Germany

Phone: +49 (0)9365-8826-0
Fax: +49 (0)9365-8826-99

Email: info@rapidbiomed.de
Web: www.rapidbiomed.de
© 2019-02-13 RAPID Biomedical GmbH

Document Version: 0.15

Technical changes reserved.

Table of Contents

1	Important Information	4
2	Introduction	9
2.1	Shipping Contents	9
2.2	Environment Storage Requirement	10
3	Coil Installation	11
3.1	Plug label and Mirror protection	11
3.2	Adapting the coil to the patient table	12
3.3	Performance / Quality Assurance (QA)	19
4	Field Replaceable Units (FRUs)	26
4.1	List of Field Replaceable Units	26
4.2	Replacement of the cushions	27
4.3	Returning for Service	28
5	Troubleshooting	29
5.1	General notes	29
5.2	Failed MCQA	29

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


2 Introduction

2.1 Shipping Contents

Description	GEHC Part #	Rapid Part No.	Illustration
3.0T 16Ch Diagnostic Breast Coil with cable	5772248-2	P-H16LE-015-01630	
16Ch Diagnostic Breast Coil positioning parts set	5772248-3	ZUB-01858	
16Ch Diagnostic Breast Coil head rest framework	5772248-4	ZUB-01859	
16Ch Diagnostic Breast Coil cushion set	5772248-5	ZUB-01860	
CD containing user manual (all languages)	N/A	RT004808	N/A
3.0 T Foreign Language Sticker Sheet	N/A	RT004812	N/A
eIFU Leaflet	N/A	RT004810	N/A

16ch Diagnostic Breast Coil service manual	N/A	RT004811	N/A
--	-----	----------	-----

2.2 Environment Storage Requirement

Store the coil and accessories in an air conditioned scan room or equipment room until the GE field engineer can perform installation. The environmental condition limits are included below.

	<i>Operating Conditions</i>	<i>Transport & Storage Conditions</i>
<i>Temperature Range</i>	+15°C to +24°C +59°F to +75.2°F	-25°C to +60°C -13°F to +140°F
<i>Relative Humidity</i>	30 % to 80 % RH	5 % to 95 % RH
<i>Air Pressure</i>	70 kPa - 107 kPa	n/a

3 Coil Installation

3.1 Plug label and Mirror protection

3.1.1.1 List of required parts

Description	Quantity	Illustration
Foreign Language Sticker Sheet	1	


3.1.1.2 Installation instruction

Step 1: Select the sticker according to the national language

Step 2: Paste the sticker to the empty field of the connector

Step 3: Peel of the protective film from the mirror of the head rest

3.2 Adapting the coil to the patient table


The housing of the 16Ch Diagnostic Breast Coil and the head rest framework are adapted to the patient table of the individual MR system by positioning parts. The following subsections give specific instructions for the necessary adjustments in correspondence to the MR system on site.

Tools are not included the package. Use appropriate screw drivers to avoid wearing out the screw heads. Do not use excessive force when tightening the screws.

3.2.1 Curved Table


This is the assembly instruction for the Curved Table. Not all of the enclosed parts are required. Select the parts shown in the list below only.

3.2.1.1 List of required parts


Description	Quantity	Illustration
Elbow bracket	4	
Short stand with inner mounting holes	1	
Short stand with outer mounting holes	1	
Long stand for head rest	2	
Long slotted flat headcheese head screw	4	
Short slotted countersunk screw	4	
knurled Thumb screw	4	

3.2.1.2 Installation instruction

The following steps show how to adapt the coil and the head rest to the Curved Table.

Step 1: Turn the coil upside down

Step 2: Unscrew the four screws marked by the screwdrivers in the image above.

Step 3: Align the two short stands in the holes where you have removed the screws. Insert the four long slotted flat head screws through the holes of the stands.


Step 4: Tighten the four screws.

Step 5: In each corner, place a short slotted countersunk screw and an elbow bracket. The bracket points towards the side of the coil. Tighten the screws.

Step 6: Continue until all four brackets are mounted as shown above.

Step 7: Take the head rest and turn it upside down.

Step 8: Take the two long stands and the four knurled thumb screws

Step 9: Tighten all four screws by hand to finish the assembly

3.2.2 GEM Table (narrow flat table)

This is the assembly instruction for the GEM Table. Not all of the enclosed parts are required. Please select the parts shown in the list below only.

3.2.2.1 List of required parts

Description	Quantity	Illustration
Elbow bracket	4	
Short slotted countersunk screw	4	

3.2.2.2 Installation instruction

Step 1: Turn the coil upside down

Step 2: In each corner, place a short slotted countersunk screw and an elbow bracket. The bracket points towards the side of the coil. Tighten the screws.


Step 3: Continue until all four brackets are mounted as shown above.

3.2.3 Fixed Table (wide flat table)

This is the assembly instruction for the wide Fixed Table. Not all of the enclosed parts are required. Please select the parts shown in the list below only.

3.2.3.1 List of required parts

Description	Quantity	Illustration
Short slotted countersunk screw	4	
Wide table bracket	2	

3.2.3.2 Installation instruction

Step 1: Turn the coil upside down

Step 2: Take the four short slotted countersunk screws and the two wide table brackets. Place them to the mounting slots and tighten the screws.


Step 3: Continue until both brackets are mounted.

3.2.4 PET/MR Table

This is the assembly instruction for the PET/MR Table. Not all of the enclosed parts are required. Please select the parts shown in the list below only.

3.2.4.1 List of required parts

Description	Quantity	Illustration
Alignment cone	2	
PET head positioning frame	2	
Mounting slot cover	2	
Short slotted countersunk screw	2	
Thumb screw	4	

3.2.4.2 Installation instruction

Step 1: Turn the coil upside down


Step 2: Screw one of the two “alignment cones” into the thread near the cable outlet..


Step 3: Continue with the second “alignment cone” on the opposite side as shown above.


Step 4: Cover the remaining two slots using “mounting slot covers” and a short slotted countersunk screw.

Step 5: Take the PET positioning frame and the thumb screw

Step 6: Put the head rest onto the frame

Step 7: Tighten the four thumb screws

3.3 Performance / Quality Assurance (QA)

We recommend regular verification of the device's proper function by performing the Coil Quality Assurance test utilizing GE Multi-Coil Quality Assurance (MCQA) Tool.

3.3.1 MCQA setup

3.3.1.1 List of required parts

Description	Quantity	Illustration
16Ch Diagnostic Breast Coil	1	
Phantom holder foam ring	1	
Phantom GEHC part no. 5342679-2 (Large Cylindrical Unified Phantom, SiOil)	1	

3.3.1.2 Phantom positioning

The patient position is prone, feet forward. This results in the designation of “left” and “right” as shown in Figure 1.

Figure 1: MCQA setup –

left: Phantom placed inside the left cavity of the coil;

right: Phantom placed inside the right cavity of the coil.

Each cavity of the coil will be tested separately. Therefore the MCQA test must be run twice to test a single coil.

First place the foam holder in one of the cavities of the coil. The foam holder is elliptic; check its orientation for an easy fit. Place the phantoms inside the holder. The phantom shall touch the surface of the patient table.

3.3.2 Description of the MCQA test procedure

3.3.2.1 MR system's calibration

All RF coil related tests have to be performed on a system that is calibrated. The command "install_in_spec" and especially EPIWP (White Pixel from install in spec) is to result in PASS.

3.3.2.2 Initiating MCQA

- 1) From Common Service Desktop (CSD), go to **Service Browser** and select **[Image Quality]** "**Multi-Coil QA Tool**" and then "**Click here to start this tool**" as shown Figure 2.


Figure 2: Multi-Coil QA Tool – "Click here to start this tool"

Note: If a "No valid MCR-V (or MCR2/3)" warning (Figure 3) pops up click **[Yes]** and proceed with test. MCR-V diagnostics must be run before turning over system to customer.

Figure 3: Warning "No valid MCR-V (or MCR2/3)"

The field “current coil” will be automatically filled in (Figure 4) with the coil name, based on the CoilID of the coil connected to the LPCA. Choose entry ending as “Left_Array” or “Right_Array” for the test with the phantom placed inside the left or the right cavity of the coil. Enter the serial number of the coil being tested in the field “Coil Serial #”.

- 2) Click **[Start]** to begin the automated test as shown in Figure 4. The test may take from 3-5 minutes depending on the MR system.


Figure 4: Selection of Current Coil for tests on the left or the right cavity of the coil.

- 3) Upon start-up, a message stating “Phantom placement and coil landmarking are critical for repeatable results” will appear. If the landmark has been set correctly and there are no air bubbles in the phantom, click **[Yes]** to continue. (Figure 5)


Figure 5: Warning Phantom placement & landmarking

Note: The Status window of the MCQA Tool GUI will continuously update to give information on what the tool is doing at any point in time. A progress bar (Figure 6) will appear, showing approximate total test time, elapsed time and percent complete.


Figure 6: Status window

Note: The test protocol might display further (error) messages. In this case the test is considered **FAIL**. See below for corrective actions.

After the test has been completed, test results are displayed on the screen (Figure 7).

The PASS/FAIL status shows **PASS** if all coil elements are functioning properly.

The MCQA Tool GUI displays **FAIL**, e.g. for one of the following reasons:

- Incorrect positioning/placement of the phantom (cf. Figure 8)
- Incorrect phantom used for the test (cf. Figure 9)
- Bad Coil Element (contact service)

Refer to section 5.2 Failed MCQA if the MCQA test fails.

More information on the MCQA test can be found on MR service methods DVD or website via the path: Troubleshooting -> System -> Multi-Coil Quality Assurance Tool


Figure 7: Tests results displayed on the screen –

In this example all tests were passed by the left and the right array half


Figure 8: Error message related to wrong selection of the cavity. Make sure the “Current coil” is selected corresponding the the placement of the phantom.


Figure 9: Error message related to wrong phantom. Place only one phantom in one of the cavities. Tests with two phantoms simultaneously are not possible.

- 4) Perform the test for the other cavity of the coil. The test for both cavities must result in **PASS**.
- 5) Click on **[Quit]** button to exit MCQA Tool.

3.3.2.3 Using MCQA Viewer

In case if the results are to be viewed at a later stage follow the below steps:

- 1) In the MCQA Tool window select File, Open Results File and select the desired coil results file select **[View Report Details]** to review the results.

Note: The Results Viewer will open as shown in Figure 10. The Results file name and Pass/Fail Results shown on the tool GUI will also be listed across the top of the viewer.


Figure 10: MCQA Viewer

- 2) Select the ISNR option and the ISNR Specs check box in the middle portion of the Results Viewer to view the results.

3.3.3 Execution of the MCQA test procedure



Run the MCQA test as described above for the left and the right cavity. Position the phantom accordingly. Each MCQA test has to result in **PASS**.


Personal Requirements	Preliminary Reqs	Procedure	Finalization
1 person	0 min	20 mins	0 min

Test ID	Parameter Description	Expected Result
"left"	"install_in_spec" passed	PASS

4 Field Replaceable Units (FRUs)


4.1 List of Field Replaceable Units

Description	GEHC Part #	Rapid Part No.	Illustration
1.5T 16Ch Diagnostic Breast Coil with cable	5772249-2	P-H16LE-015-01811	
3.0T 16Ch Diagnostic Breast Coil with cable	5772248-2	P-H16LE-030-01630	
1.5T 16Ch Diagnostic Breast Coil positioning parts set	5772249-3	ZUB-01857	
3,0T 16Ch Diagnostic Breast Coil positioning parts set	5772248-3	ZUB-01858	


16Ch Diagnostic Breast Coil head rest framework	5772248-4	ZUB-01859	
16Ch Diagnostic Breast Coil cushion set	5772248-5	ZUB-01860	

4.2 Replacement of the cushions


The pads of the head rest framework are mounted by hook-and-loop strips. Peel the individual pads off the frame and stick on the new ones. See the following photos for details:

Step 1: Peel of the forehead pad

Step 2: Peel of the cheek pads

Step 3: All pads removed.

Step 4: Attach the new pads to the head rest.

4.3 Returning for Service

In the event that a coil has to be returned for service, make sure that all parts of the FRU “16Ch Diagnostic Breast Coil positioning parts set” are disassembled according to this Service Manual. Otherwise the device will not fit into the dedicated packaging and might be damaged during transportation. For disassembling, follow the steps described in chapter 3 Coil Installation in reverse order.


Figure 11: Dedicated packaging for 16Ch Diagnostic Breast Coil with cable

5 Troubleshooting

5.1 General notes

When encountering a problem with the 16Ch Diagnostic Breast Coils, try the following troubleshooting pointers before conducting further investigation.

- There are no bent pins on the connector.
- The coil is connected and locked to the system port. The indicator on the coil connector should be showing the green, locked symbol.
- The system detects the connected coil correctly, and shows coil without error on the screen above the bore.

5.2 Failed MCQA

When the MCQA analysis has reported a failure:

- Check the procedure (phantom, positioning, landmark...).
- Redo the MCQA.
- If the MCQA still fails, return the coil (1.5T FRU 5772249-2 or 3.0T FRU 5772248-2) for repair.