

# DICOM Worklist, SW Installation Instruction

## 1. Personnel Requirements

| Required Persons | Preliminary Reqs | Procedure | Finalization   |
|------------------|------------------|-----------|----------------|
| 1                | Not Applicable   | 15 mins   | Not Applicable |

## 2. Overview

This document describes the procedure for installing the DICOM Worklist Software on to the Scanner system's host PC. This procedure assumes that the required peripheral devices, such as the CD drive have already been properly installed onto the host computer.

## 3. Preliminary Requirements

This procedure assumes that the operator has a basic understanding of personal computers and is familiar with the installation of Windows software. This installation requires a computer system that meets the Host Computer Functional Specification along with the associated software disks.

### 3.1 Consumables

| Item                               | Qty | Part#         |
|------------------------------------|-----|---------------|
| CD, SCANNER DICOM Worklist Release | 1   | 2010-4450-001 |

## 4. Procedure

Before the procedure is executed, check the software Release Notes (DOC0798893) for updates and the latest information concerning software installation.

**Execute the following installation procedure in the order it is presented.**

### 4.1 Installation of the DICOM Worklist Software

1. Boot the Scanner system Host Computer and log in as "operator" or "sdc" user when prompted by the Windows operating system.

2. With the Scanner software running in “**View**” mode, insert the Installation CD (part # 2010-4400-001) into the CD-RW drive.
3. Click the Windows **Start** Key.
4. Click on **My Computer**.
5. Double Click on the Icon for the CD drive (drive F:).
6. Double Click on the file **MWL\_setup.exe**
7. Follow the on-screen installation instructions.
8. The DICOM Worklist software is now ready for use; refer to the DICOM Worklist User’s Guide for instructions on how to use and configure the Worklist software.
9. Remove CD and store in the CD holder.

## 5. Finalization

No finalization steps.