GE Healthcare

Technical Publications

Direction 5158949-100
Revision 5

ReportCARD 3 Installation/Service Manual

Copyright© 2007 by General Electric Company Inc.
All Rights Reserved

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation “**damage in shipment**” written on **all** copies of the freight or express bill **before** delivery is accepted or “signed for” by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately** upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14** day period.

File a report with

- Call 1–800–548–3366 and use option 8.
- Contact your local service coordinator for more information on this process.

Rev. 06/09/2006

WARNING

- THIS SERVICE MANUAL IS AVAILABLE IN ENGLISH ONLY.
- IF A CUSTOMER'S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER'S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR OR PATIENT FROM ELECTRIC SHOCK, MECHANICAL OR OTHER HAZARDS.

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REVISION HISTORY

REV	DATE	PRIMARY REASON FOR CHANGE
1..	September 20, 2005	Initial Version
2...	November 21, 2005	Updated with FRU Part Number
3...	June 09, 2006.....	Added Installation Finish Tasks Section 6 and Updated Section 8 with FRU Part Number
4	October 06, 2006	Changed the FRU number
5	January 29, 2007 ...	Chnaged the FRU number

* This revision/letter corresponds to the indicated document's revision control system.

TABLE OF CONTENTS

1- INTRODUCTION	7
1-1 ReportCARD 3 Advanced Application “Clinical” Summary	7
2- PRELIMINARY REQUIREMENTS	7
2-1 License Key Number.....	7
2-2 Save Configuration	7
2-2-1 Reading Config Media → PC	8
2-2-2 Reading Config Media → AW	8
3- INSTALLATION PROCEDURE	8
4- AW AND REPORTCARD 3 ERRORS	8
5- ADVANCED SOFTWARE APPLICATION OPERATION AND DEMO EXAM INFORMATION	9
6- INSTALLATION FINISH TASKS	9
7- SPECIAL CONSIDERATIONS	9
8- FRU LIST	10

1- INTRODUCTION

The following manual gives a brief overview of the installation process for the ReportCARD 3 advanced application software option. The ReportCARD 3 application is designed to run on the AW4.2 and forward AW platforms only. This guide makes references to all the appropriate standard Advantage Workstation service manuals for more detailed information about installing an advanced application software option on the AW. The current revision of the appropriate AW Service Manual should be used when installing this software. Unique installation processes for ReportCARD 3 and NOT found in the AW service manual can be found in this document under the section entitled SPECIAL CONSIDERATIONS.

1-1 ReportCARD 3 Advanced Application “Clinical” Summary

ReportCARD 3 application is a reviewing and reporting application for Magnetic Resonance Imaging (MRI) data, in particular, Cardiac MRI.

The tool provides simultaneous reviewing and report creation of cardiac MR exams with multiple series containing DICOM images. Multi-phase cardiac images are automatically displayed as time-lapse animations with no user interaction.

This tool is capable of automated calculation of all quantitative results derived from user measurements. Also, automatic body surface area correction for appropriate cardiac analysis can be done with ReportCARD 3. Once the image analysis is completed, the physician can approve the report by using the signature feature.

2- PRELIMINARY REQUIREMENTS

2-1 License Key Number

The ReportCARD 3 option requires a license key # based upon the workstation's Host ID. This key number may be provided by:

- ❑ Paperwork shipped with the option that will include one or more 8 to 10 digit license key number(s) – created based on the workstation's/application's FDO and the workstation's Host ID/License ID.
- ❑ The software package may include a portable media (i.e. floppy disk) with the license key information on it.
- ❑ The license key information will also be available online at the eLicense website link by searching under the order number or the workstation's Host ID/License ID: <http://egems.gemedicalsystems.com/elicense/index.jsp>.

All of these elements may exist (hard copy license #(s), media, and eLicense data), and at least one of them is necessary to load the software.

2-2 Save Configuration

Create an AW 'CONFIG' media before and after the application software installation for use as needed throughout this and other AW service activities.

2-2-1 Reading Config Media → PC

Configuration information can be read on a PC by viewing the config.txt file on the media using text readers like WordPad, Internet Explorer or MS Word (MS NotePad may not work).

2-2-2 Reading Config Media → AW

To read a Save Configuration using a floppy disk insert the floppy and type the following in a c shell:

- **mount /mnt/floppy**
- **cd /mnt/floppy**
- **more config.txt**
- **cd/**
- **umount /mnt/floppy**

To read a Save Configuration using a CDROM disk insert the CD and type the following in a c shell:

- **mount /mnt/cdrom**
- **cd /mnt/cdrom**
- **more ("backup directory name")**
- **cd/**
- **umount /mnt/cdrom**

3- INSTALLATION PROCEDURE

The 'Standard' AW Software Options Installation – Easy Install process applies. For ongoing process revision viability, the details will not be printed here, but they can be found in the AW Service Manual.

The current revisions for this manual can be viewed online at: <http://aw-ib.euro.med.ge.com>.

1. Open up the appropriate Advantage Workstation (AW) Service Manual.
2. Proceed to Job Card entitled – Software Options Installation.
3. Follow the procedure using the "Easy Install" method.
4. When complete, perform a save configuration.

4- AW AND REPORTCARD 3 ERRORS

ReportCARD 3 installation logs and AW System errors can be examined and investigated in either the usual AW error logs – such as /var/adm/messages, and/or ~sdc/logfiles/EasyInstallLog .ReportCARD application errors can be found in the application's dedicated logfile: /export/home/sdc/logfiles/ReportCARD.log .

5- ADVANCED SOFTWARE APPLICATION OPERATION AND DEMO EXAM INFORMATION

To test the functionality of the software, first you must load the DEMO exam CD shipped with the product. The 'Standard' AW Software Demo Exam Installation process applies. For ongoing process revision viability, the details will not be printed here, but they can be found in the AW Service Manual.

The current revisions for this manual can be viewed online at: <http://aw-ib.euro.med.ge.com>.

1. Open up the appropriate Advantage Workstation (AW) Service Manual.
2. Proceed to Job Card entitled – Advantage Workstation Demo Exam Reinstallation.
3. Follow the procedure to load the Demo Exams from the CD provided.

Once the application is installed, select any data set in the Browser and click on ReportCARD (after ensuring that the AW software has been restarted). The ReportCARD window will open up with images from Series 1 of the selected exam. The ReportCARD window will take up the whole screen. Select File → Quit to exit the application. This ensures that the ReportCARD has been loaded correctly.

Note :

Preferences : When upgrading the ReportCARD version from RC 2.0 to RC3.0 the normal ranges for all measurements is set to default values. Please verify the normal ranges for all measurements under Preferences with the respective application personnel or specialists at the site or edit them manually as specified in the page 3-10 of the Operators' manual.

6- INSTALLATION FINISH TASKS

GIB Update - The ReportCARD 3 application must be entered into GIB (Global Installed Base) database after the initial successful installation in order to fulfill product-tracking requirements.

FAILURE TO ENTER THE SYSTEM AND ALL SOFTWARE APPLICATIONS INTO GIB will result in GE Healthcare being in violation of FDA tracking compliance, and this system will not be in the affectivity list for future upgrades.

The OLC AW Support line can direct you as how to get GIB updated correctly if you don't know. Here is the global GIB link:

GIB: http://egems.gemedicalsystems.com/gib/gib_entry.jsp

7- SPECIAL CONSIDERATIONS

There are no special considerations to deal with when installing this advanced application option.

8- FRU LIST

The replacement parts for this option are as follows:

Part # 5160890-4 – ReportCARD 3 Software and Documentation CD Rom