



## PM Report Schedule B

### LightSpeed VCT and Discovery CT750 HD

Customer Name: \_\_\_\_\_ Service Record #: \_\_\_\_\_

System ID/ Serial Number: \_\_\_\_\_ Customer Equipment ID: \_\_\_\_\_

Service Manual Document Number: \_\_\_\_\_ Software Version: \_\_\_\_\_

PM Frequency: \_\_\_\_\_ months

System Type (Select:)

#### PRELIMINARY TASKS

##### Last PM Inspection

Verify last PM Inspection

Last Inspection Date :

Last Schedule Completed:

List any Schedule B tasks that were completed earlier

PM Schedule Sys B

#### Review Error Logs

Check the tasks you completed and list any follow-up work scheduled.

| Task                                     | Completed                      | Comments/Follow-up |
|------------------------------------------|--------------------------------|--------------------|
| Table Hydraulic Errors                   | <input type="checkbox"/>       |                    |
| System Application Stop and Start Errors | <input type="checkbox"/>       |                    |
| Linux Reboot and Shutdown Errors         | <input type="checkbox"/>       |                    |
| Record Gantry Rotations                  | <input type="text" value="0"/> |                    |

#### Initial Procedures

##### Check Temperature and Humidity

Record the temperature and humidity. Indicate if these are within spec.

| Room         | Temperature | Humidity | Pass/Fail | Comments/Follow-up |
|--------------|-------------|----------|-----------|--------------------|
| Scan room    |             |          |           |                    |
| Control Room |             |          |           |                    |

#### Console

| Task                                  | Completed                | Comments/Follow-up |
|---------------------------------------|--------------------------|--------------------|
| Console Monitor Cleaning & Inspection | <input type="checkbox"/> |                    |

**General Gantry Inspection**

| Task                      | Pass/Fail | Comments/Follow-up |
|---------------------------|-----------|--------------------|
| Listen for Unusual Noises |           |                    |
| Test Display Lights       |           |                    |
| Test E-stops              |           |                    |
| Test Patient Tilt Sensors |           |                    |

**Collimator Inspection**

| Task                                                   | Pass/Fail | Comments/Follow-up |
|--------------------------------------------------------|-----------|--------------------|
| Aperture and Filter Positions and Operation            |           |                    |
| Protective Window and Opening Inspected/Cleaned Debris |           |                    |

**Tilt Speed and 30° Functional Check**

| Task                    | Pass/Fail | Comments/Follow-up |
|-------------------------|-----------|--------------------|
| Tilt Speed +10      -10 |           |                    |
| Gantry Tilts to +30°    |           |                    |
| Gantry Tilts to -30°    |           |                    |

**NGPDU Initial Checks**

| Task                                   | Pass/Fail | Comments/Follow-up |
|----------------------------------------|-----------|--------------------|
| Highest Line-to-Line Voltage           |           |                    |
| Lowest Line-to-Line Voltage            |           |                    |
| Highest Voltage <1.02 x Lowest Voltage |           |                    |

**Gantry Power Off (Covers Off)****General Tasks**

| Task                     | Pass/Fail | Comments/Follow-up |
|--------------------------|-----------|--------------------|
| Gantry Tilt Zero Check   |           |                    |
| Emergency Off Switch     |           |                    |
| UPS Power Off Inspection |           |                    |



### General Console Cleaning & Inspection

| Task                                       | Pass/Fail | Comments/Follow-up |
|--------------------------------------------|-----------|--------------------|
| Clean Component Filters, Fans & Grills     |           |                    |
| Clean, Inspect, Replace Desktop Components |           |                    |
| Visual Fan Inspection                      |           |                    |
| Seismic Anchor Check (if necessary)        |           |                    |

### Slip Ring & Brush Block Inspection and Maintenance

| Task                          | Pass/Fail | Comments/Follow-up |
|-------------------------------|-----------|--------------------|
| Remove Slip Ring Brush Debris |           |                    |
| Slip Ring Tracks Inspected    |           |                    |
| Brush Blocks Cleaning         |           |                    |
| Brush Tip Inspection          |           |                    |

### CT Belt Tension & Adjustment

| Task                            | Pass/Fail | Comments/Follow-up |
|---------------------------------|-----------|--------------------|
| Inspect and Adjust Tension Belt |           |                    |

### Gantry Filter Cleaning

| Task                                         | Pass/Fail | Comments/Follow-up |
|----------------------------------------------|-----------|--------------------|
| Clean DAS Detector Plenum Filter             |           |                    |
| Clean Detector Face Plate                    |           |                    |
| Clean Gantry Heater Filter                   |           |                    |
| Clean Top Cover Fans                         |           |                    |
| Rotating Path Components & Cables Inspection |           |                    |

### Hydraulic Tilt Pump (Inspect for Leaks)

| Task                               | Pass/Fail | Comments/Follow-up |
|------------------------------------|-----------|--------------------|
| Hydraulic Tilt Cylinder & Fittings |           |                    |
| Hydraulic Tilt Pump & Reservoir    |           |                    |

**NGPDU**

| Task                          | Pass/Fail | Comments/Follow-up |
|-------------------------------|-----------|--------------------|
| Vacuum Covers & Inside PDU    |           |                    |
| Inspect Electronic Components |           |                    |
| Torque all Connectors         |           |                    |
| Verify Tap Connections Match  |           |                    |
| Inspect Seismic Anchor Bolts  |           |                    |

**Power On (Covers Off)****Tube Heat Exchanger & Pump**

| Task                                                             | Select | Comments/Follow-up |
|------------------------------------------------------------------|--------|--------------------|
| Inspect/Clean Tube Heat Exchanger                                |        |                    |
| Inspect/Clean JEDI Inverter Fan                                  |        |                    |
| Tube Heat Exchanger Pump Inspection<br>Mfg date/last replacement |        |                    |

**General Tasks**

| Task                                                       | Select | Comments/Follow-up |
|------------------------------------------------------------|--------|--------------------|
| UPS Powered and Operational                                |        |                    |
| Plenum Fans (5) Operational                                |        |                    |
| Grease Gantry Main Bearing<br>No. Revs since last greasing |        |                    |

**Power Off (Covers On)****Tilt Speed Adjustment**

| Task                  | Select | Comments/Follow-up |
|-----------------------|--------|--------------------|
| Tilt Speed Adjustment |        |                    |

**Power On (Covers On)****General Tasks**

| Task                               | Select | Comments/Follow-up |
|------------------------------------|--------|--------------------|
| Scan Window Inspection/Replacement |        |                    |
| Inspect Rotational Clearances      |        |                    |



## Finalization

### System Scanning Test

| Task                          | QA Phantom Serial # | Location  | Gantry Serial # |
|-------------------------------|---------------------|-----------|-----------------|
|                               |                     |           |                 |
| Scans                         | Scan Failures       | Artifacts | Artifact Types  |
| Series 1 (Scout/Auto Voice)   |                     |           |                 |
| Series 2 (Axial/Tilt)         |                     |           |                 |
| Series 3 (Helical/Auto Voice) |                     |           |                 |

### General Tasks

|                                       |                          |  |
|---------------------------------------|--------------------------|--|
| X-ray On Indicators                   |                          |  |
| Scan Control Pushbuttons Test         |                          |  |
|                                       | <b>Completed</b>         |  |
| GIB Updates Completed                 | <input type="checkbox"/> |  |
| Site Log Updated                      | <input type="checkbox"/> |  |
| PM Paper work completed for this site | <input type="checkbox"/> |  |

### Is follow-up work required?

### Service Scheduled Date:

☐ This inspection passed and meets all GE Healthcare PM specifications

### Comments:

Characters Remaining:

### Tools Used:

| Description | Serial Number | Bar Code | Cal Due Date |                                            |
|-------------|---------------|----------|--------------|--------------------------------------------|
|             |               |          |              | <input type="button" value="Add Tool"/>    |
|             |               |          |              | <input type="button" value="Remove Tool"/> |



**GE Representative:**

FSEs completing this PM

1

FSE Name: \_\_\_\_\_

SSO #: \_\_\_\_\_

FSE SR #: \_\_\_\_\_

FSE Signature Data...

**Customer Signature:** ☐ Yes ☐ No

☐ Final PM Complete

Date  
Completed: \_\_\_\_\_