



PM Report Schedule D

LightSpeed VCT and Discovery CT750 HD

Customer Name: _____ Service Record #: _____

System ID/ Serial Number: _____ Customer Equipment ID: _____

Service Manual Document Number: _____ Software Version: _____

PM Frequency: _____ months

System Type (Select:)

PRELIMINARY TASKS

Last PM Inspection

Verify last PM Inspection

Last Inspection Date :

Last Schedule Completed:

List any Schedule D tasks that were completed earlier

PM Schedule Sys D

Review Error Logs

Check the tasks you completed and list any follow-up work scheduled.

Task	Completed	Comments/Follow-up
Table Hydraulic Errors	☐	
System Application Stop and Start Errors	☐	
Linux Reboot and Shutdown Errors	☐	
Record Gantry Rotations	0	

Initial Procedures

Check Temperature and Humidity

Record the temperature and humidity. Indicate if these are within spec.

Room	Temperature	Humidity	Pass/Fail	Comments/Follow-up
Scan room				
Control room				

Console

Task	Completed	Comments/Follow-up
Console Monitor Cleaning & Inspection	☐	



General Gantry Inspection

Task	Pass/Fail	Comments/Follow-up
Listen for Unusual Noises		
Test Display Lights		
Test E-stops		
Test Patient Tilt Sensors		

General Tasks

Task	Pass/Fail	Comments/Follow-up
Emergency Off Switch		
UPS Power Off Inspection		

General Console Cleaning & Inspection

Task	Pass/Fail	Comments/Follow-up
Clean Component Filters, Fans & Grills		
Clean, Inspect, Replace Desktop Components		
Visual Fan Inspection		
Seismic Anchor Check (if necessary)		

Slip Ring & Brush Block Inspection and Maintenance

Task	Pass/Fail	Comments/Follow-up
Remove Slip Ring Brush Debris		
Slip Ring Tracks Inspected		
Clean Brush Blocks		
Brush Tip Inspection		

Gantry Filter Cleaning

Task	Pass/Fail	Comments/Follow-up
Clean DAS Detector Plenum Filter		
Clean Gantry Heater Filter		
Clean Top Cover Fans		
Clean Detector Face Plate		
Rotating Path Components & Cables Inspection		

**Gantry Power On (Covers Off)****Tube Heat Exchanger & Pump**

Task	Select	Comments/Follow-up
Inspect/Clean Tube Heat Exchanger		
Inspect/Clean JEDI Inverter Fan		
Tube Heat Exchanger Pump Inspection Mfg date/last replacement		

General Tasks

Task	Select	Comments/Follow-up
UPS Powered and Operational		
Plenum Fans (5) Operational		
Grease Gantry Main Bearing No. Revs since last greasing		

Verifications (If required)

NOTE: If you are required to complete these tests, record the readings in Appendix A of this document.

Task	Select	Comments/Follow-up
HV Tank Resistor Verification		
Meter Verification		
HHS Scans		
Task (CT750 HD Only)		
GSI Meter Verification		

Power On (Covers On)**General Tasks**

Task	Select	Comments/Follow-up
Scan Window Inspection/Replacement		
Inspect Rotational Clearances		

Finalization**System Scanning Test**

Task	QA Phantom Serial #	Location	Gantry Serial #
Scans	Scan Failures	Artifacts	Artifact Types
Series 1 (Scout/Auto Voice)			
Series 2 (Axial/Tilt)			



Scans	Scan Failures	Artifacts	Artifact Types
Series 3 (Helical/Auto Voice)			

General Tasks

Task	Select	Comments/Follow-up
X-ray On Indicators		
Scan Control Pushbuttons Test		
Quality Assurance Test		
GE Phantom Serial #		
	Completed	
Site Log Updated	☐	
PM Paper work completed for this site	☐	

Appendix A: Data Forms**Meter Verification**

	DMM Reading	A/D	Delta	Limit
mA Low Scale				+ 3.0 mA
mA High Scale				+ 3.0 mA
kV				Not Applicable

HV Tank Feedback Resistor Verification

Selected kV	On Screen kV	Measured kV
80		
100		
120		
140		
Scan Timer Value		



HHS Data - Tube, HV Tank or Collimator Replacement

CT System

Small Focal Spot		On Screen Data			
kV	mA	kV	mA	5ms mA	Time
80					
100					
120					
140					

Large Focal Spot		On Screen Data			
kV	mA	kV	mA	5ms mA	Time
80					
100					
120					
140					

GSI Meter Verification - Complete for CT750 HD Only

Selected Protocol	Measured kV	Measured mA
GSI Protocol 1		
GSI Protocol 2		
GSI Protocol 3		

Is follow-up work required?

Service Scheduled Date:

☐ This inspection passed and meets all GE Healthcare PM specifications



Comments:

Characters Remaining:

Tools Used:

Description	Serial Number	Bar Code	Cal Due Date	
				Add Tool Remove Tool

GE Representative:

FSEs completing this PM

FSE Name: _____

SSO #: _____

FSE SR #: _____

FSE Signature Data...

Customer Signature: ☐ Yes ☐ No

☐ Final PM Complete

Date
Completed: _____