

ATTENTION!!!!

Note: For the 3.0T Knee coil by Invivo (MCAT M3335LB), if the system is version 14.X (HDx), the coil is plug-and-play. No additional CD is required.

This coil, if required to be installed with a HD (Excite 3) system, the following steps must be adhered to:

- The system must have 12.0_M5_0606.b software.
- If the 12.0_M5_0606.b software is present, the HD TR Knee PA, HD TRKneeServ, HD TRK1chServ Coil configuration files must be uninstalled using the coil configuration manager.
- Insert the 12.0_M5 Service Pack 1 CD in the GOC and install the service pack.
- Using Coil configuration manager, install the HD TR Knee PA, HD TRKneeServ, HD TRK1chServ Coil configuration files [the following link is useful:
<http://3.87.118.28/optec3/common/hdx/content/496626.htm>
- Reboot the system.