

Excite II and Newer GE MR PM Report

Document #: GE-0000-00-PM1

Revision: A

Keywords: pm, ge excite II

Affected Products

Excite II and Newer GE MR

Reference Material

MR Planned Maintenance Data Sheets GE Excite II and Newer

History

Date	Name	Rev	Change descriptions
3/16/2015	Document Control	A	Initial

Approval

Department	Name	Date
Dept. Manager	Ruben Beltran	3/16/2015
Engineering		
QA & Safety		
Document Control	Jhovany Hernandez	3/16/2015

Table of Contents

Preface 3

System Information 4

Measuring Equipment Used 4

Checklist 5

Data 8

System Condition 9

Preface

This document is valid proof of a Preventative Maintenance performed on the system. The Preventative Maintenance must be performed in the specified intervals, and the results are to be documented herein.

Warning: All steps must be performed in accordance with the manufacturer's service manual.

Caution: Be aware of sharps, blood, contrast, high voltage, mechanical movements, radiation exposure and cover safety latch hazards and take appropriate safety measures to prevent personal injury or equipment damage.

This document must be completely filled out by the Field Service Engineer, i.e.:

- All boxes must be filled out. Use **N/A** where appropriate.
- Enter the Serial Number and the Date in the header of each page.
- After completing the check, fill out and sign the last page.

System Information

Customer:	
Model:	
System S/N:	
System ID:	

Measuring Equipment Used

Enter the measuring equipment used in the table below (meters, scopes, phantoms, etc.).

[illegible]

Checklist

Safety & Regulatory

Passed	Cleaned	Adjusted	Replaced	Repaired	Needs Adjustment	Needs Replacement	Needs Repair	Failed	N/A	Action	Frequency
										Check O2 Monitor Operation (if applicable)	Bi-Monthly
										Check Peripheral Pulse Gating Cable	Bi-Monthly
										Check Patient Pneumatic Alert System	Bi-Monthly
										Check RF Power Output*	Bi-Monthly
										Check Emergency Release of Cradle	Bi-Monthly
										Check Emergency Release of Patient Transport	Bi-Monthly
										Check Patient Transport Casters	Bi-Monthly
										Check Patient Armboard Setscrews	Bi-Monthly
										Verify That Cradle Must Be At Home To Lower Table	Bi-Monthly
										Check Cradle Longitudinal Drive Clutch	Bi-Monthly
										Check All Cradle and Carriage Wheels	Bi-Monthly
										Check Emergency Off Circuits	Bi-Monthly
										Check Power Monitor Functional Test	6 Months
										Check Magnet Run Down Unit (ERU, MRU)	Bi-Monthly

Image Quality

Passed	Cleaned	Adjusted	Replaced	Repaired	Needs Adjustment	Needs Replacement	Needs Repair	Failed	N/A	Action	Frequency
										Check X Gradient Calibration	Bi-Monthly
										Check Y Gradient Calibration	Bi-Monthly
										Check Z Gradient Calibration	Bi-Monthly
										Check Shim *	Bi-Monthly
										Check Head SNR *	Bi-Monthly
										Check Body SNR *	Bi-Monthly
										Check All Receive Path Functionality	Bi-Monthly

Other (Available during patient scanning)

Passed	Cleaned	Adjusted	Replaced	Repaired	Needs Adjustment	Needs Replacement	Needs Repair	Failed	N/A	Action	Frequency
										Check RF Cabinet Filters (if applicable)	Bi-Monthly
										Clean RF Cabinet Filters (if applicable)	Bi-Monthly
										Check ACGD/HFD Cabinet Filters	Bi-Monthly
										Clean ACGD/HFD Cabinet Filters	Bi-Monthly
										Check TAC Cabinet Filters (if applicable)	Bi-Monthly
										Clean TAC Cabinet Filters (if applicable)	Bi-Monthly
										Check Excite/System Cabinet Fans	Bi-Monthly
										Clean Excite/System Cabinet Filters	Bi-Monthly
										Verify GOC/System Clock Accuracy	Bi-Monthly
										Check Computer Room Air Temperature	Bi-Monthly

Other (Not available during scanning)

Passed	Cleaned	Adjusted	Replaced	Repaired	Needs Adjustment	Needs Replacement	Needs Repair	Failed	N/A	Action	Frequency
										Check RF Shielding & Door Seals	Bi-monthly
										Check Patient Fans/Blowers & Clean Filter	Bi-monthly
										Check Gradient Coil Cooling Unit	Bi-monthly
										Clean Patient Table	Bi-monthly
										Check Dock Pedal Springs	6 Months
										Run AGP/APS/SCP Diagnostics	Bi-monthly
										Clean GOC and HP Computer Fans & Air Intake	Bi-monthly
										Check and Clean GOC Mouse	Bi-monthly
										Check & Delete Error Log	Bi-monthly
										Remove Excess "core" Files	Bi-monthly
										Check Gradient Cable Connections & Support	6 Months
										Check Magnet Room RF Cables & Connections	6 Months
										Check Patient Transport Hydraulic Fluid Level	6 Months
										Inspect PDU Power Connection	6 Months
										Check Patient Alignment Lights	Bi-monthly
										Check ISO Z Table Movement	Bi-monthly
										Record Magnet Pressure *	Bi-monthly
										Record Shield Cooler Hours *	Bi-monthly
										Calculate Lost Hours of Shield Cooler *	Bi-monthly
										Record Helium Level *	Bi-monthly
										Record Static Shield Cooler Compressor Pressure *	Bi-monthly
										Back-Up System Level Information	Bi-monthly
										Back-Up Customer Protocols	Bi-monthly

Data

Magnet Data

Magnet Pressure:

Lhe Level:

Compressor Pressure:

Compressor Run Hours:

Time Run Hours Logged:

Lost Hours:

Calibration Data

Inhomogeneity:

X Gradcal:

Y Gradcal:

Z Gradcal:

Head SNR:

Body SNR:

Shim Phantom TG:

DQA Isocenter TG:

**Note: If System is TwinSpeed record Whole & Zoom Values*

*System Serial #:*_____ *Date:*_____

System Condition

After performing scheduled PM, an evaluation of the system is performed with the following result (check one):

- ☐ No deficiencies. The image quality test resulted in no differences from required reference values.
- ☐ Slight deficiencies. There is no effect on operation; however they should be corrected preventively. The image quality test resulted in no differences from required reference values.
- ☐ Serious deficiencies. This must be corrected before operation of the system is permitted.

FSE

Name:	Date:
--------------	--------------

Notes (Additional work to be performed, partial PM, customer requests):

Customer

Name:	Date:
--------------	--------------
