



Service Request Number _____

Site Name: _____

Month: _____

Procedure

9x - 12x PM Checklist

OK
Cleaned
Adjusted
Needs
Correction
N/A

PM Checklist

Talk With Technologist About System Operation..... ☐ ☐ ☐ ☐ ☐

Look at Image Quality of Previous Exams..... ☐ ☐ ☐ ☐ ☐

Check and Set System Date and Time..... ☐ ☐ ☐ ☐ ☐

SAFETY AND REGULATORY

SYSTEM - Check Gating Cable/ Respiratory Bellows / PG Probe..... ☐ ☐ ☐ ☐ ☐

SYSTEM - Check Pneumatic Patient Alert System..... ☐ ☐ ☐ ☐ ☐

SYSTEM - Check Emergency Release of Cradle and Patient Transport..... ☐ ☐ ☐ ☐ ☐

SYSTEM - Check Patient Transport Casters and Arm board Set Screws..... ☐ ☐ ☐ ☐ ☐

MAGNET - Check Magnet Bore for Metal Objects and Remove..... ☐ ☐ ☐ ☐ ☐

MAGNET - Test Magnet Emergency Rundown (Test Battery, Htr LED)..... ☐ ☐ ☐ ☐ ☐

MAGNET - Check for Condensation or Leaks on Top and Around Magnet..... ☐ ☐ ☐ ☐ ☐

MAGNET - Record Cryostat Pressure and Flow Rates..... ☐ ☐ ☐ ☐ ☐

He Level 1 _____ Pressure _____

MAGNET - Compressor Pressure/Hours..... ☐ ☐ ☐ ☐ ☐

Pressure _____ Compressor -Hrs _____

MAGNET - Coldhead and Magnet Temps..... ☐ ☐ ☐ ☐ ☐

Coldhead_Ru0 _____ Recon_Si410 _____ Recon_Ru0 _____ Shield_Si410 _____ Water Temp _____ Water flow _____

Clean, Test and Check

Check Correlated Noise (RF Integrity / RF Door Seal, top hatch, clean, Check)..... ☐ ☐ ☐ ☐ ☐

Clean Light-Weight Cradle Wheels..... ☐ ☐ ☐ ☐ ☐

Check Table Functionality (Clutch, Alignment, Undock)..... ☐ ☐ ☐ ☐ ☐

Check Body Chiller Temp and Fluid..... ☐ ☐ ☐ ☐ ☐

Water Temp and flow _____

Clean Patient Blower Filter..... ☐ ☐ ☐ ☐ ☐

Clean All Cabinet Filters..... ☐ ☐ ☐ ☐ ☐

Clean Console area (Vacuum PC, Octane, or GOC)..... ☐ ☐ ☐ ☐ ☐

Semi-Annual or First PM Safety Checks

Month: _____

Test Emergency off Buttons..... ☐ ☐ ☐ ☐ ☐

Test Emergency Stop Buttons..... ☐ ☐ ☐ ☐ ☐

Inspect all System Cables and Connections for Tightness and Security..... ☐ ☐ ☐ ☐ ☐

Check Condition and Placement of Warning Labels..... ☐ ☐ ☐ ☐ ☐

Check Condition and Placement of 0.5 mT Exclusion Warnings..... ☐ ☐ ☐ ☐ ☐

Check Magnet Exhaust Vent and Vent Glass..... ☐ ☐ ☐ ☐ ☐

Clean Water Chiller Filter - If filter is really dirty then need to clean on next PM and note status below..... ☐ ☐ ☐ ☐ ☐

Tip Seal and Solenoid Valve Replacement (TwinSpeed)..... ☐ ☐ ☐ ☐ ☐

Change Out All Patient Cradle Wheels..... ☐ ☐ ☐ ☐ ☐

Inspect and Clean Erbttec Interior..... ☐ ☐ ☐ ☐ ☐

Check PDU power and cables..... ☐ ☐ ☐ ☐ ☐

Check gradient cables..... ☐ ☐ ☐ ☐ ☐

Check gradient filters and PHPS..... ☐ ☐ ☐ ☐ ☐



Service Request Number _____ Site Name: _____ Month: _____

Procedure	9x - 12x PM Checklist	OK	Cleaned	Adjusted	Needs Correction	N/A
IMAGE QUALITY						
Check ISO-Center (Check Alignment Lights).....	TG _____ AX _____ iso Z _____	☐	☐	☐	☐	☐
Check Gradient Calibration: X,Y,Z-axis.....	Cfxfull - Old _____ New _____ Cfyfull - Old _____ New _____ Cfzfull - Old _____ New _____	☐	☐	☐	☐	☐
Check Gradient Calibration zoom mode (TRM ONLY): X,Y,Z-axis.....	Cfxfull - Old _____ New _____ Cfyfull - Old _____ New _____ Cfzfull - Old _____ New _____	☐	☐	☐	☐	☐
Check SNR and RF output- Head.....	TG _____ Signal Noise SNR	☐	☐	☐	☐	☐
Check SNR and RF output- Body.....	TG _____ Signal Noise SNR	☐	☐	☐	☐	☐
Check (CTLTOP).....	Display CV (saveinter) check for grainy images	☐	☐	☐	☐	☐
Check (CTLBOT).....	Display CV (saveinter) check for grainy images	☐	☐	☐	☐	☐
Check LVShim Localizer (DC off-sets first).....	TG _____	☐	☐	☐	☐	☐

LVshim Calibrations

Shim Values Before	X	Y	Z
Shim Values After			
Inhomogeneity _____	Z1(1,0) _____	X(1,1) _____	Z2X(3,1) _____
	Z2(2,0) _____	Y(1,-1) _____	Z2Y(3,-1) _____
	Z3(3,0) _____	ZX(2,1) _____	Zxy2(3,2) _____
	Z4(4,0) _____	ZY(2,-1) _____	ZXY(3,-2) _____
	Z5(5,0) _____	X2-Y2(2,2) _____	X3(3,3) _____
	Z6(6,0) _____	XY(2,-2) _____	Y3(3,-3) _____

Zoom Mode (After)

Inhomogeneity _____

Z1(1,0) _____ X(1,1) _____ Z2X(3,1) _____

Z2(2,0) _____ Y(1,-1) _____ Z2Y(3,-1) _____

Z3(3,0) _____ ZX(2,1) _____ Zxy2(3,2) _____

Z4(4,0) _____ ZY(2,-1) _____ ZXY(3,-2) _____

Z5(5,0) _____ X2-Y2(2,2) _____ X3(3,3) _____

Z6(6,0) _____ XY(2,-2) _____ Y3(3,-3) _____

SuperCon Shim Currents

Ax1 _____	Ax2 _____	Ax3 _____	Ax4 _____	Ax5 _____	Ax6 _____
T1-1 _____	T1-2 _____	T1-3 _____	T1-4 _____	T1-5 _____	T1-6 _____
T2-1 _____	T2-2 _____	T2-3 _____	T2-4 _____	T2-5 _____	T2-6 _____

Center Frequency _____

FINAL CHECKS

Check and Delete Error Message Log.....	☐	☐	☐	☐	☐
Do Save INFO, Save GUI and Save Protocol.....	☐	☐	☐	☐	☐
RESTORATION - DQA Phantom scan, Remove Test Scans, and Check Cabinet Doors / Covers.....	☐	☐	☐	☐	☐
RESTORATION - Shut down system pre-site procedure.....	☐	☐	☐	☐	☐
RESTORATION - Update site log book, E-mail update with any follow ups an status.....	☐	☐	☐	☐	☐

Notes