



Service Request Number _____ Site Name: _____ Month: _____

Procedure PM Checklist

OK
Cleaned
Adjusted
Needs
Correction

PM Checklist

Talk With Technologist About System Operation..... ☐ ☐ ☐ ☐

Look at Image Quality of Previous Exams..... ☐ ☐ ☐ ☐

Check and Set System Date and Time..... ☐ ☐ ☐ ☐

SAFETY AND REGULATORY

CONSOLE – Check e-stops(4-gantry, 1on PDU, 1 on SCIM)..... ☐ ☐ ☐ ☐

CONSOLE- Check Audio system, from console to Gantry..... ☐ ☐ ☐ ☐

Clean, Test and Check

CONSOLE -Clean and Vacuum inside and out of console..... ☐ ☐ ☐ ☐

CONSOLE -Clean and Vacuum CPU (open it up)..... ☐ ☐ ☐ ☐

CONSOLE -Check all fans..... ☐ ☐ ☐ ☐

CONSOLE -Clean/ Vacuum MOD Tower..... ☐ ☐ ☐ ☐

CONSOLE -Clean and check Keyboard, mouse, and SCIM..... ☐ ☐ ☐ ☐

TABLE -Clean and vacuum Table ☐ ☐ ☐ ☐

TABLE -Check and verify safety switches are good and working..... ☐ ☐ ☐ ☐

TABLE -Check covers, make sure they're mounted properly..... ☐ ☐ ☐ ☐

TABLE -Check to make sure cradle is not riding on covers..... ☐ ☐ ☐ ☐

TABLE -Clean rollers on cradle drive..... ☐ ☐ ☐ ☐

TABLE -Clean bottom of cradle, make sure no tape is on..... ☐ ☐ ☐ ☐

GANTRY -Check covers are on good, locked in..... ☐ ☐ ☐ ☐

GANTRY -Check top fans..... ☐ ☐ ☐ ☐

GANTRY -Take off covers and check power switches, are they properly working..... ☐ ☐ ☐ ☐

GANTRY -Disable 120, axial, and H,V, switches- make sure gantry is off..... ☐ ☐ ☐ ☐

GANTRY -Turn on 120 switch, does gantry brake turn off?..... ☐ ☐ ☐ ☐

GANTRY -Take off front cover..... ☐ ☐ ☐ ☐

GANTRY -Take off mylar window clean and check for damage..... ☐ ☐ ☐ ☐

GANTRY -Make sure power is off and vacuum inside front half of gantry..... ☐ ☐ ☐ ☐

GANTRY -Clean and vacuum DAS filters..... ☐ ☐ ☐ ☐

GANTRY -Turn 120 on and check for fans working properly..... ☐ ☐ ☐ ☐

GANTRY -Remove slip ring covers. Clean covers and slip ring block..... ☐ ☐ ☐ ☐

GANTRY -Check brushes for wear, replace any as needed..... ☐ ☐ ☐ ☐

GANTRY -Grease gantry as needed(2,000,000 revolutions)..... ☐ ☐ ☐ ☐

Last Lubed _____

GANTRY -Replace covers and check gantry rotation..... ☐ ☐ ☐ ☐

PDU -Check fans..... ☐ ☐ ☐ ☐

PDU -Clean and vacuum..... ☐ ☐ ☐ ☐



Service Request Number _____ Site Name: _____ Month: _____

Procedure

PM Checklist

FINAL CHECKS

	OK	Cleaned	Adjusted	Needs Correction
Check Error Message Log.....	☐	☐	☐	☐
Do Save State.....	☐	☐	☐	☐
RESTORATION - Do scout, helical, & axial scan, Remove Test Scans, and Covers.....	☐	☐	☐	☐
RESTORATION - Shut down system pre-site procedure.....	☐	☐	☐	☐
RESTORATION - Update site log book, E-mail update with any follow up's and status.....	☐	☐	☐	☐

Notes:



Service Request Number _____ Site Name: _____ Month: _____

Procedure	PM Checklist	OK	Cleaned	Adjusted	Needs Correction
1st Quarter Checks					
TABLE -Vacuum power supplies inside table.....		☐	☐	☐	☐
TABLE- MEASURE TABLE POWER SUPPLY VALUE.....		☐	☐	☐	☐
TABLE -Check fans inside table.....		☐	☐	☐	☐
TABLE -Check clutch on cradle asm.....		☐	☐	☐	☐
TABLE -Grease IS.....		☐	☐	☐	☐

2nd Quarter Checks					
GANTRY -Check to make sure heat exchange pump is working.....		☐	☐	☐	☐
GANTRY -Turn axial and check belt for any damage.....		☐	☐	☐	☐
GANTRY -Check Thermistor.....		☐	☐	☐	☐
GANTRY -Clean and check for pitting or arch marks.....		☐	☐	☐	☐

3rd Quarter Checks					
GANTRY -Check Tilt functions.....		☐	☐	☐	☐
GANTRY -Check speed of Tilt.....		☐	☐	☐	☐
GANTRY -Clean Tube Heat Exchange and Heat Exchange Pump.....		☐	☐	☐	☐
GANTRY -HV Tank feedback Resistor.....		☐	☐	☐	☐
GANTRY -Meter Verification.....		☐	☐	☐	☐

4th Quarter Checks					
GANTRY -Check bolts for proper torque.....		☐	☐	☐	☐
GANTRY -Check oil inside wells, any carbon tracks, discolor.....		☐	☐	☐	☐
GANTRY -Check HV cables, pins, sleeve.....		☐	☐	☐	☐
GANTRY -Clean Collimator.....		☐	☐	☐	☐
GANTRY -Clean wells and fill with clean oil		☐	☐	☐	☐
GANTRY -Check alignment lights.....		☐	☐	☐	☐

Notes: