

Service Call Number _____

Bi-Monthly PM Checklist



Site Name: _____

Planned Maintenance Checklist

3200 NW 27th Ave Ste. 100
Pompano Beach, FL
33069
Office: 888-400-3979
aisserv.com

FE Name: _____

Month: _____

Procedure

Print Form

Reset Form

OK

Cleaned

Adjusted

Needs
Correction

PM Checklist

Talk With Technologist About System Operation.....	☐	☐	☐	☐
Look at Image Quality of Previous Exams.....	☐	☐	☐	☐
Check and Set System Date and Time.....	☐	☐	☐	☐

SAFETY AND REGULATORY

SYSTEM - Check Gating Cable/ Respiratory Bellows / PG Probe.....	☐	☐	☐	☐
SYSTEM - Check Pneumatic Patient Alert System.....	☐	☐	☐	☐
SYSTEM - Check Emergency Release of Cradle and Patient Transport.....	☐	☐	☐	☐
SYSTEM - Check Patient Transport Casters and Arm board Set Screws.....	☐	☐	☐	☐
MAGNET - Check Magnet Bore for Metal Objects and Remove.....	☐	☐	☐	☐
MAGNET - Test Magnet Emergency Rundown (Test Battery, Htr LED).....	☐	☐	☐	☐
MAGNET - Check for Condensation or Leaks on Top and Around Magnet.....	☐	☐	☐	☐
MAGNET - Record Cryostat Pressure and Flow Rates.....	☐	☐	☐	☐

He Level 1 _____

Pressure _____

MAGNET - Compressor Pressure.....	☐	☐	☐	☐
-----------------------------------	--------------------------	--------------------------	--------------------------	--------------------------

Pressure _____

MAGNET - Coldhead and Magnet Temps.....	☐	☐	☐	☐
---	--------------------------	--------------------------	--------------------------	--------------------------

Coldhead_Ru0 _____

Recon_Si410 _____

Recon_Ru0 _____

Shield_Si410 _____

Water Temp _____

Water flow _____

Clean, Test and Check

Check Correlated Noise (RF Integrity / RF Door Seal, top hatch, clean, Check).....	☐	☐	☐	☐
Clean Light-Weight Cradle Wheels.....	☐	☐	☐	☐
Check Table Functionality (Clutch, Alignment, Undock).....	☐	☐	☐	☐
Check Body Chiller Temp and Fluid.....	☐	☐	☐	☐
Clean Patient Blower Filter.....	☐	☐	☐	☐
Clean All Cabinet Filters.....	☐	☐	☐	☐
Clean Console area (Vacuum PC, Octane, or GOC).....	☐	☐	☐	☐

Water Temp and flow _____

Semi-Annual or First PM Safety Checks

Month: _____

Test Emergency off Buttons.....	☐	☐	☐	☐
Test Emergency Stop Buttons.....	☐	☐	☐	☐
Inspect all System Cables and Connections for Tightness and Security.....	☐	☐	☐	☐
Check Condition and Placement of Warning Labels.....	☐	☐	☐	☐
Check Condition and Placement of 0.5 mT Exclusion Warnings.....	☐	☐	☐	☐
Check Magnet Exhaust Vent and Vent Glass.....	☐	☐	☐	☐
Clean Water Chiller Filter - If filter is really dirty then need to clean on next PM and note status below.....	☐	☐	☐	☐
Tip Seal and Solenoid Valve Replacement (Twinspeed).....	☐	☐	☐	☐
Change Out All Patient Cradle Wheels.....	☐	☐	☐	☐
Inspect and Clean Erbttec Interior.....	☐	☐	☐	☐
Check PDU power and cables.....	☐	☐	☐	☐
Check gradient cables.....	☐	☐	☐	☐
Check gradient filters and PHPS.....	☐	☐	☐	☐

Service Call Number _____

Bi-Monthly PM Checklist



Site Name: _____

Planned Maintenance Checklist

3200 NW 27th Ave Ste. 100

Pompano Beach, FL

33069

Office: 888-400-3979

aisserv.com

FE Name: _____

Month: _____

Procedure

IMAGE QUALITY

OK
Cleaned
Adjusted
Needs
Correction

Check ISO-Center (Check Alignment Lights).....☐ ☐ ☐ ☐

TG _____ AX _____ iso Z _____

Check Gradient Calibration: X,Y,Z-axis.....☐ ☐ ☐ ☐

Cfxfull - Old _____ New _____ Cfyfull - Old _____ New _____ Cfzfull - Old _____ New _____

Check Gradient Calibration zoom mode (TRM ONLY): X,Y,Z-axis.....☐ ☐ ☐ ☐

Cfxfull - Old _____ New _____ Cfyfull - Old _____ New _____ Cfzfull - Old _____ New _____

Check SNR and RF output- Head.....☐ ☐ ☐ ☐

TG _____ Signal Noise SNR

Check SNR and RF output- Body.....☐ ☐ ☐ ☐

TG _____ Signal Noise SNR

Check (CTLTOP).....Display CV (saveinter) check for grainy images ☐ ☐ ☐ ☐

TG _____

Check (CTLBOT).....Display CV (saveinter) check for grainy images ☐ ☐ ☐ ☐

TG _____

Check LVShim Localizer (DC off-sets before centering phatom).....☐ ☐ ☐ ☐

TG _____

LVshim Calibrations

Shim Values Before

X _____ Y _____ Z _____

Shim Values After

X _____ Y _____ Z _____

Inhomogeneity _____

Z1(1,0)	X(1,1)	Z2X(3,1)
Z2(2,0)	Y(1,-1)	Z2Y(3,-1)
Z3(3,0)	ZX(2,1)	Zxy2(3,2)
Z4(4,0)	ZY(2,-1)	ZXY(3,-2)
Z5(5,0)	X2-Y2(2,2)	X3(3,3)
Z6(6,0)	XY(2,-2)	Y3(3,-3)

TRM ONLY

Zoom Mode (After)

X _____ Y _____ Z _____ Z2 _____

Inhomogeneity _____

Z1(1,0)	X(1,1)	Z2X(3,1)
Z2(2,0)	Y(1,-1)	Z2Y(3,-1)
Z3(3,0)	ZX(2,1)	Zxy2(3,2)
Z4(4,0)	ZY(2,-1)	ZXY(3,-2)
Z5(5,0)	X2-Y2(2,2)	X3(3,3)
Z6(6,0)	XY(2,-2)	Y3(3,-3)

Supercon Shim Currents

Ax1	Ax2	Ax3	Ax4	Ax5	Ax6
T1-1	T1-2	T1-3	T1-4	T1-5	T1-6
T2-1	T2-2	T2-3	T2-4	T2-5	T2-6

Center Frequency _____

FINAL CHECKS

Check and Delete Error Message Log.....☐ ☐ ☐ ☐

Do Save INFO, Save GUI and Save Protocol.....☐ ☐ ☐ ☐

RESTORATION - DQA Phantom scan, Remove Test Scans, and Check Cabinet Doors / Covers.....☐ ☐ ☐ ☐

RESTORATION - Shut down system pre site procedure.....☐ ☐ ☐ ☐

RESTORATION - Update site log book, E-mail update with any follow up's and status.....☐ ☐ ☐ ☐

Notes