

DVMR PM Checklist (450, 450W, 750, 750W)

ADVANCED
IMAGING SYSTEMS

DEDICATED
IMAGING SOLUTIONS

Service Request Number _____

Site Name: _____

FE Name: _____

Planned Maintenance Checklist

Month: _____

3200 NW 27th Ave Ste. 100
Pompano Beach, FL
33069
Office: 888-400-3979
aisserv.com

Procedure

PM Checklist

	OK	Cleaned	Adjusted	Needs Correction
Talk With Technologist About System Operation.....	☐	☐	☐	☐
Look at Image Quality of Previous Exams.....	☐	☐	☐	☐
Check and Set System Date and Time.....	☐	☐	☐	☐

SAFETY AND REGULATORY

SYSTEM - Check Gating Cable/ Respiratory Bellows / PG Probe.....	☐	☐	☐	☐
SYSTEM - Check Pneumatic Patient Alert System.....	☐	☐	☐	☐
SYSTEM - Check Emergency Release of Cradle and Patient Transport.....	☐	☐	☐	☐
SYSTEM - Check Patient Transport Casters and Arm board Set Screws.....	☐	☐	☐	☐
MAGNET - Check Magnet Bore for Metal Objects and Remove.....	☐	☐	☐	☐
MAGNET - Test Magnet Emergency Rundown (Test Battery, Htr LED).....	☐	☐	☐	☐
MAGNET - Check for Condensation or Leaks on Top and Around Magnet.....	☐	☐	☐	☐
MAGNET - Record Cryostat Pressure and Flow Rates.....	☐	☐	☐	☐

He Level _____

Pressure _____

MAGNET - Compressor Pressure (Dynamic)..... = _____ PSIG ☐ ☐ ☐ ☐

MAGNET - Coldhead and Magnet Temps..... ☐ ☐ ☐ ☐

Coldhead_Ru0 _____

Recon_Si410 _____

Recon_Ru0 _____

Shield_Si410 _____

Water Temp _____

Water flow _____

Clean, Test and Check

Check Correlated Noise (RF Integrity / RF Door Seal, top hatch, clean, Check).....	☐	☐	☐	☐
Clean Light-Weight Cradle Wheels.....	☐	☐	☐	☐
Check Table Functionality (Belt adjustment, Alignment, Undock).....	☐	☐	☐	☐
Clean Pen Wall Body Coil Filter.....	☐	☐	☐	☐
Clean All Cabinet Filters.....	☐	☐	☐	☐
Clean Console area (Vacuum PC).....	☐	☐	☐	☐

Semi-Annual or First PM Safety Checks

Month: _____

Test Emergency off Buttons.....	☐	☐	☐	☐
Test Emergency Stop Buttons.....	☐	☐	☐	☐
Inspect all System Cables and Connections for Tightness and Security.....	☐	☐	☐	☐
Check Condition and Placement of Warning Labels.....	☐	☐	☐	☐
Check Condition and Placement of 0.5 mT Exclusion Warnings.....	☐	☐	☐	☐
Check Magnet Exhaust Vent and Vent Glass.....	☐	☐	☐	☐
Cryocooler Filter check and clean in HEC cabinet	☐	☐	☐	☐
Change Out All Patient Cradle Wheels.....	☐	☐	☐	☐
Check PDU power and cables.....	☐	☐	☐	☐
Check gradient cables.....	☐	☐	☐	☐
Check gradient filters	☐	☐	☐	☐

Annual Procedure

Month: _____

HEC Cabinet Deionization..... ☐ ☐ ☐ ☐

DVMR PM Checklist (450, 450W, 750, 750W)



Service Request Number _____

Site Name: _____

FE Name: _____

Planned Maintenance Checklist

Month: _____

3200 NW 27th Ave Ste. 100
Pompano Beach, FL 33069
Office: 888-400-3979
aisserv.com

Procedure

IMAGE QUALITY

	OK	Cleaned	Adjusted	Needs Correction
Check ISO-Center (Check Alignment Lights).....	☐	☐	☐	☐
TG _____ AX _____ iso Z _____				
Check Gradient Calibration: X,Y,Z-axis.....	☐	☐	☐	☐
Cfxfull - Old _____ New _____ Cfyfull - Old _____ New _____ Cfzfull - Old _____ New _____				
Check SNR and RF output- Head.....	☐	☐	☐	☐
TG _____ Signal Noise SNR				
Check SNR and RF output- Body.....	☐	☐	☐	☐
TG _____ Signal Noise SNR				
Check PA Coil Top.....	☐	☐	☐	☐
TG _____				
Check PA Coil Bottom.....	☐	☐	☐	☐
TG _____				
Check LVShim Localizer	☐	☐	☐	☐
TG _____				

LVshim Calibrations

Shim Values Before X _____ Y _____ Z _____

Shim Values After X _____ Y _____ Z _____

Center Frequency _____

PASS

FAIL

Inhomogeneities	45x45 p2p (ppm) _____ [15.00]	45x45 rms (ppm) _____ [1.25]
	40x40 p2p (ppm) _____ [6.00]	40x40 rms (ppm) _____ [1.00]
	30x30 p2p (ppm) _____ [2.00]	30x30 rms (ppm) _____ [0.50]
	20x20 p2p (ppm) _____ [1.00]	20x20 rms (ppm) _____ [0.10]

HEC Cabinet

Facility Water

Temperature _____ Incoming Pressure _____ Outgoing Pressure _____ Pressure Differential _____

Gradient Coil

IN Temperature _____ OUT Temperature _____ Flow _____ Pressure _____ Deionization _____ kOhm/cm [4.00<x<14]

Power Electronics

IN Temperature _____ OUT Temperature _____ Flow _____ Pressure _____

Cryogen Compressor

Temperature _____ Flow _____

FINAL CHECKS

Check and Delete Error Message Log.....	☐	☐	☐	☐
Do Save INFO, Save GUI and Save Protocol.....	☐	☐	☐	☐
RESTORATION - DQA Phantom scan, Remove Test Scans, and Check Cabinet Doors / Covers.....	☐	☐	☐	☐
RESTORATION - Shut down system per site procedure.....	☐	☐	☐	☐
RESTORATION - Update site logbook, E-mail update with any follow ups and status.....	☐	☐	☐	☐

Notes