



Service Request Number _____

Site Name: _____

Month: _____

Procedure Siemens Espree PM Checklist

PM Checklist

	OK	Cleaned	Adjusted	Needs Correction	N/A
Talk With Technologist About System Operation.....	☐	☐	☐	☐	☐
Look at Image Quality of Previous Exams.....	☐	☐	☐	☐	☐
Check and Set System Date and Time.....	☐	☐	☐	☐	☐

MAINTENANCE INSTRUCTIONS

COOLING - Cooling system general checks	☐	☐	☐	☐	☐
ACC CABINET - Check the fan assembly.....	☐	☐	☐	☐	☐
GRADIENTS - Check gradient filter fans.....	☐	☐	☐	☐	☐
MAGNET - Check Magnet Bore for Metal Objects and Remove.....	☐	☐	☐	☐	☐
VACUUM PUMP (option) - Replace Vacuum pump filter..... Date last replaced	☐	☐	☐	☐	☐
SEP (option) - Clean the strainer of the primary and secondary water circuit	☐	☐	☐	☐	☐
MAGNET - Check for Condensation or Leaks on Top and Around Magnet.....	☐	☐	☐	☐	☐
COLD HEAD - Replace cold head every 24 months..... Date last replaced	☐	☐	☐	☐	☐
MAGNET - Compressor Pressure.....	☐	☐	☐	☐	☐
COMPRESSOR (CSW-71 2yrs/F70 3yrs) - Replace Absorber..... Date last replaced	☐	☐	☐	☐	☐
Check Correlated Noise (RF Integrity / RF Door Seal, top hatch, clean, Check).....	☐	☐	☐	☐	☐
COMPUTER - Clean and check the PC.....	☐	☐	☐	☐	☐
TFT MONITOR - Check and clean TFT monitor.....	☐	☐	☐	☐	☐
PATIENT TABLE - Check patient table drive unit.....	☐	☐	☐	☐	☐
SOFTWARE - Delete MR save log files.....	☐	☐	☐	☐	☐

FINAL CHECKS

Check and Delete Error Message Log.....	☐	☐	☐	☐	☐
Do Save site specific data and Save Protocol.....	☐	☐	☐	☐	☐
RESTORATION - DQA Phantom scan, Remove Test Scans, and Check Cabinet Doors / Covers.....	☐	☐	☐	☐	☐
RESTORATION - Shut down system.....	☐	☐	☐	☐	☐
RESTORATION - Update site log book, E-mail update with any follow ups and status.....	☐	☐	☐	☐	☐

Notes