

Airis Scanners Quarterly PM Report



Document #: HMS-00-00-PM1

Revision: A

Affected Products

Airis Elite, Airis 2, Airis 1

History

Date	Name	Rev	Change Descriptions
6/18/15	Document Control	A	Initial release

Approval

Department	Name	Date
Dept. Manager		
Engineering		
QA & Safety		
Document Control	Jhovany Hernandez	6/18/15

Preface

This document is valid proof of a Preventative Maintenance performed on the system. The Preventative Maintenance must be performed in the specified intervals, and the results are to be documented herein.

This section addresses planned maintenance that should be performed by service personnel. Some of it is also listed in the Operator's Manual for the user to perform on daily and weekly basis.

However, service personnel should make sure that these procedures have been followed and if not, include them on the periodic visits.

Warning: All steps must be performed in accordance with the manufacturer's service manual.

This document must be completely filled out by the Field Service Engineer, i.e.:

- All boxes must be filled out. Use N/A where appropriate.
- Enter the Serial Number and the Date in the header of each page.
- After completing the check, fill out and sign the last page.

Serial #:

Date:

System Information

Customer:

Model:

System S/N:

System ID:

OEM System ID:

Measuring Equipment Used

Enter the measuring equipment used in the table below (meters, scopes, phantoms, etc.):

Equipment:

Serial #:

Serial #:

Date:

Checklist

Console

Pass	Fail	N/A	Action	Frequency
☐	☐	☐	Inspect overall of console.	Quarterly
☐	☐	☐	Clean keyboard.	Quarterly
☐	☐	☐	Clean mouse.	Quarterly
☐	☐	☐	Clean monitor.	Quarterly
☐	☐	☐	Clean host computer.	Quarterly
☐	☐	☐	Clean MOD/DVD.	Quarterly
☐	☐	☐	Clean and inspect all fans.	Quarterly

MR Cabinet

Pass	Fail	N/A	Action	Frequency
☐	☐	☐	Clean and inspect all rear fans.	Quarterly
☐	☐	☐	Clean and inspect gradient amp fans.	Quarterly
☐	☐	☐	Clean and inspect RF amp fan(s).	Quarterly
☐	☐	☐	Clean filters.	Quarterly
☐	☐	☐	Perform general cleaning inside cabinet.	Quarterly

Scanner Gantry

Pass	Fail	N/A	Action	Frequency
☐	☐	☐	Inspect and remove any metal inside of gantry.	Quarterly
☐	☐	☐	Clean covers.	Quarterly

Serial #:

Date:

Scanner Table

Pass	Fail	N/A	Action	Frequency
☐	☐	☐	Perform general cleaning.	Quarterly
☐	☐	☐	Inspect and remove any metal.	Quarterly
☐	☐	☐	Drive table in/out/down/up.	Quarterly

Coils

Pass	Fail	N/A	Action	Frequency
☐	☐	☐	Inspect coil cables and connectors.	Quarterly
☐	☐	☐	Tighten any loose screws.	Quarterly

System Performance

Pass	Fail	N/A	Action	Frequency
☐	☐	☐	Perform a QA scan on QD head coil.	Quarterly
☐	☐	☐	Record RF amp forward/reflect: Specs: Forward=70+/-15; Reflect=5+/-10	Quarterly

Coils Performance

Pass	Fail	N/A	Action	Frequency
☐	☐	☐	QD head coil, SNR/tune/H1:	Quarterly
☐	☐	☐	QD knee coil, SNR/tune/H1:	Quarterly
☐	☐	☐	QD shoulder coil, SNR/tune/H1:	Quarterly
☐	☐	☐	QD C-spine coil, SNR/tune/H1:	Quarterly
☐	☐	☐	QD medium body coil, SNR/tune/H1:	Quarterly
☐	☐	☐	QD large body coil, SNR/tune/H1:	Quarterly
☐	☐	☐	QD X-large body coil, SNR/tune/H1:	Quarterly
☐	☐	☐	Signal channel coil, SNR/tune/H1:	Quarterly
☐	☐	☐	Any extra coils? Name of coil, SNR/tune/H1:	Quarterly

Serial #:

Date:

Final Step

Pass	Fail	N/A	Action	Frequency
☐	☐	☐	Save system backup.	Quarterly
☐	☐	☐	Restore all covers.	Quarterly

Serial #:

Date:

System Condition

After performing scheduled PM, an evaluation of the system is performed with the following result (check one):

☐

No deficiencies. The image quality test resulted in no differences from required reference values.

☐

Slight deficiencies. There is no effect on operation; however they should be corrected preventively. The image quality test resulted in no differences from required reference values.

☐

Serious deficiencies. This must be corrected before operation of the system is permitted.

Field Service Engineer

Name:

Date:

Notes (Additional work to be performed, partial PM, customer requests):

Customer

Name:

Date: