

SIEMENS MRI System PM Checklist

Site: MRI & CT Diagnostic			Date: 12/27/22									
System ID #: MR2432			Service Report #: 12008									
Service Engineer: Ray Jacobsen												
Procedure		Serial Number for tools	Pass	Cleaned	Adjusted	Replaced	Repaired	Needs Adjustment	Needs Replacement	Needs Repair	Does not Apply	Failed
Safety and Regulatory												
S.1	Table collision switch and stop button checks		☒	NA	NA	☐	☐	☐	☐	☐	☐	☐
S.2	Check patient trolley frame, wheels, brake & hydraulic		☐	☐	☐	NA	☐	☐	☐	☐	☒	☐
S.3	Check Phantoms for cracks and refill if necessary		☒	NA	NA	NA	☐	☐	☐	☐	☐	☐
S.4	Service Turret free of water & Magnet pressure in Spec.		☒	NA	NA	NA	☐	NA	NA	☐	☐	☐
S.5	Symphony annual grease cartridge replacement Date of cartridge installation		☐	NA	NA	☐	NA	NA	☐	NA	☒	☐
S.6	Patient Blower & Checking the air filter		☒	☒	NA	☐	NA	NA	☐	NA	☐	☐
S.7	Tim System annual Gradient Supervision check Date of last test 12/28/21		☒	NA	NA	NA	☐	☐	☐	☐	☐	☐
S.8	ERDU tests performed annually Date of last test 12/28/21		☐	NA	NA	NA	NA	NA	NA	☐	☐	☐
S.9	ERDU Battery 23.5Vdc or higher or replace every 2 years Voltage 24.2Vdc		☒	NA	NA	☐	NA	NA	☐	NA	☐	☐
S.10	Verify Emergency Power Off switches		☒	NA	☐	NA	NA	☐	NA	☐	☐	☐
S.11	Assurance test complete		☐	NA	☐	NA	NA	☒	NA	☐	☐	☐

fix formatting

Notes:

Body Coil failed tuning, time needs to be scheduled to re-tune.

SIEMENS MRI System PM Checklist

	Procedure	Serial Number for tools	Pass	Cleaned	Adjusted	Replaced	Repaired	Needs Adjustment	Needs Replacement	Needs Repair	Does Not Apply	Failed
Maintenance												
M.1	Vacuum dust from computer grill & graphics card		☒	☒	☐	NA	NA	NA	NA	☐	☐	☐
M.2	Vacuum Pump filter (Comfort Kit Option)		☐	NA	NA	NA	NA	NA	☐	NA	☒	☐
M.3	Record compressor hours 125554.1		☒	NA	NA	☐	NA	NA	☐	NA	☐	☐
M.4	Helium Compressor – Static Pressure 13.6 Bar		☒	NA	☐	NA	NA	☐	☐	☐	☐	☐
M.5	Check cooling fans, Cabinets, Gradient filters, clean RF air filters on 3T		☒	☐	☐	NA	NA	☐	☐	NA	☐	☐
M.6	Replace 2 AA batteries, every 2 years on Tim Trolley		☐	NA	NA	☐	NA	NA	☐	NA	☒	☐
M.7	Check primary & secondary water flow, clean filter if necessary, check water lines for leaks & condensation		☒	☐	☐	NA	NA	☐	NA	NA	☐	☐
M.8	Checking the shield temperatures		☒	NA	NA	NA	NA	NA	NA	NA	☐	☐
M.9	Cleaning up directories		☒	☐	NA	NA	NA	NA	NA	NA	☐	☐
M.10	Backup Site specific data, saving Site specific data		☒	NA	NA	NA	NA	NA	NA	NA	☐	☐

Notes:

Full CoBR/disk imaging backup was performed.

Service Engineer Signature and Date: e29b602a-0fd9-4ce6-a8a5-3c17497e0db3 Digitally signed by e29b602a-0fd9-4ce6-a8a5-3c17497e0db3
Date: 2022.12.27 21:46:39 -05'00'