GE Medical Systems

Technical Publications

**Direction 2297997
Revision 0**

Signa LX Release 8.4 ACGD M4 Patch Software Load

Copyright© 2001 by General Electric Co.

Operating Documentation

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "**damage in shipment**" written on **all** copies of the freight or express bill **before** delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately** upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14** day period.

Immediately complete a "Damage Loss Claim Form", available via MS Exchange Mail, after the damage is found.

MS Exchange Path:

Outlook/Public Folder/All Public Folders/Medical Systems/!Global Initiatives/Information Management/Forms/Common Forms/DAMAGE LOSS CLAIM FORM.

Send the completed form to the email address listed in the form.

For more information about the Transportation Claim Procedure, access the GE Medical Systems Intranet and enter the following URL address (case sensitive):

<ftp://3.87.40.2/globepro/qualsys/Docs/190016MF.PDF>

Rev. 11/15/2000

Direction 2128126 - Language Policy For Service Documentation, Rev 0

W A R N I N G

- THIS SERVICE MANUAL IS AVAILABLE IN ENGLISH ONLY.
- IF A CUSTOMER'S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER'S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR OR PATIENT FROM ELECTRIC SHOCK, MECHANICAL OR OTHER HAZARDS.

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REVISION HISTORY

<u>REV</u>	<u>DATE</u>	<u>PRIMARY REASON FOR CHANGE</u>
0	Apr. 12, 2001	Initial version of procedure.

<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>
Title Page.....	0						
Damage in Trans.	-						
Direction 2128126..	0*						
5 to 12	0						

* This revision/letter corresponds to the indicated document's revision control system.

TABLE OF CONTENTS

1- OVERVIEW	6
2- LOADING SIGNA 8.4 ACGD M4 PATCHES FROM CD-ROM	6
APPENDIX A – LX 8.4 ACGD M4 PATCH DESCRIPTIONS	11

1- OVERVIEW

The 8.4 ACGD M4 Patch CD contains fixes for the following software issues (See Appendix A for details):

- **MRlge64190** - EPI-GRE Violates dBdt 66% Limit
- **MRlge65450** - telnetd security
- **MRlge65605** - Minimum Slice thickness check wrong for spiral spectral spatial
- **MRlge62877** - Bad images in 2nd station acquisition of Bolus Chase sequence
- **MRlge63693** - Scan fails in Bolus Chase 3rd station with CERD ERROR SRV REQ
- **MRlge65449** - SmartStep Produces Non-Diagnostic Exams – Reproducible
- **MRlge63977** - Striped Artifact with Multistation/Smart Prep/62.5 kHz
- **MRlge64167** - Stage 3 of Runoff No images

Note

Signa software must be at Release 8.4 ACGD M4 to install this patch CD.

2- LOADING SIGNA 8.4 ACGD M4 PATCHES FROM CD-ROM

1. With the system running and software at the scanning level, go to the Service Desktop and select the **[C Shell]** button in the lower left toolbar.
2. Inside the C Shell window, log in as the root user by typing the following:


su <Enter>

Password: **operator <Enter>** ("operator" is the system default root password. Use whatever password is appropriate for the system you are working on if it has been changed.)

3. Place the "Signa LX 84ACGD M4 Patch" CD in the system CD drive. Wait 20 seconds for the system to see the CD, or until the amber light on the front of the CD drive stops blinking.
4. Inside the C Shell, type the following:

mediad <Enter> (This command mounts the CDRom.)


swmgr <Enter> (This command starts the software manager which is used to load the patch.)

SOFTWARE MANAGER INITIALIZATION
ILLUSTRATION 2-1

- Once the Software Manager starts, begin typing the following in the **[Software Location Edit Window]**. See Illustration 2-1.

/CDROM (The system will automatically enter **/ldist** after you type **/CDROM**.)
- Click the **[Lookup]** button found to the right of the **[Software Location Edit Window]**. See Illustration 2-1.
- Click the **[Customize Installation]** button at the top of the screen. See Illustration 2-1.


LOADING THE 8.4 ACGD M4 PATCHES
ILLUSTRATION 2-2

8. All available patches should be visible at this time in the selection window. Any patches that apply to your system and are not loaded will be selected (check marked). Make sure both "dB/dt patch for 84_ACGD" and "Security Fix for telnetd" are checked; make sure no others are checked. See Illustration 2-2. One of two things will happen next:
 - a. The **[Start]** button will activate (The letters in "Start" will change colors or shades). If the **[Start]** button activates, push it and the patch load will begin. A progress bar will open and count up from 0 to 100%. The load takes less than 1 minute. When the progress bar goes away, a confirmation window will be displayed. Select **[Exit]** to quit the Software manager.
 - b. After the patch is selected, if there is a problem, the **[Conflicts]** button will activate, (the letters in "Conflicts" change colors or shades). If this happens, this indicates that your system does not need this patch and the indicated patch will not be loaded. Push the **[Conflicts]** to view the error if desired. (See Note 1 and Note 2 .)

Note 1: If Patch_84ACGD_M4_01 "dB/dt patch for 84_ACGD" was previously installed, it will be flagged as a conflict now. If Patch_84ACGD_M4_01 was previously installed, you should not install it again.

Note 2: If Patch SG0004060 "Security Fix for telnetd" was previously installed, it will be flagged as a conflict now. If Patch SG0004060 was previously installed, you should not install it again.

9. After installation, select **Exit, Installations completed** checkbox, then select **[OK]**. See Illustration 2-3.

INSTALLATION DECISION WINDOW
ILLUSTRATION 2-3

10. In the C Shell type the following to eject the CD-ROM from the drive:

umount /CDROM <Enter> (Un-mounts the CD from the system.)
eject /CDROM <Enter> (Opens the CD drawer.)

11. If the patches were required and loaded, verify they are installed by typing in the C Shell:

getver <Enter> You should see an entry for the 84ACGD patch as follows:

Build number for MrpApps is 84_ACGD.282J_M4_0105.b
Build number for PostSdC is 84_ACGD.282J_M4_0105.b
Build number for cclass is 84_ACGD.282J_M4_0105.b
Build number for driverSupport is 84_ACGD.282J_M4_0105.b
Build number for insite is 84_ACGD.282J_M4_0105.b
Build number for install is 84_ACGD.282J_M4_0105.b
Build number for os_cd is 6.8.a
Build number for os_u_cd is 6u7.14.a
Build number for patch_84ACGD_M4_01 is 84_ACGD.282J_M4_0105.b.PA <- **84ACGD Patch entry**

Build number for pc is 84_ACGD.282J_M4_0105.b
MR Software release: 8.4._A.C.G.D.0105b

12. In the C-Shell type: **showprods | grep patch <Enter>**. The following three lines for the “telnetd” patch should be included in the output:

(Ignore other patch lines before the following.)

```
I patchSG0004060 <date> Patch SG0004060: Security Fix for telnetd (does not require a
reboot)
I patchSG0004060.eoe_sw <date> IRIX Execution Environment Software
I patchSG0004060.eoe_sw.base <date> IRIX Base Execution Environment
```

(Note: The following three lines are displayed for the 84ACGD patch.)

```
I patch_84ACGD_M4_01 <date> Patch_84ACGD_M4_01: dB/dt patch for 84_ACGD
I patch_84ACGD_M4_01.MrpApps_sw <date> MR Applications
I patch_84ACGD_M4_01.MrpApps_sw.base <date> MR Applications
```

13. Reboot Signa for the patches to take effect. (The “telnetd” patch does not require a reboot, but the “84ACGD_M4_01” patch does require a reboot.)

APPENDIX A – LX 8.4 ACGD M4 PATCH DESCRIPTIONS

This patch CD includes the fix for three problems found in 84ACGD_M4 Software. The size of the patches is approximately 19 megabytes (uncompressed).

PROBLEM 1 - dB/dt limit problems in epi

SPR - MRlge64190 EPI-GRE Violates dBdt 66% Limit
SPR - MRlge65605 Minimum Slice thickness check wrong for spiral spectral spatial

There may be an error in dB/dt calculations when the slice thickness is decreased below 2.1 mm using EPI sequences on a CRM system. In that case, the gradient strength exceeds the maximum allowable for rise time of the spatial-spectral pulse that is fixed. The fix limits the maximum gradient strength of the pulse based on the Reilly curve limit. This then limits the minimum slice thickness that can be achieved with this pulse (2.1 mm => SR150, 40 mT/m)

The fix is in: epi, epi2, epiVP, and spiral PSD files

PROBLEM 2 - Security Fix for telnetd

SPR - MRlge65450 telnetd patch

An SGI patch to the SGI software that Release 8.4 ACGD runs on is also being distributed at the same time. It is security-related.

The fix is in: telnetd

PROBLEM 3 - Multi-Station Bolus Chase Artifacts after First Station

SPRs:

MRlge62877 Bad images in 2nd station acquisition of Bolus Chase sequence
MRlge63693 Scan fails in Bolus Chase 3rd station with CERD ERROR SRV REQ
MRlge65449 SmartStep Produces Non-Diagnostic Exams - Reproducible
MRlge63977 Striped Artifact with Multistation/Smart Prep/62.5 kHz
MRlge64167 Stage 3 of Runoff No images

A striped artifact was observed to occur in the 2nd and later stations of Multi-Station EPI Bolus Chase acquisition. Problem was isolated to a filter program pointer on the CERD ASC not getting initialized when the ASC is in "EPI" mode.

The fix is in: cerd.P3.i40

Uninstall Instructions – Signa 8.4 ACGD M4 Patch

In the event that the customer is not satisfied with the new software, the patches may be uninstalled as follows.

1. With the system running and software at the scanning level, go to the Service Desktop and select **[C Shell]**.
2. Inside the C Shell window, log in as root by typing: **su <Enter>**. (Password default is: **operator <Enter>**.)

3. Inside the C Shell window, start the software manager by typing: **swmgr <Enter>**
4. Click on the **[Manage Installed Software]** button.
5. Scroll the window until you find the "Patch_84ACGD_M4_01: dB/dt patch for 84_ACGD" line. Select the checkbox on the left side of the line. Also select the checkbox for "Security Fix for telnetd". Make sure no other lines are checked.
6. Click on **[Start]**. The patch uninstall will begin. A progress bar will count up to 100%. The uninstall should take less than one minute.
7. At the completion of the uninstall, a popup panel appears. Select **"Exit, Installation Complete"**, and click **[OK]**.
8. Reboot Signa for the patches uninstall to take effect.