

TABLE OF CONTENTS

TABLE OF CONTENTS	1
1- OVERVIEW	2
2- REMOVE, INSPECT, AND REPLACE AIR FILTERS.....	2
REVISION HISTORY.....	3

Description - This section describes the air filter replacement performed on the 1.0T and 1.5T RF/PEN II cabinet.

1- OVERVIEW

The cabinet uses forced air cabinet cooling. The cabinet air passes through viscous impingement air filters that must be replaced as necessary to ensure proper cabinet cooling. **Replace filters during the initial installation of the site before turning the system over to the customer.**

When replacing air filters, be sure to use GE part number 2160200-46.

2- REMOVE, INSPECT, AND REPLACE AIR FILTERS

1. Turn the recessed screw at the top center of the front cover 1/4 turn counterclockwise to release the cover latch.
2. Tilt the top of the front cover forward, and lift up and out to remove the door for access to the filter(s).
3. Remove each filter by lifting up slightly and pulling bottom edge outward.
4. Inspect the filter for dust build-up, and replace as required.
5. Insert the filters first up and into the top bracket, then drop down into lower bracket, keeping the intake side (side with cleaning instructions) facing the front cover.
6. Replace the front cover on the cabinet by inserting the lower cover lip behind the shoulder bolts at the bottom of the cabinet and swinging the top of the cover into the spring loaded "slam-latch".

REVISION HISTORY

REV	DATE	AUTHOR	PRIMARY REASONS FOR CHANGE
0	June 12, 1998	Erbtec Engineering	Initial conversion to Word. (KK)