DV ICN Upgrade to GEN6

OPERATING DOCUMENTATION

5670087
Revision 3.0

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**WARNING
(EN)**

This service manual is available in English only.

- If a customer's service provider requires a language other than English, it is the customer's responsibility to provide translation services.
- Do not attempt to service the equipment unless this service manual has been consulted and is understood.
- Failure to heed this warning may result in injury to the service provider, operator or patient from electric shock, mechanical or other hazards.

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This page left intentionally blank.

Revision History

Revision	Date	Description
1.0	March 5, 2014	Initial Release
2.0	March 1, 2016	Updated to DV25.1_R01, GEN6, and Dell R-630. Updated material return procedures.
3.0	February 23, 2018	Chapter 1: Revised harvest information. Added MyPLCReader App instructions.

This page left intentionally blank.

Table of Contents

Chapter 1 OVERVIEW.....	15
1 ICN Upgrade Installation.....	15
1.1 Upgrade Introduction.....	15
1.2 Damage in Transportation	15
1.3 Product Locator.....	15
1.4 Required Equipment and Tools.....	16
1.5 Discarded Material Return Policy.....	17
1.5.1 Material Policy Overview	17
1.5.1.1 Returned Material	17
1.5.1.2 Disposed Material	17
1.5.2 Harvest Testing.....	18
1.5.3 Harvested Materials.....	18
1.5.4 Return Kit	18
1.5.5 Deinstallation Procedures	19
1.5.6 Product Locator Information	19
1.6 Upgrade Components.....	19
Chapter 2 POWER OFF.....	21
1 Power Off Procedure.....	21
1.1 Overview.....	21
1.2 Safety.....	21
1.3 Power Down ICNs from VRE Power Control.....	21
1.4 Power Off ICNs Manually.....	25
Chapter 3 UPGRADE PROCEDURES.....	27
1 Sun ICN Upgrade Procedures.....	27
1.1 Overview.....	27
1.2 Removal Procedures.....	27
1.2.1 Removing ICN Sun 4100.....	27
1.2.2 Removing ICN 4170.....	30

1.3 Install Dell R-630 ICN and Mounting Hardware.....	32
2 Dell ICN Upgrade Procedures.....	39
2.1 Overview.....	39
2.2 Removal Procedures.....	39
Chapter 4 POWER ON.....	43
1 Power On Procedure.....	43
Chapter 5 ICN RECONFIGURATION.....	45
1 ICN Reconfiguration.....	45

Chapter 1 Overview

1 ICN Upgrade Installation

1.1 Upgrade Introduction

This publication documents the mechanical installation information for upgrading the existing SUN ICN (GEN2, 3, 4) module or existing Dell 32 channel ICN (GEN5) in the PGR cabinet with the new Dell ICN (GEN6) module. The Dell ICN (GEN6) is only compatible with DV25.1_R01 and later software.

WARNING

FERROUS MATERIAL HAZARD!

DELIVERY EQUIPMENT AND OTHER TOOLS AND PARTS REQUIRED FOR THIS INSTALLATION CONTAIN FERROUS MATERIAL AND WILL BE STRONGLY ATTRACTED TO MAGNET AND MAY BECOME DANGEROUS PROJECTILES.

KEEP ALL FERROUS TOOLS AT LEAST 10 FEET AWAY FROM THE MAGNET.

CAUTION

Possible personal injury!

The Sun 4170 ICN weighs approximately 40 pounds.

Two people are required to lift it.

1.2 Damage in Transportation

All packages should be closely examined at time of delivery. If damage is apparent, have notation "**damage in shipment**" written on **all** copies of the freight or express bill **before** delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately** upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14** day period. File a report with:

- Call 1-800-548-3366. Select "Install Support Services for FOA and MIS"
- Contact your local service coordinator for more information on this process.

1.3 Product Locator

During installation, the Install FE shall confirm that the Model and Serial Numbers are part of the asset. If missing, use either the gib web interface (http://gib.gehealthcare.com/gib/gib_entry.jsp) or the MyPLCReader App. This information will then be loaded directly to Siebel by the IB Verification Team.

Use the following procedure to install the MyPLCReader App.

1. From the GE App Icon, search for and install the MyPLCReader App on your iPhone or iPad.

Illustration 1-1: MyPLCReader App

2. Update the Configuration Page.
 - a. The first time the App opens, the Configuration Page will open.
 - b. Enter your SSO, Country and **if in US select the edit button to change the default email address** (by default, this app will send the email to the Europe IB Team Resource.EuropeanCentralAdminGIB@ge.com). The email to be used in US is usa.ib.verification@ge.com.
3. Start your first Scan.
 - a. Go to the Menu, make the appropriate selection for the job you are performing. **New Install, Upgrade System, or Update System.**
 - b. Enter the applicable Data, Install End Date, System ID & RFS#.
 - c. Click **Scan Bar Codes** and start scanning each barcode on the PL Cards.
 - d. Click the “Email” Button.
 - e. Email is sent to the Default email address, which was set up earlier.

1.4 Required Equipment and Tools

NOTICE

Follow all required safety and PPE procedures customary for your organization when working on this product.

Table 1-1: General Safety Requirements

Item	Description
1	Wear composite safety shoes for foot protection.
2	Wear gloves for hand protection.
3	Wear safety glasses for eye protection.

FERROUS MATERIAL HAZARD!

IF MAGNET IS ENERGIZED, TOOLS, COMPONENTS, AND PARTS REQUIRED FOR THIS INSTALLATION THAT CONTAIN FERROUS MATERIAL WILL BE STRONGLY ATTRACTED TO MAGNET AND MAY BECOME DANGEROUS PROJECTILES.

KEEP ALL FERROUS TOOLS AS FAR FROM THE MAGNET AS POSSIBLE.

Tools to upgrade the system are listed in following table:

Table 1-2: Installation Equipment

Item	GE Part Number	Description
1		Standard Hand Tools
2		Screwdriver, slotted (long & narrow)
3	46-271138G1	Restricted Access Control Kit. Contains two plastic warning signs for posting at site during installation and service activity.
4		Service Methods DVD
5	46-198094P1	Wrist grounding strap
6	---	Plastic or aluminum flashlight

1.5 Discarded Material Return Policy

The equipment or material being removed by the upgrade have a potential value for service of the installed base. The disposition of such must be controlled by GE. Disposition of the removed material is per GE Medical Systems Policy and Procedures.

1.5.1 Material Policy Overview

1.5.1.1 Returned Material

Per GE Healthcare Policy 3.9: all requested material is to be returned, as directed, intact and promptly to Approved GE Recycling Center. Do not remove or swap any parts from the material to be returned.

1.5.1.2 Disposed Material

Regulations governing the disposal of this equipment or material vary from location to location and from country to country. In order to insure compliance with all local, state, or national environment, health, and safety regulations, dispose of all items per local GEHC policy and procedures and other applicable published guidance. Contact GEHC Recycling Center for assistance.

1.5.2 Harvest Testing

NOTICE

The equipment must be tested prior to the system shutdown/upgrade to allow appropriate Harvest capabilities.

Before going on site:

- Print the proper regional address labeling from DOC2100599 Harvest Return Addresses.
- Verify that you can download the Harvest test procedure from MyWorkshop using DOC1899393.
- If you have not been contacted to test the system or cannot download the document, reach out to the local Harvest team. For all other requests contact:

@HEALTH Global Harvest Team GlobalHarvestTeam@ge.com

Once testing is completed, submit the test record for review. Once these test records are reviewed, they will be processed for electronic signatures.

Package the equipment and follow your regions method of returning harvest equipment.

1.5.3 Harvested Materials

Table 1-3: List of Materials to Be Harvested

System Type for Upgrade	MR450, MR450W, MR750 and MR750W- Test Procedure DOC1899393		
Subassembly	Part Number	Description	Qty
Parts to be returned are located in PGR Cabinet	5127452-2	Image Compute Node - Gen 2 with VRF	Up to 2
	5127452-6	Image Compute Node - Gen 2 with VRF and without Infiniband	1
	5337894	Image Compute Node - DV 16CH GEN 3	1
	5337894-2	Image Compute Node - DV 32CH GEN 3	1
	5911000	Image Compute Node - Gen 4.1 DV 16 Channel	1
	5911000-2	Image Compute Node - Gen 4.1 DV 32 Channel 32GB	1
	5911000-3	Image Compute Node - Gen 4.1 DV 32 Channel 72GB	1
	5921000	Gen5 16CH ICN	1
	5921000-2	DV Gen5 32CH ICN	1
	5921000-3	DV Gen5 32CH+ ICN	1
		Mounting brackets	Up to 4
		Cables	

NOTE: Reuse the packaging that the new Dell ICN was shipped in to return the removed Sun or Dell ICNs.

1.5.4 Return Kit

If needed, packaging and/or transportation details can be obtained through the local Installation Specialist.

Reuse appropriate containers that new hardware and modules were shipped to site for return of above material.

1.5.5 Deinstallation Procedures

Refer to instructions in this manual for deinstallation procedures of site.

1.5.6 Product Locator Information

Process product locator information for all serialized components removed from the system according to policy. Refer to [Section 1.3](#).

1.6 Upgrade Components

The following components are provided for the Sun to Dell ICN upgrade:

- Dell ICN
- ICN Mounting Kit that includes:
 - Slide Bracket (2)
 - Slide Rail (1)
 - Nuts (8)
 - M5 x 12 Screw (8)
- Cable Kit that includes:
 - Ethernet Cables (3)
 - ICN Power Cord (1)

For upgrading an existing Dell 32 Channel to a Dell 32 Channel +, only the ICN will be upgraded. The provided slides, slide brackets, and cables are to be discarded.

This page left intentionally blank.

Chapter 2 Power Off

1 Power Off Procedure

1.1 Overview

This procedure documents the proper methods of shutdown for GEN2 to GEN4 ICNs. **NOTE:** The ICN modules are automatically powered on when the host computer at the GOC is powered on. There is no need to switch the ICNs on individually.

1.2 Safety

WARNING

ELECTROCUTION HAZARD!
DANGEROUS VOLTAGES ARE PRESENT THROUGHOUT THE SYSTEM CABINET WHEN POWER IS APPLIED.
REFER TO THE MR SERVICE SAFETY MANUAL, PN 5452735.

NOTICE

Follow this procedure for removing power to the ICN chassis. Failure to properly power down ICNs may result in premature ICN failure.


NOTICE

Do not turn power off to the ICNs using the power button on the front panel. The ICNs have software running on disk drives that may become corrupted if not powered down properly.

1.3 Power Down ICNs from VRE Power Control

1. At the host computer, on to the Service Browser, select *Utilities > Toolbox > VRE > Power Control*.

Illustration 2-1: Power Control

NOTE: The number configuration should match the number of VREs.


- Click [POWER OFF ICN(x)] (shown in previous illustration). Wait for the ICN to power off.

Illustration 2-2: ICN Switching Off Indication

When the Power Control link is refreshed, the screen shows the ICNs in their current state.

Illustration 2-3: ICN 1 Switched Off

NOTE: The illustration above shows images from a system with a Sun 4100 ICN. Systems with a Sun 4170 or 4170M2 ICN have only one ICN selection.

3. The Sun ICN(s) are powered down if the green disk drive LEDs on the front of the ICN(s) are no longer illuminated. The Power LEDs continue to blink even after power-down.

Illustration 2-4: ICN Front Panel (Sun 4100)


Illustration 2-5: ICN Front Panel (Sun 4170 and 4170M2)

When ICN is powered-down, Power LED will blink, but disk drive LEDs will be OFF.

4. The Dell ICN is powered down if the green power LED located on the upper left side of the front of the ICN is not illuminated.

Illustration 2-6: ICN front Panel for Dell R-630

1.4 Power Off ICNs Manually

NOTICE

If the browser or the command line are not accessible, the ICNs must be powered off manually. Try this method before powering off the GOC. This is the preferred method because it does not corrupt the files on the hard drive.

1. For Sun ICNs, momentarily press the power button. This brings down the ICN in a soft method. It takes about a minute or two for the operating system (OS) to boot down. (After power-down, the green power LED continues to slowly blink, but the disk drive LEDs go dark.)

Illustration 2-7: ICN Power LED and Power Button (Sun 4100)


Illustration 2-8: ICN Power LED and Power Button (Sun 4170 and 4170M2)

2. If the OS freezes or is not responding, then the soft method did not work properly. In this case, hold the power button down for six seconds to force the power off.
3. For Dell ICN momentarily press the power button. The green LED is also the power button.

NOTICE

Holding the power button down for several seconds will bring it down hard.

Chapter 3 Upgrade Procedures

1 Sun ICN Upgrade Procedures

1.1 Overview

This procedure provides instructions on for removing existing Image Compute Nodes (ICNs) Sun 4100, Sun 4170, and Sun 4170 M2 to be replaced with the Dell R-630.

1.2 Removal Procedures

1. Make sure that the power to the ICN(s) have been turned off. Refer to Chapter 1.
2. Perform LOTO on the reconstruction equipment. See Chapter 3, Section 2.5 of the *MR Service Safety Manual*, PN 5452735.


Refer to one of the following two procedures for removal of the existing ICN..

- [Section 1.2.1](#) for removal of Sun 4100 ICN
- [Section 1.2.2](#) for removal of Sun 4170 ICN

1.2.1 Removing ICN Sun 4100

Location of ICNs in PGR cabinet.

Illustration 3-1: ICN Location - 4100


1. On the front of the ICN, undo the thumbscrew on each side of the ICNs being removed.

Illustration 3-2: Thumbscrews on ICN Front - 4100

2. Slide out the ICN until the stop is reached at the end of the side rails.
3. Remove all the connections at the back of the ICN.

Illustration 3-3: Cable Connections on Sun 4100 ICN


CAUTION

Possible personal injury!
The Sun 4100 ICN weighs approximately 40 pounds.
Two people are required to lift it.


4. Locate the two release tabs on the slide rails. Press both release tabs toward the center of the ICN while pulling the ICN forward. Pull the ICN forward until the rails release and the ICN is free. Place the ICN to the side in a safe area.

Illustration 3-4: Slide Rails and Release Tab


5. Remove the ICN trays from the ICN mounting brackets.

Illustration 3-5: Sun 4100 ICN Trays

6. Unfasten the screws securing the ICN mounting bracket/slide assemblies from the vertical rails and remove the assemblies from each side of the cabinet.

Illustration 3-6: Side Mounting Brackets for 4100


7. Remove and save all clip nuts that were used to fasten the above slides and bracket assemblies. They will be reused for installing the new ICN.

1.2.2 Removing ICN 4170

Location of ICN in PGR cabinet.

Illustration 3-7: ICN Location - 4170

1. On the front of the ICN, undo the thumbscrew on each side of the ICN to be removed.


Illustration 3-8: Thumbscrews on ICN Front - 4170

2. Slide out the ICN until the stop is reached at the end of the side rails.

3. Remove all the connections at the back of the ICN.

Illustration 3-9: Cable Connections on Sun 4170 ICN

CAUTION


Possible personal injury!

The Sun 4170 ICN weighs approximately 40 pounds.

Two people are required to lift it.


4. Locate the two release tabs on the slide rails. Press both release tabs toward the center of the ICN while pulling the ICN forward. Pull the ICN forward until the rails release and the ICN is free. Remove the ICN and place to the side in a safe area.

Illustration 3-10: Slide Rails and Release Tab


5. Loosen the mounting screws for the network switch. Slide the switch down to give access to the mounting location for the slide rails for the 4170 ICN.

Illustration 3-11: Network Switch Screws

6. Unfasten screws and remove slide and bracket assemblies from each side of cabinet.

Illustration 3-12: ICN Slide and Bracket Removal


7. Remove and save all clip nuts that were used to fasten the above slides and bracket assemblies. They will be reused for installing the new ICN.

1.3 Install Dell R-630 ICN and Mounting Hardware

The installation of the Dell R-630 will require attaching new mounting brackets (5496069) on the front vertical cabinet rails.

Illustration 3-13: ICN Bracket Installation

1. Install the clip nuts to the locations indicated around the cut out on the front vertical rack rail as shown in the above illustration. Left side vertical rack rail is shown. Repeat for right side rail.
2. Attach brackets to each rail using the previously removed screws.

Illustration 3-14: Ethernet Switch


3. Reattach the Ethernet Switch.
4. New cables will be installed later in procedure.

Illustration 3-15: ICN Slide Rails

5. Locate the Dell mounting rails.
6. Remove the orange shipping clips.

Illustration 3-16: Align Front of Slide


7. Position the left slide rail labeled FRONT facing inward and align the end of the rail to seat in the bottom three holes of the installed bracket.

NOTICE

Verify that the new back rail bracket is aligned to the proper holes prior to seating. Use a level to make sure that the rail is aligned to the correct holes. It is difficult to remove the back rail bracket after it has been installed.

Illustration 3-17: Secure Rear of Slide

8. At the back end of the rail, align and engage the slide rail by pulling it outward, until it fully seats into the rear vertical cabinet rack rail and the latch clicks into place

Illustration 3-18: Secure Front of Slide


9. Engage the front end of the slide rail into the vertical cabinet rack rail.

Illustration 3-19: Install Right Side Slide Rail


10. Repeat previous steps for installing right side slide rail.

Illustration 3-20: Release Tab


11. The blue tab on the end of the slides must be flat with the end surface of the rail for the slide rail to be completely connected to the vertical cabinet rail or bracket. If rails need to be removed, pull the blue latch release tab on the end of each rail and unseat. The blue tab mechanism can easily be observed on the front of the rails.

Illustration 3-21: Slide Out Rails

12. Pull the inner slide rails out of the rack until they lock into place.

Illustration 3-22: Move ICN to Rails


CAUTION

Possible personal injury!
The ICN weighs approximately 40 pounds.
Two people are required to lift it.


13. Remove the new ICN from packaging. Position the rear rail standoff on each side of the ICN and lower into the rear J-slots on the slides.

Illustration 3-23: Secure ICN in Rails


14. Rotate the ICN module downward until all eight rail standoffs are seated in the J-slots. Push the ICN module inward until the two lock levers, one on each side, click into place.

Illustration 3-24: Install ICN Cables


15. Remove the existing 3 Ethernet cables that run between the ICN and Ethernet switch and install new ones that are provided.
16. Connect power cables to the power port(s) on the ICN.

Illustration 3-25: Connect Fiber Optic Cables

17. Reconnect the existing fiber optic cables to J1 and J5 as marked. Two sets of fiber optic cables will be connected if the system has a 32 Channel ICN as shown in the illustration.

Illustration 3-26: Move ICN Into Cabinet

18. Press the slide-release lock buttons on both rails and carefully slide the module into the cabinet.

2 Dell ICN Upgrade Procedures


2.1 Overview

This procedure provides instructions on for removing existing Dell R-620 or Dell R-630 32 Channel Image Compute Node (ICN) and replaced with the Dell R-630 32 Channel Plus ICN.

2.2 Removal Procedures


1. Make sure that the power to the ICN has been turned off. Refer to Chapter 1.
2. Perform LOTO on the reconstruction equipment. See the *MR Service Safety Manual*, PN 5452735.

Illustration 3-27: Move ICN out of Cabinet


3. Pull the ICN module slide rails out of the rack until they lock into place.
4. Disconnect all cables from the rear of the ICN module.

Illustration 3-28: Remove Dell ICN


5. Locate the two release tabs on the slide rails. Press both release tabs toward the center of the ICN while pulling the ICN forward. Pull the ICN forward until the rails release and the ICN is free. Pivot ICN upwards and remove.

Illustration 3-29: Move New Dell Plus ICN to Rails


6. Remove the new Dell 32 Channel Plus ICN from packaging. Position the rear rail standoff on each side of the ICN and lower into the rear J-slots on the slides.

Illustration 3-30: Secure ICN in Rails


7. Rotate the ICN module downward until all eight rail standoffs are seated in the J-slots. Push the ICN module inward until the two lock levers, one on each side, click into place.

Illustration 3-31: Connect cables to ICN

8. Reconnect the 3 ethernet cables. Connect the two power cables.

Illustration 3-32: Move ICN into Cabinet

9. Press the slide-release lock buttons on both rails and carefully slide the module into the cabinet.

This page left intentionally blank.

Chapter 4 Power On

1 Power On Procedure

1. Remove LOTO from the reconstruction equipment. See the *MR Service Safety Manual*, PN 5452735.

NOTE: The Dell ICN (GEN5) is only compatible with DV25 software. The Dell ICN (GEN6) is only compatible with DV25.1_R01 software.

2. On the Common Service Desktop, go to Utilities > Toolbox > VRE > Power Control. Click [POWER ON ICN] for the selected ICN. The new ICN is automatically powered on.

NOTE: The Dell ICN will display the boot status on the LED readout, which is located on the front panel of the ICN.

This page left intentionally blank.

Chapter 5 ICN Reconfiguration

1 ICN Reconfiguration

The Dell ICN (GEN6) is only compatible with DV25 software.

1. Gen 6 ICNs are supported only on DV25.1_R01 and later s/w. If needed perform a DV25 Host Software Load.
2. If already on DV25.1_R01 or later s/w, Perform ICN reconfiguration per Service procedures.
3. After ICN Reconfiguration completes, reboot the system.
4. Perform a Check Scan and ensure system scans successfully.
5. Perform a SAVEINFO to save new updated vre.cfg files onto Save Info Media.

This page left intentionally blank.

© 2018 General Electric Company. All Rights Reserved.

www.gehealthcare.com